COOK

MEDICAL

CE
0123

EN
3

Roadrunner® Extra Support Wire Guide

Instructions for Use

BG
5

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ROADRUNNER® EXTRA SUPPORT WIRE GUIDE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Roadrunner® Extra Support Wire Guide features a PTFE-coated nitinol (nickel-titanium alloy) shaft and a short, angled platinum coil tip. This device comes in two diameters and multiple lengths. This device is provided with an insertion tool and torque device. Refer to the product label for product specifications (e.g., wire guide diameter and length). The device is packaged in a spiral holder with a Luer lock hub.

Performance Characteristics

- Radiopaque platinum coil for positioning under fluoroscopy
- PTFE coating on main body for low-friction manipulation
- Nitinol core for flexibility
- Torque device for directional manipulation and torque control
- Insertion tool for facilitating insertion into a procedural device

Device Compatibility

- Vascular devices with a larger inner diameter than the specified outer diameter of the wire guide
- Heparinized saline
- Standard Luer syringe

Patient Population

The target patient population are adult patients with a medical condition whose treatment a physician has determined requires access to the peripheral vasculature.

Intended Users

These devices are used by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

This device contacts endoluminal surfaces in peripheral vasculature.

Operating Principle

This device is inserted (flexible end first) into the target vessel through the lumen of a catheter or other device. An insertion tool may be employed to aid in wire guide insertion (e.g., by straightening out a curved tip or advancing the wire guide through a hemostasis valve and into an access sheath or catheter). After insertion, medical imaging is used to maneuver the wire guide to the desired location within the body. A torque device may be attached to the proximal end of the wire guide to aid in rotational manipulation of the device. Once the wire guide has reached the target location, another device can be advanced over the wire guide and directed to the target anatomic site.

INTENDED USE

The Roadrunner Extra Support Wire Guide facilitates delivery of devices into the peripheral vasculature.

DEVICE AND CLINICAL BENEFITS

To assist with maintenance of vascular access for minimally invasive placement of medical devices, by providing a track across which the devices can be inserted into the vasculature (without the wire guide itself having a direct diagnostic or therapeutic function).

INDICATIONS FOR USE

This device is indicated for use in procedures that require access to the peripheral vasculature.

CONTRAINDICATIONS

This device is not intended for coronary or neurovascular use.

WARNINGS

- To prevent wire guide loss, keep at least 5 cm of the wire protruding from the device fitting at all times and maintain contact with the wire.
 - Manually altering the tip configuration may damage the wire guide.
 - The wire guide should be advanced only when visualizing its tip under imaging guidance. Do not torque wire guide without evidence of corresponding movement of distal tip.
- NOTE:** Torquing against resistance may cause vessel trauma or wire guide damage which may lead to device failure.
- To prevent possible tissue damage, care should be taken when manipulating a device over a wire guide during the device's placement and withdrawal. If resistance is felt during device placement, discontinue the procedure and determine the cause of resistance before proceeding.

If the cause of resistance cannot be determined, remove the wire guide and device as a unit to prevent possible damage and/or complications.

- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

This device contains a nickel-titanium alloy. It should not be used on patients who have allergies to any of the device components.

POTENTIAL ADVERSE EVENTS

- Bleeding
- Infection
- Injury to (or inflammation of) surrounding structures
- Wire coating shearing off
- Wire displacement, dislodgement, or loss of access
- Wire entrapment
- Wire fracture
- Wire kinking, bending, or knotting
- Wire unravelling or uncoiling
- Thrombus formation on the wire or local to the access/procedure site
- Vasospasm

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. If desired, attach torque device to wire guide.
2. Under imaging guidance, and while maintaining position of the catheter, advance wire guide to target anatomic site.
3. Confirm location of wire guide in two planes to assure placement of distal tip in desired vessel and not within a side branching vessel.
4. With the wire guide secured in place, advance any additional device over the wire guide and direct it to the target site.

Wire Guide Exchange

1. Remove torque device from wire guide.
2. Open hemostatic valve and flush line.
NOTE: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.
3. Withdraw wire guide under fluoroscopy slowly and gently to prevent trauma.
NOTE: If the wire guide will be used again for the same patient during the same procedure, after removing wire guide from catheter, wipe with sterile gauze pad moistened with heparinized saline solution and place back into holder.
4. Close hemostatic valve, flush line, and prepare a new wire guide.
NOTE: Prior to re-inserting exchanged wire guide, verify that distal tip of cardiovascular device is free within vessel. If tip is against vessel wall, vessel trauma may result upon re-introduction of wire guide.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd