EN
3

Roadrunner® Extra Support Wire Guide

Instructions for Use

CS
4

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DA
5

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Bruksanvisning

ROADRUNNER® EXTRA SUPPORT WIRE GUIDE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Roadrunner Extra Support Wire Guide features a gray PTFE coated nitinol shaft and a short, angled platinum coil tip. This device comes in two diameters and multiple lengths. To aid in gripping the device, a pin vise is included. Refer to product label for product specifications.

INTENDED USE

The Roadrunner Extra Support Wire Guide is intended for use in facilitating delivery of percutaneous catheters into the cardiovascular system.

CONTRAINDICATIONS

This device is not intended for coronary and neurovascular use.

WARNINGS

- The wire guide should be advanced only when visualizing its tip fluoroscopically. Do not torque wire guide without evidence of corresponding movement of distal tip.
- Refer to labeling of percutaneous cardiovascular access catheter for contraindications and potential complications associated with use.
- **Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.**

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in cardiovascular system access procedures. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Withdrawal or manipulation of distal spring coil portion of wire guide through needle tip may result in breakage.
- This wire guide is a delicate instrument and should be handled carefully; forceful angulation should be avoided.

INSTRUCTIONS FOR USE

1. Attach pin vise to wire guide.
2. Under fluoroscopy and while maintaining position of cardiovascular access catheter, advance wire guide to targeted site.

WARNING: Observe all wire guide movement in vessels under fluoroscopy. Do not torque wire guide without evidence of corresponding movement of distal tip. Further torquing against resistance may cause vessel trauma or wire guide damage, which may lead to device fracture.

WARNING: Wire guide should be advanced only when visualizing wire guide tip fluoroscopically. If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

3. Upon reaching desired location, advance distal tip of wire guide distal to lesion, placing a more supportive section of wire guide at targeted site.
4. Confirm location of wire guide in two planes to assure placement of distal tip in desired vessel and not within a side branching vessel.
5. With wire guide secured in place, advance cardiovascular access catheter to target site.

Wire Guide Exchange

1. Remove pin vise from wire guide.
2. Open hemostatic valve and flush line.

WARNING: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.

3. Withdraw wire guide under fluoroscopy slowly and gently to prevent trauma. **NOTE:** After removing wire guide from catheter, wipe with sterile gauze pad moistened with heparinized saline solution or sterile water and place back into holder.

4. Close hemostatic valve, flush line, and prepare a new wire guide.

WARNING: Prior to re-inserting exchanged wire guide, verify that distal tip of cardiovascular device is free within vessel. If tip is against vessel wall, vessel trauma may result upon re-introduction of wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd