EN
3

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Bruksanvisning

RING TRANSJUGULAR INTRAHEPATIC ACCESS SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Ring Transjugular Intrahepatic Access Set may consist of the following components:

- Ross modified Colapinto Needle
- PTFE catheter
- Flexor® Check-Flo® Introducer
- Dilator
- Torcon® NB Advantage Catheter
- Van Andel Dilatation Catheter
- Curved Newton Wire Guide
- Amplatz Extra Stiff Wire Guide

INTENDED USE

The Ring Transjugular Intrahepatic Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Avoid transcapsular puncture.
- Manipulation of products requires fluoroscopic guidance.
- The patient should be monitored (EKG, oximetry, blood pressure, pulse) during the procedure.
- Before introduction and periodically during the procedure, the transjugular introducer sheath should be flushed with saline.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- If the device is to remain in a vessel for an extended period of time, appropriate measures (i.e., flushing with heparinized intravenous flushing solution) should be taken to reduce risk of clot formation.

POTENTIAL ADVERSE EVENTS

- Intraperitoneal and liver hemorrhage
- Carotid artery puncture
- Pneumothorax
- Air embolism
- Puncture site hematoma
- Cardiac arrhythmia
- Arteriovenous fistula

INSTRUCTIONS FOR USE

1. Using Micropuncture or standard 18 gage access techniques, introduce an appropriate .035 inch wire guide into the inferior vena cava via jugular venous access.
2. Using the selective hepatic catheter and a wire guide of choice, catheterize the right hepatic vein (RHV) or the most adequate hepatic vein branch. Leave the wire guide in a safe, distal position, and remove the catheter.
3. Advance the transjugular introducer sheath assembly over the wire guide as distal as possible into the hepatic vein. Remove the wire and dilator.
4. Insert the Colapinto needle into the PTFE catheter, and introduce the catheter/needle assembly through the sheath. The needle should not protrude past the catheter end hole.
5. Advance the catheter/needle assembly into the distal hepatic vein.
6. Orient the catheter/needle assembly inferiorly and rotate anteriorly (arrow baseplate on needle indicates direction of needle curve).
7. Wedging the catheter/needle assembly against the vein wall, thrust the needle in one movement forward through the hepatic parenchyma and toward the portal system.
8. Connect a syringe with contrast medium to the needle, apply suction and withdraw the needle until blood is seen.
9. Inject a small amount of contrast medium in order to confirm position in the portal system.

10. Advance the catheter over the needle further into the portal vein, then remove the needle.
11. Introduce an appropriate 0.035" wire guide into the portal branch and select the main portal vein (MPV).
12. Advance the catheter and sheath into position across the parenchymal tract.
13. Advance an appropriate catheter through the PTFE catheter into the portal vein.
14. Proceed with any indicated interventional procedures.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd