**EN**  
**3**

**Rösch-Uchida Transjugular Liver  
Access Set**

Instructions for Use

**DA**  
**4**

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Bruksanvisning

## RÖSCH-UCHIDA TRANSJUGULAR LIVER ACCESS SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The Rösch-Uchida Transjugular Liver Access Set may consist of:

- Catheter with 0.038 inch flexible trocar stylet
- Stiffening cannula
- PTFE catheter
- Flexor® Check-Flo® introducer set
- Dilator

### INTENDED USE

The Rösch-Uchida Transjugular Liver Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

### CONTRAINDICATIONS

None known

### WARNINGS

- If resistance is encountered during advancement of Flexor sheath, assess cause of resistance and consider dilation of any restriction identified or consider alternate treatment strategy. If Flexor sheath is advanced through resistance, force to remove the sheath will be higher, increasing the risk of sheath material or hub separation upon withdrawal.
- Reinsertion of dilator prior to removal of Flexor sheath increases the strength of the sheath and lessens the risk of device separation. If resistance is anticipated or encountered during withdrawal of Flexor sheath, consider carefully reinserting the dilator prior to continuing removal.

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Avoid transcapsular puncture.
- Manipulation of products requires fluoroscopic guidance.
- The patient should be monitored (EKG, oximetry, blood pressure, pulse) during the procedure.
- Before introduction and periodically during the procedure, the transjugular introducer sheath should be flushed with saline.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- If the device is to remain in a vessel for an extended period of time, appropriate measures (i.e., flushing with heparinized intravenous flushing solution) should be taken to reduce risk of clot formation.

### POTENTIAL ADVERSE EVENTS

- Intraperitoneal and liver hemorrhage
- Carotid artery puncture
- Pneumothorax
- Air embolism
- Puncture site hematoma
- Cardiac arrhythmia
- Arteriovenous fistula

### INSTRUCTIONS FOR USE

1. Using Micropuncture or standard 18 gage access techniques, introduce an appropriate .035 inch wire guide into the inferior vena cava via jugular venous access.
2. Using a selective catheter and a wire guide of choice, catheterize the right hepatic vein (RHV) or the most adequate hepatic vein branch.
3. Leaving the wire guide in a safe, distal position, remove the catheter.
4. Advance the introducer sheath assembly over the wire guide into the hepatic vein. Remove the dilator from the sheath, leaving the wire guide in place.
5. Insert the stiffening cannula into the PTFE catheter, and introduce the catheter/cannula assembly over the wire into the introducer sheath, positioning the catheter/cannula assembly into the distal hepatic vein. Remove the wire.
6. Introduce the blue catheter with flexible trocar stylet assembly through the catheter/cannula assembly. **NOTE:** The stylet should not protrude past the PTFE catheter end hole.
7. Orient the catheter/cannula assembly inferiorly and rotate anteriorly (arrow baseplate on needle indicates direction of needle curve).

8. Wedging the catheter/cannula assembly against the vein wall, thrust the blue catheter/trocar stylet assembly in one movement forward through the hepatic parenchyma and toward the portal system.
9. Remove the trocar stylet from the blue catheter. Connect a syringe with contrast medium to the catheter, apply suction, and withdraw the catheter until blood is seen. Inject a small amount of contrast medium in order to confirm position in the portal system.
10. Introduce a wire guide into the portal branch and select the main portal vein (MPV).
11. Holding the stiffening cannula in stable position, advance the PTFE catheter and introducer sheath over the blue catheter and wire guide into position across the parenchymal tract.
12. Remove the stiffening cannula and blue catheter from the assembly. The wire guide and PTFE catheter may be removed or utilized as required for further diagnostic or interventional procedures.
13. Following completion of diagnostic and interventional procedures, remove the introducer sheath. Avoid applying traction to hub during removal. If resistance is anticipated or encountered during withdrawal of the introducer sheath, consider reinserting the dilator and removing the sheath and dilator as a unit.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd