COOK
MEDICAL

CE
0123

EN
3

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Bruksanvisning

T _ S A A _ R E V 1

SIDE-ARM ADAPTER MALE LUER LOCK AND TUOHY-BORST ADAPTER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Side-Arm Adapter is an endoscopic adapter designed to facilitate the passage of instruments ranging in size from 3 to 6 French in diameter during a procedure. The adapter can be connected to any standard female Luer lock fitting.

INTENDED USE

This device is intended to resist the backflow of fluid around an instrument inserted through the working channel of a ureteroscope, cystoscope, or other endoscope.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in urological procedures. Standard techniques should be employed.
- The potential effects of phthalates on children or pregnant/nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Secure the adapter onto the working channel of an endoscope via a standard female Luer lock fitting.
2. Introduce catheter or other desired instrument through the port of the adapter as necessary to complete the intended procedure.
3. After removal of the inserted device and completion of the procedure, carefully remove the adapter from the endoscope.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd