EN
3

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T _ S A F W G _ R E V 1

SAFETY WIRE GUIDE INTRODUCER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

- Radiopaque obturator
- Radiopaque sheath

NOTE: Placement of the obturator should be performed over a wire guide not to exceed 0.038" in diameter.

INTENDED USE

The Safety Wire Guide Introducer is used for antegrade placement of a safety wire guide prior to endourological procedures.

CONTRAINDICATIONS

This product is contraindicated in the presence of conditions which create unacceptable risk during endourological procedures.

WARNINGS

Do not force the device during the procedure; proceed with caution if excessive resistance is encountered.

PRECAUTIONS

The product is intended for use by physicians trained and experienced in urological procedures. Standard techniques for placement of ureteral introducers should be employed.

POTENTIAL ADVERSE EVENTS

Complications associated with endourological procedures may include, but are not limited to:

- Bleeding
- Urinary tract perforation
- Stricture
- Urinary retention.

INSTRUCTIONS FOR USE

After gaining percutaneous access to the renal pelvis with a wire guide:

1. Over a wire guide, dilate the nephrostomy tract to 12.0 French with the included obturator and sheath.
2. Leaving the sheath and wire guide in place, remove the obturator.
3. Insert an additional wire guide through the sheath, advancing it to an appropriate distance down the ureter.
4. Remove the sheath leaving both wire guides.
5. Additional devices compatible with endourological wire guides may now be placed as necessary to complete the intended procedure.
6. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not

use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd