COOK

MEDICAL

CE
0123

EN
3

Entuit™ Secure Adjustable Gastrointestinal Suture Anchor Set

Instructions for Use

DA
4

Entuit™ Secure justerbart gastrointestinalt suturankersæt

Brugsanvisning

DE
6

Entuit™ Secure verstellbares gastrointestinales Fadenanker-Set

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EL
9

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
1

ENTUIT™ SECURE ADJUSTABLE GASTROINTESTINAL SUTURE ANCHOR SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Entuit Secure Adjustable Gastrointestinal Suture Anchor Set consists of two or three introducer needles with preloaded anchors and a .035 inch wire guide.

INTENDED USE

The Entuit Secure Adjustable Gastrointestinal Suture Anchor Set is intended for anchoring the anterior wall of the stomach to the abdominal wall prior to introduction of interventional catheters.

CONTRAINDICATIONS

None known

WARNINGS

- Do not leave the red plunger in place after the procedure is finished. Leaving the plunger in place may allow the anchor to become repositioned, resulting in leakage into the peritoneal cavity.
- Gastropexy site should be checked daily for early signs of pressure ulcers, including discoloration and tenderness around gastropexy site.
- Patients with ascites may be at an increased risk of complications, including pressure ulcers. Additional site monitoring should be implemented for these patients.
- If signs of pressure ulcer development are detected at gastropexy site, use the extra red plunger provided in the set to adjust the tension of the retention bolster.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of gastrointestinal suture anchors should be employed.
- Manipulation of the product requires imaging guidance.
- Withdrawal or manipulation of the wire guide through the needle tip may result in breakage.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Entuit Secure Adjustable Gastrointestinal Suture Anchor Set is **MR Conditional** according to ASTM F2503. A patient with this metallic anchor and suture retention mechanism may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient of 1600 Gauss/cm (16.0 T/m) or less
- Normal Operating Mode: Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence)

MRI-Related Heating

Under the scan conditions defined above, the Entuit Secure Adjustable Gastrointestinal Suture Anchor Set is not expected to result in a temperature rise of more than 1.8 °C after 15 minutes of continuous scanning.

Image Artifact

MR image quality might be compromised if the area of interest is within approximately 78 mm of the position of the suture retention mechanism of the Entuit Secure Adjustable Gastrointestinal Suture Anchor Set as found during nonclinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a MR 3.0 Tesla system. Therefore, it might be necessary to optimize MR imaging parameters for the presence of this metallic anchor and suture retention mechanism.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

INSTRUCTIONS FOR USE

1. Insufflate the stomach with approximately 500 mL of air through a nasogastric tube.
2. Localize the puncture site and administer local anesthetic to the peritoneum. **NOTE:** The preferred puncture site into the gastric lumen is at the junction of the left anterior rib cage and the lateral border of the rectus abdominis muscle, provided there is no interposition of the liver or colon.
3. Attach a 10 mL syringe, half filled with contrast medium or sterile water, to the preloaded introducer needle.
4. Using endoscopic or fluoroscopic visualization, advance the introducer needle through the anterior abdominal wall and then rapidly thrust the needle another 2-5 cm into the air-filled stomach. **NOTE: Air should be freely aspirated at this point. Injection of contrast medium under fluoroscopy will aid in confirming proper intragastric needle position.**
5. Remove the 10 mL syringe from the needle hub.
6. While maintaining slight tension on the trailing suture, introduce the distal spring coil of the wire guide into the needle and use it to push the anchor out of the needle into the stomach cavity. **Maintain slight tension on the suture during anchor introduction.**
7. Remove the introducer needle over the wire guide.
8. With the wire guide still in position, apply traction to the suture to pull the anterior wall of the stomach against the abdominal wall.
9. While maintaining traction on the suture anchor, remove the wire guide.
10. While compressing the red plunger against the external gastropexy bolster, slide the external gastropexy bolster down the suture to achieve the desired tension. Bolster position may be adjusted to relieve or increase tension as needed, as long as the red plunger is compressed against the bolster. The gastropexy bolster should contact the skin but not create any visible skin indentation when at optimal tension.
11. Once the optimal tension is achieved, stabilize the external bolster manually and pull the red plunger away from the main body of the bolster. **NOTE: Do not** leave the plunger attached to the bolster after the procedure is complete. If further adjustment of gastropexy tension is required, reattach the plunger or use the spare plunger included with the device to make the adjustment. **(Fig. 1)**
12. Using fluoroscopic or endoscopic visualization, place the gastropexy anchors in a 2-, 3-, or 4-point configuration. Be sure the anchors are positioned such that they will not come into contact with the inflated gastrostomy catheter balloon. **NOTE:** Use of a single-point gastropexy anchor is not recommended.
13. Cut the suture, leaving enough length attached to the gastropexy bolster to allow for future adjustments.
14. Introduce an .035 inch wire guide-compatible access needle through the skin and into the gastric lumen centered within the gastropexy configuration.
15. Confirm gastric position with either endoscopic or fluoroscopic visualization. Syringe aspiration of air and/or injection of contrast medium may be used to help confirm gastric positioning.
16. Advance an .035 inch wire guide through the access needle into the gastric lumen. Confirm position within the gastric lumen with endoscopic or fluoroscopic visualization.
17. Maintaining wire guide position, remove the access needle.
18. The wire guide may now be used to facilitate introduction of fascial dilators, peel-away sheaths, and the gastrostomy catheter. **NOTE:** If desired, a scalpel may be used for a dermatotomy at the puncture site to ease placement of dilators, sheaths, and catheters.

NOTE: The sutures may be left in place for 10 to 14 days while tract formation is completed, or they may be cut after gastrostomy catheter placement. This releases the anchors into the stomach, allowing their passage via the gastrointestinal system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd