EN  
4

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

1

2

## S~CURVE® URETHRAL DILATORS AND S~CURVE® URETHRAL DILATOR SET

**Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The S~Curve® Urethral Dilators with AQ® Hydrophilic Coating are manufactured from radiopaque tubing with a length of 37 cm. The dilators are shaped in an S-Curve formation and have a tapered tip. They are available as individual dilators or as a set. The individual dilators are available in 22, 24, and 26 Fr sizes. The set includes 7 dilators in sizes 8, 10, 12, 14, 16, 18, 20 Fr. AQ® Hydrophilic coating is applied to the surface of the dilator for ease of entry to the access site. A single sideport is placed near the tapered tip. This sideport is angled toward the tapered tip and placed on the inside of the curve. The device French size is indicated on the end of the dilator without a taper. The French size label aligns with the sideport to act as an indicator of the sideport's location during insertion.

### Performance Characteristics

- S-curve shape helps traverse and support the natural curvature of the male urethra.
- Radiopaque material to aid in positioning under imaging techniques.
- Tapered tip to facilitate insertion into anatomy.
- Hydrophilic coating to enhance lubricity of device.

### Device Compatibility

These devices are compatible with up to a 0.038 inch diameter wire guide.

### Patient Population

The target patient population for these devices are men with urological conditions who require dilation of the urethra.

### Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

### Contact with Body Tissue

These devices interact with the urethral tract including the urethra and bladder.

### Operating Principle

These devices are intended to be introduced into the urethra over a previously placed wire guide. To accommodate the patient's anatomy, the dilator should be inserted with the penis in an upward position and the dilator oriented with the French size label on the end without the tapered tip facing upward. Alternatively, the dilator can be inserted with the penis in a downward position and the dilator oriented with the French size label on the end without a tapered tip facing downward. In the latter placement method, when the tapered tip is at the bulbar urethra, the dilator should be rotated 180° before advancing into the bladder.

The dilators are intended to be introduced sequentially from smallest to largest to gradually dilate the urethra. The sideport allows for urine to flow into the dilator to confirm placement in the bladder.

### INTENDED USE

These devices are S-Curve dilators which are used for dilation of the male urethra in adult patients by retrograde insertion of the devices.

### INDICATIONS FOR USE

These devices are used for dilating urethral strictures in adult patients.

## CLINICAL BENEFITS

Urethral dilation to widen the diameter of the urethra.

## CONTRAINDICATIONS

None

## WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

## PRECAUTIONS

No known Precautions

## POTENTIAL ADVERSE EVENTS

- Bleeding
- False passage
- Fistula
- Hematoma
- Hematuria
- Infection
- Pain
- Perforation, urethral
- Rectal injury
- Urethral injury
- Lower urinary tract symptoms (poor/diminished urinary stream, dysuria, incomplete emptying of the bladder, (recurrent) urinary tract infections, overflow incontinence, frequency, hesitancy, urgency, urinary retention, increasing post void residual, and/or pain)

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

## INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

## INSTRUCTIONS FOR USE

**NOTE:** Prior to use, immerse the AQ® dilator in sterile water or isotonic saline to allow the hydrophilic surface to absorb liquid and become lubricious. This will ease placement under standard conditions.

1. Under direct vision or fluoroscopic control, position a wire guide through the urethra and into the bladder.
2. Introduce the smallest appropriately sized dilator, using one of the two methods below:
  - a. Hold the penis in an upward position. (**Fig 1**)
  - b. Orient the dilator so the French Size label on the end without a tapered tip faces upward, and introduce the dilator over the previously placed wire guide.

OR

- a. Hold the penis in a downward position. (**Fig 2**)
- b. Orient the dilator so the French Size label on the end without a tapered tip faces downward, and introduce the dilator over the previously placed wire guide.

- c. Advance the dilator to the level of the bulbar urethra.
- d. With the dilator positioned at the level of the bulbar urethra, rotate the dilator 180 degrees so the dilator is now oriented with the French Size label facing upward.
3. Continue to advance the dilator into the bladder. Urine drainage will confirm proper placement.
4. While maintaining the wire guide position, remove the dilator currently in place and replace it with the next appropriately sized dilator. Repeat steps 2-3 as needed, progressing from smallest to largest appropriately sized dilator until desired urethral dilation is achieved.

## **DISPOSAL OF DEVICES**

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

## **REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

## **PATIENT COUNSELING INFORMATION**

Please inform the patient as necessary of the relevant warnings, precautions contraindications, measures to be taken, and limitations of use that the patient should be aware of.

## **SERIOUS INCIDENT REPORTING**

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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**Cook Medical Europe Ltd.**
**William A. Cook Australia Pty. Ltd.**