COOK

MEDICAL

CE
0088

EN
3

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Bruksanvisning

T _ S D A S _ R E V 2

KOLENDA SALIVARY ACCESS INTRODUCER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Salivary Access Introducers are designed to perform as a guiding sheath and/or introducer sheath.

INTENDED USE

Salivary Access Introducers are intended for introduction of diagnostic or interventional devices in the submandibular salivary duct.

CONTRAINDICATIONS

None known

WARNINGS

- This product is not intended for use in the vasculature.
- This product is not intended for use in the parotid or sublingual salivary glands.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of salivary access sheaths should be employed.
- Ensure compatibility of instruments to be introduced through the introducer. All instruments or catheters used with this product should move freely through the hub and sheath. Damage to the hub/introducer may result when the fit is tight.
- When inserting, manipulating or withdrawing a device through an introducer, always maintain introducer position.
- When puncturing, suturing or incising the tissue near the introducer, use caution to avoid damaging the introducer.
- Do not attempt to insert or withdraw the wire guide and/or introducer if resistance is felt.
- Do not attempt if proximal duct stones have been previously identified.
- Do not attempt in the presence or suspected presence of unfavorable patient anatomy.

INSTRUCTIONS FOR USE

1. Upon removal from package, ensure the inner diameter (ID) of the introducer is appropriate for the maximum diameter of the instrument or device to be introduced.
2. Remove dilator from the distal end of the sheath and insert it completely into the proximal fitting.
3. If desired, lubricate the outside diameter of the sheath with water or other suitable lubricant, to ease introduction.
4. Introduce the dilator and sheath into the duct over an appropriately sized wire guide. If resistance is encountered, do not force; damage to the sheath and/or duct may occur.

NOTE: The duct may be pre-dilated to increase ease of introduction.

5. Remove the dilator from the sheath by pulling the dilator fitting to the side to release the dilator from the sheath.

NOTE: Wire guide may be removed if it is not being used for additional guidance of interventional devices.

6. Once the sheath is inserted, the side-arm connection may be connected to suction to remove any excessive fluids, stone fragments, or irrigation.

NOTE: If stone fragments occlude the suction path, close or shut off suction and continue with irrigation to flush out the stone fragments. Once cleared, resume suction as desired.

NOTE: If resistance is felt when advancing instruments through the sheath or withdrawing the sheath, do not force; damage to the sheath and/or duct may occur.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-pouch packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2015
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.