- EN**
3 **Sof-Flex® Double Pigtail Stent Set**
Instructions for Use
- CS**
6 **Souprava double pigtailového stentu Sof-Flex®**
Návod k použití
- DA**
9 **Sof-Flex® stentsæt med dobbelt grisehale**
Brugsanvisning
- DE**
12 **Sof-Flex® Doppelpigtail-Stentset**
Gebrauchsanweisung
- EL**
15 **Σετ ενδοπρόσθεσης διπλού σπειροειδούς άκρου Sof-Flex®**
Οδηγίες χρήσης
- ES**
19 **Equipo de stent de doble pigtail Sof-Flex®**
Instrucciones de uso
- FR**
22 **Set d'endoprothèse double pigtail Sof-Flex®**
Mode d'emploi
- HU**
25 **Sof-Flex® kettős pigtail sztentkészlet**
Használati utasítás
- IT**
29 **Set con stent a doppio pigtail Sof-Flex®**
Istruzioni per l'uso
- NL**
32 **Sof-Flex® stentset met dubbele pigtail**
Gebruiksaanwijzing
- NO**
35 **Sof-Flex® stentsett med dobbel pigtail**
Bruksanvisning
- PL**
38 **Zestaw stentu Sof-Flex® z dwiema końcówkami pigtail**
Instrukcja użycia
- PT**
42 **Conjunto de stent de duplo pigtail Sof-Flex®**
Instruções de utilização
- SV**
45 **Sof-Flex® stentset med dubbla pigtails**
Bruksanvisning

SOF-FLEX® DOUBLE PIGTAIL STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

Set Contains:

- Radiopaque stent
- Stent positioner
- Wire guide (*size and configuration may vary*)

NOTE: Set components may vary.

NOTE: If no wire guide is provided with this set, the following is recommended:

3.0 Fr. accepts .018"	8.2 Fr. accepts .038"
4.7 Fr. accepts .035"	10 Fr. accepts .038"
6.0 Fr. accepts .038"	12 Fr. accepts .038"
7.0 Fr. accepts .038"	

INTENDED USE

The Sof-Flex Double Pigtail Stent Set is used for temporary internal drainage from the ureteropelvic junction to the bladder. Sof-Flex Double Pigtail Stent Sets have been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques.

3 French stents are indicated for pediatric patients.

CONTRAINDICATIONS

- Do not attempt stent placement in a patient with suspected ureteral avulsion,
- Do not use, when, in the judgement of the physician, such a procedure would be contrary to the best interest of the patient.

WARNINGS

None known

PRECAUTIONS

- Manipulation of the wire guide requires appropriate imaging control. Use caution not to force or overmanipulate the wire guide when gaining access.
- When using a wire guide through a metal cannula/needle, use caution as damage may occur to the outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- These wire guides are not intended for PTCA use.
- Remove tether if stent is to remain indwelling for more than 14 days.
- Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to each patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- These stents must not remain indwelling more than six months. If the patient's status permits, the stent may be replaced with a new stent.
- These stents are not intended as permanent indwelling devices.

- Do not force components during removal or replacement. If resistance is encountered, stop. Determine the cause of the resistance before proceeding.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- Hydrophilically coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Service department at the address/phone listed at www.cookmedical.com.
- The potential effects of phthalates on children or pregnant or nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.
- The patients should be instructed in terms that they understand to inform the physician if they are experiencing any pain, cloudy urine, bladder irritation or any sign or symptoms that they are having difficulty with urination.
- If the stent is placed in a child the parents should be given these instructions.

POTENTIAL ADVERSE EVENTS

- Extravasation
- Occlusion
- Migration
- Hemorrhage
- Sepsis
- Perforation of the urinary tract
- Peritonitis
- Encrustation
- Urinary tract infection
- Loss of renal function

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Sof-Flex Double Pigtail Stents are **MR Conditional** according to ASTM F2503. A patient with this device can be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or 1.5 Tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 4.0 W/kg (Normal Operating or First Operating Mode)

Nonclinical testing demonstrated that the Sof-Flex Double Pigtail Stents are not expected to induce RF heating greater than that of biological tissue.

The image artifact extends approximately 6.8 mm from the Sof-Flex Double Pigtail Stents as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

INSTRUCTIONS FOR USE

Wire Guide Preparation

NOTE: If your device did not include a wire guide, then the instructions below are not applicable. If your device did include a wire guide, please refer to the instructions below for proper use.

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. If using a hydrophilically coated wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment or, after extended use, replace it with a new hydrophilically coated wire guide.
3. Remove the wire guide from the holder.
4. If using a movable core wire guide, adjust mandril position by advancing or withdrawing the mandril within the wire guide to achieve the required degree of flexibility.

Endoscopic Stent Placement

1. Pass the flexible wire guide tip beyond the obstruction to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using the wire guide and an open-ended ureteral catheter in combination.

NOTE: Endhole size and length of the device must be taken into consideration to ensure proper fit between the wire guide and the device.

2. Using a baseline pyelogram, estimate the length of stent needed by taking the ureteral length and adding 1 cm. Accurate measurement enhances drainage efficiency and patient comfort.
3. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
4. Watch for the distal end of the stent at the ureterovesical junction. At that point, halt advancement of the stent. As an assistant removes the wire guide, hold the stent in position with the positioner. The stent pigtail will form spontaneously. Carefully remove the positioner from the cystoscope.

NOTE: If necessary, final adjustment can be made with endoscopic forceps or the tether. Fluoroscopy facilitates stent placement; however, standard radiography may be used.

Stent Removal

- The stent may be removed by gentle withdrawal traction using endoscopic forceps or tether.
- Retrieve cystoscopically by gently pulling on the stent, and if resistance is encountered during removal, fluoroscopically determine stent position and cause of resistance.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M. F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single End Ureteral Stent," *Urology*, 13 (1979), 516-520.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd