EN  
3

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
Figure A


Figure B

## SOF-FLEX® MULTI-LENGTH URETERAL STENT SET

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

### DEVICE DESCRIPTION

Set Contains:

- Radiopaque stent
- Stent positioner
- Appropriate diameter stainless steel wire guide

**NOTE:** The multi-length stents will accommodate ureters from approximately 4 - 32 cm long.

**NOTE:** Set components may vary.

### RECOMMENDED WIRE GUIDE

Stent Size	Accepts Wire Guide
3.0 French	.018 inch
4.7 French	.035 inch
6.0 French	.038 inch
7.0 French	.038 inch

### INTENDED USE

The Sof-Flex Multi-Length Ureteral Stent Set is intended for temporary internal drainage from the ureteropelvic junction to the bladder. Sof-Flex Multi-Length Ureteral Stent Sets have been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques.

3 French stents are indicated for pediatric patients.

### CONTRAINDICATIONS

- Do not attempt stent placement in a patient with suspected ureteral avulsion,
- Do not use, when, in the judgement of the physician, such a procedure would be contrary to the best interest of the patient.

### WARNINGS

Formation of knots in multi-length stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.

### PRECAUTIONS

- Manipulation of the wire guide requires appropriate imaging control. Use caution not to force or overmanipulate the wire guide when gaining access.
- When using a wire guide through a metal cannula/needle, use caution as damage may occur to the outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- These wire guides are not intended for PTCA use.
- Complications of ureteral stent placement are documented in literature. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.

- These stents must not remain indwelling more than six (6) months. If the patient's status permits, the stent may be replaced with a new stent.
- The tether should be removed if the stent is to remain indwelling longer than 14 days.
- These stents are not intended as permanent indwelling devices.
- Do not force components during removal or replacement. If resistance is encountered, stop. Determine the cause of the resistance before proceeding.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- The potential effects of phthalates on children or pregnant or nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.
- Hydrophilically coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- The patients should be instructed in terms that they understand to inform the physician if they are experiencing any pain, cloudy urine, bladder irritation or any sign or symptoms that they are having difficulty with urination.
- If the stent is placed in a child the parents should be given these instructions.

#### POTENTIAL ADVERSE EVENTS

- Extravasation
- Occlusion
- Migration
- Hemorrhage
- Sepsis
- Perforation of the urinary tract
- Peritonitis
- Encrustation
- Urinary tract infection
- Loss of renal function

#### MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Sof-Flex Multi-Length Ureteral Stents are **MR Conditional** according to ASTM F2503. A patient with this device can be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or 1.5 Tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 4.0 W/kg (Normal Operating or First Operating Mode)

Nonclinical testing demonstrated that the Sof-Flex Multi-Length Ureteral Stents are not expected to induce RF heating greater than that of biological tissue.

The image artifact extends approximately 6.8 mm from the Sof-Flex Multi-Length Ureteral Stents found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

## INSTRUCTIONS FOR USE

### Wire Guide Preparation

**NOTE:** If your device did not include a wire guide, then the instructions below are not applicable. If your device did include a wire guide, please refer to the instructions below for proper use.

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. If using a hydrophilically coated wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment or, after extended use, replace it with a new hydrophilically coated wire guide.
3. Remove the wire guide from the holder.
4. If using a movable core wire guide, adjust mandril position by advancing or withdrawing the mandril within the wire guide to achieve the required degree of flexibility.

### Endoscopic Sof-Flex Multi-Length Stent Placement

1. Pass the flexible wire guide tip beyond the obstruction to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using the wire guide and an open-end ureteral catheter in combination. Remove catheter before attempting to place stent.
2. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
3. When the renal pelvis has been entered (**Figure A**) (this can be confirmed fluoroscopically), gently remove the wire guide while firmly holding the positioning catheter. The stent coil will form spontaneously in the renal pelvis. Carefully remove the positioning catheter (**Figure B**).
  - For ureters longer than 22 cm, endoscopic forceps can be used to adjust the stent position, leaving approximately equal coils in the renal pelvis and bladder.

**NOTE:** If problems occur using this device, please call your Cook Incorporated sales representative or contact our Customer Quality Assurance department at the address/phone number listed within this booklet.

### Stent Removal

- The stent may be removed by gentle withdrawal traction using endoscopic forceps or tether.
- Retrieve cystoscopically by gently pulling on the stent, and if resistance is encountered during removal, fluoroscopically determine stent position and cause of resistance.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd