- EN**
3
Flexor® Shuttle® - SL Introducer
Flexor® Shuttle Select® - Introducer
Instructions for Use
- BG**
5
Въвеждащо устройство Flexor® Shuttle® - SL
Въвеждащо устройство Flexor® Shuttle Select®
Инструкции за употреба
- CS**
9
Zavaděč Flexor® Shuttle® - SL
Zavaděč Flexor® Shuttle Select®
Návod k použití
- DA**
12
Flexor® Shuttle® - SL Indfører
Flexor® Shuttle Select® Indfører
Brugsanvisning
- DE**
15
Flexor® Shuttle® - SL Einführschleuse
Flexor® Shuttle Select® Einführschleuse
Gebrauchsanweisung
- EL**
18
Εισαγωγέας Flexor® Shuttle® - SL
Εισαγωγέας Flexor® Shuttle Select®
Οδηγίες χρήσης
- ES**
22
Introductor Flexor® Shuttle® - SL
Introductor Flexor® Shuttle Select®
Instrucciones de uso
- FR**
25
Introduceur Flexor® Shuttle® - SL
Introduceur Flexor® Shuttle Select®
Mode d'emploi
- HU**
28
Flexor® Shuttle® - SL bevezetőeszköz
Flexor® Shuttle Select® bevezetőeszköz
Használati utasítás
- IT**
31
Introduttore Flexor® Shuttle® - SL
Introduttore Flexor® Shuttle Select®
Istruzioni per l'uso
- LT**
34
Intubatorius Flexor® Shuttle® - SL
Intubatorius Flexor® Shuttle Select®
Naudojimo instrukcija
- LV**
37
Flexor® Shuttle® - SL ievadītājs
Flexor® Shuttle Select® ievadītājs
Lietošanas instrukcija
- NL**
40
Flexor® Shuttle® - SL Introducer
Flexor® Shuttle Select® - Introducer
Gebruiksaanwijzing
- NO**
43
Flexor® Shuttle® - SL innføringsenhet
Flexor® Shuttle Select® - innføringsenhet
Bruksanvisning

- PL**
46 **Introduktor Flexor® Shuttle® - SL**
Introduktor Flexor® Shuttle Select®
Instrukcja użycia
- PT**
49 **Introdutor Flexor® Shuttle® - SL**
Introdutor Flexor® Shuttle Select®
Instruções de utilização
- SK**
53 **Zavádzač Flexor® Shuttle® - SL**
Zavádzač Flexor® Shuttle Select®
Návod na použitie
- SL**
56 **Uvajalo Flexor® Shuttle® - SL**
Uvajalo Flexor® Shuttle Select®
Navodila za uporabo
- SV**
58 **Flexor® Shuttle® - SL införare**
Flexor® Shuttle Select® införare
Bruksanvisning
- TR**
61 **Flexor® Shuttle® - SL İntrodüser**
Flexor® Shuttle Select® İntrodüser
Kullanma Talimatı

FLEXOR® SHUTTLE® - SL INTRODUCER**FLEXOR® SHUTTLE SELECT® - INTRODUCER**

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Flexor® Shuttle® - SL Introducers and Flexor® Shuttle Select® - Introducers are designed to perform as a guiding sheath and/or introducer sheath. These devices incorporate a Flexor shaft with a hemostasis valve and are provided with one or more dilators. These devices are available in various sizes, lengths, and configurations. Configurations include differences in shaft construction, such as varying shaft stiffness and distal tip material, dilator material, and hydrophilic coating. Refer to the product label for product specifications (i.e., sheath inner diameter, dilator inner diameter).

Performance Characteristics

- Stainless steel coiling incorporated into sheath material for resistance to kinking.
- Radiopaque marker band incorporated into sheath material for distal tip positioning under imaging techniques.
- Hemostasis valve to minimize leakage.
- Rounded sheath tip and tapered dilator tip to facilitate insertion into patient anatomy.
- Multiple durometers of nylon sheath shaft material for pushability on proximal end and flexibility on distal end.
- If applicable, hydrophilic coating to enhance lubricity of device.

Device Compatibility

- Standard luer syringe
- Heparinized saline
- Contrast media
- If applicable, 0.018" Wire Guide
- If applicable, 0.035" Wire Guide
- If applicable, 0.038" Wire Guide
- Vascular devices with a smaller outer diameter than the specified inner diameter of the sheath

Patient Population

The intended patient population are adult patients with a medical condition whose treatment a physician has determined requires access to the vasculature. The introducer is not intended to be placed in the coronary or neuro vasculature.

Intended Users

These devices are intended for use by physicians trained and experienced in diagnostic and interventional techniques. Endovascular access procedures that require the use of the Flexor Introducers and Guiding Sheaths should only be performed by licensed and trained clinical practitioners, such as vascular surgeons, interventional radiologists, and interventional cardiologists.

Contact with Body Tissue

These devices are for single-use, categorized as external communicating devices in contact with blood and limited (< 24 hours) human contact.

Operating Principle

These devices include a single lumen sheath ranging from 4 Fr to 8 Fr, and depending on the device specification, one or two dilators. The sheaths are available in lengths from 80 cm to 110 cm. The sheath is made of a nylon material which is supported by coils to provide flexibility while minimizing kinking or compression. Additionally, hydrophilic coating is applied to the distal portion of the sheath and dilator. There is a radiopaque marker band located approximately 2-3 mm from the end of the sheath to provide precise positioning. The dilator tip extends past the sheath to facilitate positioning of the introducer assembly. The proximal end of the sheath is comprised of a Tuohy-Borst haemostasis valve which accepts a range of devices while minimizing blood reflux, as well as a stopcock on the side-arm.

INTENDED USE

The intended purpose of these devices is to serve as a pathway for other devices to be inserted into the vascular system. The introducer is not intended to be placed in the coronary or neuro vasculature.

INDICATIONS FOR USE

These devices are indicated when therapeutic or diagnostic devices are to be introduced into the vasculature. The introducer is not intended to be placed in the coronary or neuro vasculature.

CLINICAL BENEFITS

These devices provide a minimally-invasive pathway to access the vascular system for diagnostic or therapeutic device placement.

CONTRAINDICATIONS

- Use in pregnant women.
- Use in pediatric patients.
- Presence of pseudoaneurysm, hematoma, or active infection in the area of the access vessel.

WARNINGS

- If resistance is encountered during advancement of Flexor sheath, assess cause of resistance and consider dilation of any restriction identified or consider alternate treatment strategy. If Flexor sheath is advanced through resistance, force to remove the sheath will be higher, increasing the risk of sheath material or hub separation upon withdrawal.
- Reinsertion of dilator prior to removal of Flexor sheath increases the strength of the sheath and lessens the risk of device separation. If resistance is anticipated or encountered during withdrawal of Flexor sheath, consider carefully reinserting the dilator prior to continuing removal.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and or to re-use may lead to chemical contamination, device failure and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.
- In order to ensure device compatibility, choose a sheath size large enough to accommodate the maximum outer diameter of any devices that will be placed through the sheath.
- All interventional or diagnostic instruments used with this product should move freely through the valve and sheath to avoid damage.
- Do not attempt to heat or reshape the device.
- When inserting, manipulating or withdrawing a device through the sheath, always maintain sheath position.
- When puncturing, suturing or incising the tissue near the sheath, use caution to avoid damaging the sheath.

POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of an introducer set:

- Access site complications (including hematoma, infection, local inflammation, local pain, pseudoaneurysm)
- Bleeding
- Death
- Distal embolization
- Extravasation
- Hematoma
- Infection
- Injury to adjacent structures
- Ischemia
- Nerve injury
- Stroke or cerebrovascular event
- Vascular thrombosis or occlusion
- Vessel injury (including dissection, laceration, perforation, vasospasm)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

Sheath Introduction:

1. Ensure that the inner diameter (ID) of the sheath is appropriate for the maximum diameter of the instruments to be introduced.
2. Using the side-arm of the valve, flush the sheath by filling the sheath assembly completely with heparinized saline.
3. Flush the dilator with heparinized saline.
4. Insert the dilator completely into the sheath. If the sheath has a Tuohy-Borst valve, tighten the valve around the dilator.
5. If the device has hydrophilic coating, activate the coating by wetting the outer surface of the device with heparinized saline. **NOTE:** For best results, maintain wetted condition of device during placement.
6. Using standard Seldinger technique, access the target vessel with the appropriate needle.
7. Insert a wire guide into the vessel through the needle.
8. Leave the wire guide in place and remove the needle.
9. Insert the dilator/sheath assembly over the wire guide.
10. Remove the wire guide and dilator, then aspirate and flush through the sheath side-arm.
11. Insert appropriately sized device(s) as needed.

Sheath Removal

1. Insert a wire guide until its tip extends at least 10 cm past the tip of the sheath.
2. Remove the sheath. Avoid applying traction to hub during removal. If resistance is anticipated or encountered during withdrawal of the Flexor sheath, consider reinserting the dilator and removing the sheath and dilator as a unit.
3. Remove the wire guide.

DISPOSAL OF DEVICES

Upon completion of procedure, dispose of all used products and associated packaging per local medical waste disposal procedures.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device - (0827002CINC050-032_AHS) - to find the Summary of Safety and Clinical Performance (SSCP) on the Eudamed website. When the Eudamed is available, use the following link: ec.europa.eu/tools/eudamed

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2022
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd