EN 3 **Salle Intraoperative Pyeloplasty Stent Set**

Instructions for Use

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Bruksanvisning

T _ S I P S _ R E V 1

SALLE INTRAOPERATIVE PYELOPLASTY STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Multi-length Sof-Flex® stent and nephrostomy catheter with positioning stylet
- TFE-coated wire guide
- Connecting tube
- Tuohy-Borst adapter with cap
- Silicone retention disc with pull tie

NOTE: Set components may vary.

INTENDED USE

Used to establish and maintain nephrostomy drainage following pyeloplasty.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The maximum indwelling time for this device is four weeks.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Pass the positioning stylet through the open renal pelvis; push it through the desired calyx.
2. If the stylet “tents” the renal capsule, make a small incision over the stylet to allow for smooth passage of the nephrostomy catheter through the parenchyma.
3. Utilizing the wire guide, place the stent down the ureter to the bladder in a typical intraoperative fashion.
4. The posterior aspect of the anastomosis should be completed and the stent introduced and positioned in the bladder. If difficulties are encountered while passing the distal end through the ureterovesical junction, the lower coil can be excised, leaving the straight portion of the stent in the ureter below the anastomosis.
5. Once this part of the procedure is complete, check for full coil in the renal pelvis, then bring the nephrostomy catheter portion of the stent out through a separate external incision.
6. Remove the positioning stylet by transecting the catheter at the most proximal ink mark.
7. Attach the Tuohy-Borst adapter to the catheter.
8. Secure the catheter to the skin with the retention disc and pull tie.
9. Use the connecting tube to connect the catheter to a drainage bag.

NOTE: The correct position of the stent may be confirmed by radiographic means, utilizing the radiopaque marks at the renal and bladder coils.

NOTE: Before discharge, the exterior part of the stent can be partially excised, leaving a 3-4 cm segment closed and secured to the skin.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.