EN
3 **Salle Intraoperative Pyeloplasty Stent Set**

Instructions for Use

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Bruksanvisning

T _ S I P S _ R E V 3

SALLE INTRAOPERATIVE PYELOPLASTY STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Multi-length Sof-Flex® stent with nephrostomy catheter and positioning stylet
- PTFE-coated, stainless steel wire guide
- Connecting tube
- Tuohy-Borst adapter with cap
- Silicone retention disc with pull tie

INTENDED USE

This device is intended to establish and maintain nephrostomy and internal drainage following pyeloplasty for pediatric patients. This device is intended for patients 6 months and older.

CONTRAINDICATIONS

- Absence or severe reduction of kidney function
- Patients with a ureteral/bladder fistula to the bowel or vagina
- Current untreated urinary tract infection
- Oliguria or anuria

WARNINGS

Formation of knots in multi-length stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.

PRECAUTION

The maximum indwelling time for this device is four weeks.

POTENTIAL ADVERSE EVENTS

- Infection
- Migration
- Encrustation
- Breakage/Erosion
- Perforation
- Urgency/frequency
- Mucosal irritation or inflammation
- Obstruction

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated the Salle Intraoperative Pyeloplasty Stent is **MR conditional**. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0T/m) or less

Nonclinical testing demonstrated that the Salle Intraoperative Pyeloplasty Stent is not expected to induce RF heating greater than that of biological tissue.

In nonclinical testing, the image artifact caused by the device extends approximately 2 mm from the Salle Intraoperative Pyeloplasty Stent when imaged with a gradient echo pulse sequence and a 3.0 Tesla MR system.

INSTRUCTIONS FOR USE

1. The posterior aspect of the anastomosis should be completed prior to the stent being introduced.
2. Insert the wire guide through a sideport of the stent to straighten the distal pigtail. **NOTE:** Utilize the most proximal sideport in the shaft of the stent and not a sideport in the proximal pigtail.
3. In a typical intraoperative fashion, place the stent down the ureter and into the bladder. **NOTE:** If difficulties are encountered while passing the distal pigtail through the ureterovesical junction or to accommodate patient anatomy, the distal pigtail can be transected, leaving the straight portion of the stent in the ureter below the anastomosis.
4. With the distal pigtail positioned in the bladder, slightly remove the wire guide to allow the retentive coil of the stent to form. Once you have confirmed appropriate positioning of the bladder pigtail, the wire guide can be removed.
5. Pass the positioning stylet of the stent through the open renal pelvis; push it through the desired calyx.
6. If the stylet "tents" the renal capsule, make a small incision over the stylet to allow for smooth passage of the nephrostomy catheter portion of the stent through the parenchyma.
7. Verify full coil in the renal pelvis and confirm appropriate positioning of the pigtail.
8. Remove the positioning stylet by transecting the catheter at the most proximal ink mark.
9. Attach the Tuohy-Borst adapter to the catheter.
10. Secure the catheter to the skin with the retention disc and pull tie.
11. Use the connecting tube to connect the catheter to a drainage bag.

NOTE: The correct position of the stent may be confirmed fluoroscopically, utilizing the renal and bladder coils.

NOTE: Before discharge, the exterior part of the stent can be transected, leaving a 3-4 cm segment. This should be closed and secured to the skin.

NOTE: The stent may be removed by gentle withdrawal traction of the extracorporeal portion of the stent through the puncture site.

Troubleshooting:

For Customer support please contact Cook Medical at:

Cook Medical Inc.
P.O. Box 4195
Bloomington, IN 47402
U.S.A.
Phone: 812.339.2235
Alternative Phone: 800.457.4500 (toll free)
Fax: 800.554.8335 (toll free)
Email: customersupport@cookmedical.com

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd