**COOK**

**MEDICAL**

**EN**  
**3**

**Torq-Flex® Wire Guide TFE Coated**

Instructions for Use

**DA**  
**6**

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T \_ T H T F \_ R E V 2

## TORQ-FLEX® WIRE GUIDE TFE COATED

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The TFE coated Torq-Flex Wire Guide is etch marked in five centimeter increments from the distal tip.

### INTENDED USE

The TFE coated Torq-Flex Wire Guide is intended for placement of Turbo-Flo® Over-The-Wire Peripherally Inserted Central Venous Catheter.

### CONTRAINDICATIONS

None known

### WARNINGS

None known

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in venous access techniques. Standard techniques for placement of peripherally inserted central venous catheters should be employed.
- When exchanging a Twin Turbo-Flo® double lumen catheter, ensure that the wire guide is inserted into the hub marked "LARGE" (round lumen).

### PRODUCT RECOMMENDATIONS

This wire guide is recommended for use only with Turbo-Flo® peripherally inserted central venous catheters.

### INSTRUCTIONS FOR USE

Turbo-Flo Over-The-Wire Peripherally Inserted Central Venous Catheter placement using the Mircopuncture® Peel-Away® Introducer technique with the graduated TFE coated Torq-Flex Wire Guide.


Fig. 1

1. After prepping the access site, introduce a thinwall percutaneous entry needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to locate the tip of the needle during access into the vessel. (Fig. 1)


Fig. 2

2. Using fluoroscopic guidance, introduce the wire guide through the needle; advance the wire guide 15 to 20 cm into the vessel. (Fig. 2)


Fig. 3

3. Leaving the wire guide in place, withdraw the needle. If necessary, enlarge the puncture site with a #11 scalpel blade. (Fig. 3)


Fig. 4

4. Introduce a sheath-introducer assembly which has an internal diameter of .018 inch over the wire guide. With a twisting motion, advance the assembly into the vessel. (Fig. 4)


Fig. 5

- Using fluoroscopic control, determine the correct catheter length by advancing the wire guide to the desired catheter tip location. Once the wire guide tip is in proper position, and using the proximal portion of the wire that is external to the patient, measure the distance on the wire guide from the triple etch marks to the puncture site. (The distal 60 cm of the wire guide is marked in 5 cm increments, with the triple etch marks positioned proximally at 60 cm.) (**Fig. 5**) Trim the Turbo-Flo catheter to the appropriate (60 cm minus X cm) length. A Turbo-Flo catheter's original length is 60 cm.
- Continue catheter placement as designated by institutional protocol.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.