EN
3

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TRANSJUGULAR LIVER ACCESS AND BIOPSY SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Transjugular Liver Access and Biopsy Set consists of a transjugular biopsy needle and a catheter.

INTENDED USE

The Transjugular Liver Access and Biopsy Set is intended for transjugular liver access and biopsy in diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Manipulation of products requires fluoroscopic and/or ultrasound guidance.
- When introducing the needle through the catheter, the needle should be advanced over the wire guide in order to prevent catheter damage.

INSTRUCTIONS FOR USE

1. Using standard venous access technique, introduce a wire guide into the jugular vein and advance it into the inferior vena cava.
2. Advance an appropriately sized introducer over the wire guide and remove the dilator, leaving the sheath in place.
3. Insert the preassembled needle/catheter as a unit over the wire guide into the sheath; advance the catheter to protect the needle tip.
4. Introduce the needle/catheter assembly into the middle hepatic vein, then remove the wire guide and inject contrast to confirm position.

For Liver Biopsy

5. Attach a syringe with saline solution to the needle and advance the needle/catheter assembly to a wedged position while the patient holds his/her breath in inspiration.
6. Holding the catheter steady, advance the needle with a thrust approximately 2 cm into the liver parenchyma.
7. Immediately withdraw the needle, while applying suction to the syringe. Discontinue suction just before the needle tip re-enters the hepatic vein, ensuring that the biopsy specimen remains in the needle. Leaving the catheter in position, withdraw the needle from the catheter.
8. Gently flush out the specimen with saline solution.
9. If further specimens are required, reinsert the biopsy needle into the catheter over a wire guide and repeat steps 4-8.

For Portal Branch Access

10. Direct the needle/catheter assembly toward the selected portal branch and gently advance the needle.
11. Aspirate to confirm portal vein blood.
12. Reinsert the wire guide and advance it into the portal system.
13. Access to the portal vein has now been established, and the patient is ready for the planned procedure.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd