Transseptal Needle Set

Instructions for Use

EN
3

Transseptalt nålesæt

Brugsanvisning

DA
5

Transseptalkanülenet

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DE
7

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T - T S N C - R E V . 1

TRANSSEPTAL NEEDLE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Transseptal Needle Set consists of a transseptal needle and obturator.

INTENDED USE

The Transseptal Needle is intended for transseptal left heart access in both diagnostic and interventional procedures.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

Absolute

- Recent systemic embolization
- Left atrial myxoma
- Prothrombin time more than twice the normal value or INR > 2.2 in patients with anticoagulation treatment

Relative

- Severe deformation of the chest
- Congenital defects with markedly changed position of heart cavities
- Absence of inferior vena cava (superior vena cava approach)
- Extreme atrial dilation

WARNINGS

- The transseptal needle must only be rotated counterclockwise to find the puncture site, in order to minimize the risk of puncturing the ascending aorta.
- Never puncture the septum without monitoring the pressure at the needle tip.
- Do not heparinize the patient until the catheter has been safely inserted into the left atrium.

PRECAUTIONS

Manipulation of product requires biplane fluoroscopic control.

PRODUCT RECOMMENDATIONS

For transseptal left heart catheterization, a Mullins introducer sheath assembly is recommended. If a Mullins introducer is used, it is important to measure the distance between the hub of the introducer and the metallic arrow of the transseptal needle before use. The distance should be measured with the needle tip positioned just proximal to the tip of the introducer. This distance, which is usually one to two finger breadths, should be kept in mind to ensure proper fit of the assembly during the procedure.

INSTRUCTIONS FOR USE

Needle Introduction

1. Perform a standard vein puncture of the right femoral vein and advance a wire guide to the superior vena cava.
2. When using a Mullins introducer sheath assembly, advance it over the wire guide into the superior vena cava. Remove the wire guide.
3. Aspirate and irrigate the introducer to remove remaining air.
4. Flush the transseptal needle and introduce the obturator, which will prevent damage to or puncture of the introducer during introduction of the needle.

5. Under fluoroscopic guidance, advance the needle through the introducer.
NOTE: The needle must move freely within the introducer.
6. With the needle tip and hub of the introducer correctly aligned, remove the obturator, leaving the needle tip within the introducer. Under continuous pressure monitoring, withdraw the unit caudally under anteroposterior fluoroscopic guidance.

Right Atrium Landmarks

1. As the distal tip enters the right atrium, the first of three landmarks of descent is apparent as the tip is deflected leftward when visualized in the anteroposterior plane. Constant right atrial pressure should be monitored.
2. Free mobility of the tip within the right atrium should be confirmed by gentle rotation. When central septal bulge and adjacent gutters have been created by left atrial enlargement, the relation of the introducer tip to the central septum can sometimes be appreciated by this maneuver.
3. Using the needle arrow and anteroposterior, then lateral, fluoroscopy, direct the needle tip posteromedially at 45-65 degrees to the horizontal plane.
4. Under anteroposterior fluoroscopy and continued pressure monitoring, again move the needle and the introducer assembly inferiorly. A second gradual leftward deflection of the tip is apparent as the introducer traverses the bulge of the ascending aorta, which may transmit a perceptible systolic pulsation. As the introducer tip is further withdrawn, a distinct third leftward deflection occurs as the tip clears the limbic ledge and engages the fossa ovalis. This is variably displaced depending upon the underlying valvular disease and its effect on septal size, orientation and relations.
5. The introducer tip position should then be checked in the lateral plane. Slight posterior angulation should be observed at the level of the fossa, but this again varies dependent upon the existing valvular disease.

Transseptal Puncture

1. Move the introducer/needle assembly superiorly against the septum with 5-10 degrees rotation of 4 o'clock position. This may permit the assembly to pass directly to the left atrium without the need for needle advancement. Probing will frequently detect a patent foramen ovale.
2. If the septum has not been crossed, gently advance the needle tip through the introducer tip and into the interatrial septum, at which time the right atrial pressure trace may disappear.
3. Left atrial pressure must be recorded before proceeding. Left atrial entry can be further confirmed by contrast injection and oximetry of left atrial blood.
4. Gently advance the needle and introducer as a unit.
5. Slightly withdraw the needle into the introducer and advance the sheath about 4 cm over the introducer, into the left atrium. Once the sheath is in the left atrium, remove the needle and the introducer and allow the sheath to backbleed before flushing.
6. Left atrial pressure is then recorded through the side-arm, and once satisfactory positioning within the left atrium is confirmed, the patient may receive systemic heparin and the procedure can continue.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.