Günther Tulip Vena Cava Filter Retrieval Set for Jugular Vein Approach

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GÜNTHER TULIP VENA CAVA FILTER RETRIEVAL SET

For Jugular Vein Approach

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Günther Tulip Vena Cava Filter Retrieval Set for jugular vein approach consists of a retrieval loop system with a braided platinum retrieval loop, an 11 French coaxial retrieval sheath system, an 18 gauge entry needle, a 0.035 inch wire guide, a 12 French pre-dilator, and a three-way stopcock. The coaxial retrieval sheath system has a radiopaque band on the outer sheath to identify the precise location of the distal tip of the sheath for positioning accuracy.

The Günther Tulip Vena Cava Filter Retrieval Set and Cook Celect Vena Cava Filters (including both the Cook Celect and Cook Celect Platinum Vena Cava Filters) are designed to act as permanent or retrievable filters; the filters may be retrieved if clinically indicated.

SET COMPONENTS


Fig. 1

INTENDED USE

The product has been designed for retrieval of implanted Günther Tulip Vena Cava Filter Retrieval Set and Cook Celect Vena Cava Filters in patients who no longer require a filter. Retrieval of the filter can be performed only by jugular approach.

CONTRAINDICATIONS

- Retrieval of the filter with significant amounts of trapped thrombus (greater than 25% of the volume of the cone).
- Retrieval of the filter for patients with an ongoing high risk for pulmonary embolism.

WARNINGS

- An inferior vena caval imaging evaluation for residual captured thrombus should be performed prior to an attempted retrieval.
- Excessive force should not be exerted to retrieve the filter.
- Never attempt to re-deploy a retrieved filter.
- Specific for the Gunther Tulip Vena Cava Filter Retrieval Set: do not advance the inner sheath of the coaxial retrieval sheath over the anchors of the filter; doing so may cause particles to scratch off the sheath.
- Specific for the Cook Celect Vena Cava Filters (i.e., Celect and Celect Platinum): the filter can be completely collapsed in the inner sheath before subsequently advancing the outer sheath.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in interventional endovascular techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Product modification or alteration is not recommended, as the product's safety and effectiveness have not been established following any modifications.
- Manipulation of products requires imaging control.
- Possible allergic reactions (e.g., to natural rubber latex in the side-arm adapter, or to cobalt chromium, nickel and platinum in the wire loop) should be considered.
- Physician practice guidelines and published guidance from regulatory agencies recommend that patients with indwelling filters undergo routine follow-up. The risks/benefits of filter retrieval should be considered for each patient during follow-up. Refer to the "**REFERENCES**" section of the Instructions for Use for citations that include recommendations related to filter follow-up and retrieval.
- Once protection from pulmonary embolism (PE) is no longer necessary, filter retrieval should be considered. Filter retrieval should be attempted when feasible and clinically indicated. Filter retrieval is a patient-specific, clinically complex decision; the decision to remove a filter should be based on each patient's individual risk/benefit profile (e.g., a patient's continued need for protection from PE compared to their experience with and (or) ongoing risk of experiencing filter-related complications). For all retrievable IVC filters, retrieval becomes more challenging with time, and this is commonly due to encapsulation of the filter legs or hook (in a tilted filter) by tissue ingrowth. Please refer to the "**CLINICAL STUDIES**" section of the respective instructions for Use for the Gunther Tulip Vena Cava Filter Retrieval Set or Cook Celect Vena Cava Filters, which include filter retrieval data from clinical studies and (or) published clinical literature.
- For filter retrieval, the right jugular vein is usually preferred due to its straighter route to the vena cava.
- Do not pull on the filter beyond what is required to keep tension on the loop. Doing so may cause damage to the caval wall.
- Do not retract the retrieval loop system before the tip of the coaxial retrieval sheath is at the filter anchors. Doing so may cause damage to the caval wall.
- Before the final retrieval step, ensure that the entire filter is covered by the outer sheath to prevent the filter anchors damaging the caval wall.
- The published clinical literature includes descriptions of alternative techniques for filter retrieval (e.g.,

techniques that vary from that described in this Instructions for Use); use of these techniques varies according to physician experience, patient anatomy, and filter position. **The safety or effectiveness of these alternative retrieval techniques has not been established.** The “REFERENCES” section of this booklet includes citations that describe alternative retrieval techniques; this information is provided as reference.

- Filter tilt, vena caval wall penetration/perforation, filter fracture, and filter or filter fragment migration and (or) embolization have been reported. These observations may be symptomatic or asymptomatic and could result from a retrieval attempt and (or) contribute to difficult or failed retrieval. These observations may contribute to a physician's decision to attempt filter retrieval and (or) the outcome of a filter retrieval attempt.

POTENTIAL ADVERSE EVENTS

Potential adverse events may be observed at the time of filter retrieval attempt and (or) in association with a retrieval attempt. Potential adverse events that may occur include, but are not limited to, the following:

- Access site-related complications
- Air embolism
- Arrhythmia
- Blood loss
- Blood vessel injury (including vena cava, jugular and innominate veins)
- Cardiac injury
- Cardiac tamponade
- Edema
- Device or filter fragment embolization
- Hematoma
- Hemorrhage
- Pain (including back and abdominal)
- Pulmonary embolism
- Injury to adjacent structures
- Vena cava stenosis
- Vessel occlusion/thrombosis

CLINICAL STUDIES

Please refer to the “CLINICAL STUDIES” section of the Instructions for Use for Gunther Tulip Vena Cava Filter Retrieval Set or Cook Celect Vena Cava Filters for further data on filter retrieval; clinical outcomes are described for the Gunther Tulip Vena Cava Filter Retrieval Set and Cook Celect Vena Cava Filters.

INSTRUCTIONS FOR USE

The product is intended for use by physicians trained and experienced in diagnostic and interventional endovascular techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed. Local anesthesia, sedation and analgesia will be used as required according to institution standard of care.

PREPARATION

1. Flush the pre-dilator, coaxial retrieval sheath system, and retrieval loop system.
2. On the retrieval loop system, ensure that the retrieval loop is fully expanded and tighten the pin vise. (Fig. 2)


Fig. 2

3. Hold the Y-fitting and pull back the plastic pin vise to retract the retrieval loop into the catheter. Tighten the screw of the Y-fitting to keep the retrieval loop inside the catheter. **(Fig. 3)**


Fig. 3

FILTER RETRIEVAL

4. Access the chosen jugular vein using the Seldinger technique.
5. Place a supportive 0.035 inch wire guide in the IVC.
6. If necessary, dilate the puncture site with the 12 French pre-dilator over the wire guide.
7. Remove the pre-dilator and advance a flush catheter (not included in set).
8. Remove the wire guide.
9. Perform diagnostic imaging to evaluate the filter position and check for thrombus.
10. After confirming the plan for filter retrieval on diagnostic imaging, insert the wire guide.
11. Exchange the flush catheter for the coaxial retrieval sheath system, advancing it over the wire guide to a point approximately 3-5 cm above the hook of the filter.
12. Remove the introducer dilator and the wire guide from the coaxial retrieval sheath system. **(Fig. 4)** Verify the position of the coaxial retrieval sheath by diagnostic imaging. The 3-Way-Stopcock should be attached to the inner sheath hub during injection of contrast media.

Fig. 4

13. Introduce the retrieval loop system through the coaxial retrieval sheath, then advance and connect the side-arm fitting of the loop system to the sheath. Tighten the side-arm fitting adapter around the catheter to prevent loss of blood. **(Fig. 5)**

Fig. 5

14. Loosen the screw of the Y-fitting. Hold the Y-fitting and advance the pin vise until the retrieval loop has fully expanded inside the vena cava and surrounds the filter. **(Fig. 6)**


Fig. 6

15. Pull back the retrieval loop until it engages the hook of the filter. **(Fig. 7)**


Fig. 7

CAUTION: Do not pull on the filter beyond what is required to keep tension on the retrieval loop. **Doing so may cause damage to the caval wall.**

16. Hold the retrieval loop steady with the pin vise, then push the Y-fitting with the catheter forward, still ensuring that the filter hook is caught in the retrieval loop. Lock the Y-fitting to secure the snare around the filter hook. (Fig. 8)

NOTE: If the retrieval loop loses its shape during the retrieval attempt, it can be removed and gently reshaped. After reshaping, clean the retrieval loop and proceed from step 13.


Fig. 8

17. While holding the retrieval loop and Y-fitting steady, advance the side-arm fitting to collapse the filter into the coaxial retrieval sheath and disengage the anchors of the filter from the caval wall. (Fig. 9)

CAUTION: Do not retract the retrieval loop system at this step. **Doing so may cause damage to the caval wall.**


Fig. 9

18. When the tip of the coaxial retrieval sheath is at the filter anchors, loosen the hub of the outer sheath and advance the outer sheath forward to cover the entire filter. (Fig. 10)

WARNING: Specific for the Gunther Tulip Vena Cava Filter Retrieval Set: do not advance the inner sheath over the anchors of the filter; doing so may cause particles to scratch off the sheath.

CAUTION: Ensure that the entire filter is covered by the outer sheath to prevent the filter anchors from damaging the caval wall.


Fig. 10

19. Retrieve the retrieval loop system and Coaxial Retrieval Sheath assembly.
20. Perform diagnostic imaging to assess the inferior vena cava for potential injury following filter retrieval.

NOTE: Hospital standard of care should be followed for removing the complete assembly and providing hemostasis to prevent bleeding at the vascular access site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. IVC filter guidelines, ISO 25539-3, and regulatory safety communications regarding IVC filters. Refer to your local Cook sales representative for information on available literature.

RECOMMENDATIONS RELATED TO FILTER FOLLOW-UP AND RETRIEVAL:

- Caplin et al.; Society of interventional Radiology Standards for Practice Committee. Quality improvement guidelines for the performance of inferior vena cava filter placement for the prevention of pulmonary embolism. J Vasc Interv Radiol, 2011;22:1499-1506.
- ISO 25539-3:2011. Cardiovascular implants – Endovascular devices – Part 3: Vena cava filters.
- Removing retrievable inferior vena cava filters: FDA Safety Communication; Issued May 6, 2014.

- Retrievable inferior vena cava (IVC) filters - serious complications associated with attempted IVC filter retrieval. MHRA Medical Device Alert; Issued May 2, 2013.

Filter retrieval is a patient specific, clinically complex decision; the decision to remove a filter should be based on each patient's individual risk/benefit profile (e.g., a patient's continued need for protection from PE compared to their experience with and (or) ongoing risk of experiencing filter-related complications).

For all retrievable IVC filters, retrieval becomes more challenging with time, and this is commonly due to encapsulation of the filter legs or hook (in a tilted filter) by tissue ingrowth.

The following references include descriptions of alternative techniques for filter retrieval. **The safety or effectiveness of these alternative retrieval techniques has not been established.** Use of these techniques varies according to physician experience, patient anatomy, and filter position.

- Al-Hakim et al. The hangman technique: a modified loop snare technique for the retrieval of inferior vena cava filters with embedded hooks. *J Vasc Interv Radiol.* 2015;26(1):107-10.
- Cho et al. Failed inferior vena cava filter retrieval by conventional method: Analysis of its causes and retrieval of it by modified double-loop technique. *Phlebology.* 2015;30(8):549-56.
- Foley et al. A "fall-back" technique for difficult inferior vena cava filter retrieval. *J Vasc Surg.* 2012;56(6):1629-33.
- Kuo et al. Excimer laser-assisted removal of embedded inferior vena cava filters: A single-center prospective study. *Circ Cardiovasc Interv.* 2013;6(5):560-6.
- Stavropoulos et al. Retrieval of tip-embedded inferior vena cava filters by using the endobronchial forceps technique: Experience at a single institution. *Radiology.* 2015;275(3):900-7.

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- Al-Hakim et coll. The hangman technique: a modified loop snare technique for the retrieval of inferior vena cava filters with embedded hooks. *J Vasc Interv Radiol*. 2015; 26(1) : 107-110.
- Cho et coll. Failed inferior vena cava filter retrieval by conventional method: Analysis of its causes and retrieval of it by modified double-loop technique. *Phlebology*. 2015; 30(8) : 549-556.
- Foley et coll. A "fall-back" technique for difficult inferior vena cava filter retrieval. *J Vasc Surg*. 2012; 56(6) : 1629-1633.
- Kuo et coll. Excimer laser-assisted removal of embedded inferior vena cava filters: A single-center prospective study. *Circ Cardiovasc Interv*. 2013; 6(5) : 560-566.
- Stavropoulos et coll. Retrieval of tip-embedded inferior vena cava filters by using the endobronchial forceps technique: Experience at a single institution. *Radiology*. 2015; 275(3) : 900-907.

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Bloomington, IN 47404 U.S.A.