EN
3

Universal Check-Flo® Adapter
Instructions for Use

CS
3

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T _ U C F A _ R E V 2

UNIVERSAL CHECK-FLO® ADAPTER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Check-Flo Adapter is an endoscopic adapter designed to facilitate passage of instruments ranging in size up to 9 French in diameter during a procedure. The self-sealing adapter can be connected to any standard female Luer lock fitting or bulb-type endoscopic fitting.

INTENDED USE

This device is intended to resist the backflow of fluid around an instrument inserted through the working channel of a ureteroscope, cystoscope, or other endoscope.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

This device is intended for use by physicians trained and experienced in urological procedures. Standard techniques should be employed.

INSTRUCTIONS FOR USE

1. Introduce the adapter onto the working channel of an endoscope via a standard female Luer lock fitting or bulb-type endoscopic fitting.
2. Introduce catheter or other desired instrument through the port of the adapter as necessary to complete the intended procedure.
3. After removal of the inserted device and completion of the procedure, carefully remove the adapter from the endoscope.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.