EN 3	Urethral Dilator Set Instructions for Use
BG 5	Комплект уретрални дилатори Инструкции за употреба
CS 7	Souprava uretrálních dilatátorů Návod k použití
DA 9	Dilatatorsæt til urethra Brugsanvisning
DE 11	Urethra-Dilatator-Set Gebrauchsanweisung
EL 13	Σετ διαστολέων ουρήθρας Οδηγίες χρήσης
ES 15	Equipo de dilatadores uretrales Instrucciones de uso
ET 17	Ureetra dilataatori komplekt Kasutusjuhend
FI 19	Virtsaputken laajenninsetti Käyttöohjeet
FR 21	Set de dilateurs urétraux Mode d'emploi
HR 23	Komplet dilatatora uretre Upute za uporabu
HU 25	Urethrális dilatátorkészlet Használati utasítás
IT 27	Set di dilatatori uretrali Istruzioni per l'uso
LT 30	Šlaplės plėtiklių rinkinys Naudojimo instrukcija
LV 32	Urīnizvadkanāla dilatatoru komplekts Lietošanas instrukcija
NL 34	Urethrale dilatatorset Gebruiksaanwijzing
NO 36	Urinrørdilatatorsett Bruksanvisning

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50 Kullanma Talimatı

URETHRAL DILATOR SET

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Urethral Dilator Set includes 9 dilators in 8, 10, 12, 14, 16, 18, 20, 22, and 24 Fr sizes and in a length of 37 cm. Dilators manufactured in 8 - 14 Fr are designed with a single sideport angled toward the tapered tip. The sideport is aligned with the device French size indicated on the end of the dilator without a tapered tip.

Performance Characteristics

- Tapered tip facilitates insertion into anatomy.
- Radiopaque material aids in positioning under imaging techniques.

Device Compatibility

These dilators are compatible with up to a 0.038 inch diameter wire guide.

Patient Population

The target patient population for these devices are adult patients who require dilation of the urethra.

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the urethral tract including the urethra and bladder.

Operating Principle

These devices are intended to be introduced into the urethra over a wire guide. Each dilator features a tapered tip to provide a gradual transition to the maximum outer diameter. The dilators are intended to be introduced sequentially from smallest to largest to gradually dilate the urethra. The sideport allows for urine to flow into the dilator to confirm placement in the bladder.

INTENDED USE

These devices are dilators intended for mechanical dilation of the urethra by the sequential introduction of dilators of gradually increasing French sizes.

INDICATIONS FOR USE

These devices are used for adult patients that require dilating and traversing of their urethra for cystoscopy.

CLINICAL BENEFITS

Urethral dilation to widen the diameter of the urethra.

CONTRAINDICATIONS

None

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

No known Precautions

POTENTIAL ADVERSE EVENTS

- Bleeding
- False passage
- Fistula
- Hematoma
- Hematuria
- Infection
- Pain
- Perforation, urethral
- Rectal injury
- Urethral injury
- Lower urinary tract symptoms (poor/diminished urinary stream, dysuria, incomplete emptying of the bladder, (recurrent) urinary tract infections, overflow incontinence, frequency, hesitancy, urgency, urinary retention, increasing post void residual, and/or pain)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Advance the wire guide directly from the enclosed holder by introducing the flexible tip of the wire guide into the urethral meatus, and gently manipulate it beyond the obstruction and into the bladder.

NOTE: The wire guide can also be placed through a cystoscope and maintained in position as the scope is removed.

2. Progressing from the smallest to the largest appropriate size, pass the radiopaque dilators over the wire guide while maintaining wire guide position.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.