EN
3

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URETERAL DILATOR SET AND AQ® HYDROPHILIC DILATOR

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The ureteral dilators are available as an individual dilator and in three sets.

The single dilator (part number AQ-076012-TP) has an outer diameter that tapers from 12 to 6 Fr.

The devices identified by part numbers 076000 and AQ-076000 each include 9 dilators in 6, 8, 9, 10, 11, 12, 14, 16, and 18 Fr sizes. The device identified by part number 076600 includes 7 dilators in 6, 8, 10, 12, 14, 16, and 18 Fr sizes.

All dilators have a length of 60 cm. The end of the dilators are tapered. The device French size is indicated on the end of the dilator without a taper. The devices identified by part numbers AQ-076012-TP and AQ-076000 have a hydrophilic coating that is applied to the surface of the dilator for ease of entry to the access site.

Performance Characteristics

- Tapered tip to facilitate insertion into anatomy.
- Radiopaque material to aid in positioning under imaging techniques.
- Hydrophilic Coating (Only applicable to part numbers: AQ-076012-TP and AQ-076000) to enhance lubricity of device.

Device Compatibility

These devices are compatible with up to a 0.038 inch diameter wire guide.

Patient Population

The target patient population for these devices are adult patients who require dilation of the ureter.

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the urinary tract including the ureter, bladder and urethra.

Operating Principle

These devices are intended to be introduced into the ureter over a previously placed wire guide. Each dilator features a tapered tip to provide a gradual transition to the maximum outer diameter. Additionally, some dilators feature a hydrophilic coating to reduce friction during insertion. Hydrophilic coating is activated by immersing the dilator in sterile water or saline prior to insertion. The dilators are intended to be introduced sequentially from smallest to largest to gradually dilate the ureter.

INTENDED USE

These devices are dilators intended for mechanical dilation of the ureter by the sequential use of dilators of gradually increasing French sizes, prior to ureteroscopy and/or stone manipulation.

INDICATIONS FOR USE

These devices are used for adult patients that require dilating and traversing of their ureter for ureteroscopy procedures.

CLINICAL BENEFITS

Ureteral dilation to permit ureteroscope insertion.

CONTRAINDICATIONS

None

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

No known Precautions

POTENTIAL ADVERSE EVENTS

- Ureteral perforations
- Hematomas
- Mucosal splitting/tear
- Failed access to the stone
- Ischemic damage and stricture formation (ureteral stricture)
- Failure to dilate and the need for a second procedure: ureteric stenting or loss of access with emergent percutaneous nephrostomy tube
- Transmural extrusion of a stone

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

NOTE: Prior to use, immerse AQ® dilators in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

1. Under x-ray or fluoroscopic control, cystoscopically pass the wire guide the desired distance into the ureter.
2. Progressing from the smallest to the largest appropriate size, pass the radiopaque dilators over the wire guide while maintaining the wire guide's position. **NOTE:** Dilators up to approximately 12 Fr may be passed through the working channel of the cystoscope. However, the telescope and inner scope elements must be removed—leaving the sheath—to accommodate the larger sizes.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.