- EN**
3 **Ultraxx™ Nephrostomy Balloon and Set**
Instructions for Use
- CS**
5 **Nefrostomický balónek Ultraxx™ se soupravou**
Návod k použití
- DA**
8 **Ultraxx™ nefrostomiballon og sæt**
Brugsanvisning
- DE**
10 **Ultraxx™-Nephrostomie-Ballon und -Set**
Gebrauchsanweisung
- EL**
13 **Μπαλόνι και σετ νεφροστομίας Ultraxx™**
Οδηγίες χρήσης
- ES**
16 **Balón y equipo de nefrostomía Ultraxx™**
Instrucciones de uso
- FR**
19 **Set et ballonnet de néphrostomie Ultraxx™**
Mode d'emploi
- HU**
21 **Ultraxx™ nephrostomiás ballon és készlet**
Használati utasítás
- IT**
24 **Palloncino e set per nefrostomia Ultraxx™**
Istruzioni per l'uso
- NL**
27 **Ultraxx™ nefrostomieballon en -set**
Gebruiksaanwijzing
- NO**
30 **Ultraxx™ nefrostomiballong og -sett**
Bruksanvisning
- PL**
32 **Balon i zestaw Ultraxx™ do nefrostomii**
Instrukcja użycia
- PT**
35 **Balão e conjunto de nefrostomia Ultraxx™**
Instruções de utilização
- SV**
38 **Ultraxx™ nefrostomiballong och -set**
Bruksanvisning

ULTRAXX™ NEPHROSTOMY BALLOON AND SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Ultraxx™ Nephrostomy Balloon and Set

Set contains (*set components may vary*)

- Radiopaque balloon catheter
- Radiopaque sheath
- Inflation device (*optional*)
 - Balloon catheters accept .038 inch diameter wire guides

Sphere Inflation Device

The Cook Sphere Inflation Device is a one-piece, plastic, disposable inflation device with a lock lever design that controls the piston, a manometer, and a connecting tube with a male rotating adapter. An optional 3-way stopcock may also be enclosed for use during preparation of the device. The manometer measures pressures ranging from vacuum to gauge capacity; the gauge is marked in 1 atm increments. The gauge also has an inner scale of comparable PSI measurements. The accuracy of the manometer has been determined to be within 1 atm over the range.

INTENDED USE

Ultraxx™ Nephrostomy Balloon and Set

Used to dilate the musculofascia, renal capsule, and parenchyma to establish and maintain a percutaneous tract.

Sphere Inflation Device

The inflation device is recommended for use while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon.

CONTRAINDICATIONS

None known

WARNINGS

- Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.
- Do not overinflate. Using excessive pressure to inflate the balloon on this device can cause the balloon to rupture.

PRECAUTIONS

Ultraxx™ Nephrostomy Balloon and Set

- Balloon dilation catheters are intended for use by physicians trained and experienced in techniques for balloon-catheter dilation.
- Do not exceed the maximum rated burst pressure (listed on label) for this balloon device.
- Do not preinflate the balloon.
- Refer to product label or the inflation check valve on the balloon device for appropriate balloon volume
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

Sphere Inflation Device

- Before use, ensure the connector tubing is completely free of air.
- Do not reuse or resterilize the Sphere Inflation Device. Read instructions prior to use.

POTENTIAL ADVERSE EVENTS

Complications that might occur during this procedure include overinflation of the balloon, which could result in trauma to the surrounding tissue.

INSTRUCTIONS FOR USE

Balloon Catheter Preparation

1. Using an aseptic technique, remove the Ultraxx™ balloon tray from its outer package and place in the sterile field.
2. Upon removal from package, inspect the balloon catheter to ensure no damage has occurred during shipment. Defective devices should be returned to Cook.
3. Lock a 10 ml syringe onto the stopcock and draw back the syringe to apply negative pressure. Close the stopcock and detach the syringe. Remove the protective sleeve from the balloon and discard.
4. You may lubricate the deflated balloon with a water-soluble lubricant.

Balloon Catheter Introduction

1. Prior to use, remove the inner protective tubing from the clear Amplatz sheath.
2. Backload the Amplatz sheath over the balloon and secure it to the sheath lock at the proximal end of the catheter.
3. Verify that the sheath is properly oriented, with the beveled tip distal on the balloon catheter.
4. Carefully advance the balloon catheter over a wire guide and into final position in the kidney. Utilize the radiopaque marker at the distal portion of the balloon to ensure proper positioning of the balloon catheter within the kidney.

NOTE: If your device did not include a Sphere Inflation Device, the following instructions may not apply. Follow the instructions for use provided by the manufacturer of the inflation device you are using.

Inflation Device Preparation and Balloon Inflation

NOTE: Make all aspiration and injection maneuvers with the lock lever pushed left, i.e., unlocked. In this position, you can freely pull the piston back for aspiration, or push it forward for injection. To lock the piston in position, slide the lever right to the straight up position.

1. Prepare a solution of contrast medium and normal saline in a small sterile bowl. Check contrast medium instructions for specific contrast mixture recommendations.
2. Orient the tubing downward into the contrast medium.
3. Push the release lever left and aspirate enough solution to fill the syringe. (Attach stopcock, if applicable.)
4. Hold the device upright to purge the air from the syringe and connecting tube. Tap the syringe lightly, if necessary, to remove all the air bubbles and fill the connecting tube completely.
5. Inspect the syringe and tubing (and stopcock, if applicable) to ensure that the device has been completely purged of air bubbles.
6. Adjust the syringe volume to the desired amount. If more contrast is needed, submerge the syringe tip into the basin of solution and aspirate. (Close stopcock, if applicable.)
7. Attach the inflation device to the balloon catheter. If a separate syringe was used to prepare the balloon catheter, remove it. When a stopcock is installed on the end of the inflation device connecting tube, it should be opened and purged with contrast media from the inflation device to eliminate air. Create a fluid-fluid connection between the balloon and the stopcock or connecting tube (male rotating adapter) of the inflation device by placing a drop of contrast solution from the syringe into each hub.
8. Hand-tighten the hubs securely.
9. Release the lock lever and allow the piston to move forward into neutral position (0 atm).
10. To inflate the balloon, engage the lock lever, turn the palm grip on the piston clockwise slowly until the desired inflation pressure is reached. Close the stopcock, (lock lever) on the inflation device to maintain pressure.

11. Advance the sheath over the inflated balloon and into the proper position. **NOTE:** The tolerances of the sheath and inflated balloon are very close. If resistance is encountered when advancing the sheath, it may be necessary to decrease balloon inflation to facilitate sheath passage.

Balloon Deflation and Withdrawal

1. Deflate the balloon by syringe aspiration. Aspirate the balloon completely before gently removing the device. Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal.
 - a. To gradually deflate the balloon, turn the palm grip on the piston counterclockwise slowly until the desired deflation pressure is reached.
 - b. To rapidly deflate the balloon, push the lock lever left, releasing the piston, and pull back. Slide the lock lever back to lock, if desired.
2. While removing the balloon catheter, carefully hold the sheath in place, ensuring its position has not changed.
3. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd