EN
3

Ultraxx™ Nephrostomy Balloon

Instructions for Use

CS
4

Nefrostomický balóněk Ultraxx™

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DA
6

Ultraxx™ nefrostomiballon

Brugsanvisning

DE
8

Ultraxx™ Nephrostomie-Ballon

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EL
10

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T _ U N B _ R E V O

ULTRAXX™ NEPHROSTOMY BALLOON

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

Set contains:

- A radiopaque balloon catheter (accepts .038 inch diameter wire guide)
- A radiopaque sheath

INTENDED USE

This product is intended to dilate the musculofascial tract, renal capsule, and parenchyma to establish and maintain a percutaneous tract.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during percutaneous nephrostomy.

WARNINGS

- Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.
- Do not overinflate the balloon. Using excessive pressure to inflate the balloon on this device can cause the balloon to rupture.
- If resistance is felt when removing the balloon from the working sheath, stop and consider removing them as a single unit to prevent damage to the product.

PRECAUTIONS

- Balloon dilation catheters are intended for use by physicians trained and experienced in techniques for balloon catheter dilation. Standard techniques should be employed.
- Do not exceed the maximum rated burst pressure (listed on label) for this balloon catheter.
- Do not preinflate the balloon.
- Refer to the product label or the inflation check valve on the balloon catheter for appropriate balloon volume.
- If loss of pressure in the balloon occurs or the balloon ruptures, immediately stop the procedure, deflate the balloon and remove carefully. Do not attempt to reinflate the balloon.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- The balloon catheter is not intended for multiple inflations.

POTENTIAL ADVERSE EVENTS

- Tissue trauma
- Bleeding
- Perforation
- Infection
- Urine leakage
- Fistula
- Pain

INSTRUCTIONS FOR USE

Balloon Catheter Preparation

1. Using aseptic technique, remove the balloon tray from its outer package and place it in the sterile field.
2. Inspect the balloon catheter to ensure no damage has occurred during shipment. Defective devices should be returned to Cook.
3. Lock a 10 mL syringe onto the stopcock and draw back the syringe to apply negative pressure. Close the stopcock and detach the syringe. Remove the protective sleeve from the balloon and discard the sleeve.
4. You may lubricate the deflated balloon with a water-soluble lubricant.

Balloon Catheter Introduction

1. Prior to use, remove the inner protective tubing from the clear Amplatz sheath.
2. Backload the Amplatz sheath over the balloon and secure it to the sheath lock at the proximal end of the catheter.
3. Verify that the sheath is properly oriented, with the beveled tip distal on the balloon catheter.
4. Carefully advance the balloon catheter over a wire guide and into final position in the kidney. Utilize the radiopaque marker at the distal portion of the balloon to ensure proper positioning of the balloon catheter within the kidney.

Balloon Inflation

1. Prepare a solution of contrast medium and normal saline in a small sterile bowl. Check the contrast medium instructions for specific contrast mixture recommendations.
2. Follow instructions provided by the manufacturer of the inflation device and fill the balloon to the desired pressure (≤ 20 atm).
3. Advance the sheath over the inflated balloon and into the proper position. **NOTE:** The tolerances of the sheath and inflated balloon are very close. If you encounter resistance when advancing the sheath, it may be necessary to decrease balloon inflation to facilitate sheath passage.

Balloon Deflation and Withdrawal

1. Deflate the balloon by syringe aspiration. Aspirate the balloon completely before gently removing the device. Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal.
2. While removing the balloon catheter, carefully hold the sheath in place, ensuring its position has not changed.
3. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.