EN 4	Turbo-Flo™ Peripherally Inserted Central Venous Catheters Instructions for Use
DA 9	Turbo-Flo™ perifert indførte centralvenekatetre Brugsanvisning
DE 14	Turbo-Flo™ peripher eingeführte zentrale Venenkatheter Gebrauchsanweisung
EL 20	Περιφερικά εισαγομενοι κεντρικοί φλεβικοί καθετηρες Turbo-Flo™ Οδηγίες χρήσης
ES 27	Catéteres venosos centrales de inserción periférica Turbo-Flo™ Instrucciones de uso
FR 32	Cathéters veineux centraux insérés en périphérie Turbo-Flo™ Mode d'emploi
IT 38	Cateteri venosi centrali ad inserimento periferico Turbo-Flo™ Istruzioni per l'uso
NL 44	Turbo-Flo™ perifeer ingebrachte centraalveneuze katheters Gebruiksaanwijzing
PT 50	Cateteres venosos centrais de introdução periférica Turbo-Flo™ Instruções de utilização
SV 56	Turbo-Flo™ perifert insatta centrala venkatetrar Bruksanvisning


Fig. 1


Fig. 2


Fig. 3

TURBO-FLO™ PERIPHERALLY INSERTED CENTRAL VENOUS CATHETERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Turbo-Flo Peripherally Inserted Central Venous Catheter (PICC) is a polyurethane catheter with a polyurethane suture wing, reinforced extensions and plastic clamps. Catheters can be trimmed to fit the patient's anatomy. PICC sets and trays also contain a Peel-Away® introducer, appropriately sized access needle, wire guides, positive pressure injection caps, and other accessories for percutaneous vascular placement.

Single-Lumen Catheters

French Size	Lumen Volume (ml)	Lumen Diameter (inches)	Length (cm)
4.0	0.9	.033	60
5.0	1.3	.048	60

Double-Lumen Catheters

French Size	Small Lumen Diameter (inches)	Small Lumen Volume (ml)	Large Lumen Diameter (inches)	Large Lumen Volume (ml)	Length (cm)
5.0	.015	0.6	.025	0.62	60

INTENDED USE

Turbo-Flo Peripherally Inserted Central Venous Catheter (PICC) Sets and Trays are intended for venous pressure monitoring, blood sampling and administration of drugs and fluids.

CONTRAINDICATIONS

None known

WARNINGS

- Do not power inject contrast medium through catheter. Catheter rupture may result. Use of a 10 ml or larger syringe will reduce the risk of catheter rupture.
- **Peripherally Inserted Central Venous Catheters play an important role in treatment of patients. However, catheter tips can erode or perforate vascular walls. Extreme caution must be used in placement and monitoring of catheters.**

- Catheter tip position should be verified by X-ray and monitored on a routine basis. Periodic lateral view X-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall. (**Reference 1**)

PRECAUTIONS

- The product is intended for use by healthcare practitioners trained and experienced in proper positioning of catheters in the central venous system using percutaneous entry (Seldinger) technique. Standard technique for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Patient movement can cause catheter tip displacement. Catheters placed via an antecubital vein have shown tip movement of up to 10 cm with motion of the extremity.
- Catheter size should be as small as the use will allow.
- When exchanging a Twin Turbo-Flo Double-Lumen Catheter, ensure that the wire guide is inserted into the hub marked "LARGE."

PRODUCT RECOMMENDATIONS

Catheter Size and Puncture Site


Preliminary reports indicate that catheter size can influence clotting. Larger diameter catheters have more tendency to promote clots. As reported by Amplatz, Gianturco and others, clot formation has less relation to type of catheter material than to size of catheter.

The angle of the catheter tip to the vessel wall should be checked carefully. Blackshear reviewed the medical literature of catheter perforations, which have confirming X-rays, and found that an incident angle of the catheter to vessel wall greater than 40 degrees was more likely to perforate. (**Reference 2**)

The following variables must also be considered in selecting appropriate catheter and length:

1. Patient history
2. Patient age and size
3. Access site available
4. Unusual anatomical variables
5. Proposed use and duration of treatment plan

Catheter Tip Positioning

Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.

Catheter Maintenance

Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization. After catheter placement and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, catheter tip position should be immediately reevaluated.** If catheter is not to be used immediately, its lumen should be maintained by continuous saline or heparinized saline drip or locked with heparinized saline solution. Catheter heparinization should be determined by institutional protocol and clinical judgement. Heparin concentrations of 10 Units/ml to 100 Units/ml have been reported adequate to maintain lumen patency. **(Reference 3)** Heparin lock should be reestablished after every use or at least every 24 hours if unused. Before using catheter lumen already locked with heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumen should be flushed with normal saline between administration of different infusates. After use, lumen should again be

flushed with twice the indicated lumen volume using normal saline before reestablishing heparin lock. Strict aseptic technique must be adhered to while using and maintaining catheter.

INSTRUCTIONS FOR USE

Catheter Obturator Preparation

1. Flush catheter with saline solution, heparinized saline solution or sterile water. **NOTE:** The catheter may be trimmed if a shorter length is required.
2. Attach a syringe with saline solution, heparinized saline solution or sterile water to the Luer lock fitting of the catheter obturator holder.
3. Inject enough solution to wet the obturator surface entirely. This will activate the AQ® hydrophilic coating, making the obturator surface very lubricious.
4. Remove the obturator from its holder and insert it into the catheter, locking it into position. **NOTE:** When placing a Twin Turbo-Flo Double-Lumen Catheter, insert obturator into the hub marked "LARGE." If the catheter has been trimmed, advance the obturator only to the distal end of the catheter. **NOTE:** If the surface of the obturator becomes dry after removal from the holder, wetting with additional saline or sterile water will renew the hydrophilic effect.
5. The catheter/obturator assembly can now be introduced as described in the following section.

Catheter Placement (Fluoroscopic Method)

6. After prepping the access site, introduce the access needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to visualize the tip of the needle during vessel access.
7. Using fluoroscopic guidance, introduce the wire guide through the needle and advance it 15-20 cm into the vessel.
8. Withdraw the needle, leaving wire guide in place. If necessary, enlarge the puncture site with scalpel blade.
9. Introduce the Peel-Away introducer assembly (sheath and dilator) over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 1**)
10. Using fluoroscopic control, determine the correct catheter length by advancing the wire guide to the desired catheter tip location. Once the wire guide tip is in proper position, mark the length by clamping forceps onto the wire guide at the skin site.
11. Withdraw the wire guide and measure it from the forceps to the distal tip to determine correct catheter length. Trim catheter if necessary. **NOTE:** Remove obturator prior to trimming the catheter; reinsert for catheter introduction.
12. Leaving the sheath in place, remove the dilator. (**Fig. 2**) **NOTE:** To prevent inadvertent air aspiration after removing wire guide and dilator, place thumb or finger over the cuffed proximal end of the sheath.
13. Introduce the catheter/obturator assembly into the sheath as far as possible. **NOTE:** Resistance may be felt approximately 7 cm distal to the catheter suture wing when introducing the assembly into the sheath due to an increase in outer diameter. (**Fig. 3**)

14. Peel the sheath away from the catheter by grasping the two wings of the sheath and pulling outward and upward. **NOTE:** Final catheter position is accomplished by alternately advancing the catheter into the sheath and then further pulling on the two wings. Once the sheath is removed, a slight advancement of the catheter may be needed for final positioning.
15. After catheter is in final position, remove obturator, secure catheter to skin and dress in standard fashion.
16. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC.

Catheter Placement (Non-Fluoroscopic Method)

17. After prepping the access site, introduce the access needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to locate the tip of the needle during vessel access.
18. Introduce the wire guide through the needle and advance it 15-20 cm into the vessel.
19. Withdraw the needle, leaving wire guide in place. If necessary, enlarge the puncture site with scalpel blade.
20. Introduce the Peel-Away introducer assembly (sheath and dilator) over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 1**)
21. Using a Tyvek® tape measure, clinical judgment, or other institutional protocol, determine the correct catheter length and trim catheter as needed. **NOTE:** Remove obturator prior to trimming catheter; reinsert for catheter introduction.
22. Leaving the sheath in place, remove the dilator. (**Fig. 2**) **NOTE:** To prevent inadvertent air aspiration after removing wire guide and dilator, place thumb or finger over the proximal end of the sheath.
23. Introduce the catheter/obturator assembly into the sheath as far as possible. **NOTE:** Resistance may be felt approximately 7 cm distal to the catheter suture wing when introducing the assembly into the sheath due to an increase in outer diameter. (**Fig. 3**)
24. Peel the sheath away from the catheter by grasping the two wings of the sheath and pulling outward and upward. **NOTE:** Final catheter position is accomplished by alternately advancing the catheter into the sheath and then further pulling on the two wings. Once the sheath is removed, a slight advancement of the catheter may be needed for final positioning.
25. After catheter is in final position, remove obturator, secure catheter to skin and dress in standard fashion.
26. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

1. Tocino IM, Watanabe A: "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition," *American Journal of Roentgenology*, 1986; 146:487-490.
2. Welch RH, Gravenstein N, Blackshear RH: "Multilumen Central Venous Catheters in Children: Relative Potential to Perforate Vessels. An In Vitro Study," *Journal of Clinical Monitoring*, 1997; 13:75-59.
3. Kelly C, Dumenko L, McGregor SE, McHutchion ME: "Change in Flushing Protocols of Central Venous Catheters," *Oncology Nursing Forum*, 1992; 19:599-605.

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Bloomington, IN 47404 U.S.A.