
- EN**
3 **UPJ Occlusion Balloon Catheter Set**
Instruction for Use
- CS**
5 **Souprava balónkového katetru pro okluzi UPJ**
Návod k použití
- DA**
7 **UPJ okklusjonsballonkatetersæt**
Brugsanvisning
- DE**
9 **Ballonkatheter-Set zur Okklusion des pyeloureteralen Übergangs**
Gebrauchsanweisung
- EL**
11 **Σετ καθετήρα με μπαλόني για την απόφραξη της πυελοουρητηρικής συμβολής**
Οδηγίες χρήσης
- ES**
13 **Equipo de catéter balón para oclusión de la unión ureteropélvica**
Instrucciones de uso
- FR**
15 **Set de sonde à ballonnet pour occlusion de la jonction pyélo-urétérale**
Mode d'emploi
- HU**
17 **UPJ occlusió ballonkatéterkészlet**
Használati utasítás
- IT**
19 **Set con catetere a palloncino per l'occlusione della giunzione ureteropelvica**
Istruzioni per l'uso
- NL**
21 **Ballonkatheterset voor occlusie van de pyelo-ureterale overgang**
Gebruiksaanwijzing
- NO**
23 **UPJ-okklusjonsballongkatetersett**
Bruksanvisning
- PL**
25 **Zestaw okluzyjnego cewnika balonowego do połączenia miedniczkowo-moczowodowego**
Instrukcja użycia
- PT**
27 **Conjunto de cateter com balão para oclusão da junção ureteropélvica**
Instruções de utilização
- SV**
29 **Ballongkateterset för UPJ-okklusion**
Bruksanvisning


Tuohy-Borst Adapter
 Adaptér Tuohy-Borst
 Tuohy-Borst adapter
 Tuohy-Borst-Adapter
 Προσαρμογέας Tuohy-Borst
 Adaptador Tuohy-Borst
 Adaptateur Tuohy-Borst
 Tuohy-Borst-adapter
 Adattatore Tuohy-Borst
 Tuohy-Borst-adapter
 Tuohy-Borst-adapter
 Łącznik Tuohy-Borst
 Adaptador Tuohy-Borst
 Tuohy-Borst-adapter

Occlusion Balloon Catheter
 Okluzní balónkový katetr
 Okklusionsballonkateter
 Okklusions-Ballonkatether
 Καθετήρας απόφραξης με μπαλόνι
 Catéter balón para oclusión
 Sonde à ballonnet pour occlusion de la jonction pyélo-urétérale
 Occlusiós ballonkatéter
 Catetere di occlusione a palloncino
 Occlusieballonkatheter
 Okklusjonsballongkateter
 Okluzyjny cewnik balonowy
 Cateter de balão para oclusão
 Ballongkateter för ocklusion

Wire Guide
 Vodící drát
 Kateterfeder
 Führungsdraht
 Συρμάτινος οδηγός
 Guia
 Guide
 Vezetődrót
 Guida
 Voerdraad
 Ledevaier
 Prowadnik
 Fio guida
 Ledare

Inflation Check Valve
 Pojistný plnicí ventil
 Inflationskontraventil
 Inflationsrückschlagventil
 Ανεπίστροφη βαλβίδα πλήρωσης
 Válvula de hinchado
 Valve anti-retour de gonflage
 Feltöltési visszacsapó szelep
 Valvola anti-reflusso per il gonfiaggio
 Vulkeerklep
 Antirefluksventil
 Zastawka zwrotna do napełniania
 Válvula anti-refluxo
 Fyllningsspärventil

Inflation/Injection Adapter
 Plnicí/injekční adaptér
 Inflations-/injektionsadapter
 Inflations-/injektionsadapter
 Προσαρμογέας πλήρωσης/έγχυσης
 Adaptador de hinchado/inyección
 Adaptateur de gonflage/injection
 Feltöltési/befecskendezési adapter
 Adattatore di gonfiaggio/iniezione
 Vul-/injectieadapter
 Fyllings-/injiseringsadapter
 Adapter do napełniania/wstrzykiwania
 Adaptador para enchimento/injecção
 Fyllnings-/injektionsadapter

Occlusion Balloon Catheter
 Okluzní balónkový katetr
 Okklusionsballonkateter
 Okklusions-Ballonkatether
 Καθετήρας απόφραξης με μπαλόνι
 Catéter balón para oclusión
 Sonde à ballonnet pour occlusion de la jonction pyélo-urétérale
 Occlusiós ballonkatéter
 Catetere di occlusione a palloncino
 Occlusieballonkatheter
 Okklusjonsballongkateter
 Okluzyjny cewnik balonowy
 Cateter de balão para oclusão
 Ballongkateter för ocklusion

UPJ OCCLUSION BALLOON CATHETER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

Set Contains:

- Balloon catheter
- Appropriate diameter PTFE coated stainless steel wire guide 80 cm long
- Removable T-shaped inflation/injection adapter
 - The removable T-shaped adapter allows the cystoscope to be removed after the catheter is in place.
- Syringe

INTENDED USE

Used to temporarily occlude the ureteropelvic junction to prevent stone fragments from entering the ureter during percutaneous lithotripsy, and for injection of contrast medium.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions, which create unacceptable risk during catheterization.

WARNINGS

- Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.
- Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. This device is intended for one-time use.

PRECAUTIONS

- Over inflation of the balloon on this device can cause the balloon to rupture.
- Do not overfill greater than 1.0 mL or risk device bursting.
- Only water soluble lubricants are recommended for use with the balloon catheter. Failure to use water soluble lubricants may cause the inflated balloon to burst.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Tissue trauma
- Perforation
- Infection
- Hematuria due to local ureteral bleeding

INSTRUCTIONS FOR USE

Catheter and Balloon Preparation

1. Remove the PTFE sheath which protects the balloon catheter; discard the sheath.
2. (See Figure A) Loosen the Tuohy-Borst adapter and remove the wire guide and hub.
3. (See Figure B) Replace the Tuohy-Borst adapter with the T-shaped inflation/injection adapter. Advance the T-shaped adapter until you feel resistance. Tighten the caps on both ends of the T-shaped adapter.
4. Fill the syringe with dilute contrast medium. Eliminate air from the syringe.

5. Press the syringe tip into the inflation check valve. **NOTE:** If the syringe is not inserted far enough, the balloon will not inflate or deflate when needed. The syringe will fit securely in the inflation check valve.
6. Inflate the balloon, ensuring not to inflate above 1.0 mL. Inflating and deflating the balloon for one cycle prior to patient contact will verify the functionality of the balloon.
7. After inflating and deflating the balloon, remove the syringe from the T-shaped adapter.
8. Loosen the caps on the T-shaped adapter, but do not remove them.
9. Remove the T-shaped adapter and replace it with the provided wire guide. Secure the wire guide to the balloon catheter by tightening the Tuohy-Borst adapter. The wire should protrude approximately 1-2 cm from the balloon catheter tip to act as a filiform during placement.

Access and Inflation

1. (See **Figure A**) Verify that the wire guide is secured to the balloon catheter by tightening the Tuohy-Borst adapter. The wire should protrude approximately 1-2 cm from the balloon catheter tip to act as a filiform during placement.
2. Under fluoroscopic control, cystoscopically pass the balloon catheter with wire guide into the ureter and up to the UPJ.
3. After positioning and before inflating the balloon, remove the wire guide and then remove the cystoscope.
4. (See **Figure B**) Replace the Tuohy-Borst adapter with the T-shaped adapter. Advance the T-shaped adapter until you feel resistance. Tighten the caps on both ends of the T-shaped adapter.
5. Fill the syringe with dilute contrast medium. Eliminate air from the syringe.
6. Press the syringe tip into the inflation check valve. **NOTE:** If the syringe is not inserted far enough, the balloon will not inflate or deflate when needed. The syringe will fit securely in the inflation check valve.
7. Inject 0.4 mL of contrast medium to start expanding the balloon and occlude the ureter, ensuring not to inflate above 1.0 mL. Confirm the amount of balloon inflation with fluoroscopy.
8. Remove the syringe from the inflation check valve to lock the injected contrast medium in the balloon for the duration of the procedure.

Deflation and Removal

1. Press the syringe tip into the inflation check valve until back pressure begins to fill the syringe with contrast medium.
2. Pull back on syringe plunger until all fluid is removed. **NOTE:** If unable to deflate the balloon with the syringe, loosen the caps on both ends of the T-shaped adapter and contrast medium will drain from the adapter.
3. Under fluoroscopic control confirm balloon is completely deflated.
4. Once deflated, slowly remove the balloon catheter.
5. If you feel resistance during removal, stop and use fluoroscopy to determine what is causing the resistance.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd