**COOK**  
MEDICAL

**CE**  
0123

**EN**  
3

## Ureteral Catheters

Instructions for Use

**CS**  
4

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**DA**  
6

## Ureterkatetre

Brugsanvisning

**DE**  
8

## Ureterkatheter

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**EL**  
10

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T \_ U R E C A T \_ R E V 1

## URETERAL CATHETERS

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

Ureteral catheters are available in a variety of French sizes, lengths, and distal tip configurations. The choice of catheter should be based on the physician's preference and clinical situation.

### INTENDED USE

Ureteral catheters are indicated for access and catheterization of the urinary tract, including the following applications:

- Delivery of contrast media
- Drainage of fluids from the urinary tract
- Delivery of irrigation fluids to the urinary tract
- Navigation of a tortuous ureter
- Access, advancement, or exchange of wire guides (open-ended catheters only)

### CONTRAINDICATIONS

These devices are contraindicated for use in patients unsuitable for catheterization.

### WARNINGS

None known

### PRECAUTIONS

- These devices are intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.
- Avoid bending or kinking the catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.
- If you encounter resistance while advancing or withdrawing the catheter, stop. Determine the cause of the resistance before proceeding.
- The potential effects of phthalates on children or pregnant or nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.
- Do not withdraw catheter while deflected in cystoscope/endoscope.
- Do not over-tighten the catheter adapter, as it may occlude the lumen.

### POTENTIAL ADVERSE EVENTS

Adverse events that may result from ureteral catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Peritonitis
- Urinary tract infection
- Abrasion

## INSTRUCTIONS FOR USE

### Placement with a Wire Guide (Open-Ended Catheters)

1. Place the cystoscope within the bladder to locate the ureteral orifice.
2. Flush the ureteral catheter of choice.
3. Advance the catheter through the cystoscope and into the bladder.
4. Advance the wire guide through the catheter. Slightly advance and retract the wire guide to assist in gaining access to the ureteral orifice and into the ureter.
5. Once the wire is advanced to the desired location within the ureter, advance the catheter over the wire and to the desired location.
6. Remove the wire guide.
7. Under fluoroscopy, inject contrast media through the catheter to visualize the ureter and to determine required intervention.

### Placement without a Wire Guide (Closed-Ended Catheters)

1. With a cystoscope in the bladder, locate the ureteral orifice to determine choice of ureteral catheter.
2. Flush the ureteral catheter.
3. Under direct vision, advance the catheter through the cystoscope and position the catheter tip in the desired location within the ureter.
4. Under fluoroscopy, inject contrast media through the catheter to visualize the ureter and to determine required intervention.

**NOTE:** Ensure sufficient length of catheter is inserted into the adapter to prevent inadvertent removal of adapter from catheter tubing. The recommended location of the proximal end of the catheter tubing is approximately the middle of the adapter.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd