Universa® Stent Set

Instruction for Use

Souprava stentu Universa®

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Universa® stentsæt

Brugsanvisning

Universa® Stent-Set

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T _ U S S _ R E V 3

Universa® STENT SOFT OR FIRM

Used for temporary internal drainage from the ureteropelvic junction to the bladder. Ureteral stents have been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques. The multi-length stents will accommodate ureters from approximately 22-32cm in length.

CAUTION: Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.

NOTE: Prior to use, immerse stent in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

SET CONTAINS *(set components may vary)*

- Radiopaque double pigtail stent
- Stainless steel wire guide 145 cm long with 3 cm flexible tip
- Radiopaque stent positioner with radiopaque tip
- Pigtail straightener
- Lockable pusher (-LP stents only)

CONTRAINDICATIONS

There are no known contraindications.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

NOTE: If no wire guide is provided with this set, the following is recommended:

5.0 French stents accept .038"

6.0 French stents accept .038"

7.0 French stents accept .038"

8.0 French stents accepts .038"

WARNINGS

Formation of knots in multi-length stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.

SUGGESTED INSTRUCTIONS FOR USING Universa® STENT SETS

ENDOSCOPIC PLACEMENT

1. Pass a flexible wire guide tip to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using a wire guide and an open-end ureteral catheter in combination.
2. Using a baseline pyelogram, estimate the proper stent length; add 1 cm to that estimated ureteral measurement. Accurate measurement enhances drainage efficiency and patient comfort.
3. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wireguide into the renal parenchyma.
4. Watch for the distal end of the stent at the ureterovesical junction. At that point, halt advancement of the stent. As an assistant removes the wire guide, hold the stent in position with the positioner. The stent pigtail will form spontaneously. Carefully remove the positioner from the cystoscope.

NOTE: If necessary, final adjustment can be made with endoscopic forceps. The stent may be removed easily by gentle withdrawal traction using endoscopic forceps. Fluoroscopy facilitates stent placement; however, standard radiography may be used.

CAUTIONS

Intended for one-time use only.

The Universa® Firm stents must not remain indwelling more than twelve (12) months. The Universa® Soft stents must not remain indwelling more than six (6) months. If the patient's status permits, the stent may be replaced with a new stent.

These stents are not intended as permanent indwelling devices.

Tether should be removed if stent is to remain indwelling longer than 14 days.

Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.

A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.

Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period.

Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended.

Scopes larger than 21.0 French are suggested.

Individual variations of interaction between stents and the urinary system are unpredictable.

Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.

If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Quality Assurance department at the address/phone listed within this booklet.

REFERENCES

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," The Journal of Urology, 119 (1978), 731.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," The Journal of Urology, 121 (1979), 17.

M. F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single End Ureteral Stent," Urology, 13 (1979), 516-520.

M. A. Koyle, A. D. Smith: "Ureteral Stents," Chapter 31 in A. D. Smith, W. R. Castañeda-Zuñiga, J. G. Bronson: Endourology: Principles and Practice. Thieme Inc., New York, New York, 1986.

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REFERANSER

T.W. Hepperlen, H.K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," The Journal of Urology, 119 (1978), 731.

T.W. Hepperlen, H.K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," The Journal of Urology, 121 (1979), 17.

M.F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V.A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification to Single End Ureteral Stent," Urology, 13 (1979), 516-520.

M.A. Koyle, A.D. Smith: "Ureteral Stents," Chapter 31 in A.D. Smith, W.R. Castaneda-Zuniga, J.G. Bronson: Endourology: Principles and Practice. Thieme Inc., New York, New York, 1986.

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© COOK 2016
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.