COOK

MEDICAL

- EN**
3
Universa® Soft Ureteral Stents and Stent Sets
Instructions for Use
- CS**
5
Měkké ureterální stenty a sady stentů
Universa®
Návod k použití
- DA**
7
Universa® bløde ureterstents og stentsæt
Brugsanvisning
- DE**
9
Weiche Universa® Ureterstents und Stent-Sets
Gebrauchsanweisung
- EL**
11
Μαλακές ενδοπροσθέσεις ουρητήρα και
σετ μαλακών ενδοπροσθέσεων ουρητήρα
Universa®
Οδηγίες χρήσης
- ES**
14
Stents y equipos de stents ureterales blandos
Universa®
Instrucciones de uso
- FR**
16
Endoprothèses urétérales et sets
d'endoprothèses urétérales souples Universa®
Mode d'emploi
- HU**
18
Universa® lágy ureterális sztentek és
sztentkészletek
Használati utasítás
- IT**
20
Stent e set con stent ureterale morbido
Universa®
Istruzioni per l'uso
- NL**
23
Universa® zachte ureterale stents en stentsets
Gebruiksaanwijzing
- NO**
25
Universa® myke ureterstenter og stentsett
Bruksanvisning
- PL**
27
Miękkie stenty moczowodowe Universa®,
pojedyncze i w zestawach
Instrukcja użycia
- PT**
29
Stents ureterais maleáveis e Conjuntos de
stents Universa®
Instruções de utilização
- SV**
32
Universa® mjuka uretärstentar och stentset
Bruksanvisning

T _ U S U S N O C E _ R E V O

UNIVERSA® SOFT URETERAL STENTS AND STENT SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Universa Soft Ureteral Stents and Stent Sets are flexible tubular devices made from radiopaque polyurethane. The stents are available with or without hydrophilic coating. Each end of the Universa Soft Ureteral Stent consists of a curled, tapered pigtail. Sideports are placed on each end of the stent, including on the pigtails of the stent. A braided or monofilament tether for repositioning and removal of the device is located on the proximal pigtail of the stent. Radiopaque graduation marks are located on the stent to increase fluoroscopic visualization during stent advancement and placement.

The Universa Soft Ureteral Stent Set includes (set components may vary):

- Radiopaque double pigtail stent
- Stainless steel wire guide
- Radiopaque stent positioner with radiopaque tip
- Pigtail straightener

INTENDED USE

The Universa Soft Ureteral Stents and Stent Sets are used for temporary internal drainage from the ureteropelvic junction to the bladder. Ureteral stents have been used to relieve obstruction in a variety of benign, malignant and post-traumatic conditions. The stents may be placed using endoscopic, percutaneous or open surgical techniques.

CONTRAINDICATIONS

None known

WARNINGS

- Formation of knots in multi-length stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.
- After use this device may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and with applicable local, state and federal laws and regulations.

PRECAUTIONS

- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for the reproductive and developmental effects.
- The Universa Soft Stents must not remain indwelling more than six months. If the patient's status permits, the stent may be replaced with a new stent.
- These stents are not intended as permanent indwelling devices.
- **The tether should be removed if the stent is to remain indwelling longer than 14 days.**
- Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period.
- Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Ureteral stents should be checked periodically for signs of encrustation and proper function. Periodic checks of the stent by cystoscopic and / or radiographic procedures are recommended at intervals deemed to be appropriate by the physician in consideration of the individual patient's condition and other patient specific factors. The stent is not intended as a permanent indwelling device which should not exceed 180 days.
- Use of the device should be based upon consideration of risk-benefit factors as they apply to your patient.

POTENTIAL ADVERSE EVENTS

Complications of ureteral stent placement are documented. These complications include, but are not limited to:

- Extravasations
- Occlusion
- Migration
- Hemorrhage
- Sepsis
- Perforation of the urinary tract

- Peritonitis
- Encrustation
- Urinary tract infection
- Loss of renal function.

PRODUCT RECOMMENDATIONS

If no wire guide is provided with this set, the following is recommended:

- 5.0 French stents accept .035"
- 6.0 French stents accept .038"
- 7.0 French stents accept .038"
- 8.0 French stents accept .038"

MRI INFORMATION

Nonclinical testing has demonstrated that the Universa Ureteral Soft Stents are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 tesla or 3.0 tesla only
- Maximum spatial gradient magnetic field of 1600 gauss/cm (16.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 2.0 W/kg (Normal Operating Mode) for 15 minutes of continuous scanning

Under the scan conditions defined above, Universa Soft Ureteral Stents are expected to produce a maximum temperature rise of 1.6°C after 15 minutes of continuous scanning.

The image artifact extends approximately 6 mm for the Universa Soft Ureteral Stents as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 tesla MR system. The image artifact obscures the device lumen.

INSTRUCTIONS FOR USE

Endoscopic Placement

NOTE: Prior to use, immerse stent in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

1. Pass a flexible wire guide tip to the renal pelvis. Tortuosity in the obstructed ureter can be resolved using a wire guide and an open-end ureteral catheter in combination.
2. Using a baseline pyelogram, estimate the proper stent length; add 1 cm to that estimated ureteral measurement. Accurate measurement enhances drainage efficiency and patient comfort.
3. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
4. Watch for the distal end of the stent to appear at the ureterovesical junction. At that point, halt advancement of the stent.
5. As an assistant removes the wire guide, hold the stent in position with the positioner. The stent pigtail will form spontaneously.
6. Carefully remove the positioner from the cystoscope.

NOTE: If necessary, final adjustment can be made with endoscopic forceps.

NOTE: The stent may be removed easily by gentle withdrawal traction using endoscopic forceps.

NOTE: Fluoroscopy facilitates stent placement; however, standard radiography may be used.

If problems occur using this device, please call your Cook sales representative or contact our Customer Relations Department.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

T.W. Hepperlen, H.K.Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T.W. Hepperlen, H.K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M.F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V.A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification to Single End Ureteral Stent," *Urology*, 13 (1979), 516-520.

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Bloomington, IN 47404 U.S.A.