
VueLite™
LED LIGHT SOURCE

www.cookmedical.com

Language Index

English	1-1
VueLite™ LED LIGHT SOURCE USER MANUAL	
Česky	2-1
VueLite™ SVETELNÝ ZDROJ LED UŽIVATELSKÁ PŘÍRUČKA	
Dansk	3-1
VueLite™ LED-LYSKILDE BRUGERVEJLEDNING	
Deutsch	4-1
VueLite™ LED-LICHTQUELLE BENUTZERHANDBUCH	
Ελληνικά	5-1
ΠΗΓΗ ΦΩΤΟΣ LED VueLite™ ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ	
Español	6-1
FUENTE DE LUZ LED VueLite™ MANUAL DEL USUARIO	
Suomi	7-1
VueLite™ LED -VALONLÄHDE KÄYTTÖOPAS	
Français	8-1
SOURCE LUMINEUSE À LED VueLite™ MANUEL D'UTILISATION	
Magyar	9-1
VueLite™ LED FÉNYFORRÁS FELHASZNÁLÓI KÉZIKÖNYV	
Italiano	10-1
SORGENTE LUMINOSA LED VueLite™ MANUALE D'USO	
Nederlands	11-1
VueLite™ LED LICHTBRON GEBRUIKERSHANDLEIDING	
Norsk	12-1
VueLite™ LED-LYSKILDE BRUKERHÅNDBOK	
Polski	13-1
VueLite™ ŹRÓDŁO ŚWIATŁA LED PODRĘCZNIK UŻYTKOWNIKA	
Português	14-1
VueLite™ FONTE DE ILUMINAÇÃO LED MANUAL DO UTILIZADOR	
Svenska	15-1
VueLite™ LED-LJUSKÄLLA ANVÄNDARHANDBOK	

Introduction

This manual describes the proper procedures for using the VueLite LED Light Source. The manual also contains pertinent information on the proper care and handling of the device during use, reprocessing and storage. Please read the entire manual carefully before using the VueLite. If you have any questions about the material in this manual or the operation or safety of the equipment, please contact Cook Medical’s customer service.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Table of Contents

VueLite™ Description, Specifications and Accessories	1-2	Agency Compliance Statements and Guidance Tables	1-11
Body Diagram and Descriptions	1-2	Guidance and Manufacturer’s Declaration – Emissions	1-11
Specifications	1-3	Guidance and Manufacturer’s Declaration – Immunity	1-11 - 1-12
Intended Use	1-3	Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the VueLite™	1-12
Warnings and Precautions	1-4		
VueLite™ Instructions for Use	1-5	Repair Service	1-13
Preparation	1-5	Customer Service	1-13
Operation	1-5	Returning Goods to Cook Medical	1-13
Battery and Charger	1-6 - 1-7		
Care and Maintenance	1-8	Warranty	1-14
Reprocessing	1-8		
Reprocessing Protocol	1-8 - 1-9	Symbol Descriptions	1-15
Precleaning, Enzymatic Cleaning and Rinsing, High-Level Disinfection			
Drying	1-9	Regulatory	1-16
Troubleshooting	1-10	Vendor Information	1-16

This device is a **type BF** applied part. It should be connected only to medically approved type BF endoscopes that conform to IEC 60601-1 and IEC 60601-2-18.

IPX7 This equipment is rated IPX7.

Body Diagram and Descriptions

Fig. 1

Adapter Nut: The adapter nut allows the VueLite to attach to the user's endoscope. Rotate the nut clockwise to attach the device to an endoscope and counterclockwise to remove the device. The standard configuration of the VueLite is with an adapter nut that threads onto the standard universal light guide post's STORZ adapter thread, found on most endoscopes.

Light Intensity Switch: The light intensity switch controls the on/off and intensity functions of the VueLite. Moving the lever toward the wider end of the intensity indicator turns the LED on and increases the light intensity; moving the lever toward the narrower end of the intensity indicator decreases the light intensity and turns the LED off.

Battery Cap: The battery cap seals the battery and electronics within the control body.

Specifications

Parameter	Specification
Body Diameter	27 mm (1.06 inch)
Overall Length	89 mm (3.51 inch)
Weight (with Battery)	90 g (3.2 oz)
Battery Life at Full LED Power	45 minutes
Battery Charge Time	2 hours
Battery Type	Lithium-ion NP-700 3.6 V, 700 mAH minimum
Power Source	Internally powered
Mode of Operation	Continuous
Safe Operating Ambient Temperature Range	15-33°C (59-91°F)
Safe Storage and Transport Temperature Range	-25-50°C (-13-122°F)
Safe Operating, Storage and Transport Relative Humidity Range	0%-95% RH

Intended Use

The VueLite LED Light Source is intended to provide illumination for examination, diagnostic and therapeutic applications.

Please be sure to handle this optical instrument with care at all times. The electrical and optical components can be damaged by physical trauma, extreme temperatures or fluid invasion into the device.

Warnings

- This device is not supplied sterile. Clean and disinfect the VueLite prior to use, employing the procedure described in this manual.
- This equipment is not suitable for use in the presence of flammable mixtures.
- The following methods of sterilization and disinfection will result in damage to the VueLite and should not be used: ethylene oxide, autoclave, ultrasonic, chlorides, formaldehyde and hydrogen peroxide.
- To avoid potential safety hazards, the user of this device should consult the manuals of all medical electrical equipment used in conjunction with this device. If this device is used with high-frequency equipment, the attached endoscope should be medically approved for use with high-frequency equipment, and the user should consult the manuals of all devices for safe operating voltage levels.
- Avoid direct viewing of the high-intensity light at the front of the VueLite. Always connect and disconnect the VueLite with the light turned off.

Precautions

- This instrument is intended for use by persons thoroughly trained in the techniques that employ this light source.
- The VueLite may be damaged if used improperly. Read this User's Manual thoroughly before attempting to use the device.
- Check all items upon receipt to ensure damage has not occurred during shipment.
- Avoid storing or using the VueLite in areas of heavy traffic where the device may sustain physical damage.
- Avoid immersing the VueLite for periods of time greater than recommended by the disinfectant manufacturer. Prolonged immersion may damage the exterior of the device and could result in fluid invasion into the light source.
- Do not attempt to disassemble the VueLite in any way. There are no user-serviceable parts, and disassembly will void all warranties.
- Avoid looking directly into the LED when it is on.
- Prior to use, check the outer surface of the VueLite, endoscope and any endoscopic accessories for rough surfaces, sharp edges or protrusions that may cause a safety hazard.

Preparation

Before use, the VueLite should be reprocessed. Please refer to the Care and Maintenance section of this manual for proper reprocessing protocol.

Operation

Fig. 2. Attaching the VueLite to an endoscope

Fig. 3. Intensity control switch

- A.** Hold the VueLite so that the main body fits comfortably in one hand, allowing easy manipulation of the adapter nut with the thumb and forefinger.
- B.** Tighten the adapter nut onto the endoscope light post by turning the nut clockwise (Fig. 2). Do not thread the VueLite by turning the device body—rotate only the adapter nut.
- C.** Turn on the LED by rotating the light intensity switch toward the wider end of the intensity indicator (Fig. 3). Adjust the brightness to the desired level by rotating the intensity switch further.
- D.** At the completion of use, turn the LED off by rotating the light intensity switch toward the narrower end of the intensity indicator.
- E.** Remove the VueLite from the endoscope by rotating the adapter nut counterclockwise.
- F.** Avoid leaving the LED on for extended periods of time when not in use because this will drain the battery unnecessarily.

Battery and Charger

Fig. 4. Removing the battery cap and inserting/removing the Li-ion battery

The device's internal LED operates on a single lithium-ion battery and provides approximately 45 minutes to two hours of continuous operation (depending on intensity) when a new battery is fully charged. To insert a new battery, align the battery with the cutout in the battery compartment, making sure the end of the battery with the two gold electrical contacts is inserted into the battery compartment first. Once a charged battery is inserted, replace the battery cap by screwing the cap clockwise until a tight seal is made. To remove the battery, simply unscrew the battery cap and turn the device upright until the battery slides out. Please take care not to let the battery fall to the floor or onto a hard object, which could damage the battery.

Caution

- Do not incinerate the battery.
- Do not expose the battery to high temperatures above 122°F/50°C.
- Do not disassemble the battery.
- Do not short-circuit the battery.

Warning: The battery should be removed if the equipment is not likely to be used for some time.

To Charge Battery


Fig. 5. Battery charger and end views of battery

- A. Plug the AC adapter into a wall outlet.
- B. Connect the AC power adapter to the DC IN jack located at the back of the charger.
The **status** indicator will be solid **red** when the charger is ready, and the **charge** indicator will be flashing **green** when the charger is waiting for a battery to charge.
- C. Slide the battery into the charger, making sure the end of the battery with the two gold electrical contacts is inserted into the charger first. The **charge** indicator will turn **red** when the battery is charging.
- D. The **charge** indicator will alternately flash **red** and **green** when the battery has been charged to over 90%.
- E. When the battery is fully charged, the **charge** indicator will turn solid green.
- F. After the battery has been fully charged, slide the battery out of the charger and unplug the AC adapter for the charger.

Caution

- Do not expose the charger to fire, water or moisture.
- To avoid electrical shock, do not disassemble the charger.
- No user-serviceable parts are inside the charger.

Warning: The battery charger is not a medical device and should not be used in a room where a patient may be present.

When the end user wishes to discard this product, it must be sent to a separate collection facility for recovery and recycling.

Reprocessing

Warning: Failure to comply with the following may result in damage to the VueLite and will void the product warranty.

The VueLite LED is rated IPX7 for protection against ingress of water. The entire device is submersible in water and disinfectant solutions. However, the VueLite should not remain immersed in water or cleaning solutions longer than the manufacturer's specifications for the particular cleaning solution. After the disinfection cycle has concluded, remove the device from the disinfection solution, and dry and store the device in a clean, dry environment.

Thorough cleaning and rinsing are the first and most important steps in the reprocessing of the VueLite. Without thorough cleaning and rinsing, it might not be possible to achieve high-level disinfection. Cleaning is the removal of all adherent visible soil (i.e., blood, protein substances and other debris) from the surface, crevices and joints of the device.

Use caution when cleaning and disinfecting the VueLite; some methods may be harmful to the device and could result in extensive damage. Manual cleaning is the recommended method for cleaning. To be effective, cleaning agents must assist in the removal of residual organic debris without damaging the device. No single cleaning agent removes all types of debris or is safe to use with every type of reusable device. Certain cleaning agents may damage the device's materials. Cleaning agents should be used in the correct dilution/concentration and at the correct temperature in accordance with the cleaning agent manufacturer's directions. Contact your local Cook Medical representative to verify the compatibility of a cleaning method not listed in this manual.

Warning: Ensure that the battery cap is fully tightened on the O-ring seal to prevent liquid intrusion into the battery compartment.

Reprocessing Protocol

1. Precleaning

Cook Medical recommends the VueLite be precleaned immediately after the device has become contaminated. This will prevent patient material from adhering to the device. Dried material is difficult to remove and can render the disinfection or sterilization process ineffective. Gently wipe all debris from the exterior surfaces using a soft, lint-free cloth and water to which you have added a low-sudsing detergent, diluted in accordance with the detergent manufacturer's instructions.

2. Enzymatic Cleaning and Rinsing

- A. Fill a basin with a freshly made solution of water and a low-sudsing detergent, diluted in accordance with the detergent manufacturer's instructions.
- B. Immerse the VueLite. Wash all debris from the exterior of the device by gently brushing and wiping the instrument while it is submerged in the detergent solution.
- C. Rinse the VueLite in clean water.
- D. Wipe dry the exterior of the device with a soft, clean cloth.

3. High-Level Disinfection with .55% Ortho-Phthalaldehyde Solution:

- E. Prepare the disinfectant solution in accordance with the solution manufacturer's label.
- F. Completely immerse the VueLite in the disinfectant solution for the recommended time and temperature as indicated on the label of the disinfectant.
- G. Thoroughly rinse the exterior of the device with large amounts of clean water.
- H. A final rinse using 70% isopropyl alcohol solution can be used to speed the drying process and reduce the numbers of any organisms present as a result of rinsing with potable water.
- I. Dry the exterior of the device with a soft, clean towel.

Drying

If the VueLite becomes wet, either from reprocessing or cleaning fluids, thoroughly dry the exterior immediately. Rinsing with isopropyl alcohol followed by drying with a compressed gas, such as compressed air or commercial compressed gas canisters, is highly recommended to reduce contamination of the glass light guide and reduce the possibility of oxidation of the materials used in the construction of the device.

Warning: Failure to comply with the following may result in damage to the VueLite and will void the product warranty.

In addition, if the VueLite becomes wet, thoroughly dry the battery cap seal around the O-ring, with compressed gas or by thoroughly wiping any excess liquid from this area, prior to removing the battery cap. This will prevent any trapped fluid from inadvertently entering the battery compartment and exposing the electronics to fluids.

Warning: Ensure there is no trapped liquid between the battery cap and VueLite body prior to removing the battery cap.

Troubleshooting

Problem	Possible Cause	Action
Light goes out or dims.	Battery power is low.	Replace battery with a freshly charged battery.
	Battery cannot supply sufficient current at this power setting.	Change battery for a freshly charged battery. Reduce the light intensity and continue using the current battery. Recharge the battery as soon as possible.
Adapter nut does not attach properly to the endoscope.	The endoscope has an adapter on its light post.	Refer to the guide below and remove the adapter.
	The endoscope and device do not have a compatible coupling system.	Use an adapter sleeve or switch to device with a different type of adapter nut.
Light intensity is low.	Foreign material or film is on LED output taper.	Clean the taper with a cotton swab wetted with isopropyl alcohol until the debris or film is cleared off. Dry with compressed air or gas.
Device becomes hot after prolonged use at full power.	Device has been left on for more than 30 minutes at full power.	Turn the LED to a lower intensity setting if the examination must continue past 30 minutes. Or turn the device off for a few minutes to permit adequate cooldown of the device.

The standard light guide post found on most endoscopes consists of a snap-ring groove for an ACMI fitting connection, a Wolf adapter thread and a STORZ adapter thread. The VueLite attaches to the STORZ adapter thread, providing a secure connection between the VueLite and light guide within the scope. If an adapter sleeve (see Fig. 7 below) is on the light guide post, it must be removed prior to connecting the VueLite to the STORZ adapter thread.


Fig. 6. Standard universal light guide post


Fig. 7. Standard universal light guide post (left), post with Wolf adapter (center), and post with STORZ adapter (right)

Agency Compliance Statements and Guidance Tables

Medical electrical equipment needs special precautions regarding EMC and must be installed and put into service according to the EMC information provided in this section.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Warning: The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

The VueLite is intended for use in the electromagnetic environments specified below. The customer or user of the device should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Emissions

Emission Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The VueLite uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The VueLite is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	N/A	
Flicker IEC 61000-3-3	N/A	

Guidance and Manufacturer's Declaration - Immunity

Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD EN/IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the RH should be at least 30%.
EFT EN/IEC 61000-4-4	± 2 kV mains ± 1 kV I/Os	N/A – VueLite is powered by internal battery and has no I/Os.	N/A
Surge EN/IEC 61000-4-5	± 1 kV differential ± 2 kV common	N/A – VueLite is powered by internal battery and has no I/Os.	N/A
Voltage Dips/Dropout EN/IEC 61000-4-11	> 95% dip for 0.5 cycle 60% dip for 5 cycles 30% dip for 25 cycles > 95% dip for 5 seconds	N/A – VueLite is powered by internal battery and has no I/Os.	N/A
Power Frequency 50/60 Hz Magnetic Field EN/IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Guidance and Manufacturer's Declaration - Immunity (cont.)

Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A – VueLite powered by internal battery; no I/Os.	Portable and mobile communications equipment should be separated from the VueLite LED by no less than the distances calculated or listed below:
Radiated RF EN/IEC 61000-4-3	3 Vrms 80 MHz to 2.5 GHz	(E1) 3 V/m	$D = (3.5/E1)(\sqrt{P})$ 80-800 MHz $D = (7/E1)(\sqrt{P})$ 800 MHz to 2.5 GHz Where P is the max power in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (E1). Interference may occur in the vicinity of equipment containing a transmitter.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the VueLite: Equipment and Systems that Are Not Life-Supporting

The VueLite is intended for use in an electromagnetic environment in which radiated disturbances are controlled. The customer or user of the VueLite can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the VueLite as recommended below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150 kHz to 80 MHz $D = (3.5/\sqrt{P})(\sqrt{P})$	Separation (m) 80-800 MHz $D = (3.5/\sqrt{P})(\sqrt{P})$	Separation (m) 800 MHz to 2.5 GHz $D = (3.5/\sqrt{P})(\sqrt{P})$
.01	N/A*	.1166	.2333
.1	N/A*	.3689	.7378
1	N/A*	1.1666	2.3333
10	N/A*	3.6893	7.3786
100	N/A*	11.666	23.3333

*The VueLite is internally battery powered and has no I/Os. Conducted RF immunity testing does not apply, resulting in no separation data from 150 kHz to 80 MHz.

Customer Service

When it is determined that a VueLite LED Light Source needs servicing, the customer should fill out the Cook Complaint Communication Form in its entirety for submittal to Cook for service. We ask that all service requests go through Cook Medical Customer Relations via the following contact information:

Global Location	Contact Information
U.S. and Canada South/Latin America APAC	customerrelationsna@cookmedical.com
	Call 800.468.1379 x102146 to relay the customer complaint information over the phone.
Europe Middle East Africa	european.complaints@cookmedical.com
	Call the Shared Service Centre using a local country phone number, and a representative will take the details and log the complaint.

Returning Goods to Cook Medical

Ship the light source and any associated accessories in the original corrugated box to prevent damage during shipment. Please ship to the following addresses based on the global location:

Global Location	Ship To
U.S. South/Latin America APAC	Cook Medical Attn: Q.E. Complaints 1100 W. Morgan Street Spencer, IN 47460 USA
Canada	Cook Medical 65 Mostar Street Stouffville, Ontario, L4A 0Y2 Canada
Europe Middle East Africa	Cook Medical Attn: Complaints O'Halloran Road National Technology Park Limerick, Ireland

Attention: If the VueLite has been used in a clinical setting, the device must be reprocessed according to procedures outlined in the Care and Maintenance section of this manual before shipment to Cook Medical.

Upon evaluation, the customer will be contacted and advised of the findings and estimated repair cost. Repairs will not begin until authorization or a purchase order is issued indicating the approval of charges.

Warranty

The VueLite LED Light Source is warranted to be free from defects in materials and workmanship for a period of one year from the date of purchase.

All nonwarranty repairs will be warranted to be free from defects in materials and workmanship for a period of ninety days from the date of the invoice.

Accidental damage and damage resulting from misuse, abuse, excessive sterilization, disinfection and sterilization methods not approved by Cook Medical, as well as fluid invasion and normal wear and tear, will be subject to prevailing repair charges. Disassembly, alteration or repair performed by any person not authorized by Cook Medical will result in immediate loss of warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF ALL OTHER WARRANTIES EITHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

Suitability of use of the medical device for any surgical procedure shall be determined by the user. Cook Medical shall not be liable for incidental or consequential damages of any kind.

All shipping charges to Cook Medical's facility are the responsibility of the customer.

Symbol Descriptions

Attention: Read user's manual for warnings, cautions and instructions for use.

IPX7

This equipment is rated IPX7 and is protected against the effect of immersion between 15 cm and 1 m.

S/N

This symbol indicates the product's serial number.

This symbol indicates that this product is a type BF applied part.

Light intensity indicator.

Light intensity switch.

In the European Union, this symbol indicates that when the end user wishes to discard this product, it must be sent to a separate collection facility for recovery and recycling. By separating this product from other household waste, the volume of waste sent to incinerators or landfills will be reduced and natural resources will thus be conserved.

Regulatory

The VueLite LED has been tested to, and found in conformance with, the following standards:

IEC 60601-1
IEC 60601-1-1
IEC 60601-2-18
CSA C22.2 #601.1.1

MANUFACTURER
COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

www.cookmedical.com
© COOK 2012

EC REPRESENTATIVE
COOK IRELAND LTD.
O'Halloran Road
National Technology Park
Limerick, Ireland

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