Vascular Retrieval Forceps

**EN
3**

Instructions for Use

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**DA
4**

Brugsanvisning

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**DE
5**

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T - V R F - R E V 4

VASCULAR RETRIEVAL FORCEPS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).


Fig. 1

DEVICE DESCRIPTION

The Vascular Retrieval Forceps consist of a handle, a TFE-sheathed, braided stainless steel shaft, forceps with grasping jaws (**Fig. 1**) and an inserter for use during insertion through a hemostasis valve.

INTENDED USE

The Vascular Retrieval Forceps are intended for retrieval of fragments of catheter tubing, wire guides and other foreign objects from the vascular system.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Manipulation of the product requires fluoroscopic control.
- Visually inspect the product before use to make sure the shaft is free of kinks or severe angulation.
- The recommended introducer sheath should be one French size larger than the forceps size. Always check fit of forceps through intended guiding catheter or sheath prior to use.

INSTRUCTIONS FOR USE

1. Carefully remove the plastic protector which covers the forceps and distal spring coil tip.
2. Introduce the retrieval forceps into an appropriate access sheath or guiding catheter and, utilizing fluoroscopic control, advance them to the target site.
NOTE: If an introducer sheath with a hemostasis valve is utilized for access, use the supplied inserter for protection of the spring coil tip during insertion through the valve.
3. Push the button forward on the handle assembly to open the jaw of the forceps and engage the foreign body.
4. Close the jaw of the forceps by pulling the button of the handle assembly backward, which should grasp the foreign body for retrieval. Maintain pressure to keep jaws closed during retrieval.
5. With the foreign body firmly engaged, pull back the retrieval forceps into the access sheath and remove as one unit. **NOTE:** If multiple fragments are being retrieved, rinse the forceps with heparinized saline between retrievals.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.