COOK

MEDICAL

CE
0123

EN 2	CloverSnare® 4-Loop Vascular Retriever Instructions for Use
CS 4	Cévní extraktor CloverSnare® se 4 smyčkami Návod k použití
DA 7	CloverSnare® vaskulær udtrækker med 4 løkker Brugsanvisning
DE 9	CloverSnare® intravaskulärer Rückholer mit 4 Schlaufen Gebrauchsanweisung
EL 12	Αγγειακή συσκευή ανάκτησης 4 βρόχων CloverSnare® Οδηγίες χρήσης
ES 15	Recuperador vascular de 4 lazos CloverSnare® Instrucciones de uso
FI 17	4-silmukkainen vaskulaarinen CloverSnare®-poistolaite Käyttöohjeet
FR 20	Dispositif de récupération intravasculaire à 4 boucles CloverSnare® Mode d'emploi
HR 23	CloverSnare® vaskularni uređaj za uklanjanje s 4 petlje Upute za uporabu
HU 25	CloverSnare® négyhurkú vaskuláris eltávolító Használati utasítás
IT 28	Recuperatore vascolare a 4 anse CloverSnare® Istruzioni per l'uso
NL 30	CloverSnare® vasculair verwijderingsinstrument met 4 lussen Gebruiksaanwijzing
NO 33	CloverSnare® vaskulær uttrekker med 4 løkker Bruksanvisning
PL 36	4-pętlowe urządzenie do wydobywania w obrębie układu naczyniowego CloverSnare® Instrukcja użycia
PT 38	Recuperador vascular de 4 ansas CloverSnare® Instruções de utilização
SK 41	Štvorslučkový vaskulárny extraktor CloverSnare® Návod na použitie
SL 43	Vaskularna naprava za vračanje s 4 zankami CloverSnare® Navodila za uporabo
SV 46	CloverSnare® vaskulär fragmentfångare med 4 öglor Bruksanvisning
TR 48	CloverSnare® 4 Halkalı Vasküler Geri Alıcı Kullanma Talimatı

T _ V R S _ M _ R E V O

CLOVERSNARE® 4-LOOP VASCULAR RETRIEVER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The CloverSnare® 4-Loop Vascular Retriever consists of four petals composed of nitinol, for shape memory, and tantalum, for visibility. It is packaged with a coaxial sheath system with dilator. A radiopaque tip is incorporated into each sheath of the coaxial system to identify the distal ends of each sheath.

Performance Characteristics

- Snare expands when extended from and collapses when retracted into catheter to facilitate foreign body capture.
- Radiopaque wire incorporated into snare material and radiopaque material incorporated into sheath tips for positioning under imaging techniques.
- Tuohy-Borst valve aids in manipulation or retrieval of foreign body by maintaining tension on snare.
- Hemostasis valves on inner sheath and snare catheter minimize leakage.
- Tapered dilator tip and rounded sheath tips to facilitate atraumatic insertion into patient vasculature.
- Hubs of introducer dilator, inner sheath, and outer sheath connect to facilitate insertion into patient vasculature.

Device Compatibility

- Standard luer syringe.
- 0.035 inch wire guide.
- Diagnostic catheter with an outer diameter less than or equal to 2.50 mm (7.5 Fr).
- Heparinized saline.

NOTE: All components are compatible with a standard luer except the outer sheath. The outer sheath may be flushed with a standard luer syringe, but it will not form a leak-free connection with a luer-lock connector.

Patient Population

The target population for this device is adult patients with an intravascular foreign body who require and are suitable for endovascular intervention based on the assessment of their treating physician/healthcare provider. The underlying patient condition, the coagulation status of the patient, and the location of the foreign body should be considered.

Intended User

This device is intended for use by physicians trained and experienced in intravascular diagnostic and interventional techniques.

Contact with Body Tissue

This device contacts circulating blood for limited duration (≤ 24 hours).

Operating Principle

This device manipulates or retrieves a foreign body (including inferior vena cava filters) from the intended vasculature by capturing the body in wire loops and either repositions the object or retrieves the object through an intravascular sheath. The shape memory properties of the wire loops and ability of the wire loops to be collapsed by the extension of the sheath facilitate the capture of the foreign body. Fluoroscopic imaging identifies and facilitates control of both the wire loops and distal tip of the inner and outer sheaths throughout the procedure.

INTENDED USE

This device is intended to manipulate or retrieve a foreign body within the central and peripheral circulatory systems by capturing the foreign body in wire loops.

INDICATIONS FOR USE

This device is indicated for use in adult patients with medical conditions related to the presence of a vascular foreign body and for whom manipulation or retrieval of the foreign body would provide a clinical benefit.

CLINICAL BENEFITS

This device provides a minimally-invasive means to manipulate or retrieve an intravascular foreign body.

CONTRAINDICATIONS

- This device is not intended to be placed in the coronary or neuro vasculature.

For optional inferior vena cava (IVC) filter retrieval:

- Retrieval of a Cook IVC filter with significant amounts of trapped thrombus (greater than 25% of the volume of the cone).
- Retrieval of the filter for patients with an on-going high risk for pulmonary embolism.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- This device contains nitinol, a nickel titanium alloy. Allergic reactions to nickel should be considered.
- This device contains cobalt as a component of a stainless-steel alloy. Cobalt is a carcinogen and a substance toxic to reproduction. However, this stainless-steel alloy does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- This device should be used in patients who are suitable for endovascular interventions with particular attention paid to the patient's coagulation status, presence of excessive thrombus formation on the foreign body, and the risk of vessel injury posed by retrieval.
- This device should be used in patients with vessel diameters compatible with associated device components.
- Excessive force should not be used to manipulate or retrieve foreign objects.
- Do not attempt to shape the catheter tip or snare, as doing so may damage the device. Damage may include, but not be limited to, separation of the nitinol snare(s) from the snare catheter.
- Manipulation of the product requires fluoroscopic control.
- When performing inferior vena cava imaging, an appropriately sized diagnostic catheter should be used.

For optional inferior vena cava (IVC) filter retrieval:

- During filter retrieval, an inferior vena cava imaging evaluation for residual captured thrombus should be performed prior to attempted retrieval.
- Never attempt to re-deploy a retrieved filter.

POTENTIAL ADVERSE EVENTS

- Arrhythmia
- Bleeding / Hemorrhage
- Device and foreign body entrapment
- Embolization
- Heart injury
- Infection
- Myocardial infarction
- Pulmonary embolism
- Stroke
- Thrombosis
- Vessel injury

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. After removing product from the package, the sheaths, dilator, and snare should be individually flushed with heparinized saline.
2. Gain vascular access using Seldinger technique.
3. Introduce the coaxial sheaths with dilator over the wire guide and advance to the desired location.
4. Remove the wire guide and dilator, then insert the snare catheter and advance it to the desired location.
5. Loosen the screw of the clear Y-fitting to enable advancement of the snare loops. Holding the clear Y-fitting steady, advance the white pin vise. The white pin vise can be used to advance, retract, and manipulate the loops to snare the foreign body.

6. Once the foreign body is snared by the loops, hold the pin vise steady and tighten the screw of the clear Y-fitting.
NOTE: If at any time during the procedure, the distance between the pin vise and the clear Y-fitting changes, the screw at the top of the clear Y-fitting should be further tightened.
7. While holding the clear Y-fitting steady, advance the coaxial sheath system over the foreign body.
NOTE: The outer sheath of the coaxial system may be advanced over the tip of the inner sheath to cover any portion of the foreign body not contained inside the distal tip of the inner sheath.
8. Separate the hub of the inner and outer sheaths and remove inner sheath, snare catheter, and foreign body.
9. After removing the sheath, hospital standard of care should be followed for providing hemostasis to prevent bleeding at vascular access site.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device – (0827002CINC050-211_A-2J7) – to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: ec.europa.eu/tools/eudamed

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2022
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd