- EN**
3 **CloverSnare™ 4-Loop Vascular Retriever**
Instructions for Use
- DA**
4 **CloverSnare™ vaskulær udtrækker med 4 løkker**
Brugsanvisning
- DE**
5 **CloverSnare™ intravaskulärer Rückholer mit 4 Schlaufen**
Gebrauchsanweisung
- EL**
6 **Αγγειακή συσκευή ανάκτησης 4 βρόχων CloverSnare™**
Οδηγίες χρήσης
- ES**
7 **Recuperador vascular de 4 lazos CloverSnare™**
Instrucciones de uso
- FR**
9 **Anse d'extraction vasculaire à 4 boucles CloverSnare™**
Mode d'emploi
- IT**
10 **Recuperatore vascolare a 4 anse CloverSnare™**
Istruzioni per l'uso
- NL**
11 **CloverSnare™ vasculair verwijderingsinstrument met 4 lussen**
Gebruiksaanwijzing
- PT**
12 **Recuperador vascular de 4 ansas CloverSnare™**
Instruções de utilização
- SV**
14 **CloverSnare™ vaskulär hämtnara med 4 öglor**
Bruksanvisning

CLOVERSNARE™ 4-LOOP VASCULAR RETRIEVER

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The CloverSnare 4-Loop Vascular Retriever consists of four petals composed of nitinol, for shape memory, and tantalum, for visibility. It is packaged with a coaxial sheath system with dilator. Each of the sheaths of the coaxial system has radiopaque bands incorporated to identify the location of the distal tip of each sheath.

INTENDED USE

The CloverSnare 4-Loop Vascular Retriever is intended for use in the cardiovascular system to manipulate and retrieve foreign objects, including, but not limited to, wire guides, coils, balloons, catheters, and filters.

CONTRAINDICATIONS

None known

WARNINGS

This product contains nitinol, a nickel titanium alloy. Allergic reactions to nickel should be considered.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Manipulation of the product requires fluoroscopic control.
- Excessive force should not be used to manipulate or retrieve foreign objects.
- Do not attempt to shape the catheter tip or snare, as doing so may damage the device. Damage may include, but not be limited to, separation of the nitinol snare(s) from the snare catheter.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

Potential adverse procedural events related to foreign body manipulation and retrieval in the vasculature can include, but are not limited to:

- Embolization
- Stroke
- Pulmonary embolism
- Myocardial infarction
- Vessel perforation
- Device entrapment

INSTRUCTIONS FOR USE

1. After removing product from the package, the sheaths, dilator, and snare should be flushed with heparinized saline.
2. Gain vascular access using Seldinger technique.
3. Introduce the coaxial sheaths with dilator over the wire guide and advance to the desired location.
4. Remove the wire guide and dilator, then insert the snare catheter and advance it to the desired location.
5. Loosen the screw of the clear Y-fitting to enable advancement of the snare loops. Holding the clear Y-fitting steady, advance the white pin vise. The pin vise can be used to advance, retract, and manipulate loops to capture/surround foreign object.
6. Once the foreign object is captured/surrounded by the loops of the snare, hold the pin vise steady and tighten the screw of the clear Y-fitting.

NOTE: If at any time during the procedure, the distance between the pin vise and the clear Y-fitting changes, the screw at the top of the clear Y-fitting should be further tightened.
7. While holding the clear Y-fitting steady, advance the coaxial sheath system over the foreign object.

NOTE: The outer sheath of the coaxial system may be advanced over the tip of the inner sheath to cover any portion of the foreign object not contained inside the distal tip of the inner sheath.
8. Separate the hub of the inner and outer sheaths and remove inner sheath, snare catheter, and foreign object.
9. After retrieval, hospital standard of care should be followed for removing the sheath and providing hemostasis to prevent bleeding at vascular access site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a

dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd