**EN**  
**3**

# **Willscher Vas Reapproximator Set**

Instructions for Use

**CS**  
**4**

# **Souprava Willscher pro reapproximaci chámovodu**

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**DA**  
**5**

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Fig. 2

## WILLSCHER VAS REAPPROXIMATOR SET

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

### DEVICE DESCRIPTION

- One vas reapproximator plate, 3.5 cm x 5.5 cm
- Two vas skewers, .017 inch (0.43 mm) diameter, 4.5 cm long with-pencil point tip

### INTENDED USE

Used to stabilize the vas ends during microscopic vasovasostomy.

### CONTRAINDICATIONS

None known

### WARNINGS

None known

### PRECAUTIONS

The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### INSTRUCTIONS FOR USE

1. Make a 3 cm incision on each side of the scrotum and locate the previously cut vas ends.
2. Once the vas ends have been located, place a suture through the adventitia of the vas.
3. Draw the ends of the vas up through the window of the vas reapproximator plate. (**Fig.1**)
4. Pass the skewer through the adventitia of each vas a few millimeters from the end and rest the skewers in a V-notch of the plate to stabilize and approximate the vas for anastomosis. (**Fig. 2**)  
**NOTE:** If the vas ends are loose and redundant, the skewers can be widely separated on the plate. If the vas is quite taut due to extensive vas resection during vasectomy, the skewers will have to be placed close to each other in adjacent grooves.
5. Anastomose the vas using appropriately sized sutures.
6. Following anastomosis, withdraw the skewers. The newly anastomosed vas will fall back into the scrotum. Close the wound site.

### STERILIZATION INFORMATION

- Prior to reuse, the items in the set may be sterilized using ethylene oxide (ETO) gas not exceeding 120°F.
- Routine ETO gas sterilization is the sterilization method of choice.
- ETO gas sterilization and aeration should be performed in accordance with AAMI Standards and Recommended Practices with a validated sterilization cycle, obtaining a 10<sup>-6</sup> sterility assurance level.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Reusable. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

M. K. Willscher, D. E. Novick: "Simplified Technique for Microscopic Vaso-vasostomy," Urology, 15 (1980), 147.

AAMI/ISO 11135-1994-Medical Devices-Validation and Routine Control of Ethylene Oxide Sterilization.

ANSI/AAMI ST41:1999-Ethylene Oxide Sterilization in Health Care Facilities: Safety and Effectiveness.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.