EN 16 **Zenith® Spiral-Z AAA Iliac Leg with the Z-Trak™ Introduction System**

Instructions for Use

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Bruksanvisning

TABLE OF CONTENTS

Illustrations.....	11-15
1 DEVICE DESCRIPTION.....	16
1.1 Zenith Spiral-Z AAA Iliac Leg.....	16
1.2 Delivery System.....	16
2 INDICATIONS FOR USE.....	16
3 CONTRAINDICATIONS.....	16
4 WARNINGS AND PRECAUTIONS.....	16
4.1 General.....	16
4.2 Patient Selection, Treatment and Follow-Up.....	16
4.3 Pre-Procedure Measurement Techniques and Imaging.....	16
4.4 Device Selection.....	17
4.5 Implant Procedure.....	17
4.6 Molding Balloon Use.....	17
4.7 MRI Information.....	17
Static Magnetic Field.....	17
MRI-Related Heating.....	17
Image Artifact.....	17
5 ADVERSE EVENTS.....	17
5.1 Observed Adverse Events.....	17
5.2 Potential Adverse Events.....	17
Device-Related Adverse Event Reporting.....	18
6 SUMMARY OF CLINICAL STUDIES.....	18
7 PATIENT SELECTION AND TREATMENT.....	18
7.1 Individualization of Treatment.....	18
8 PATIENT COUNSELING INFORMATION.....	18
9 HOW SUPPLIED.....	18
10 CLINICAL USE INFORMATION.....	18
10.1 Physician Training.....	18
10.2 Inspection Prior to Use.....	18
10.3 Materials Required.....	18
10.4 Materials Recommended.....	18
10.5 Device Diameter Sizing Guidelines.....	19
Table 10.5.1 Spiral-Z AAA Iliac Leg Graft Sizing Guide.....	19
11 DIRECTIONS FOR USE.....	19
Anatomical Requirements.....	19
General Use Information.....	19
Pre-Implant Determinants.....	19
Patient Preparation.....	19
11.1 Zenith Spiral-Z AAA Iliac Leg System.....	19
11.1.1 Contralateral Iliac Leg Preparation/Flush.....	19
11.1.2 Ipsilateral Iliac Leg Preparation/Flush.....	19
11.1.3 Vascular Access and Angiography.....	19
11.1.4 Contralateral Iliac Leg Placement and Deployment.....	19
11.1.5 Ipsilateral Iliac Leg Placement and Deployment.....	19
11.1.6 Ipsilateral Iliac Leg Placement and Development in Conjunction with Zenith Alpha Abdominal Endovascular Graft or Zenith Low Profile Endovascular Graft.....	20
11.1.7 Molding Balloon Insertion.....	20
Final Angiogram.....	20
12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.....	20
12.1 General.....	20
12.2 Additional Surveillance and Treatment.....	20
12.3 MRI Information.....	20
Static Magnetic Field.....	20
MRI-Related Heating.....	20
Image Artifact.....	20
13 PATIENT TRACKING INFORMATION.....	21

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









11

12

ZENITH® SPIRAL-Z AAA ILIAC LEG WITH THE Z-TRAK™ INTRODUCTION SYSTEM

Suggested Instructions for Use

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CAUTION: All contents of the outer pouch (including the introduction system and endovascular grafts) are supplied sterile, for single use only.

For the Zenith product line there are several applicable Suggested Instructions for Use (IFU). This IFU describes the Suggested Instructions for Use for the Zenith Spiral-Z AAA Iliac Leg. For information on other applicable Zenith components, please refer to the following Instructions for Use:

- Zenith AAA Endovascular Graft (Zenith AAA Endovascular Graft main body);
- Zenith Flex® AAA Endovascular Graft Main Body Component;
- Zenith Alpha™ Abdominal Endovascular Graft;
- Zenith Low Profile AAA Endovascular Graft;
- Zenith Renu® AAA Ancillary Graft (main body extension and converter configurations);
- Zenith Fenestrated AAA Endovascular Graft;
- Zenith Universal Distal Body Endovascular Graft;
- Zenith Flex® AUI AAA Endovascular Graft with the Z-Trak™ Introduction System
- Zenith Branch Iliac Endovascular Graft;
- Zenith AAA Endovascular Graft Ancillary Components (Zenith AAA Endovascular Graft main body extensions, iliac leg extensions, converters and iliac plugs); and
- Coda® Balloon Catheter.

1 DEVICE DESCRIPTION

1.1 Zenith Spiral-Z AAA Iliac Leg

The Zenith Spiral-Z AAA Iliac Leg should be used in conjunction with the Zenith AAA main bodies (Flex, Fenestrated, Low Profile, Alpha, Universal Distal Body, and Flex AUI), Branch, or Renu and is a part of a modular system consisting of multiple components, most typically a bifurcated main body and two iliac legs. (**Fig. 1**) The iliac legs are constructed of full-thickness woven polyester fabric sewn to two self-expanding stainless steel Cook-Z® stents and a continuous nitinol spiral stent with braided polyester and monofilament polypropylene suture. The graft is fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents placed at the ends of the graft provide the necessary attachment and seal of the graft to the vessel wall.

1.2 Delivery System

The Zenith Spiral-Z AAA Iliac Leg is shipped pre-loaded onto a 14 French (5.3 mm OD) or 16 French (6.0 mm OD) Z-Trak Introduction System. (**Fig. 2**) The delivery system is designed for ease of use with minimal preparation. All systems are compatible with a .035 inch wire guide.

For added hemostasis, the Captor™ Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. In addition, the delivery system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated.

Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

2 INDICATIONS FOR USE

The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is indicated for use with the Zenith AAA Endovascular Graft family of products, including the Zenith Flex AAA Endovascular Graft, Zenith Alpha Abdominal Endovascular Graft, Zenith Low Profile AAA Endovascular Graft, Zenith Renu Ancillary Graft, Zenith Fenestrated AAA Endovascular Graft, Zenith Universal Distal Body Endovascular Graft, Zenith Flex AUI, or Zenith Branch Iliac Endovascular Graft, during either a primary or a secondary procedure in patients who have adequate iliac/femoral access compatible with the required introduction systems. The graft is used in combination with these products for the endovascular treatment of abdominal aortic and aorto-iliac aneurysms.

3 CONTRAINDICATIONS

The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), polypropylene, nitinol, polytetrafluoroethylene (PTFE), or gold.
- Patients with a systemic or local infection that may increase the risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 10.1, Physician Training**.
- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.

4.2 Patient Selection, Treatment and Follow-Up

- Zenith Spiral-Z iliac artery distal fixation site greater than 10 mm in length and 7.5-20 mm in diameter (measured outer wall to outer wall) is required. These sizing measurements are critical to the performance of the endovascular repair.
- For sizing requirements and a list of key anatomical elements that may affect successful aneurysm exclusion using a device from the Zenith AAA Endovascular Graft family of products, refer to the appropriate Instructions for Use.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and delivery systems of a 14 French to 16 French vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization. A vascular conduit technique may be necessary to achieve success in some patients.
- Pre-existing regions of stenosis/narrowing (less than approximately 20 mm ID in the aorta or 7 to 8 mm ID in the iliacs) have been shown to increase the risk of a thromboembolic event (e.g., graft limb occlusion). The potential for this increased risk in these patients may preclude placement of an endovascular graft. Dilatation of these regions with a noncompliant balloon and/or stent placement may be necessary to help assure maintained graft patency and to reduce the risk of a thromboembolic event. Additionally, the completion angiogram (with stiff wire guides removed) should be reviewed carefully to determine if further treatment in these regions is necessary (e.g., adjunctive ballooning or stenting). Failure to remove the stiff wire guide prior to the angiogram could mask any limb kinking or narrowing that might occur when the wire guide is removed.
- Follow-up imaging should be carefully reviewed for narrowing within the graft leg. Patients with a graft leg lumen of less than approximately 5 mm ID may be at increased risk of a thromboembolic event (e.g., graft limb occlusion). Reintervention (e.g., noncompliant ballooning or stenting in these regions) should be considered to help assure maintained graft patency and to reduce the risk of a thromboembolic event.
- Patients with poor outflow or a hypercoagulable state (e.g., cancer) may be at an increased risk of a thromboembolic event.
- The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensible inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The Zenith Spiral-Z AAA Iliac Leg has not been explicitly evaluated clinically; however, its performance is represented by the Zenith AAA Endovascular Graft Iliac Leg (a previous version of the device), which has not been evaluated in the following patient populations:
 - traumatic aortic injury
 - leaking, pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed endovascular grafts
 - uncorrectable coagulopathy
 - indispensible mesenteric artery
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - morbidly obese patients
 - less than 18 years of age
 - key anatomical elements that fall outside the sizing requirements specified in the appropriate main body or Renu Instructions for Use.
- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3, Pre-Procedure Measurement Techniques and Imaging**.
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith Spiral-Z AAA Iliac Leg. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend that angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are well defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg components.

Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with

proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

Lengths

Utilizing CT, length measurements should be determined to accurately assess length as well as planning Zenith Spiral-Z AAA Iliac Leg components. These reconstructions should be performed in sagittal, coronal, and 3-D.

- All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components or stent fracture) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

4.4 Device Selection

- Strict adherence to the Zenith Spiral-Z AAA Iliac Leg IFU sizing guide is strongly recommended when selecting the appropriate device size (**Table 10.5.1**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

4.5 Implant Procedure

(Refer to **Section 11, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith Spiral-Z AAA Iliac Leg and assure accurate apposition to the vessel wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Spiral-Z AAA Iliac Leg Graft.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis, or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith Spiral-Z AAA Iliac Leg in a location that will occlude arteries necessary to supply blood flow to organs or extremities.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Inaccurate placement and/or incomplete sealing of the Zenith Spiral-Z AAA Iliac Leg within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the internal iliac arteries.
- Inadequate overlap of the Zenith Spiral-Z AAA Iliac Leg may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implant procedure based on hospital- and physician-preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise (e.g., adequate hydration).
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization, or may rupture the aneurysm.
- Avoid damaging the graft or disturbing graft positioning after placement in the event reinstrumentation (secondary intervention) of the graft is necessary.
- Verify that the appropriate iliac leg is selected for insertion on the contralateral side of the patient before implantation.
- Excessive overlap 10 mm above the main body bifurcation may increase the risk of limb thrombosis.

4.6 Molding Balloon Use

- Do not inflate the balloon in the vessel outside of the graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care in inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.

- Confirm complete deflation of the balloon prior to repositioning.
- For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

4.7 MRI Information

Nonclinical testing has demonstrated that the Zenith Spiral-Z AAA Iliac Leg Graft is MR Conditional. It can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less.
- Highest spatial magnetic gradient field of 720 Gauss/cm

Nonclinical evaluation was conducted in a 3.0 Tesla MR system (General Electric Excite) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith Spiral-Z AAA Iliac Leg produced a temperature rise of less than or equal to 2.1 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg, for 15 minutes of MR scanning in a 1.5 Tesla, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.9 W/kg, which corresponds to a calorimetry measured value of 2.1 W/kg.

3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith Spiral-Z AAA Iliac Leg produced a temperature rise of less than or equal to 2.6 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Healthcare, 14X.M5 Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.9 W/kg, which corresponds to a calorimetry measured value of 2.7 W/kg.

Image Artifact

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 5 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region, depending on distance from the device to the area of interest.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners.

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382
Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US
Fax: 209-669-2450
Web: www.medicalert.org

5 ADVERSE EVENTS

5.1 Observed Adverse Events

For information on observed adverse events in patients receiving Zenith AAA Endovascular Grafts, please refer to the appropriate IFU from the Zenith AAA Endovascular Graft family of products. A copy is available online at www.cookmedical.com.

5.2 Potential Adverse Events

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; component separation from another graft component; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; and corrosion
- Fever and localized inflammation

- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)

Device-Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith Spiral-Z AAA Iliac Leg should be reported to Cook immediately. For customers inside the United States, to report an incident, call the Customer Relations Department at 1-800-457-4500 (24 hour) or 1-812-339-2235. For customers outside the United States, please call your distributor.

6 SUMMARY OF CLINICAL STUDIES

For clinical study information regarding patients receiving Zenith AAA Endovascular Grafts, please refer to the appropriate IFU from the Zenith AAA Endovascular Graft family of products. A copy is available online at www.cookmedical.com.

7 PATIENT SELECTION AND TREATMENT

(See **Section 4, WARNINGS AND PRECAUTIONS**)

7.1 Individualization of Treatment

Cook recommends that the Zenith Spiral-Z AAA Iliac Leg diameters be selected as described in **Table 10.5.1**. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The risks and benefits should be carefully considered for each patient before use of the Zenith Spiral-Z AAA Iliac Leg. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Spiral-Z AAA Iliac Leg
- Patient's ability to tolerate general, regional or local anesthesia
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 14 French to 16 French vascular introducer sheath
- Zenith Spiral-Z iliac artery distal fixation site greater than 10 mm in length and 7.5-20 mm in diameter (measured outer wall to outer wall)
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft

The final treatment decision is at the discretion of the physician and patient.

8 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.
- Physicians must advise all patients that it is important to seek prompt medical attention if they experience signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.

- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant. Men who undergo endovascular or open surgical repair may experience impotence.

Physicians should refer patients to the Patient Guide regarding risks occurring during or after implantation of the device. Procedure-related risks include cardiac, pulmonary, neurologic, bowel and bleeding complications. Device-related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See **Section 5.1, Observed Adverse Events** and **Section 5.2, Potential Adverse Events**). The physician should complete the Patient I.D. Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

9 HOW SUPPLIED

- The Zenith Spiral-Z AAA Iliac Legs are sterilized by ethylene oxide gas, are pre-loaded into the Z-Trak Introduction System, and are supplied in peel-open packages.
- The devices are intended for single use only. Do not re-sterilize the devices.
- The product is sterile if the package is unopened and undamaged. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- The Zenith Spiral-Z AAA Iliac Leg is loaded onto a 14 or 16 French Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a sterile gauze pad soaked in saline solution under sterile conditions.
- Do not use after the "USE BY" (expiration) date printed on the label.
- Store in a cool, dry place.

10 CLINICAL USE INFORMATION

10.1 Physician Training

CAUTION: Always have a vascular surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System are outlined below.

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection, planning and sizing.

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

10.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

10.3 Materials Required

- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

10.4 Materials Recommended

The following products are recommended for implantation of any component in the Zenith product line. For information on the use of these products, refer to the individual product's suggested instructions for use.

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm;
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide;
 - Cook .035 inch wire guides
 - Cook Nimble™ Wire Guides
- Molding Balloons;
 - Cook Coda Balloon Catheter
- Introducer sets;
 - Cook Check-Flo Introducer Sets
 - Cook Extra Large Check-Flo Introducer Sets
 - Cook Flexor Balkin Up & Over Contralateral Introducers
- Sizing catheter;
 - Cook Aorous Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters;
 - Cook Beacon Tip Angiographic Catheters

- Cook Beacon Tip Royal Flush Catheters
- Entry needles;
 - Cook single wall entry needles
- Endovascular Dilators;
 - Cook Endovascular Dilator Sets

10.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Table 10.5.1 Spiral-Z AAA Iliac Leg Graft Sizing Guide*

Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Leg Diameter ³ (mm)	Iliac Leg Label Length ⁴ (mm)	Introducer Sheath (Fr)
≤8	9	39, 56, 74, 90, 107, 122	14
9	11	39, 56, 74, 90, 107, 122	14
10 - 12	13	39, 56, 74, 90, 107, 122	14
13 - 15	16	39, 56, 74, 90	14
16 - 18	20	39, 56, 74, 90	16
19 - 20	24	39, 56, 74, 90	16

¹Maximum diameter along the proximal fixation site.

²Round measured aortic diameter to nearest mm.

³Additional considerations may affect choice of diameter.

⁴Overall leg length = label length + 22 mm docking stent.

*All dimensions are nominal.

11 DIRECTIONS FOR USE

Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.
- For additional anatomical requirements, please refer to the appropriate device IFU from the Zenith AAA Endovascular Graft family of products. A copy is available online at www.cookmedical.com.

Prior to use of the Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System, review this Suggested Instructions for Use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

General Use Information

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System. The Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is compatible with .035 inch diameter wire guides.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected.

Determinants include:

1. Femoral artery selection for introduction of the delivery system (i.e., define respective contralateral and ipsilateral iliac arteries).
2. Angulation of aortic neck, aneurysm and iliac arteries.
3. Diameters of infrarenal aortic neck and distal iliac arteries.
4. Length from the aortic bifurcation of a previously placed main body or Renu from the Zenith AAA Endovascular Graft family of products to the internal iliac arteries/attachment site(s).
5. Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
6. Degree of vascular calcification.

Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Expose selected common femoral artery using standard surgical technique.
4. Establish adequate proximal and distal vascular control of selected femoral vessel.

11.1 Zenith Spiral-Z AAA Iliac Leg System (Fig 2)

NOTE: For directions on how to place a compatible device from the Zenith AAA Endovascular Graft family of products, refer to the Instructions for Use included with that device.

11.1.1 Contralateral Iliac Leg Preparation/Flush

1. If applicable, remove gray-hubbed inner stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (Fig. 3) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (Fig. 4) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on the connecting tube.
 - NOTE:** Graft flushing solution of heparinized saline is often used.
2. Attach syringe with heparinized saline to the hub on the distal inner cannula. Flush until fluid exits the distal dilator tip. (Fig. 5)
 - NOTE:** When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

11.1.2 Ipsilateral Iliac Leg Preparation/Flush

Follow the instructions in the previous section, Contralateral Iliac Leg Preparation/Flush, to ensure proper flushing of the ipsilateral iliac leg graft and activation of the hydrophilic coating.

11.1.3 Vascular Access and Angiography

1. Puncture the selected common femoral arteries using standard technique with an 18 or 19 UT gage arterial needle. Upon vessel entry, insert:
 - Wire guides – standard .035 inch diameter, 145 cm long, J tip or Bentson Wire Guide
 - Appropriate size sheaths (e.g., 6 or 8 French)
 - Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)
2. Perform angiography to identify level(s) of aortic bifurcation and iliac bifurcations.
 - NOTE:** If fluoroscope angulation is used with an angulated neck, it may be necessary to perform angiograms using various projections.
 - NOTE:** Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

11.1.4 Contralateral Iliac Leg Placement and Deployment

1. Position the image intensifier to show both the contralateral internal iliac artery and contralateral common iliac artery.
2. Prior to the introduction of the contralateral iliac leg delivery system, inject contrast through the contralateral femoral sheath to locate the contralateral internal iliac artery.
3. Introduce the contralateral iliac leg delivery system into the artery. Advance slowly until at least one stent of the iliac leg graft overlaps within the main body and not past the radiopaque marker band positioned 30 mm from the proximal end of the iliac leg graft inside the contralateral limb of the main body. (Fig. 6) If there is any tendency for the main body graft to move during this maneuver, hold it in position by stabilizing the gray positioner on the ipsilateral side.
 - NOTE:** A radiopaque marker band is positioned 30 mm from the proximal end of the iliac leg graft to identify the maximum amount of overlap.
 - NOTE:** If difficulty is encountered advancing the iliac leg delivery system, exchange to a more supportive wire guide. In tortuous vessels the anatomy may alter significantly with the introduction of the rigid wires and sheath systems.
4. Confirm position of distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure internal iliac patency, a minimum overlap of one stent, and a maximum overlap of 30 mm within the main body endovascular graft.
 - NOTE:** Ensure the Captor Hemostatic Valve on the iliac leg introducer sheath is turned to the open position. (Fig. 7)
5. To deploy, hold the iliac leg graft in position with the gripper on the gray positioner while withdrawing the sheath. (Figs. 8 and 9) Ensure overlap is maintained.
6. Stop withdrawing the sheath as soon as the distal end of the iliac leg graft is released.
7. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise and retract inner cannula to dock tapered dilator to gray positioner. Tighten pin vise. Maintain sheath position while withdrawing gray positioner with secured inner cannula.
8. Close the Captor Hemostatic Valve on the iliac leg introducer sheath by turning it in a clockwise direction until it stops. (Fig. 10)
9. Re-check the position of the wire guide.

11.1.5 Ipsilateral Iliac Leg Placement and Deployment

NOTE: If using this device in conjunction with a Zenith Alpha Abdominal Endovascular Graft or Zenith Low Profile Endovascular Graft, proceed to Section 11.1.6, Ipsilateral Iliac Leg Placement and Deployment in Conjunction with Zenith Alpha Abdominal Endovascular Graft or Zenith Low Profile Endovascular Graft. For all other systems, continue with steps 1-7 below.

NOTE: Ensure the Captor Hemostatic Valve on the main body introducer sheath is turned to the open position. (Fig. 11)

1. Utilize the main body graft wire and sheath assembly to introduce the ipsilateral iliac leg graft. Advance dilator and sheath assembly into the main body sheath.
 - NOTE:** In tortuous vessels, the position of the internal iliac arteries may alter significantly with the introduction of the rigid wires and sheath systems.
2. Advance slowly until the ipsilateral iliac leg graft overlaps a minimum of one stent inside the ipsilateral limb of the main body. (Fig. 12)
 - NOTE:** If an overlap of greater than 55 mm is required, it may be necessary to consider use of a leg extension in the bifurcation area of the opposite side.
 - NOTE:** For use with Renu converter or Flex AUI, ensure iliac leg overlaps a minimum of one full iliac leg stent (i.e., proximal stent of iliac leg graft) inside the Renu converter.
3. Confirm position of distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure internal iliac patency.
 - NOTE:** Ensure the Captor Hemostatic Valve on the iliac leg introducer sheath is turned to the open position. (Fig. 7)
4. To deploy, stabilize the iliac leg graft with the gripper on the gray positioner while withdrawing the iliac leg sheath. If necessary, withdraw the main body sheath. (Figs. 8 and 13)
5. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to gray positioner. Tighten pin vise. Maintain main body sheath position while withdrawing the iliac leg sheath and gray positioner with secured inner cannula.
6. Close the Captor Hemostatic Valve on the main body introducer sheath by turning it in a clockwise direction until it stops.
7. Re-check the position of the wire guides. Leave sheath and wire guide in place.

11.1.6 Ipsilateral Iliac Leg Placement and Deployment in Conjunction with Zenith Alpha Abdominal Endovascular Graft or Zenith Low Profile Endovascular Graft

1. Position the image intensifier to show both the ipsilateral internal iliac artery and ipsilateral common iliac artery.
2. Prior to removing the main body sheath, inject contrast through the sheath to locate the ipsilateral internal iliac artery.
3. Remove the main body sheath.
4. Introduce the ipsilateral iliac leg delivery system, and continue advancing slowly until the proximal edge of the ipsilateral leg graft aligns with the proximal edge of the previously-placed contralateral leg graft. (Fig. 14)
5. Confirm position of the distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure internal iliac patency.
NOTE: Ensure the Captor Hemostatic Valve on the iliac leg introducer sheath is turned to the open position.
6. To deploy, stabilize the iliac leg graft with the gripper on the gray positioner while withdrawing the iliac leg sheath.
7. Under fluoroscopy and after verification of iliac leg graft position, loosen the pin vise, and retract the inner cannula to dock the tapered dilator to the gray positioner. Tighten the pin vise. Maintain sheath position while withdrawing the gray positioner with secured inner cannula.
8. Close the Captor Hemostatic Valve on the iliac leg introducer sheath by turning it in a clockwise direction until it stops.
9. Re-check the position of the wire guides. Leave the sheath and wire guide in place.

11.1.7 Molding Balloon Insertion

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counter-clockwise. (Fig. 11)
3. Advance the molding balloon over the wire guide and through the Captor Hemostatic Valve of the main body introduction system to the level of the renal arteries. Maintain proper sheath position.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise. (Fig. 10)
CAUTION: Do not inflate the balloon in the vessel outside of the graft.
5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the most proximal covered stent and the infrarenal neck, starting proximally and working in the distal direction. (Fig. 15)
CAUTION: Confirm complete deflation of balloon prior to repositioning.
CAUTION: Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
6. Withdraw the molding balloon to the ipsilateral limb overlap region and expand.
CAUTION: Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
7. Withdraw the molding balloon to the ipsilateral distal fixation site and expand.
CAUTION: Do not inflate the balloon in the vessel outside of the graft.
CAUTION: Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
8. Deflate and remove molding balloon. Transfer molding balloon onto the contralateral wire guide and into the contralateral iliac leg introduction system. Advance molding balloon to the contralateral limb overlap and expand.
CAUTION: Confirm complete deflation of balloon prior to repositioning.
CAUTION: Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
9. Withdraw the molding balloon to the contralateral iliac leg/vessel distal fixation site and expand. (Fig. 15)
CAUTION: Do not inflate the balloon in the vessel outside of the graft.
10. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
11. Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

Final Angiogram

1. Position angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of internal iliac arteries.
2. Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.
NOTE: If endoleaks or other problems are observed, refer to the Suggested Instructions for Use for the Zenith AAA Endovascular Graft Ancillary Components.
3. Repair vessels and close in standard surgical fashion.

12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

For information on imaging guidelines and postoperative follow-up, please refer to the Instructions for Use for the Zenith AAA device that was used. A copy is available online at www.cookmedical.com.

12.1 General

- The long-term performance of endovascular grafts with secondary endovascular intervention using additional components has not yet been established.
- All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.

- Physicians should evaluate patients on an individual basis and prescribe follow-up relative to the needs and circumstances of each individual patient. The minimum requirement for patient follow-up (described in the Instructions for Use for the Zenith AAA device that was used) should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (e.g., separation between components, stent fracture).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT.
- The minimum imaging follow-up for patients with Zenith AAA stent grafts is described in the Instructions for Use for the Zenith AAA device that was used, which can be found at www.cookmedical.com. Patients requiring enhanced follow-up should have interim evaluations.

12.2 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement ≥ 5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counseled that the need for subsequent reinterventions, including catheter-based and open surgical conversion, are possible following endograft placement.

12.3 MRI Information

Nonclinical testing has demonstrated that the Zenith Spiral-Z AAA Iliac Leg Graft is MR Conditional. It can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less
- Highest spatial magnetic gradient field of 720 Gauss/cm

Nonclinical evaluation was conducted in a 3.0 Tesla MR system (General Electric Excite) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith Spiral-Z AAA Iliac Leg produced a temperature rise of less than or equal to 2.1 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg, for 15 minutes of MR scanning in a 1.5 Tesla, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.9 W/kg, which corresponds to a calorimetry measured value of 2.1 W/kg.

3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In nonclinical testing, the Zenith Spiral-Z AAA Iliac Leg produced a temperature rise of less than or equal to 2.6 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, 14X.M5 Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.9 W/kg, which corresponds to a calorimetry measured value of 2.7 W/kg.

Image Artifact

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 5 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region, depending on distance from the device to the area of interest.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners.

Mail: MedAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382
Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US
Fax: 209-669-2450
Web: www.medicalert.org

13 PATIENT TRACKING INFORMATION

In addition to these Instructions for Use, the Zenith Spiral-Z AAA Iliac Leg with the Z-Trak Introduction System is packaged with a Device Tracking Form, which the hospital staff is required to complete and forward to Cook for the purposes of tracking all patients who receive the Zenith Spiral-Z AAA Iliac Leg (as required by U.S. Federal Regulation).

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd