EN
17

Zenith Alpha® Spiral-Z® Endovascular Leg

Instructions for Use

CS
22

Endovaskulární rameno Zenith Alpha® Spiral-Z®

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DA
27

Zenith Alpha® Spiral-Z® endovaskulært ben

Brugsanvisning

DE
32

Zenith Alpha® Spiral-Z® endovaskulärer Schenkel

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EL
38

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Zenith Alpha® Spiral-Z® endovaskulärt graftben

Bruksanvisning

ENGLISH	
TABLE OF CONTENTS	
Illustrations	11-16
1 DEVICE DESCRIPTION	17
1.1 Zenith Alpha Spiral-Z Endovascular Leg	17
1.2 Delivery System.....	17
2 INDICATIONS FOR USE	17
3 CONTRAINDICATIONS.....	17
4 WARNINGS AND PRECAUTIONS	17
4.1 General.....	17
4.2 Patient Selection, Treatment and Follow-Up	17
4.3 Pre-Procedure Measurement Techniques and Imaging.....	17
4.4 Device Selection	18
4.5 Implant Procedure	18
4.6 Molding Balloon Use	18
4.7 MRI Information	18
Static Magnetic Field.....	18
MRI-Related Heating	18
Image Artifact.....	18
5 ADVERSE EVENTS.....	18
5.1 Potential Adverse Events	18
Device-Related Adverse Event Reporting	18
6 PATIENT SELECTION AND TREATMENT	19
6.1 Individualization of Treatment	19
7 PATIENT COUNSELING INFORMATION.....	19
8 HOW SUPPLIED.....	19
9 CLINICAL USE INFORMATION.....	19
9.1 Physician Training.....	19
9.2 Inspection Prior to Use.....	19
9.3 Materials Required	19
9.4 Materials Recommended	19
9.5 Device Diameter Sizing Guidelines	19
Table 9.5.1 Zenith Alpha Spiral-Z Endovascular Leg	
Graft Sizing Guide	20
10 DIRECTIONS FOR USE.....	20
Anatomical Requirements	20
General Use Information	20
Pre-Implant Determinants	20
Patient Preparation	20
10.1 Zenith Alpha Spiral-Z Endovascular Leg System	20
10.1.1 Contralateral Iliac Leg Preparation/Flush	20
10.1.2 Ipsilateral Iliac Leg Preparation/Flush.....	20
10.1.3 Vascular Access and Angiography	20
10.1.4 Contralateral Iliac Leg Placement and Deployment	20
10.1.5 Ipsilateral Iliac Leg Placement and Deployment	20
10.1.6 Molding Balloon Insertion	21
Final Angiogram	21
11 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP	21
11.1 General	21
11.2 Additional Surveillance and Treatment.....	21
11.3 MRI Information.....	21
Static Magnetic Field	
MRI-Related Heating	21
Image Artifact.....	21
12 PATIENT TRACKING INFORMATION.....	21
13 REFERENCES	21

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














ZENITH ALPHA® SPIRAL-Z® ENDOVASCULAR LEG

Suggested Instructions for Use

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

CAUTION: All contents of the outer pouch (including the introduction system and endovascular graft) are supplied sterile, for single use only.

For the Zenith product line there are several applicable Suggested Instructions for Use (IFU). This IFU describes the Suggested Instructions for Use for the Zenith Alpha Spiral-Z Endovascular Leg. For information on other applicable Zenith components, please refer to the following Instructions for Use:

- Zenith® Low Profile AAA/Zenith Alpha™ Abdominal Ancillary Components;
- Zenith Alpha® Abdominal Endovascular Graft;
- Zenith Flex® AAA Endovascular Graft;
- Zenith® Fenestrated AAA Endovascular Graft;
- Zenith® Branch Endovascular Graft-Iliac Bifurcation;
- Zenith® Universal Distal Endovascular Graft;
- Zenith® Renu™ AAA Ancillary Graft; and
- Coda® Balloon Catheter.

1 DEVICE DESCRIPTION

1.1 Zenith Alpha Spiral-Z Endovascular Leg

The Zenith Alpha Spiral-Z Endovascular Leg is part of a modular system consisting of multiple components, most typically a bifurcated main body and two iliac legs. (Fig. 1) It can be used in conjunction with systems listed above. The iliac legs are constructed of woven polyester fabric sewn to five self-expanding nitinol Cook-Z® stents and a continuous nitinol spiral stent with braided polyester and monofilament polypropylene suture, providing a conduit that is intended to exclude the aneurysm from blood flow. The graft is fully stented to provide stability and the expansive force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents placed at the ends of the graft provide the necessary attachment and seal of the graft to the vessel wall.

On each iliac leg graft there are three gold markers, as explained below:

Location	
Proximal Marker 1	Proximal edge
Proximal Marker 2	16 mm from proximal edge
Distal Marker	Distal edge

1.2 Delivery System

The Zenith Alpha Spiral-Z Endovascular Leg is shipped pre-loaded into a 12 Fr (4.7 mm OD) or 14 Fr (5.3 mm OD) introduction system. (Fig. 2) The introduction system is designed for ease of use with minimal preparation. Both systems are compatible with a .035 inch wire guide.

For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. In addition, the delivery system features a Flexor introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

2 INDICATIONS FOR USE

The Zenith Alpha Spiral-Z Endovascular Leg is indicated for use with the Zenith Alpha Abdominal Endovascular Graft, Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components, Zenith Flex AAA Endovascular Graft, Zenith Renu AAA Ancillary Graft, Zenith Fenestrated AAA Endovascular Graft and Zenith Branch Endovascular Graft-Iliac Bifurcation, during either a primary or a secondary procedure in patients who have adequate iliac/femoral access compatible with the required introduction systems. The graft is used in combination with these products for the endovascular treatment of abdominal aortic and aorto-iliac aneurysms.

3 CONTRAINDICATIONS

The Zenith Alpha Spiral-Z Endovascular Leg is contraindicated in:

- Patients with known sensitivities or allergies to polyester, polypropylene, nitinol, urethane, polytetrafluoroethylene (PTFE), nylon, or gold.
- Patients with a systemic or local infection that may increase the risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith Alpha Spiral-Z Endovascular Leg should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 9.1, Physician Training.**
- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- When a 42 or 59 mm leg graft is used on the ipsilateral side, contralateral leg overlap into the contralateral main body limb should be limited to 16 mm. Failure to do so may result in occlusion of the ipsilateral limb.

4.2 Patient Selection, Treatment and Follow-Up

- Iliac artery distal fixation site greater than 10 mm in length and 8-20 mm in diameter (measured outer wall to outer wall) is required. These sizing measurements are critical to the performance of the endovascular repair.
- For sizing requirements and a list of key anatomical elements that may affect successful aneurysm exclusion, refer to the appropriate Zenith Instructions for Use.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and 12 Fr (4.7 mm OD) or 14 Fr (5.3 mm OD) introduction systems. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization or thrombosis. A vascular conduit technique may be necessary to achieve success in some patients.
- Pre-existing regions of stenosis/narrowing (less than approximately 20 mm ID in the aorta or 7 to 8 mm ID in the iliacs) have been shown to increase the risk of a thromboembolic event (e.g., graft limb occlusion). The potential for this increased risk in these patients may preclude placement of an endovascular graft. Dilatation of these regions with a noncompliant balloon and/or stent placement may be necessary to help assure maintained graft patency and to reduce the risk of a thromboembolic event. Additionally, the completion angiogram (with stiff wire guides removed) should be reviewed carefully to determine if further treatment in these regions is necessary (e.g., adjunctive ballooning or stenting). Failure to remove the stiff wire guide prior to the angiogram could mask any limb kinking or narrowing that might occur when the wire guide is removed.
- Follow-up imaging should be carefully reviewed for narrowing within the graft leg. Patients with a graft leg lumen of less than approximately 5 mm ID may be at increased risk of a thromboembolic event (e.g., graft limb occlusion). Reintervention (e.g., noncompliant ballooning or stenting in these regions) should be considered to help assure maintained graft patency and to reduce the risk of a thromboembolic event.
- Patients with poor outflow or a hypercoagulable state (e.g., cancer) may be at an increased risk of a thromboembolic event.
- The Zenith Alpha Spiral-Z Endovascular Leg is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith Alpha Spiral-Z Endovascular Leg is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensible inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The Zenith AAA family of grafts has not been formally tested in the following patient populations:
 - traumatic aortic injury
 - leaking, pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed endovascular grafts
 - uncorrectable coagulopathy
 - indispensible mesenteric artery
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - morbidly obese patients
 - patients less than 18 years of age
 - key anatomical elements that fall outside the sizing requirements specified in the appropriate main body or Renu Instructions for Use.
- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3, Pre-Procedure Measurement Techniques and Imaging.**
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in suboptimal device sizing or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith Alpha Spiral-Z Endovascular Leg. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend that angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are well defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg components.

Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

Lengths

Utilizing CT, length measurements should be determined to accurately assess length as well as planning Zenith Alpha Spiral-Z Endovascular Leg

components. These reconstructions should be performed in sagittal, coronal, and 3-D.

- **The long-term performance of this endovascular graft has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.**
- The Zenith Alpha Spiral-Z Endovascular Leg is not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.**
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components or stent fracture) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

4.4 Device Selection

- Strict adherence to the Zenith Alpha Spiral-Z Endovascular Leg IFU sizing guide is strongly recommended when selecting the appropriate device size (**Table 9.5.1**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

4.5 Implant Procedure

(Refer to **Section 10, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith Alpha Spiral-Z Endovascular Leg and assure accurate apposition to the vessel wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Alpha Spiral-Z Endovascular Leg Graft.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all components of the system together (from outer sheath to inner cannula).
- To avoid damage to the sheath, be careful to advance all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis, or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith Alpha Spiral-Z Endovascular Leg in a location that will occlude arteries necessary to supply blood flow to organs or extremities.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Inaccurate placement and/or incomplete sealing of the Zenith Alpha Spiral-Z Endovascular Leg within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the internal iliac arteries.
- Inadequate overlap of the Zenith Alpha Spiral-Z Endovascular Leg with the main body graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implant procedure based on hospital- and physician-preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Alpha Spiral-Z Endovascular Leg requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise.
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization, or may rupture the aneurysm.
- Avoid damaging the graft or disturbing graft positioning after placement in the event reinstrumentation (secondary intervention) of the graft is necessary.
- Excessive overlap of 12 mm above the main body graft bifurcation may increase the risk of limb thrombosis.

4.6 Molding Balloon Use

- Do not inflate the balloon in the vessel outside of the graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care in inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.
- For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

4.7 MRI Information

NOTE: If using this device in conjunction with another endovascular graft from the Zenith family, refer to the appropriate device's IFU for additional MRI information.

Nonclinical testing has demonstrated that the Zenith Alpha Abdominal Endovascular Graft used with Zenith Alpha Spiral-Z Endovascular Legs is **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned under the following conditions.

Static Magnetic Field

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum spatial magnetic gradient of 1600 gauss/cm (16.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode) for 15 minutes of continuous scanning.

MRI-Related Heating

Under the scan conditions provided above, the Zenith Alpha Abdominal Endovascular Graft used with Zenith Alpha Spiral-Z Endovascular Legs is not expected to result in a temperature rise of more than 1.3°C after 15 minutes of continuous scanning.

Image Artifact

The image artifact extends approximately 9.5 mm from the Zenith Alpha Abdominal Endovascular Graft used with Zenith Alpha Spiral-Z Endovascular Legs as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. The image artifact obscures the device lumen.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail:	MedicAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 (from outside the US)
Fax:	209-669-2450
Web:	www.medicalert.org

5 ADVERSE EVENTS

5.1 Potential Adverse Events

Adverse events that may occur and/or require intervention include:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Aortic fistulation to adjacent organs/anatomic structures
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; component separation from another graft component; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Limb stenosis
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Radiation injury and/or late malignancy
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)

Device-Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith Alpha Spiral-Z Endovascular Leg should be reported to Cook immediately. For customers inside the United States, to report an incident, call the Customer Relations Department at 1-800-457-4500 (24 hour) or 1-812-339-2235. For customers outside the United States, please call your distributor.

6 PATIENT SELECTION AND TREATMENT

(See Section 4, WARNINGS AND PRECAUTIONS)

6.1 Individualization of Treatment

Cook recommends that the Zenith Alpha Spiral-Z Endovascular Leg diameters be selected as described in Table 9.5.1. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The risks and benefits should be carefully considered for each patient before use of the Zenith Alpha Spiral-Z Endovascular Leg. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Alpha Spiral-Z Endovascular Leg
- Patient's ability to tolerate general, regional or local anesthesia
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories with a delivery profile of 12 Fr to 14 Fr
- Iliac artery distal fixation site greater than 10 mm in length and 8-20 mm in diameter (measured outer wall to outer wall)
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft

The final treatment decision is at the discretion of the physician and patient.

7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of this endovascular graft has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.**
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise every patient that it is important to seek prompt medical attention if they experience signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant.
- Men who undergo endovascular or open surgical repair may experience impotence.

Physicians should refer patients to the the appropriate patient guide for the associated main body or Renu graft that is being used in conjunction with this iliac leg graft regarding risks occurring during or after implantation of the device. Procedure-related risks include cardiac, pulmonary, neurologic, bowel and bleeding complications. Device-related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See **Section 5.1, Potential Adverse Events**). The physician should complete the Patient I.D. Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

8 HOW SUPPLIED

- **STERILE - DO NOT RESTERILIZE - SINGLE USE ONLY**
- The Zenith Alpha Spiral-Z Endovascular Leg is sterilized by ethylene oxide gas, is pre-loaded into an introduction system, and is supplied in peel-open packages.
- The devices are intended for single use only. Do not re-sterilize the devices.
- The product is sterile if the package is unopened and undamaged. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

- The Zenith Alpha Spiral-Z Endovascular Leg is loaded into a 12 Fr (4.7 mm OD) or 14 Fr (5.3 mm OD) Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a sterile gauze pad soaked in saline solution under sterile conditions.
- Do not use after the "USE BY" (expiration) date printed on the label.
- Store in a cool, dry place.

9 CLINICAL USE INFORMATION

9.1 Physician Training

CAUTION: Always have a vascular surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Alpha Spiral-Z Endovascular Leg should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Alpha Spiral-Z Endovascular Leg are outlined below.

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection, planning and sizing.

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

9.3 Materials Required

- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

9.4 Materials Recommended

The following products are recommended for implantation of any component in the Zenith product line. For information on the use of these products, refer to the individual product's suggested instructions for use.

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm;
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide;
 - Cook .035 inch wire guides
 - Cook Nimble™ Wire Guides
- Molding Balloons;
 - Cook Coda Balloon Catheter
- Introducer sets;
 - Cook Check-Flo Introducer Sets
 - Cook Extra Large Check-Flo Introducer Sets
 - Cook Flexor Balkin Up & Over Contralateral Introducers
- Sizing catheter;
 - Cook Auros Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters;
 - Cook Beacon Tip Angiographic Catheters
 - Cook Beacon Tip Royal Flush Catheters
- Entry needles;
 - Cook single wall entry needles
- Endovascular Dilators;
 - Cook Endovascular Dilator Sets

9.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Table 9.5.1 Zenith Alpha Spiral-Z Endovascular Leg Graft Sizing Guide*

Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Leg Diameter ³ (mm)	Iliac Leg Label Length ⁴ (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
8	9	42, 59, 77, 93, 110, 125	12	4.0, 4.7
9	11	42, 59, 77, 93, 110, 125	12	4.0, 4.7
10 - 12	13	42, 59, 77, 93, 110, 125	12	4.0, 4.7
13 - 15	16	42, 59, 77, 93	12	4.0, 4.7
16 - 18	20	42, 59, 77, 93	14	4.7, 5.3
19 - 20	24	42, 59, 77, 93	14	4.7, 5.3

¹Maximum diameter along the distal fixation site.

²Round measured aortic diameter to nearest mm.

³Additional considerations may affect choice of diameter.

⁴Approximate overall leg length = label length + 28 mm.

*All dimensions are nominal.

10 DIRECTIONS FOR USE

Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.
- Iliac artery distal fixation site should be greater than 10 mm in length and 8 to 20 mm in diameter (measured outer wall to outer wall).
- For additional anatomical requirements, please refer to the appropriate main body or Renu IFU from the Zenith AAA family of products. A copy is available online at www.cookmedical.com.

Prior to use of the Zenith Alpha Spiral-Z Endovascular Leg, review this Suggested Instructions for Use booklet. The following instructions embody a guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

General Use Information

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Alpha Spiral-Z Endovascular Leg. The Zenith Alpha Spiral-Z Endovascular Leg is compatible with .035 inch diameter wire guides.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected.

Determinants include:

1. Femoral artery selection for introduction of the delivery system (i.e., define respective contralateral and ipsilateral iliac arteries).
2. Angulation of aortic neck, aneurysm and iliac arteries.
3. Diameters of infrarenal aortic neck and distal iliac arteries.
4. Length from the aortic bifurcation of a previously placed main body or Renu from the Zenith AAA Endovascular Graft family of products to the internal iliac arteries/attachment site(s).
5. Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
6. Consider the degree of vascular calcification, stenosis and narrowing.

Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position the patient on the imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Expose the selected common femoral artery using standard surgical technique.
4. Establish adequate proximal and distal vascular control of the selected femoral vessel.

10.1 Zenith Alpha Spiral-Z Endovascular Leg System (Fig 2)

NOTE: For directions on how to place a main body or other Zenith devices refer to the appropriate Instructions for Use included with that Zenith device.

10.1.1 Contralateral Iliac Leg Preparation/Flush

1. Remove the Peel-Away sheath from the back of the hemostatic valve. (Fig. 3) Elevate the distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (Fig. 4) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close the stopcock on the connecting tube.
NOTE: Graft flushing solution of heparinized saline is often used.
2. Attach a syringe with heparinized saline to the black hub of the inner cannula. Flush until fluid exits the distal dilator tip. (Fig. 5)
NOTE: When flushing the system, elevate the distal end of system to facilitate removal of air.
3. Soak sterile gauze pads in saline solution and use them to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both the sheath and dilator liberally.

10.1.2 Ipsilateral Iliac Leg Preparation/Flush

Follow the instructions in the previous section, Contralateral Iliac Leg Preparation/Flush, to ensure proper flushing of the ipsilateral iliac leg graft and activation of the hydrophilic coating.

10.1.3 Vascular Access and Angiography

1. Puncture the selected common femoral arteries using standard technique with an 18 or 19 UT gage arterial needle. Upon vessel entry, insert:

- Wire guides – standard .035 inch diameter, 145 cm long
 - Appropriate size sheaths (e.g., 6 or 8 Fr)
 - Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)
2. Perform angiography to identify level(s) of aortic bifurcation and iliac bifurcations.
NOTE: If fluoroscope angulation is used with an angulated neck, it may be necessary to perform angiograms using various projections.

10.1.4 Contralateral Iliac Leg Placement and Deployment

CAUTION: Verify that the contralateral iliac leg is selected.

NOTE: When using a 42 or 59 mm ipsilateral leg graft, contralateral leg overlap should be limited to 16 mm.

1. Position the image intensifier to show both the contralateral internal iliac artery and contralateral common iliac artery.
2. Prior to introduction of the contralateral iliac leg delivery system, inject contrast through the contralateral femoral sheath to locate the contralateral internal iliac artery.
3. Remove the femoral sheath and introduce the contralateral iliac leg delivery system into the artery. Advance slowly until the zone between the two proximal gold markers aligns with the gold radiopaque marker at the bifurcation of the main body graft. This will allow between 16 mm and 32 mm of overlap between components. (Fig. 6) If there is any tendency for the main body graft to move during this maneuver, hold it in position by stabilizing the positioner on the ipsilateral side.

NOTE: If used in conjunction with a Zenith Flex main body, Zenith Renu, Zenith Fenestrated distal body, or Zenith Universal Distal body, advance the introduction system slowly until the proximal gold marker of the leg graft aligns with the proximal edge of the first limb stent. This provides the minimum one stent overlap. Advancing the proximal gold marker an additional one-half stent length into the limb, provides the maximum one and one-half stent overlap.

NOTE: On each iliac leg there are three gold markers, as explained below:

	Location	Graft Overlap
Proximal Marker 1	Proximal edge	When marker is aligned with the main body bifurcation marker, there is an overlap of 16 mm on the contralateral side (minimum recommended overlap)
Proximal Marker 2	16 mm from proximal edge	When this marker is aligned with the main body bifurcation marker, there is an overlap of 32 mm on the contralateral side (maximum recommended overlap)
Distal Marker	Distal edge	N/A

- NOTE:** If difficulty is encountered advancing the iliac leg delivery system, exchange to a more supportive wire guide. In tortuous vessels, the anatomy may alter significantly with the introduction of the rigid wires and sheath systems.
4. Confirm position of the distal end of the contralateral iliac leg graft. Reposition the contralateral iliac leg graft if necessary to ensure both internal iliac patency and minimum overlap of 2 stents (16 mm) within the main body endovascular graft.
 5. To deploy, hold the contralateral iliac leg graft in position with the gray positioner while withdrawing the sheath approximately 10 mm. (Figs. 8 and 9)
 6. Check the graft position and reposition if necessary.
 7. Continue to deploy the graft by withdrawing the sheath while continuously checking the position of the graft.
 8. Stop withdrawing the sheath as soon as the distal end of the contralateral iliac leg graft is released.
 9. Under fluoroscopy and after verification of iliac leg graft position, loosen the pin vise and retract the inner cannula to dock the tapered dilator to the positioner. Tighten the pin vise. Maintain sheath position while withdrawing the gray positioner with secured inner cannula. (Fig. 10)
 10. Re-check the position of the wire guide.

10.1.5 Ipsilateral Iliac Leg Placement and Deployment

NOTE: Ensure the Captor Hemostatic Valve on the introducer sheath is turned to the open position. (Fig. 11)

1. Position the image intensifier to show both the ipsilateral internal iliac artery and ipsilateral common iliac artery.
2. Prior to the introduction of the ipsilateral iliac leg delivery system, inject contrast through the main body femoral sheath to locate the ipsilateral internal iliac artery.

- Utilize the main body graft wire and sheath assembly to introduce the ipsilateral iliac leg graft. Advance the dilator and sheath assembly into the main body sheath.
NOTE: In tortuous vessels, the position of the internal iliac arteries may alter significantly with the introduction of the rigid wires and sheath systems.
- Continue advancing slowly until the zone between the two proximal gold markers aligns with the gold radiopaque marker at the bifurcation of the main body graft. Closely align the proximal edge of the ipsilateral leg graft with the proximal edge of the previously placed contralateral leg graft. (**Fig. 12 and 13**)
NOTE: If this device is used in conjunction with a Zenith Flex main body or Zenith Universal Distal body, advance the ipsilateral leg graft until it is closely aligned with the previously placed contralateral leg graft.
NOTE: For use with a Renu converter, ensure that the iliac leg overlaps a minimum of two full iliac leg stents (i.e., proximal stents of the iliac leg graft) inside the Renu converter.
- Confirm position of the distal end of the iliac leg graft. Using the gold markers as reference points, reposition the iliac leg graft if necessary to ensure internal iliac patency.
- To deploy, hold the iliac leg graft in position with the gray positioner while withdrawing the sheath approximately 10 mm. (**Figs. 14 and 15**)
- Check the graft position and reposition if necessary.
- Continue to deploy the graft by withdrawing the sheath while continuously checking the position of the graft.
- Under fluoroscopy and after verification of iliac leg graft position, loosen the pin vise, and retract the inner cannula to dock the tapered dilator to the positioner. Tighten the pin vise. Maintain sheath position while withdrawing the positioner with secured inner cannula. (**Fig. 16**)
- Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops. (**Fig. 7**)
- Re-check the position of the wire guides. Leave the sheath and wire guides in place.

10.1.6 Molding Balloon Insertion

- Prepare the molding balloon as follows:
 - Flush the wire lumen with heparinized saline.
 - Remove all air from the balloon.
- In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counter-clockwise.
- Advance the molding balloon over the wire guide and through the Captor Hemostatic Valve of the main body introduction system to the level of the renal arteries. Maintain proper sheath position.
- Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.
CAUTION: Do not inflate the balloon in the vessel outside of the graft.
- Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the most proximal covered stent and the infrarenal neck, starting proximally and working in the distal direction. (**Fig. 17**)
CAUTION: Confirm complete deflation of the balloon prior to repositioning.
CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
- Withdraw the molding balloon to the ipsilateral limb overlap region and expand.
CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
- Withdraw the molding balloon to the ipsilateral distal fixation site and expand.
CAUTION: Do not inflate the balloon in the vessel outside of the graft.
CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
- Deflate and remove the molding balloon. Transfer the molding balloon onto the contralateral wire guide and into the contralateral iliac leg introduction system. Advance the molding balloon to the contralateral limb overlap and expand.
CAUTION: Confirm complete deflation of the balloon prior to repositioning.
CAUTION: Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.
- Withdraw the molding balloon to the contralateral iliac leg/vessel distal fixation site and expand. (**Fig. 17**)
CAUTION: Do not inflate the balloon in the vessel outside of the graft.
- Remove the molding balloon and replace it with an angiographic catheter to perform completion angiograms.
- Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

Final Angiogram

- Position an angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of internal iliac arteries.
- Confirm there are no endoleaks or kinks and verify position of the proximal gold radiopaque markers. Remove the sheaths, wires and catheters.
NOTE: If endoleaks or other problems are observed, refer to the Suggested Instructions for Use for the appropriate main body or Renu IFU.
- Repair vessels and close in standard surgical fashion.

11 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

For information on imaging guidelines and postoperative follow-up, please refer to the Instructions for Use for the Zenith devices that were used. A copy is available online at www.cookmedical.com.

11.1 General

- The long-term performance of endovascular grafts with secondary endovascular intervention using additional components has not yet been established.

- All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.
- Physicians should evaluate patients on an individual basis and prescribe follow-up relative to the needs and circumstances of each individual patient. The minimum requirement for patient follow-up (described in the Instructions for Use for the Zenith AAA device that was used) should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (e.g., separation between components, stent fracture).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound.
- The minimum imaging follow-up for patients with Zenith AAA stent grafts is described in the Instructions for Use for the Zenith AAA main body or Renu that was used, which can be found at www.cookmedical.com. Patients requiring enhanced follow-up should have interim evaluations.

11.2 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement $\geq$ 5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counseled that the need for subsequent reinterventions, including catheter-based and open surgical conversion, are possible following endograft placement.

11.3 MRI Information

NOTE: If using this device in conjunction with another endovascular graft from the Zenith family, refer to the appropriate device's IFU for additional MRI information.

Nonclinical testing has demonstrated that the Zenith Alpha Abdominal Endovascular Graft used with Zenith Alpha Spiral-Z Endovascular Legs is **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned under the following conditions.

Static Magnetic Field

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum spatial magnetic gradient of 1600 gauss/cm (16.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of $\leq$ 2.0 W/kg (Normal Operating Mode) for 15 minutes of continuous scanning.

MRI-Related Heating

Under the scan conditions provided above, the Zenith Alpha Abdominal Endovascular Graft used with Zenith Alpha Spiral-Z Endovascular Legs is not expected to result in a temperature rise of more than 1.3°C after 15 minutes of continuous scanning.

Image Artifact

The image artifact extends approximately 9.5 mm from the Zenith Alpha Abdominal Endovascular Graft used with Zenith Alpha Spiral-Z Endovascular Legs as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. The image artifact obscures the device lumen.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail:	MediAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medicalert.org

12 PATIENT TRACKING INFORMATION

In addition to these Instructions for Use, the Zenith Alpha Spiral-Z Endovascular Leg is packaged with a Device Tracking Form, which the hospital staff is required to complete and forward to Cook for the purposes of tracking all patients who receive the Zenith Alpha Spiral-Z Endovascular Leg (as required by U.S. Federal Regulation).

13 REFERENCES

These instructions for use are based on experience from physicians (and/or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020