EN 3	Zaontz Urethral Stent Instructions for Use
CS 4	Uretrální stent Zaontz Návod k použití
DA 5	Zaontz urethrastent Brugsanvisning
DE 6	Zaontz Urethrastent Gebrauchsanweisung
EL 7	Ουρηθρική ενδοπρόσθεση Zaontz Οδηγίες χρήσης
ES 8	Stent uretral Zaontz Instrucciones de uso
FR 9	Endoprothèse urétrale de Zaontz Mode d'emploi
HU 10	Zaontz urethrális sztent Használati utasítás
IT 11	Stent uretrale Zaontz Istruzioni per l'uso
NL 12	Zaontz urethrastent Gebruiksaanwijzing
NO 13	Zaontz uretrastent Bruksanvisning
PL 14	Stent Zaontz do cewki moczowej Instrukcja użycia
PT 15	Stent uretral Zaontz Instruções de utilização
SV 16	Zaontz uretral stent Bruksanvisning

ZAONTZ URETHRAL STENT

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

This device is constructed from C-Flex tubing. The stent is available in 6, 8, and 10 French diameter sizes and the length is 12 cm. The stent has two drainage sideports at the distal tip. The stent's proximal tip is flared. The 6 French stent has 2 suture sideports and the 8 and 10 French stents have four suture sideports at the proximal flared tip.

INTENDED USE

This device is used for stenting the urethra during hypospadias or epispadias repair and to allow postoperative drainage of the bladder in infants (29 days to less than 2 years old) and children (2 years old to less than 12 years old).

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions that create unacceptable risk during catheterization.

WARNINGS

None known

PRECAUTIONS

This device is intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.

POTENTIAL ADVERSE EVENTS

- Urinary tract infection
- Perforation or laceration of the urethra
- Perforation or laceration of the bladder
- Abrasions

INSTRUCTIONS FOR USE

1. Following the procedure, advance the stent into the urethra and to the bladder.
2. Using absorbable suture, anchor the stent to the glans meatus. **NOTE:** The excess portion of the stent can be trimmed if necessary.

NOTE: The stent must not remain indwelling for more than 21 days. During the healing process, the sutures will be absorbed, and the stent should be expelled spontaneously during micturition.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd