EN
3**Holmium Laser Fiber (Single-Use)**

Instructions for Use

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
Fig. 2

HOLMIUM LASER FIBER (SINGLE-USE)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTENDED USE

Used with holmium laser systems for the treatment of kidney, ureter, and bladder stones, as well as soft tissue ablation.

CONTRAINDICATIONS

None known

WARNINGS

All personnel present in the laser hazard zone must wear all suggested protective devices. Wear laser safety eyewear per the laser manufacturer's specifications.

PRECAUTIONS

- A number of different factors affect the life of any particular fiber, including:
 - Extended lasing at high power
 - Continuous lasing with the fiber tip in contact with tissue
 - Lasing with a contaminated or damaged proximal end
 - Improper handling
 - Poor laser beam alignment or focus
- Never subject fiberoptics to sharp bends in handling, use, or storage.
- Always keep connector-end dry and free from contaminants.
- Discard any fiberoptic assembly that is cracked or broken, or does not meet minimum transmission standards.
- Do not exceed recommended power limits.

Fiber Core Size	Power Limit (Single Cavity Lasers)
150 μm	6 W
200 μm	10 W
273 μm	15 W
365 μm	20 W

NOTE: In addition to above recommended fiber handling protocol and precautions, all other applicable hospital protocols required to ensure patient and hospital personnel safety should be followed.

NOTE: Due to the extremely small fiber core size, the output power from 150 or 200 μm core fibers may be more than 20% lower than set power.

INSTRUCTIONS FOR USE

1. Open and remove the packaging tray from the sterile pouch.
2. Remove the Fiber Clip from the packaging tray.
3. Peel and remove the backing from the Fiber Clip. (**Fig. 1**)
4. Adhere the clip to the desired location.

- Remove the proximal connector **and tubing** from the packaging tray.
- Hand the proximal connector to a non-sterile attendant.
NOTE: The laser system is NOT sterile. Anyone who touches the laser system to insert the fiber will no longer be "scrubbed in."
- Open any laser aperture cover or door and insert the proximal connector into the laser system; screw in finger-tight.
- Pull the coiled tubing into the sterile field and allow the fiber to slide out of the tubing.
NOTE: Keep output of tubing in line with fiber to minimize drag as the fiber is exiting the coil.
- Slide the fiber into the lower clip to hold the fiber in place.
NOTE: Use the upper clip to hold the fiber tip when not in use. (Fig. 2)
- Follow the laser manufacturer's directions for use.
- Remove the fiber from the Fiber Clip and discard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404, U.S.A
COOK IRELAND LTD.