
- EN**  
2  
**Entuit® Gastrostomy BR** Balloon Retention Feeding Tube with **ENFit®** Connection  
Instructions for Use
- CS**  
6  
**Balónková retenční výživová sonda Entuit® Gastrostomy BR** s konektorem **ENFit®**  
Návod k použití
- DA**  
10  
**Entuit® Gastrostomy BR** ernæringssonde med retentionsballon og **ENFit®** konnektor  
Brugsanvisning
- DE**  
14  
**Entuit® Gastrostomy BR** Ernährungssonde mit Retentionsballon und **ENFit®** Verbindung  
Gebrauchsanweisung
- EL**  
18  
Σωληνας σιτισης με μπαλονι συγκρατησης **Entuit® Gastrostomy BR** με συνδεση **ENFit®**  
Οδηγιες χρησης
- ES**  
22  
Sonda de alimentación de retención con balón **Entuit® Gastrostomy BR** con conexión **ENFit®**  
Instrucciones de uso
- FR**  
26  
Sonde d'alimentation avec ballonnet de rétention **Entuit® Gastrostomy BR** et connexion **ENFit®**  
Mode d'emploi
- HU**  
30  
**Entuit® Gastrostomy BR** Ballonnal megtartott táplálószonda **ENFit®** csatlakozóval  
Használati utasítás
- IT**  
34  
Sondino con sistema di ritenzione a palloncino **Entuit® Gastrostomy BR** con connettori **ENFit®**  
Istruzioni per l'uso
- KO**  
38  
**ENFit®** 연결 장치 통합형 **Entuit® Gastrostomy BR** 풍선 유지 피딩 튜브  
사용 설명서
- NL**  
42  
**Entuit® Gastrostomy BR**-voedingssonde met fixatieballon en **ENFit®**-aansluiting  
Gebruiksaanwijzing
- NO**  
46  
**Entuit® Gastrostomy BR** ernæringssonde med ballongretensjon og **ENFit®**-tilkobling  
Bruksanvisning
- PL**  
50  
**Entuit® Gastrostomy BR** - cewnik żywieniowy utrzymywany balonem z łącznikiem **ENFit®**  
Instrukcja użycia
- PT**  
54  
Tubo de alimentação com balão de retenção **Entuit® Gastrostomy BR** com conexão **ENFit®**  
Instruções de utilização
- SV**  
58  
**Entuit® Gastrostomy BR** matnings slang med retentionsballong och **ENFit®**-anslutning  
Bruksanvisning

## ENTUIT® GASTROSTOMY BR BALLOON RETENTION FEEDING TUBE WITH ENFit® CONNECTION

### Instructions For Use:

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

### Contents of package:

1 - Entuit® Gastrostomy BR Balloon Retention Feeding Tube with ENFit® Connection

### DEVICE DESCRIPTION

The feeding tube is a sterile device consisting of a silicone balloon, bi-lumen shaft, three-port funnel and bolster that allows for proper retention during enteral feeding, medication administration and decompression. The device contains one valve that allows for inflation and deflation of the silicone balloon.

The device is intended to connect to the following enteral devices:

- Enteral giving sets with ENFit connector for feeding administration
- Syringes with ENFit connector for feeding / medication administration
- Inflation syringe with Luer connector for balloon inflation and deflation

### INTENDED USE/INDICATIONS FOR USE

The feeding tube is intended to provide gastric access for enteral feeding, medication administration and decompression through an established gastrointestinal stoma tract.

The feeding tube is indicated for use in percutaneous placement of an enteral feeding tube in adult and pediatric patients that require enteral feeding, medication administration or decompression through an established gastrointestinal stoma tract.

### CONTRAINDICATIONS

Device placement is contraindicated for patients with evidence of granulation tissue, infection, and/or irritation in the stoma site.

### WARNING

- Inspect package integrity before use. Do not use device if package is damaged or if the sterile barrier has been compromised. Do not use if labeling is incomplete or illegible.
- Do not resterilize or reprocess this medical device as this may have an adverse effect on the known characteristics of the structural integrity, performance and biocompatibility of the device.
- For single patient use only. Do not reuse this medical device as this may increase the risk of contamination leading to transmission of infectious diseases which has the potential of resulting in patient injury, illness or death.
- **FOR ENTERAL USE ONLY.** The device is intended to connect to enteral giving sets for enteral feeding and syringes for feeding, medication and inflation of balloon.
- This device incorporates ENFit connectors which are designed to reduce the likelihood of tubing misconnections; however, the potential to misconnect this device with medical devices of other healthcare applications still exists.
- This device contains small parts that are a potential choking hazard.
- Do not attempt to force non-enteral connectors to this device.
- Do not overtighten ENFit connectors. Full engagement of connector is achieved once resistance is felt, typically after half to one full turn of the connector, as illustrated below.

- Clean device/connectors with **ONLY** mild soap and water
- Do **NOT** use this product in the vasculature.
- Once used, dispose of packaging and package contents in accordance with healthcare institution guidelines and/or local government policy.

### PRECAUTIONS

- This device should only be used for its intended purpose by or under the supervision of trained healthcare professionals with a comprehensive understanding in clinical principles, procedures and risks associated with percutaneous placement of enteral feeding devices. It is recommended to adhere to the instructions for use provided with this device, the enteral feeding instructions for the healthcare facility and instructions recommended by physicians.

#### COMPLICATIONS/ADVERSE EVENTS

Complications and adverse events associated with placement and use of a balloon gastrostomy tube include, but are not limited to:

- Aspiration, reflux, sepsis, ascites, bleeding, peritonitis and perforation
- Granulation tissue, pressure necrosis and ulcers
- Irritation and infection such as redness, edema or purulent drainage
- Severe gastroesophageal reflux or diffuse inflammatory, infectious, or neoplastic disease involving the walls of the abdomen or anterior stomach
- Gastrointestinal obstruction and proximal small bowel fistula
- Tube clogging, kinking, malposition, migration, leakage and unintended tube dislodgement

#### HOW SUPPLIED

The feeding tube is provided sterile with one (1) Instructions for Use for each device.

#### HANDLING AND STORAGE

Store in a cool, dry, and secure environment to prevent damage.

#### DEVICE PREPARATION

**Warning:** *An inappropriately sized tube may cause leakage of gastric contents, necrosis, granulation tissue, sepsis, associated sequelae and/or buried bumper syndrome.*

1. Select the appropriately sized feeding tube based on the diameter of patient's stoma.
2. Inspect package integrity before use. Do not use device if package is damaged or if the sterile barrier has been compromised.
3. Remove device from package and inspect contents for damage. Do not use contents if damaged.
4. Slide external bolster up and down shaft several times to ensure bolster can be adjusted to properly secure device when placed.
5. Inflate the feeding tube balloon to recommended fill volume (*printed on the device*) with sterile or distilled water using a sterile syringe.

**Caution:** *Only use sterile or distilled water for balloon inflation. Do not use air, saline, feeding formula, medication, or radiopaque contrast for balloon inflation. Inflate balloon slowly to prevent buildup of excess pressure in inflation system.*

6. Verify balloon integrity by visually inspecting inflated balloon for symmetry and by gently squeezing the balloon to determine if the balloon is leaking. Balloon symmetry may be improved by gently manipulating the inflated balloon up and down. Do not use device if any suspected characteristics that may affect the performance of the balloon are noted.

#### DEVICE PREPARATION (continued)

7. Deflate balloon by reinserting syringe and removing all water from balloon. Do not attempt to use device if balloon cannot be deflated.
8. If necessary, lubricate the tip of the device with a water soluble lubricant for ease of insertion. Sterile gauze pads may be used for ease of tip lubrication.

**Note:** *Do not use petroleum jelly or mineral oil for tip lubrication.*

#### DEVICE PLACEMENT

**Caution:** *Prior to device placement, a Gastrostomy must have been performed to affix the stomach wall to the abdominal wall, an incision site must be identified for creation of stoma tract, then the stoma must be dilated and finally measured, in order to ensure patient safety and comfort.*

**Caution:** *Confirm that the balloon of the feeding tube will not occlude the Pylorus when inflated, specifically, when device will be inserted into infants or children.*

**Warning:** *Do not attempt to affix the stomach to the abdominal wall using the feeding tube. Excessive tension, tissue necrosis and inadvertent tube removal may result in separation of the gastric and abdominal wall.*

1. **For Device Replacement:** Clean stoma site with mild soap and water. Allow area to dry.
2. Select and prepare the appropriately sized feeding tube in accordance with the "DEVICE PREPARATION" instructions in the section above.
3. Carefully insert device through the established stoma and into the stomach lumen.  
**Note:** *Maintain insertion angle perpendicular to the surface of the skin.*
4. Verify that device is inside the stomach lumen using fluoroscopic procedures to determine position of device using the radiopaque mark located on the tip of the device or by following the instructions provided in the "VERIFICATION OF TUBE PLACEMENT" section listed below.
5. Inflate the balloon to recommended fill volume (printed on the device) with sterile or distilled water using a sterile syringe.  
**Caution:** *Only use sterile or distilled water for balloon inflation. Do not use air, saline, feeding formula, medication, or radiopaque contrast for balloon inflation. Inflate balloon slowly to prevent buildup of excess pressure in inflation system.*
6. Gently pull device toward abdomen until balloon is up against the inner stomach wall.
7. Slide the external bolster down to the skin.
8. Verify that there is no leakage of gastric contents around the stoma site and that there is not excessive pressure between the balloon and external bolster.

**Warning:** *Inflating balloon with less volume than the recommended volume may result in leakage of gastric contents and/or device migration. Inflating balloon with more volume than the recommended volume may decrease balloon performance.*


9. Clean any remaining contents around the stoma site with mild soap and water.

#### VERIFICATION OF TUBE PLACEMENT

**Warning:** Verify proper tube placement before feeding administration.

**Note:** This device incorporates ENFit connectors from the standards ISO 80369-1 and AAMI/CN3-2014 Part 3 for improved patient safety. The ENFit connectors are designed to reduce the risk of tubing misconnections between unrelated medical devices.

1. Remove the cap from the device by rotating cap counter-clockwise.

2. Draw 10mL of water into an ENFit compatible syringe and connect to the feeding port of the funnel by rotating the syringe clockwise (**do not overtighten**).
3. Gently pull back on syringe plunger until stomach contents (typically yellow or clear unless food is in patient's stomach) are aspirated and visible in the feeding lumen.

**Note:** The presence of gastric contents in the feeding lumen confirms correct positioning and that feeding tube is inside the stomach lumen. Flush contents with the 10mL of water in the syringe.

4. Flush contents with the 10mL of water in the syringe.
5. Disconnect syringe from device by rotating syringe counter-clockwise and reconnect cap by rotating cap clockwise (**do not overtighten**).

#### FEEDING ADMINISTRATION

**Caution:** Only begin feeding administration after confirmation of proper feeding tube placement and patency. Administration of feeding should be performed according to instructions provided by the physician.

**Caution:** Do not attempt to administer solid nutrients through feeding lumen.

**Warning:** This device is intended **FOR ENTERAL USE ONLY**.

**Note:** Verify the ENFit connector is connected to an enteral port only and not to an I.V. set prior to feeding.

1. Remove the cap from the "FEED" port by rotating the cap counter-clockwise.
2. If bolus feeding with an ENFit compliant syringe:
  - 2.1 Fill ENFit compliant syringe with water and connect to the feeding port by rotating syringe clockwise (**do not overtighten**).
  - 2.2 Prime the feeding lumen with water using the water-filled syringe and inspect for leakage.  
**Note:** The device may need to be replaced if leakage is present.
- 2.3 Disconnect syringe by rotating counter-clockwise. Remove remaining water, if any, and plunger from syringe.
- 2.4 Reconnect the syringe without the plunger to the feeding port by rotating the syringe clockwise (**do not overtighten**).
- 2.5 Pour feeding formula into syringe.
- 2.6 Raise or lower syringe to increase or decrease feeding flow rate.
- 2.7 Flush the feeding tube with the prescribed amount of water when feeding is complete.
- 2.8 Disconnect syringe by rotating counter-clockwise and reconnect the cap by rotating cap clockwise (**do not overtighten**).

#### FEEDING ADMINISTRATION (continued)

3. If feeding continuously:
  - 3.1 Attach enteral giving set connector to the feeding port by rotating connector clockwise (**do not overtighten**).
  - 3.2 Fill feeding container with feeding formula.
  - 3.3 Allow for filling of feeding lumen. Inspect for leakage.  
**Note:** The device may need to be replaced if leakage is present.
  - 3.4 Connect enteral giving set to pump.
  - 3.5 Set desired flow rate and begin feeding.
  - 3.6 If feeding continuously for long periods of time, flush feeding tube with prescribed amount of water and rinse connectors (on device and giving set) with warm water and/or mild soap every 4-6 hours.
  - 3.7 Disconnect enteral giving set from pump and feeding port (by rotating counter-clockwise) when feeding is complete.
  - 3.8 Flush the feeding tube with prescribed amount of water, rinse connectors (on device and giving set) with warm water and/or mild soap as needed, and reconnect cap to the device by rotating cap clockwise (**do not overtighten**).

#### MEDICATION ADMINISTRATION

**Caution:** Only administer medication after confirmation of proper feeding tube placement and patency. Administration of medication should be performed according to instructions provided by physician.

**Caution:** Use medication in liquid form whenever possible. If only solid medication is available, seek advice from a physician to determine if it is safe to crush solid medication and mix with water.

**Warning:** Do not crush medication that has an enteric coating or attempt to mix medication with feeding formula.

**Warning:** This device is intended **FOR ENTERAL USE ONLY**.

**Note:** Verify the ENFit connector is connected to an enteral port only and not to an I.V. set prior to administering medication.

1. Remove cap from the "MED" port by rotating cap counter-clockwise.
2. Before administering medication, flush feeding tube with prescribed amount of water using an ENFit compliant syringe.
3. Medication may be administered by using an ENFit compliant syringe and connecting to the medication port of the device (by rotating syringe clockwise - **do not overtighten**).  
**Note:** If safe, solid medication may be pulverized into a fine powder and dissolved with water before being administered.
4. Disconnect syringe by rotating counter-clockwise.
5. After medication administration is complete, flush the route used for medication with prescribed amount of water and rinse connectors (on device and giving set) with warm water and/or mild soap as needed.
6. Reconnect cap by rotating clockwise (**do not overtighten**).

## DECOMPRESSION

1. Remove cap from either the "FEED" or "MED" port of the device by rotating cap counter-clockwise.
2. Perform decompression according to physician instructions.
3. After decompression is performed, flush device with prescribed amount of water.
4. Reconnect the cap to the device by rotating cap clockwise (**do not overtighten**).

## TUBE REMOVAL

**Caution:** Never attempt to remove feeding tube unless trained by physician or healthcare provider. Before removing feeding tube, verify that the feeding tube may be replaced at bedside.

1. Obtain equipment and supplies required for removal of device. Clean hands and use powder-free gloves before touching stoma site.
2. Gently rotate the feeding tube to ensure the tube moves freely and easily.  
**Warning:** Do not apply excessive force when removing device. If resistance is felt, apply water soluble lubricant and manipulate tube by pushing and rotating until tube moves freely and easily. Contact physician if not able to remove feeding tube.
3. Remove any contents inside the balloon using a syringe.
4. Pull device out of stoma while simultaneously applying counter-pressure to the patient's abdomen.
5. Clean stoma site with mild soap and water. Allow area to dry.
6. Dispose of device following the guidelines for biohazardous waste provided by the healthcare institution.
7. If a replacement device is required, follow instructions provided in "DEVICE PREPARATION" and "DEVICE PLACEMENT" instructions above.

**Caution:** Stoma tract may begin to close in the first hour following device removal.

## TUBE MAINTENANCE

1. The stoma site should be inspected and cleaned regularly. Cleaning should be performed with mild soap and water, remove moisture afterwards. Notify physician if signs of infection or irritation are observed.
2. Clean residual contents from the external area of the feeding tube with mild soap and water as needed; remove moisture afterwards.
3. Rotate device daily to allow for adequate air circulation. Consult with physician if device does not rotate freely.
4. In order to prevent feeding tube clogging, always flush feeding lumen with prescribed amount of water after feeding administration, every 4-6 hours during continuous feeding, before administering medication and after medication administration is complete.  
**Note:** The amount of water used for flushing will depend on patient's need, clinical condition and if patient is an infant or adult.
5. Assess the feeding tube daily for damage, clogging or abnormal discoloration.
6. Check balloon volume weekly, or as recommended by the physician, and compare to initial volume used to inflate balloon in order to determine balloon integrity.
7. A loss of volume over time may be compensated for by refilling the balloon. A large drop in fill volume may indicate the need for device replacement.
8. Replace device as or if necessary following the "DEVICE PREPARATION" and "DEVICE PLACEMENT" instructions in the section above.

## MRI SAFETY INFORMATION

This device is MR Conditional.

Non-clinical testing demonstrated that this device is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or 1.5 Tesla
- Maximum spatial gradient magnetic field of 1,500-Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode of operation for the MR system

### MRI-RELATED HEATING

Under the scan conditions above, this device is expected to produce a maximum temperature rise of less than 2.3°C after 15 minutes of continuous scanning.

### ARTIFACT INFORMATION

The maximum artifact size as seen on the gradient echo pulse sequence at 3-Tesla extends approximately 30-mm relative to the size of the shape of this device.

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Solutions Pty Ltd
GmbH