COOK

MEDICAL

CE
0088

EN
3

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Bruksanvisning

C _ T _ C P D _ R E V 3

TENCKHOFF CHRONIC PERITONEAL DIALYSIS CATHETERS AND SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Tenckhoff Chronic Peritoneal Dialysis Catheters are sideported silicone catheters with dual retention cuffs available in a range of lengths and French sizes and in either a straight or spiral tip configuration. Sets also contain an appropriately sized access needle, a wire guide and accessories for use in catheter placement.

INTENDED USE

Tenckhoff Chronic Peritoneal Dialysis Catheters are intended for chronic access to the peritoneal cavity.

MRI INFORMATION

Nonclinical testing has demonstrated that the Tenckhoff Chronic Peritoneal Dialysis Catheter is **MR Conditional**. It can be safely scanned after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient field of 1600 Gauss/cm or less
- Normal operating mode: Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence)

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

In nonclinical testing, the Tenckhoff Chronic Peritoneal Dialysis Catheter produced a maximum temperature rise of 1.7 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (Siemens Magnetom Trio, A Tim System, Software Numaris/4) at an MR system reported whole-body-averaged SAR of 2.1 W/kg (associated with a calorimetry measured whole-body-averaged value of 1.7 W/kg).

Image Artifact

The image artifact does not extend from the device, but the device lumen is obscured when scanned in nonclinical testing using: GRE sequence in a 3.0 Tesla Siemens Magnetom Trio, A Tim System (Software Numaris/4) MR system with the body coil.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic, interventional and surgical techniques. Standard techniques for placement of percutaneous catheters should be employed.
- This catheter is not intended for use in patients who are not suitable for peritoneal dialysis.

POTENTIAL ADVERSE EVENTS

- Bleeding
- Catheter migration
- Catheter obstruction
- Cuff erosion
- Cuff separation
- Dialysate leakage

- Dislodgement
- Fluid flow pain
- Fracture
- Hernia
- Ileus
- Infection (exit site and tunnel)
- Intestinal perforation
- Malposition
- Omental wrapping around catheter
- Peritonitis
- Sepsis
- Subcutaneous cuff extrusion

INSTRUCTIONS FOR USE

1. Select a placement site that avoids both the inferior epigastric artery and any abdominal scars, so as to minimize the risk of entering the bowel adherent to the anterior abdominal wall. Both midline and paramedian abdominal sites below the umbilicus and low lateral sites have been used successfully.
2. After prepping the skin as for any sterile surgical procedure, infiltrate local anesthetic.
3. Make a 1 cm skin incision, followed by minimal blunt dissection down through the subcutaneous tissue.
4. Fill a 20 ml syringe with sterile saline and attach it to the access needle. Advance the needle through the incision until resistance is encountered. **NOTE:** At this point, the needle tip should be less than 1 cm from entering the peritoneal cavity.
5. Instruct the patient to tense his or her abdominal muscles, then advance the needle tip into the peritoneal cavity while simultaneously depressing the plunger on the syringe. **NOTE:** Depressing the plunger on the syringe while advancing the needle should result in a jet of saline forcefully exiting the needle tip; this will make damage to any intraperitoneal soft tissues less likely.
6. Leaving the needle in place, remove the syringe and advance the wire guide through the needle and into the peritoneal cavity.
7. Remove the needle.
8. Introduce the Peel-Away® sheath-introducer assembly over the wire guide well into the peritoneal cavity.
9. Remove the wire guide and introducer, leaving the Peel-Away sheath in place.
10. Fill or at least flush the peritoneal dialysis catheter with a heparin solution. This will reduce the likelihood of fibrin formation in the catheter during the insertion procedure.
11. Insert the catheter through the sheath and into the peritoneal cavity, then peel away the sheath. Position the distal catheter cuff just above the linea alba (for a median insertion) or just above the rectus abdominus muscle (for a paramedian insertion). **NOTE:** Some blunt dissection may be necessary to create a space for the cuff.
12. Select an appropriate catheter exit site above or below, but not at, the belt line. The proximal catheter cuff should be positioned 1-2 cm from the point of exit through the skin, without dislodging the distal cuff.
13. After infiltrating sufficient local anesthetic, introduce a tunneling tool through the initial skin incision. Tunnel subcutaneously out to the selected site and puncture through the skin to create an exit.
14. Pressfit the peritoneal dialysis catheter to the other end of the tunneling tool; pull the catheter through the tunnel and out through the exit site. Some blunt dissection may again be necessary to create space for the proximal cuff. Note that the tunneling tool has been designed to produce a snug fit of the skin around the exit site so as to decrease the risk of leakage and exit site infection.
15. If necessary, one or two resorbable subcutaneous sutures may be placed before closing the initial surgical skin incision.
16. Test the catheter for inflow and outflow dynamics.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from packages, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2015
Bloomington, IN 47404 U.S.A.