COOK

MEDICAL

**EN
3**

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

1

2

FROVA INTUBATING INTRODUCER WITH RAPI-FIT® ADAPTERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Frova Intubating Introducer is an 8.0 French, 35 cm long radiopaque catheter introducer with centimeter markings and a blunt, curved tip that can be passed blindly into the trachea to facilitate intubation. Contains no natural rubber latex.

The device consists of the following components:

- Radiopaque catheter
- Rapi-Fit adapter – 15 mm connector
- Rapi-Fit adapter – Luer lock connector
- Stiffening stylet

INTENDED USE

The Frova Intubating Introducer is intended to facilitate endotracheal intubation in adult and pediatric patients where visualization of the glottis is inadequate.

The 8 French catheter is recommended for placement of an endotracheal tube with an inner diameter of 3 mm or larger.

When used for high-pressure oxygenation with a Luer Lock connector, the 8 French catheter is recommended for patients older than 1 month of age.

CONTRAINDICATIONS

Do not use this device if the epiglottis cannot be visualized when performing direct laryngoscopy (i.e., Grade IV Cormack & Lehane laryngoscopy classification).

WARNINGS

- Advancing the catheter beyond the carina may introduce additional procedural risks and trauma.
- Attention should be paid to the insertion depth of the Frova Intubating Introducer into the patient's airway and correct tracheal position of replacement endotracheal tube (ETT). Markers on the Frova Intubating Introducer refer to distance from the distal tip of catheter.
- Take care to avoid injuring the epiglottis, glottis, trachea, bronchi, or lung parenchyma, and to avoid perforating the sinus pyriformis while using this device.
- Ensure proper sizing of the Frova Intubating Introducer within the endotracheal tube.
- To avoid barotrauma, ensure that the tip of the Frova is always above the carina, preferably 2-3 cm.
- Use of the Rapi-Fit adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should be considered only if the patient has sufficient egression of the insufflated gas volume.
- If a high-pressure oxygen source is used for insufflation (e.g., jet ventilator), begin at a lower pressure (i.e., 5 psi) and work up gradually. Rising chest wall, pulse oximetry, and oral air flow should be carefully monitored.
- Ensure that the Rapi-Fit Adapter is securely connected to the Frova Intubating Introducer prior to oxygen delivery. Failure to properly secure the adapter to the introducer may result in hypoxia, hypoxemia, and serious adverse events.

PRECAUTIONS

- This product is intended for use by clinicians trained and experienced in airway management techniques.
- This product is not intended for intravascular use.
- This product should not be used in replacement of endotracheal tubes whose inner diameter is smaller than 3 mm.
- The use of capnography is highly recommended to confirm intubation.
- It is recommended that a sterile lubricant be applied to the Frova and new endotracheal tube prior to their introduction into the patient.
- This device is not recommended for definitive, long-term airway management.

POTENTIAL ADVERSE EVENTS

- Barotrauma
- Pneumothorax
- Hypoxia
- Perforation of the trachea, bronchi, or lung parenchyma
- Injury to the epiglottis
- Catheter migration
- Failed endotracheal tube placement

INSTRUCTIONS FOR USE

1. Lubricate an appropriately sized single lumen endotracheal tube and the distal tip of the catheter introducer.
2. If desired, prior to placement, the endotracheal tube may be preloaded onto the proximal portion of the catheter introducer.
3. Using laryngoscopy, introduce the tip of the introducer beyond the epiglottis and advance it in a straight line toward the glottis.
4. Advance the introducer into the trachea, approximately 2-3 cm. If resistance is encountered, **do not force the introducer**; instead, it should be gently rotated and advanced.
5. Remove the stiffening stylet and continue to advance the catheter introducer into the trachea for no more than 10 cm, the distance being determined by the size of the patient. **NOTE:** During introduction, you may feel tracheal rings.
6. Confirm catheter position using standard methods (e.g., capnography, breath sounds, x-ray).
7. Maintaining position of the introducer, advance the endotracheal tube into the trachea to the appropriate distance.
8. Maintaining position of the endotracheal tube, remove the introducer and laryngoscope.
9. Confirm correct position of the endotracheal tube using standard methods, (e.g., breath sounds, capnography, X-ray, and endotracheal tube device labeling).

NOTE: If necessary, oxygen may be given during the procedure, by utilizing the Rapi-Fit Adapters

Use of the Rapi-Fit® Adapter

Rapi-Fit adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient’s lungs may require more time.

1. To attach the Rapi-Fit Adapter, position the adapter on the Frova, then push the white collar forward and lock into position. **(Fig. 1)**
2. To remove the adapter, pull the white collar back to release, and then remove from the Frova. **(Fig. 2)**

CATHETER OXYGENATION

The Frova Intubation Catheter is designed for both positive airway pressure ventilation (Rapi-Fit Adapter - 15 mm connector) and jet ventilation (Rapi-Fit Adapter - Luer lock connector) for both adult and pediatric patients. In the table below, delivered minute volume, and measured average maximum airway pressure are given for jet ventilation in adult and pediatric patients with healthy lung tissue.

Catheter Oxygenation

Real Part Number (RPN)	Patient Subgroup and Age Range	Delivered Minute Volume ¹ (L/min)	Measured Airway Pressure ¹ (cm H ₂ O)	
			Mean	Maximum
C-CAE-8.0-35-FII-01	Infant > 1 month to 2 years	0.7	10.3	15.0
	Child > 2 years to 12 years	0.3	3.7	13.4
	Adolescent > 12 years through 21 years	0.2	3.4	4.0
	Adult > 21 years	0.8	4.0	5.2

¹ The test conditions were used in an active model. See *Testing Conditions* table for additional details.

The following testing conditions were used in an active model mode, tested with ASL 5000, Ingmar Medical, Ltd. for each patient subgroup.

Testing Conditions

Input pressure was set to 50 psi for both patient type models.

Patient type	Model Body Weight (kg)	Inspira-tory time (s)	Expira-tory time (s)	Breaths per minute	Resis-tance (cm H ₂ O/ L/s)	Lung Compli-ance (mL/cm H ₂ O)
Infant (> 1 month to 2 years)	10	0.46	1.54	30	25	10
Child (> 2 years to 12 years)	20	0.55	1.85	25	20	20
Adolescent (> 12 years to 21 years)	50	0.9	2.1	20	5	40
Adult (> 21 years)	80	1.0	4.0	12	3	100

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.