EN
3

**Biodesign®
Surgisis® Tissue Graft**

Instructions for Use

CS
4

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Bruksanvisning

PL
18

**Biodesign®
Wszczep tkanki Surgisis®**

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PT
19

**Biodesign®
Enxerto de pele Surgisis®**

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SV
21

**Biodesign®
Vävnadstransplantat Surgisis®**

Bruksanvisning

FP0005-4D

BIODESIGN® SURGISIS® TISSUE GRAFT

INTENDED USE

The Cook® Biodesign® Surgisis® Tissue Graft is used for implantation to reinforce soft tissue. The graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

SURGISIS® 1-LAYER TISSUE GRAFT This symbol means the following: Surgisis® 1-Layer Tissue Graft

SURGISIS® 4-LAYER TISSUE GRAFT This symbol means the following: Surgisis® 4-Layer Tissue Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

PRECAUTIONS

- **Do not resterilize.** Discard all open and unused portions.
- Graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Single-layer graft should not be used in applications requiring high strength.
- The graft may not have sufficient strength to support stresses encountered in some ventral hernias or large-area, body-wall repairs.
- Ensure that the graft is rehydrated prior to suturing or stapling.
- **Graft performance has not been evaluated with suture spacing greater than 2 mm.**
- **Ensure that all layers of the graft are secured when suturing or stapling.**

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials. If any of these conditions occur, the graft should be removed:

- Infection
- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

STORAGE

This graft should be stored in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide.

USE OF ANTIMICROBIALS

Because the graft is at times used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.¹ In these cases both antibiotic prophylaxis of the patient and antimicrobial soaking of the graft have been used. Typical flora can be expected to include a variety of aerobic and facultative anaerobic organisms, including, but not limited to, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Escherichia coli*. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.
- Antibacterial prophylaxis, if chosen, should be started prior to surgery and continued post-operatively.¹

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the graft. For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,²⁻⁴ while povidone iodine,⁵ bacitracin,^{3,5} polymyxin B,⁶ and vancomycin⁷ have all been reported to slow or inhibit wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with the graft with respect to graft placement.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede the institutional protocols or professional clinical judgment concerning patient care.

NOTE: Handle the graft using aseptic technique. Minimize contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
 - Sterile forceps
 - Rehydration fluid: at least 50 ml of room temperature sterile saline or sterile lactated Ringer's solution for each graft sheet
1. Using aseptic technique, remove the graft inner pouch from its outer bag, and place the inner pouch in the sterile field.
 2. Open the inner pouch carefully, and aseptically remove the graft sheet with the sterile forceps.
 3. Place the graft sheet into the sterile dish in the sterile field. (Multiple graft sheets may be rehydrated simultaneously in the same dish.)
 4. Add to the dish at least 50 mL of the rehydration fluid for each graft sheet.
 5. Allow a 1-layer graft sheet to rehydrate for at least three minutes. Allow a 4-layer graft sheet to rehydrate for at least ten minutes.
 6. Prepare the graft site using standard surgical techniques.
 7. Using aseptic technique, trim the graft sheet to fit the site, providing a small allowance for overlap. **NOTE:** An alternative method is to cut the graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the graft sheet prior to suturing or stapling it into place. See step 5.
 8. Using aseptic technique, transfer the graft sheet to the graft site and suture or staple into place, avoiding excess tension. **NOTE:** Surgical experience indicates that suturing or stapling graft sheets with close tissue approximation produces better outcomes.
 9. Complete the standard surgical procedure.
 10. Discard any unused portions of the graft sheet according to institutional guidelines for medical waste.

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.