**EN
3**

**Biodesign®
Surgisis® Anterior Pelvic Floor Graft
Surgisis® Posterior Pelvic Floor Graft**

Instructions for Use

**CS
4**

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BIODESIGN® SURGISIS® ANTERIOR AND POSTERIOR PELVIC FLOOR GRAFTS

INTENDED USE

Cook® Biodesign® Surgisis® Anterior and Posterior Pelvic Floor Grafts are used for soft tissue repair of pelvic floor defects such as: cystocele, rectocele, enterocele, sacrocolpopexy, and/or intra-operative bladder neck suspension. The graft is supplied sterile and intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

SURGISIS® ANTERIOR PELVIC FLOOR GRAFT This symbol means the following: Surgisis® Anterior Pelvic Floor Graft

SURGISIS® POSTERIOR PELVIC FLOOR GRAFT This symbol means the following: Surgisis® Posterior Pelvic Floor Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

PRECAUTIONS

- **Do not resterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling and possible damage has occurred, or if the graft exceeds the labeled expiration date.
- Ensure that the graft is rehydrated prior to cutting, suturing, stapling, or loading the graft laparoscopically.
- Ensure that all layers of the graft are secured when suturing or stapling.
- Suturing or stapling more than one graft together may decrease graft performance.
- The graft performance has not been evaluated with suture spacing greater than 2 mm.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the graft.
- Users should exercise good surgical practice for the management of contaminated or infected wounds.
- The potential for infection of the graft material following implantation may be reduced by the use of prophylactic antibiotics.
- Users should exercise caution to avoid vessel, bowel, and bladder perforation.
- **IMPORTANT:** Users should counsel patients on abstaining from heavy lifting, strenuous exercise, and sexual intercourse for a period of four (4) to six (6) weeks after surgery.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials:

- Bleeding
- Infection
- Adhesions
- Sterile effusion
- Chronic inflammation
- Allergic reaction
- Delayed or failed incorporation of graft

If conditions of infection, inflammation, or allergic reaction cannot be resolved, consider removal of the graft.

STORAGE

This graft should be stored in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide.

USE OF ANTIMICROBIALS

Because the graft is at times used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.¹ In these cases both antibiotic prophylaxis of the patient and antimicrobial soaking of the graft have been used. Typical flora can be expected to include a variety of aerobic and facultative anaerobic organisms, including, but not limited to, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Escherichia coli*. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.

- Antibacterial prophylaxis, if chosen, should be started prior to surgery and continued post-operatively.¹

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the graft. For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,²⁻⁴ while povidone iodine,⁵ bacitracin,^{3,5} polymyxin B,⁶ and vancomycin⁷ have all been reported to slow or inhibit wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with the graft with respect to graft placement.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Always handle the graft using aseptic technique. Minimize contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
 - Sterile forceps
 - Rehydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution.
1. Using aseptic techniques, remove the graft inner pouch from its outer pouch, and place the inner pouch in the sterile field.
 2. Open the inner pouch carefully, and aseptically remove the graft with sterile forceps.
 3. Place the graft into the sterile dish in the sterile field. Multiple grafts may be rehydrated simultaneously in the same dish.
 4. Add to the dish at least 50 mL of the rehydration fluid for each graft.
 5. Allow the graft to rehydrate for at least ten (10) minutes.
 6. Prepare the surgical site using standard surgical techniques.
 7. Using aseptic technique, trim the graft to the desired shape and size, providing an allowance for overlap great enough to prevent excess tension from being placed on the suture line.

IMPORTANT: If the graft is cut too small for the defect, excess tension may be placed on the suture line as the graft incorporates into the surrounding tissues. This can result in recurrence of the original tissue defect or development of a defect in the surrounding tissues.

8. Using aseptic technique, transfer the graft to the surgical site and suture or staple into place, avoiding excess tension.

NOTE: Surgical experience indicates that suturing or stapling the graft with close tissue approximation produces better outcomes. Fundamental surgical principles suggest suture spacing approximately equal to suture bite depth.

9. Complete the standard surgical procedure.
10. Discard any unused portions following institution guidelines for disposal of medical waste.

NOTE: Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.