EN
3

**Biodesign®
Surgisis® Fistula Plug**

Instructions for Use

CS
7

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
BIODESIGN® SURGISIS® FISTULA PLUG

INTENDED USE

The Cook® Biodesign® Surgisis® Fistula Plug is for implantation to reinforce soft tissue for repair of recto-vaginal or anorectal fistulas. The plug is supplied sterile and is intended for one time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

SURGISIS® FISTULA PLUG This symbol means the following: Surgisis® Fistula Plug

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

- This plug is derived from a porcine source and should not be used for patients sensitive to porcine materials.
- **Not for vascular use.**

PRECAUTIONS

- **Do not resterilize.** Discard all open and unused portions.
- Plug is sterile if the package is dry, unopened and undamaged.
- Do not use if the package seal is broken.
- Discard plug if mishandling has caused possible damage or contamination or if the plug is past its expiration date.
- Do not implant the plug in a grossly infected or abscessed fistula tract.
- **In fistula cases involving evidence of acute inflammation, purulence, or excessive drainage, a seton should be used to allow the tract to mature and stabilize for six to eight weeks before placing the plug.**
- Ensure that the plug is rehydrated prior to placement, cutting, or suturing.
- Placement of the plug in tracts less than 1 cm in length can result in incomplete incorporation and/or expulsion of the plug.

GENERAL

- Users should be familiar with surgical technique for recto-vaginal or anorectal fistula repair.
- Users should exercise good surgical practice for the management of clean-contaminated, contaminated or infected fields.
- The potential for infection of the graft material following implantation may be reduced by the use of prophylactic antibiotics and cleaning of the fistula tract. **(See Use of Antimicrobials)**
- The plug should always be introduced through the rectal (primary) fistula opening.
- The plug should be drawn into the fistula tract completely, until the button is set against the rectal wall.
- The external/vaginal (secondary) fistula opening should remain open to allow drainage to occur.
- **IMPORTANT:** A rectal exam is recommended at the end of the 4th post-operative week to confirm separation of the non-resorbing button from the rectal wall. It is also recommended that users counsel patients on abstaining from strenuous physical activity for at least 2 weeks following anorectal fistula repair, or 6 weeks following recto-vaginal fistula repair. Refer to the section on Post-Operative Care.

POTENTIAL COMPLICATIONS

Complications that can occur with the plug include, but are not limited to:

- Inflammation
- Induration
- Migration
- Extrusion
- Seroma formation
- Infection
- Abscess
- Fistula recurrence
- Delayed or failed incorporation of the plug

If any of the following conditions occur and cannot be resolved, plug removal should be considered:

- Infection
- Abscess
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

STORAGE

This plug should be stored in a clean, dry location at room temperature.

STERILIZATION

This plug has been sterilized with ethylene oxide and should not be resterilized.

USE OF ANTIMICROBIALS

Because the plug is used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.¹⁻³ Both mechanical bowel cleansing and antibiotic prophylaxis of the patient have been used successfully, and the U.S. Centers for Disease Control recommends bowel cleansing and administration of prophylactic enteric and intravenous antimicrobial agents before elective colorectal operations.³ (See reference for specific dose, timing, and drug choices.)

Typical gastrointestinal flora can be expected to include a variety of aerobic and anaerobic organisms. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.⁴
- Mechanically prepare the colon by appropriate use of enemas and cathartic agents.
- A dose of prophylactic antimicrobial agent should be given intravenously and timed such that a bactericidal concentration of the drug is established in serum and tissues when the plug is implanted.
- Maintain therapeutic levels of the agent in serum and tissues throughout the operation.³

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the repairing tissue.⁵⁻⁷ For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,⁶ while povidone iodine,⁸ bacitracin,^{5,8} polymyxin B,⁹ and vancomycin¹⁰ have all been reported to slow or inhibit wound healing. Careful consideration is required before using any antimicrobial or antiseptic (topical or systemic) that has not been proven compatible with surgical implantation and wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with plug placement.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede the institutional protocols or professional clinical judgment concerning patient care.

NOTE: Handle plug using aseptic technique. Minimize contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
- Rehydration fluid: at least 200 ml of room temperature sterile saline or sterile lactated Ringer's solution
- Suitable resorbable suture, such as: 0 chromic, 2-0 or 0 coated polyglycolic acid suture (coated PGA)
- A short 14ga or 16ga catheter
- Saline or hydrogen peroxide for flushing

PREPARATION

1. Select the proper plug size according to Procedural step 4.
2. Remove the packaging containing the plug from the box.
3. Remove the inner pouch containing the plug from the outer package using aseptic technique. Place the inner pouch into the sterile field.
4. Using sterile gloved hands carefully open the inner pouch and remove the card containing the plug. Place the plug and card combination into the sterile dish in the sterile field.
5. Add enough rehydration fluid to the dish to fully submerge the plug. Allow the plug to rehydrate, fully submerged, until the desired handling characteristics are achieved but no longer than 5 minutes.
6. Fold the card at the indicated line and carefully remove the hydrated plug from the card.
7. Attach a suitable resorbable suture (approximately 30 cm in length) through the tail (narrow end) of the plug for pulling it into the fistula tract. Obtain at least a 4mm suture bite into the end of plug.
8. Prepare the patient and surgical site using standard surgical techniques appropriate for the fistula repair.

NOTE: The recommended practice for preoperative bowel preparation in elective colorectal surgery includes mechanical bowel cleansing through the use of enemas and cathartic agents, and administration of prophylactic antimicrobial agents, oral or intravenous. Insufficient cleansing or inadequate antibacterial prophylaxis can predispose the patient to infections.¹⁻³ (**See Use of Antimicrobials**)

PROCEDURAL

1. Perform under local, regional, or general anesthesia.
2. If a seton is not already in place, identify the rectal (primary) fistula opening by inserting a sterile probe into the secondary opening and navigating it through the fistula tract. Alternatively, injection of appropriate sterile fluids (saline or hydrogen peroxide) into the secondary opening of the fistula tract, and identifying the site of emergence at the primary opening, can also be used to assist in identifying the location of the rectal fistula opening.


NOTE: Failure to locate the rectal opening may lead to persistence of the fistula. If the rectal opening cannot be reliably identified, an alternative method of treatment should be considered.


3. Gently clean/debride the fistula tract with a probe or small brush to remove non-vascularized tissue from the fistula tract. Flush the tract as thoroughly as possible with sterile saline, or comparable solution. Do not enlarge the tract.
4. Use a sterile probe, or appropriate measuring instrument, to approximate the diameter of the fistula tract opening in the rectum. Select the proper size of plug according to the table below. Refer to the section on Preparation.

Rectal Fistula Diameter (mm)	Fistula Plug size (cm)
>1 to 2	0.2 (2 mm)
>2 to 4	0.4 (4 mm)
>4 to 7	0.7 (7 mm)


- 5a. If a seton is not in place, insert a fistula probe, seton, or suitable instrument through the secondary opening and exiting via the rectal (primary) opening. Attach the suture connected to the plug to the instrument at the primary opening end.
- 5b. If a seton is in place, cut the seton and attach the suture to the seton at the rectal (primary) opening end.

6. Draw the suture tie attached to the plug into the primary fistula opening and through the fistula tract until the plug is drawn completely into the fistula tract and the plug button is against the rectal wall.


Recto-vaginal fistula repair

Anorectal fistula repair

IMPORTANT: The rectal opening is the high-pressure zone of the fistula, as well as the site of ingress of fecal debris. The flange end of the plug must therefore be securely snug against the rectal opening to prevent ingress of fecal debris. In addition, the higher pressures within the rectum and anal canal assist in maintaining the plug in the fistula tract by simple mechanical force.

7. When the plug is properly positioned, grip the tether used to draw the plug into the fistula tract and apply gentle tension to it. This will hold the plug in place during fixation of the plug on the rectal side. Refer to the illustration below.

IMPORTANT: The button of the plug should closely approximate the rectal wall in order to minimize contact of the plug material with bowel contents.



8. Holes are provided in the button to allow fixation to the rectal wall. Using a 2-0 coated PGA or comparable suture, affix the button to the rectal wall by stitching through the large button holes and directly through the rectal wall. Obtain an adequate bite of tissue with the suture and secure the button to the tissue by tying the two suture ends together.

9. If the secondary opening is external (skin level), place a sterile dressing over it.


OPTIONAL FIXATION FOR RECTO-VAGINAL FISTULAS ONLY

10. Using a 2-0 coated PGA or comparable suture, place a minimum of two stitches through the middle of the plug at the vaginal fistula opening. The suture should be placed along the length of the plug. Pass the suture through the plug and into tissue obtaining adequate bites of both plug and tissue with the suture. Allow a 2 mm suture spacing between each stitch.

NOTE: Do not close the secondary opening completely when suturing the plug to tissue. This is to allow continued drainage of the fistula tract. Complete obstruction of the fistula opening may result in accumulation of fluid, infection, or abscess.

11. When the plug is properly secured, trim away and discard any remaining portion of the plug that is not implanted within the fistula tract. Leave at least 2 mm of material between the end of the trimmed plug and suture to provide the appropriate suture retention strength.

12. Discard any unused portions of the plug following standard technique for disposal of medical waste.

POST-OPERATIVE CARE

To provide the best environment for tissue integration into the plug, patient activity should be minimized. Provide patients with a list of post-procedural care recommendations. The following patient guidelines should be considered.

1. No strenuous physical activity beyond a gentle walk for at least 2 weeks following anorectal fistula repair, or 6 weeks following recto-vaginal fistula repair.
2. No lifting items over 10 lbs (5kg) for at least two weeks following anorectal fistula repair, or at least 6 weeks following recto-vaginal fistula repair.
3. Use of stool softeners for 2 weeks after surgery.
4. Observing a liquid diet for the first 2 days after surgery.
5. Beginning a high fiber diet 2 days after surgery.
6. Abstaining from sexual intercourse and other forms of vaginal/rectal insertion (e.g. tampons) for at least 2 weeks following anorectal fistula repair, or 6 weeks following recto-vaginal fistula repair.

REFERENCES

1. Nichols RL, Smith JW, Garcia RY, et al. *Current practices of preoperative bowel preparation among North American colorectal surgeons*. Clin Infect Dis 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. *A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations*. J Med Dent Sci 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. *Guideline for Prevention of Surgical Site Infection*, 1999. Centers for Disease Control and Prevention (CDC).
4. Aldridge KE, et al. *Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of Bacteroides fragilis group, Prevotella, Fusobacterium, Porphyromonas, and Peptostreptococcus species*. Antimicrob Agents Chemother 2001;45(4):1238-43.
5. Petroutsos G, et al. *Antibiotics and corneal epithelial wound healing*. Arch Ophthalmol 1983;101(11):1775-8.
6. Bang K, et al. *Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice*. Scand J Plast Reconstr Surg Hand Surg 1998;32(2):129-33.
7. Nelson JD, et al. *Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics*. Curr Eye Res 1990;9(3):277-85.
8. Kjolseth D, et al. *Comparison of the effects of commonly used wound agents on epithelialization and neovascularization*. J Am Coll Surg 1994;179(3):305-12.
9. Nakamura M, et al. *Effects of antimicrobials on corneal epithelial migration*. Curr Eye Res 1993;12(8):733-40.
10. Petroutsos G, Guimaraes R, and Pouliquen Y. *The effect of concentrated antibiotics on the rabbit's corneal epithelium*. Int Ophthalmol 1984;7(2):65-9.

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© COOK BIOTECH 2010
West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.