**EN
3**

**Biodesign®
Surgisis® Staple Line Reinforcement**

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
Bruksanvisning

FP0059-01C

BIODESIGN® SURGISIS® STAPLE LINE REINFORCEMENT

Staple Line Reinforcement with Self-Adhesive

INTENDED USE

The Biodesign® Surgisis® Staple Line Reinforcement is intended for use as a prosthesis for the surgical repair of soft tissue deficiencies using surgical staplers. The Staple Line Reinforcement may be used for buttressing and reinforcing staple lines during lung resection (e.g., wedge resection, blebectomy, lobectomy, bullectomy, bronchial resection, segmentectomy, pneumonectomy/pneumectomy, pneumoreduction) and other incisions and excisions of the lung and bronchus. The Staple Line Reinforcement can be used for the reinforcement of the gastric staple line during bariatric surgical procedures of gastric bypass, and for reinforcement of staple lines during small bowel, mesentery, colon and colorectal procedures. The Staple Line Reinforcement may be used with anastomotic staplers or with non-anastomotic staplers.

For Use With This symbol means the following: The surgical stapler for which this device is indicated for use.

QTY This symbol means the following: Quantity of devices per box.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

SURGISIS® STAPLE LINE REINFORCEMENT This symbol means the following: Surgisis® Staple Line Reinforcement

This product is intended for use by trained medical professionals.

This device is supplied sterile and is intended for one-time use.

CONTRAINDICATIONS

This Staple Line Reinforcement is derived from a porcine source and should not be used in patients with known sensitivity to porcine materials.

This Staple Line Reinforcement contains fructose and should not be used in patients with known sensitivity to fructose.

PRECAUTIONS

- Do not resterilize. Discard open and unused portions.
- Staple Line Reinforcement is sterile if the package is dry, unopened, and undamaged. Do not use if package seal is broken.
- Discard Staple Line Reinforcement if mishandling has caused possible damage or contamination or if Staple Line Reinforcement is past its expiration date.

NOTE: Staple Line Reinforcement may get damaged if foam applicator is folded.

- Follow the manufacturer's instructions regarding the use of the stapler, including precautions regarding use with staple reinforcement materials.
- Use stapler loads that allow for combined thickness of tissue and the Staple Line Reinforcement.*
- Use only specified staplers with appropriate sized Staple Line Reinforcement.
- Ensure that prior to firing the surgical stapler the entire active portions of stapler jaws are covered by the reinforcement device.
- Fire the stapler within 20 minutes after applying the Staple Line Reinforcement. Otherwise, the device may be difficult to dislodge from the stapler jaws after firing.
- Avoid dropping the tray-device assembly from the pouch onto the sterile field.
- Avoid using the tail to remove the Staple Line Reinforcement from the tray as it could lead to premature detachment from the foam applicator.
- Avoid trying to re-adhere the Staple Line Reinforcement to the stapler if it comes off.

- Avoid excessive manipulation of tissue with a stapler loaded with the Staple Line Reinforcement.
- To prevent loss of adhesive, avoid dipping the Staple Line Reinforcement in water/saline. Dip the stapler jaws instead.

**The average total nominal thickness (anvil+cartridge side including adhesive) of the Staple Line Reinforcement is indicated on the label as dimension "C".*

POTENTIAL COMPLICATIONS

- Infection
- Acute or chronic inflammation (Initial application of graft materials may be associated with transient, mild, localized inflammation)
- Allergic Reaction

STORAGE

This Staple Line Reinforcement should be stored in a clean, dry location at room temperature.

STERILIZATION

The Staple Line Reinforcement has been sterilized using E-beam radiation.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

REQUIRED MATERIALS

- Compatible surgical stapler
- Warm sterile water or saline
- Sterile scissors
- Sterile basin
- Sterile gauze or sponge

NOTE: Always handle the Staple Line Reinforcement using aseptic technique, minimizing contact with latex gloves.


1. Aseptically remove the tray from inside the pouch using sterile gloved hands and place into the sterile field.
2. Wet the jaws of the surgical stapler or re-load by compressing the jaws across a piece of gauze (or a sponge) that has been soaked with sterile water or saline.


NOTE: Liquid can be warmed up to 38°C (100°F) to enhance adhesion.

An alternate wetting method is to quickly dip the jaws of the surgical stapler or re-load in sterile water or saline and allow any excess water to drip off.


3. Pull the Staple Line Reinforcement out of the tray by the finger grips on wings of the foam applicator. **Do Not Pull out of the Tray by the Tail.**
4. Place the applicator inside **the wet jaws of the surgical stapler**. Visually or palpably align the parallel edges of the applicator with the edges of the stapler jaws. Gently insert the applicator to the back of the stapler jaws until the applicator can not go any further. Notches in the applicator provide for a stopping point to indicate how far back the applicator should be positioned in the stapler jaws.


5. Close the stapler on the applicator. **Keep tightly closed for at least 10 seconds.**

6. Grip the tail of the foam applicator and detach it to allow the device to be released from the foam applicator. **Discard the tail of the applicator.**

7. Open the stapler jaws, **remove and discard the foam applicator**. Inspect device alignment to ensure complete coverage of both cartridge and anvil stapler jaws. Inspect the adhered Staple Line Reinforcement to ensure that none of the foam applicator parts are stuck. If required, using sterile scissors, cut off any excess length of the Staple Line Reinforcement strip that protrudes from the tip of the stapler jaws, being certain not to dislodge the strip from the stapler jaws.


NOTE: If the stapler is a two-part stapler, use sterile scissors to cut the staple line reinforcement at the hinged fold between the two stapler jaws.

8. Move the stapler or re-load to the operative site and continue with the procedure. Use caution to avoid displacing the Staple Line Reinforcement from the stapler jaws. **Use the stapler within 20 minutes after applying the Staple Line Reinforcement to the stapler.**

NOTE: Before applying a new staple line reinforcement, between fires, wipe the anvil side of the stapler jaws with a dry piece of gauze or sterile cloth.

9. Discard any unused portions of the Staple Line Reinforcement following standard technique for disposal of medical waste.

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SURGISIS® STAPLE LINE REINFORCEMENT

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.