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
Fig. 16

Fig. 17

# Zenith® TX2® TAA Endovascular Graft with Pro-Form™ and the Z-Trak® Plus Introduction System

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgement concerning patient care.

**CAUTION:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

## 1. DEVICE DESCRIPTION

### 1.1 Zenith TX2 TAA Endovascular Graft with Pro-Form and the Z-Trak Plus Introduction System

The Zenith TX2 TAA Endovascular Graft with Pro-Form is a two-piece cylindrical endovascular graft consisting of proximal and distal components. The proximal components can be either non-tapered or tapered. The stent grafts are constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z stents with braided polyester and monofilament polypropylene suture. (Fig. 1) The Zenith TX2 TAA Endovascular Graft with Pro-Form is fully stented to provide stability and the expansile force to open the lumen of the graft during deployment. Additionally, the Cook-Z stents provide the attachment and seal of the graft to the vessel wall.

For added fixation, the covered stent at the proximal end of the proximal component contains barbs placed at a 2 mm stagger, which protrude through the graft material. In addition, the bare stent at the distal end of the distal component also contains barbs. To facilitate fluoroscopic visualization of the stent graft, four radiopaque gold markers are positioned on each end of the proximal and distal components. These markers are placed in a circumferential orientation within 1 mm of the most proximal aspect of the graft material and within 1 mm of the most distal aspect of the graft material.

The Zenith TX2 TAA Endovascular Graft with Pro-Form is shipped preloaded onto either a 20 French or 22 French Z-Trak Plus Introduction System. (Fig. 2) It has a sequential deployment method with built-in features to provide continuous control of the endovascular graft throughout the deployment procedure. The Z-Trak Plus Introduction System is designed for precise positioning before deployment of the proximal and/or distal components. The proximal component uses a single trigger-wire release mechanism. The distal component uses a dual trigger-wire release mechanism. The trigger-wires secure the endovascular graft onto the delivery system until released by the physician. (Fig. 3) All delivery systems feature Flexor® introducer sheaths, which are designed to resist kinking and are hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and thoracic aorta.

### 1.2 Zenith TX2 TAA Endovascular Graft with Pro-Form Ancillary Components

Ancillary endovascular components (proximal and distal body extensions) are available. (Fig. 1) The Zenith TX2 TAA Endovascular Graft with Pro-Form Ancillary Components are cylindrical components constructed from the same polyester fabric, self-expanding stainless steel Cook-Z stents, and polypropylene suture used in constructing the principal graft components. At the distal and proximal graft margins, the stents are attached to the inner surface. Elsewhere the stents are sutured on the external surface. The proximal extension contains proximal attachment barbs and the distal extension does not have barbs. Both the proximal and distal main body extensions can be used to provide additional length to their respective portions of the endovascular graft. Additionally, the distal main body extension can be used to increase the overlap length between components.

#### 1.2.1 Zenith TX2 TAA Endovascular Graft with Pro-Form Proximal Extensions

The Zenith TX2 TAA Endovascular Graft with Pro-Form Proximal Extension is deployed from either a 20 French or 22 French Z-Trak Plus Introduction System. (Fig. 2) A single trigger-wire release mechanism locks the endovascular graft onto the delivery system until released by the physician. All systems are compatible with a .035 inch wire guide.

The covered stent at the proximal end of the proximal extension contains barbs placed at a 2 mm stagger, which protrude through the graft material. To facilitate fluoroscopic visualization of the proximal extension, four radiopaque markers are positioned on the ends of the graft in a circumferential orientation within 1 mm of the most proximal and distal aspects of the graft material.

#### 1.2.2 Zenith TX2 TAA Endovascular Graft with Pro-Form Distal Extensions

The Zenith TX2 TAA Endovascular Graft with Pro-Form Distal Extension is deployed from either a 20 French or 22 French Z-Trak Plus Introduction System. (Fig. 2) A single trigger-wire release mechanism locks the endovascular graft onto the delivery system until released by the physician. All systems are compatible with a .035 inch wire guide.

To facilitate fluoroscopic visualization of the distal extension, four radiopaque markers are positioned on the ends of the graft in a circumferential orientation within 1 mm of the most proximal and distal aspects of the graft material.

## 2. INTENDED USE

The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System is intended for treatment of patients with atherosclerotic aneurysms, symptomatic acute or chronic dissections, contained ruptures, growing aneurysms and/or resulting in distal ischemia, in the descending thoracic aorta having vascular morphology suitable for endovascular repair (Fig. 5), including:

- Adequate iliac/femoral access compatible with the required introduction systems,
- Radius of curvature greater than 35 mm along the entire length of aorta intended to be treated,
- Non-aneurysmal aortic segments (fixation sites) proximal and distal to the aneurysm:
  - with a length of at least 20 mm
  - with a diameter measured outer wall to outer wall of no greater than 38 mm and no less than 20 mm, and
  - with an angle less than 45 degrees.

## 3. CONTRAINDICATIONS

The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), polypropylene, nitinol, or gold.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

## 4. WARNINGS AND PRECAUTIONS

### 4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 10.1, Physician Training**.
- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft and/or leaks may be required to undergo secondary endovascular interventions or surgical procedures.

#### 4.2 Patient Selection, Treatment and Follow-Up

- The Zenith TX2 TAA Endovascular Graft with Pro-Form is designed to treat aortic neck diameters no smaller than 20 mm and no larger than 38 mm. The Zenith TX2 TAA Endovascular Graft with Pro-Form is designed to treat proximal aortic necks (distal to either the left subclavian or left common carotid artery) of at least 20 mm in length. Additional proximal aortic neck length may be gained by covering the left subclavian artery (with or without discretionary transposition) when necessary to optimize device fixation and maximize aortic neck length. Distal aortic neck length of at least 20 mm proximal to the celiac axis is required. These sizing measurements are critical to the performance of the endovascular repair.
- Key anatomic elements that may affect successful exclusion of the aneurysm include a radius of curvature <35 mm; localized aortic neck angulation >45 degrees; short proximal or distal fixation sites (<20 mm); an inverted funnel shape at the proximal fixation site or a funnel shape at the distal fixation site (greater than 10% change in diameter over 20 mm of fixation site length); and circumferential thrombus and/or calcification at the arterial fixation sites. In the presence of anatomical limitations, a longer neck length may be required to obtain adequate sealing and fixation. Irregular calcification and/or plaque may compromise the attachment and sealing at the fixation sites. Necks exhibiting these key anatomic elements may be more conducive to graft migration or endoleak.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Careful evaluation of vessel size, anatomy and disease state is required to assure successful sheath introduction and subsequent withdrawal, as vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude introduction of the endovascular graft and/or may increase the risk of embolization. A vascular conduit technique may be needed to achieve access in some patients.
- The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System is not recommended for patients who cannot tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System is not recommended for patients whose weight or size would compromise or prevent the necessary imaging requirements.
- Graft implantation may increase the risk of paraplegia or paraparesis where graft exclusion covers the origins of dominant spinal cord or intercostal arteries.
- The safety and effectiveness of the Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System has not been evaluated in the following patient populations:
  - aortobronchial and aorto-esophageal fistulas
  - aortitis or inflammatory aneurysms
  - diagnosed or suspected congenital degenerative collagen disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
  - females that are pregnant, breast-feeding, or planning on becoming pregnant within 24 months
  - leaking, pending rupture or ruptured aneurysm
  - less than 18 years of age
  - mycotic aneurysms
  - pseudoaneurysms resulting from previous graft placement
  - systemic infection (e.g., sepsis)
  - traumatic aortic injury
- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3 Pre-Procedure Measurement Techniques and Imaging**.
- If occlusion of the left subclavian artery ostium is required to obtain adequate neck length for fixation and sealing, transposition or bypass of the left subclavian artery may be warranted.
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intra-operative flexibility to achieve optimal procedural outcomes.

#### 4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in suboptimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith TX2 TAA Endovascular Graft with Pro-Form. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning of the image intensifier (C-arm) so that it is perpendicular to the aneurysm neck, typically 45-75 degrees left anterior oblique (LAO) for the arch.

**Diameter:** A contrast-enhanced spiral CTA is strongly recommended for aortic diameter measurements. Diameter measurements should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. The spiral CTA scan must include the great vessels through the femoral heads at an axial slice thickness of 3 mm or less.

**Length:** Clinical experience indicates that 3-D CTA reconstruction is the strongly recommended imaging modality to accurately assess proximal and distal neck lengths for the Zenith TX2 TAA Endovascular Graft with Pro-Form. These reconstructions should be performed in sagittal, coronal and varying oblique views depending upon individual patient anatomy. If 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.**
- The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System is not recommended in patients unable to undergo, or who will not be compliant with, the necessary pre-operative and post-operative imaging and implantation studies as described in **Section 12, IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.**
- After endovascular graft placement, patients should be regularly monitored for endoleak flow, aneurysm growth, or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) chest radiographs to examine device integrity (separation between components, stent fracture, device position, and/or barb separation); and 2) contrast and non-contrast CT to examine aneurysm changes, endoleak flow, patency, tortuosity, device position and progressive disease. If renal complications or other factors preclude the use of image contrast media, chest radiographs and non-contrast CT may be used in combination with transesophageal echocardiography (for endoleak assessment) to provide similar, although suboptimal, information.

#### 4.4 Device Selection

- The recommended amount of overlap between devices is 3-4 stents. However, the proximal sealing stent of the proximal component or distal sealing stent of the distal component should not be overlapped, as doing so may cause malapposition to the vessel wall. The minimum required amount of overlap between devices is 2 stents (~50 mm) – less than 2 stents may result in endoleak (with or without component separation). Device lengths should be selected accordingly.
- Strict adherence to the Zenith TX2 TAA Endovascular Graft with Pro-Form IFU sizing guide is strongly recommended when selecting the appropriate device size (**Tables 10.1 and 10.2**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding, or compression.

#### 4.5 Implant Procedure

(Refer to **Section 10, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith TX2 TAA Endovascular Graft with Pro-Form in the neck and to assure appropriate apposition to the aortic wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith TX2 TAA Endovascular Graft with Pro-Form.
- To avoid twisting the endovascular graft, never rotate the delivery system during the procedure. Allow the device to conform naturally to the curves and tortuosity of the vessels.
- Do not continue advancing the wire guide or any portion of the delivery system if resistance is felt. Stop and assess the cause of resistance; vessel, catheter, or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis, or calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith TX2 TAA Endovascular Graft with Pro-Form in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant arch or mesenteric arteries (exception may be the left subclavian artery) with the endoprosthesis. Vessel occlusion may occur. If a left subclavian artery is to be covered with the device, the clinician should be aware of the possibility of compromise to cerebral and upper-limb circulation and collateral circulation to the spinal cord.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Repositioning the stent graft distally after partial deployment of the covered proximal stent may result in damage to the stent graft and/or vessel injury.
- During sheath withdrawal, the proximal barbs are exposed and are in contact with the vessel wall. At this stage it may be possible to advance the device, but retraction may cause aortic wall damage.
- Landing the proximal and distal ends of the device in parallel aortic neck segments without acute angulation (>45°) or circumferential thrombus/calcification is important to ensuring fixation and seal.
- Landing the proximal or distal ends of the device in an aortic neck segment with a diameter that differs from that to which the graft was sized initially may potentially result in inadequate sizing (<10% or >25%) and therefore migration, endoleak, aneurysm growth, or increased risk of thrombosis.
- Inaccurate placement and/or incomplete sealing of the Zenith TX2 TAA Endovascular Graft with Pro-Form within the vessel may result in increased risk of endoleak, migration, or inadvertent occlusion of the left subclavian, left common carotid, and/or celiac arteries.
- Inadequate fixation of the Zenith TX2 TAA Endovascular Graft with Pro-Form may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician-preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be used.
- To activate the hydrophilic coating on the outside of the sheath, the surface must be wiped with 4X4 gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Always use fluoroscopy for guidance, delivery, and observation of the Zenith TX2 TAA Endovascular Graft with Pro-Form within the vasculature.
- The use of the Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast medium used during the procedure and to observe preventative methods of treatment to decrease renal compromise (e.g., adequate hydration).
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- The Zenith TX2 TAA Endovascular Graft with Pro-Form incorporates a covered proximal stent (on the proximal component) with fixation barbs and an uncovered distal stent (on the distal component) with fixation barbs. Exercise extreme caution when manipulating interventional and angiographic devices in the region of the covered proximal stent and uncovered distal stent. Do not place the barbs and the uncovered distal stent in an area with dissection.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus or plaque, which can cause distal or cerebral embolization, or cause rupture of the aneurysm.

- Avoid damaging the graft or disturbing graft position after placement in the event reinstrumentation (secondary intervention) of the graft is necessary.
- Care should be taken not to advance the sheath while the stent graft is still within it. Advancing the sheath at this stage may cause the barbs to perforate the introducer sheath.
- To avoid entangling any catheters left *in situ*, rotate the delivery system during withdrawal.

#### 4.6 Molding Balloon Use (Optional)

- Do not inflate the balloon in the aorta outside of the graft, as doing so may cause damage to the aorta. Use the balloon in accordance with its labeling.
- Use care in inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the aorta.
- Confirm complete deflation of balloon prior to repositioning.
- For added hemostasis, the Captor® Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

#### 4.7 MRI Safety and Compatibility

**Non-clinical testing has demonstrated that the Zenith TX2 TAA Endovascular Graft with Pro-Form is MR Conditional. It can be scanned safely under the following conditions:**

##### 1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Spatial gradient field of 450 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

In non-clinical testing, the Zenith TX2 TAA Endovascular Graft with Pro-Form produced a temperature rise of less than 1.4 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom MR scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

##### 3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Spatial gradient field of 720 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith TX2 TAA Endovascular Graft with Pro-Form produced a temperature rise of less than 1.9 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of MR scanning in a 3.0 Tesla, Excite, GE Electric Healthcare MR scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in non-clinical testing using the sequence: Fast spin echo in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal, upper extremity, and head and neck region, depending on distance from the device to the area of interest.

Clinical information is available on six patients who received MRI scans during the course of the clinical trial. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been approximately 3,000 patients implanted with Zenith TAA Endovascular Grafts worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

## 5. POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Aorto-bronchial fistula
- Aorto-esophageal fistula
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma, or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, tamponade, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Compartment Syndrome
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration and/or separation; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Femoral neuropathy
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, urinary incontinence, hematuria, infection)
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula, lymphocele)

- Local or systemic neurologic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis/spinal cord shock, paralysis)
- Occlusion of device or native vessel
- Pulmonary Embolism
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vascular spasm or vascular trauma (e.g., ilio-femoral vessel dissection, bleeding, rupture, death)
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

## 6. PATIENT SELECTION AND TREATMENT

(See Section 4, WARNINGS AND PRECAUTIONS)

### 6.1 Individualization of Treatment

Cook recommends that the Zenith TX2 TAA Endovascular Graft with Pro-Form component diameters be selected as described in **Tables 10.1** and **10.2**. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intra-operative flexibility to achieve optimal procedural outcomes. The risks and benefits should be carefully considered for each patient before use of the Zenith TX2 TAA Endovascular Graft with Pro-Form. Additional considerations for patient selection include but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith TX2 TAA Endovascular Graft with Pro-Form
- Ability to tolerate general, regional, or local anesthesia
- Ability and willingness to undergo and comply with the required follow-up
- Ilio-femoral access vessel size and morphology (thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 20 French to 22 French vascular introducer sheath
- Vascular morphology suitable for endovascular repair, including:
  - Adequate iliac/femoral access compatible with the required introduction systems,
  - Radius of curvature greater than 35 mm along the entire length of aorta intended to be treated.
- Non-aneurysmal aortic segments (fixation sites) proximal and distal to the aneurysm:
  - with a length of at least 20 mm,
  - with a diameter measured outer wall to outer wall of no greater than 38 mm and no less than 20 mm, and
  - with an angle less than 45 degrees.

The final treatment decision is at the discretion of the physician and patient.

## 7. PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and open surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment to and compliance with post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.**
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of thoracic aortic aneurysms. At a minimum, annual imaging and adherence to routine post-operative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise every patient that it is important to seek prompt medical attention if he/she experiences signs of graft occlusion, aneurysm enlargement or rupture. Symptoms of graft occlusion include, but may not be limited to, pulse-less legs, pain, ischemia of intestines, and cold extremities. Aneurysm rupture may be asymptomatic, but usually presents as back or chest pain, persistent cough, dizziness, fainting, rapid heartbeat, or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant. Men who undergo endovascular or open surgical repair may experience impotence.

The physician should complete the Patient Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI). For additional information, please refer to the Zenith TX2 TAA Endovascular Graft with Pro-Form Patient Guide.

## 8. HOW SUPPLIED

- The Zenith TX2 TAA Endovascular Graft with Pro-Form is supplied ETO-sterilized and preloaded in peel-open packages. The device is designed for single-use only. Attempts to reprocess, re-sterilize and/or re-use the device may lead to device failure and/or transmission of disease.
- Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to COOK.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- The device is loaded into a 20 French or 22 French Flexor® Introducer Sheath. Its surface is treated with a hydrophilic coating that, when hydrated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a 4X4 gauze pad soaked in saline solution.
- Do not use after the expiration date printed on the label.
- Store in a dark, cool, dry place.

## 9. CLINICAL USE INFORMATION

### 9.1 Physician Training

**CAUTION: Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.**

**CAUTION: The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System should only be used by physicians and teams trained in vascular interventional techniques (endovascular and surgical) and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System are outlined below:**

#### Patient selection:

- Knowledge of the natural history of thoracic aortic aneurysms (TAA), and co-morbidities associated with TAA repair.
- Knowledge of radiographic image interpretation, patient selection, device selection, planning and sizing.

#### A multidisciplinary team that has combined procedural experience with:

- Femoral and brachial cutdown, arteriotomy, and repair or conduit technique
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

### 9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to COOK.

Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

### 9.3 Materials Required

(Not included in two-piece modular system)

- A selection of Zenith TX2 TAA Endovascular Graft with Pro-Form proximal and distal ancillary components in diameters compatible with the two-piece system are available.
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Power injector
- Syringe
- Heparinized saline solution
- Sterile 4X4 gauze pads

### 9.4 Materials Recommended

(Not included in two-piece modular system)

The following products are recommended for implantation of any component in the Zenith product line. For information on these products, refer to the individual product's Suggested Instructions For Use.

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
  - Cook Amplatz Ultra Stiff Wire Guides (AUS)
  - Cook Lunderquist Extra Stiff Wire Guides (LESDC)
- .035 inch (0.89 mm) standard wire guide; for example:
  - Cook .035 inch wire guides
  - Cook .035 inch Bentson Wire Guide
  - Cook Nimble® Wire Guides
- Molding Balloons; for example:
  - Cook CODA® Balloon Catheter
- Introducer sets; for example:
  - Cook Check-Flo® Introducer Sets
- Sizing catheter; for example:
  - Cook Aurous® Centimeter Sizing Catheters
- Angiographic radiopaque marker catheters; for example:
  - Cook Beacon® Tip Angiographic Catheters
  - Cook Beacon® Tip Royal Flush Catheters
- Entry needles; for example:
  - Cook single wall entry needles

### 9.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and **not** the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow. In order to ensure accurate diameter measurements for the purpose of graft sizing, particularly when in curved segments of the aorta, measuring the aortic diameter using 3D reconstructed views perpendicular to the aortic centerline of flow may be important.

**Table 9.1 Main Body Graft Diameter Sizing Guide\***

Intended Aortic Vessel Diameter <sup>1,2</sup> (mm)	Graft Diameter <sup>3</sup> (mm)	Overall Length of Proximal Component (mm)	Overall Length of Distal Component (mm)	Overall length of Proximal Tapered Component (mm)	Introducer Sheath (Fr)	Introducer Sheath Outer Diameter (OD) (mm)
20	22	115	NA	NA	20	7.7
21	24	115	NA	NA	20	7.7
22/23	26	115	NA	NA	20	7.7
24	28	120/140/200	147/207	N/A	20	7.7
25	30	120/140/200	147/207	N/A	20	7.7
26	30	120/140/200	147/207	N/A	20	7.7
27	30	120/140/200	147/207	N/A	20	7.7
28	32	120/140/200	147/207	160/200	20	7.7
29	32	120/140/200	147/207	160/200	20	7.7
30	34	127/152/202	136/186	157/197	20	7.7
31	36	127/152/202	136/186	157/197	22	8.6
32	36	127/152/202	136/186	157/197	22	8.6
33	38	127/152/202	136/186	152/202	22	8.6
34	38	127/152/202	136/186	152/202	22	8.6
35	40	135/162/216	144/198	158/208	22	8.6
36	40	135/162/216	144/198	158/208	22	8.6
37	42	135/162/216	144/198	158/208	22	8.6
38	42	135/162/216	144/198	158/208	22	8.6

\*All dimensions are nominal.

<sup>1</sup>Maximum diameter along the fixation site, measured outer wall to outer wall.

<sup>2</sup>Round measured aortic diameter to nearest mm.

<sup>3</sup>Additional considerations may affect choice of diameter.

**Table 9.2 Proximal and Distal Extension Graft Diameter Sizing Guide\***

Intended Aortic Vessel Diameter <sup>1,2</sup> (mm)	Graft Diameter <sup>3</sup> (mm)	Overall Length of Component (mm)	Introducer Sheath (Fr)	Introducer Sheath Outer Diameter (OD) (mm)
20	22	80	20	7.7
21	24	80	20	7.7
22/23	26	80	20	7.7
24	28	80	20	7.7
25	30	80	20	7.7
26	30	80	20	7.7
27	30	80	20	7.7
28	32	80	20	7.7
29	32	80	20	7.7
30	34	77	20	7.7
31	36	77	22	8.6
32	36	77	22	8.6
33	38	77	22	8.6
34	38	77	22	8.6
35	40	81	22	8.6
36	40	81	22	8.6
37	42	81	22	8.6
38	42	81	22	8.6

\*All dimensions are nominal.

<sup>1</sup>Maximum diameter along the fixation site, measured outer wall to outer wall.

<sup>2</sup>Round measured aortic diameter to nearest mm.

<sup>3</sup>Additional considerations may affect choice of diameter.

## 10. DIRECTIONS FOR USE

### Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit technique may be required.
- Proximal and distal aortic neck lengths should be a minimum of 20 mm.
- Aortic neck diameters measured outer wall to outer wall between 20-38 mm.
- Proximal neck diameter of at least 4 mm larger than the distal neck diameter requires the use of a proximal tapered component.
- Measurements to be taken during the pre-treatment assessment are described in **Fig. 5** and **6**.

## Proximal and Distal Component Overlap

A minimum overlap of 2 stents (~50 mm) is required; a 3-4 stent (~75-100 mm) overlap is recommended, however, the proximal sealing stent of the proximal component or distal sealing stent of the distal component should not be overlapped.

The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

## General Use Information

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System. The Zenith TX2 TAA Endovascular Graft with Pro-Form with the Z-Trak Plus Introduction System is compatible with .035 inch diameter wire guides.

## Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the delivery system(s)
2. Angulation of aorta, aneurysm and iliac arteries
3. Quality of the proximal and distal fixation sites
4. Diameters of proximal and distal fixation sites and distal iliac arteries
5. Length of proximal and distal fixation sites

## Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation, and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Expose femoral artery using standard surgical technique.
4. Establish adequate proximal and distal vascular control of femoral artery.

## 10.1 Zenith TX2 TAA Endovascular Graft with Pro-Form System Component

### Preparation/Flush/Placement – Proximal and Distal Components

1. Remove yellow-hubbed shipping stylet. Remove cannula protector tube. Remove Peel-Away<sup>®</sup> sheath from back of valve assembly. (**Fig. 7**)
2. Elevate distal tip of system and flush through the hemostatic valve until fluid emerges from the tip of the introducer sheath. (**Fig. 8**) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.  
**NOTE:** Ensure that the side-arm adapter is securely connected to the side of the valve body.  
**NOTE:** Graft flushing solution of heparinized saline is often used.
3. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal sideports and dilator tip. (**Fig. 9**)
4. Soak 4X4 gauze pads in saline solution and use to wipe the Flexor Introducer Sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

### 10.1.1 Placement of Proximal Component

1. Puncture the selected artery using standard technique with an 18 gage access needle. Upon vessel entry, insert:
  - Wire guide – standard .035 inch, 260 cm, 15 mm J tip or Bentson wire guide
  - Appropriate size sheath (e.g., 5.0 French)
  - Pigtail flush catheter (often radiopaque-banded sizing catheters; i.e., Cook Centimeter Sizing CSC-20 catheter)
2. Perform angiography at the appropriate level. Using radiopaque markers, adjust position as necessary and repeat angiography.
3. Ensure graft system has been flushed and primed with heparinized saline (appropriate flush solution), and all air has been removed.
4. Give systemic heparin. Flush all catheters and wet all wire guides with heparinized saline. This should be repeated following each exchange.
5. Replace the standard wire guide with a stiff .035 inch, 300 cm –LESDC wire guide and advance through the catheter and up to the aortic arch.
6. Remove pigtail flush catheter and sheath.  
**NOTE:** At this stage, the second femoral artery can be accessed for angiographic catheter placement. Alternatively, a brachial approach may be considered.
7. Introduce the freshly hydrated delivery system over the wire guide and advance until the desired graft position is reached.  
**CAUTION: To avoid twisting the endovascular graft, never rotate the delivery system during the procedure. Allow the device to conform naturally to the curves and tortuosity of the vessels.**  
**NOTE:** The dilator tip will soften at body temperature.  
**NOTE:** To facilitate introduction of the wire guide into the delivery system, it may be necessary to slightly straighten the delivery system dilator tip.
8. Verify wire guide position in the aortic arch. Ensure correct graft position.  
**CAUTION: Care should be taken not to advance the sheath while the stent graft is still within it. Advancing the sheath at this stage may cause the barbs to perforate the introducer sheath.**
9. Ensure that the Captor Hemostatic Valve on the Flexor Introducer Sheath is turned to the open position. (**Fig. 10**)

10. Stabilize the grey positioner (delivery system shaft) and withdraw the sheath until the graft is fully expanded and the valve assembly docks with the control handle.  
**CAUTION: As the sheath is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.**  
**CAUTION: During sheath withdrawal, the proximal barbs are exposed and are in contact with the vessel wall. At this stage it may be possible to advance the device, but retraction may cause aortic wall damage.**  
**NOTE:** If extreme difficulty is encountered when attempting to withdraw the sheath, place the device in a less tortuous position which enables the sheath to be retracted. Very carefully withdraw the sheath until it begins to retract, and stop. Move device back to original position and continue deployment.
11. Verify graft position and adjust it forward, if necessary. Recheck graft position with angiography.  
**NOTE:** If an angiographic catheter is placed parallel to the stent graft, use this to perform position angiography.
12. Loosen the safety lock from the green trigger-wire release mechanism. Withdraw the trigger-wire slowly until the proximal end of the graft opens. (Fig. 11) Withdraw the trigger-wire completely to release the distal attachment to the introducer.  
**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.
13. Remove the introduction system, leaving the wire guide in the graft.

### 10.1.2 Placement of Distal Component

1. If an angiographic catheter in the femoral artery is in use, it should be withdrawn to a position to demonstrate the aortic anatomy where the distal component is to be deployed.
2. Introduce the freshly hydrated delivery system over the wire guide until the desired graft position is reached, with a recommended 3-4 stent overlap (75-100 mm), but no less than a 2 stent overlap (50 mm) with the proximal component. Do not overlap proximal and distal sealing stents.  
**NOTE:** To facilitate introduction of the wire guide into the delivery system, it may be necessary to slightly straighten the delivery system dilator tip.
3. Check the position by angiography and adjust if necessary.
4. Ensure that the Captor Hemostatic Valve on the Flexor Introducer Sheath is turned counter-clockwise to the open position. (Fig. 10)
5. Stabilize the grey positioner (delivery system shaft) and begin withdrawing the sheath.  
**CAUTION: As the sheath or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.**  
**NOTE:** If extreme difficulty is encountered when attempting to withdraw the sheath, place the device in a less tortuous position which enables the sheath to be retracted. Very carefully withdraw the sheath until it just begins to retract, and stop instantly. Move back to original position and continue deployment.
6. Withdraw the sheath until the graft is fully expanded. Continue sheath withdrawal until the valve assembly docks with the control handle.
7. Release the distal attachment by first unscrewing the trigger-wire safety lock on the white trigger-wire release mechanism (labeled number "1"). (Figs. 12 and 13)
8. Unscrew and remove the safety lock on the telescoping handle (labeled number "2"). (Figs. 14 and 15)
9. Stabilize the delivery system and slide the telescoping handle together with the grey tube and the outer sheath in a distal direction until the distal attachment stent is released. The telescoping handle should be retracted as far as it will travel distally until it locks automatically into position. (Fig. 16)
10. Loosen the safety lock from the green trigger-wire release mechanism. Withdraw the trigger-wire slowly until the proximal end of the graft opens, then withdraw and remove the trigger-wire and release mechanism (labeled number "3"). (Fig. 17)  
**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.
11. Remove the inner introduction system entirely, leaving the sheath and wire guide in the graft.
12. Close the Captor Hemostatic Valve on the Flexor Introducer Sheath by turning it in a clockwise direction until it stops.  
**CAUTION: To avoid entangling any catheters left *in situ*, rotate the delivery system during withdrawal.**

### 10.1.3 Main Body Molding Balloon Insertion (Optional)

1. Prepare molding balloon as follows and/or per the manufacturer's instructions.
  - Flush wire lumen with heparinized saline.
  - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning it counter-clockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve of the main body introduction system to the level of the proximal fixation/seal site. Maintain proper sheath positioning.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.  
**CAUTION: Do not inflate balloon in the aorta outside of graft.**
5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the proximal covered stent, starting proximally and working in the distal direction.  
**CAUTION: Confirm complete deflation of balloon prior to repositioning.**
6. If applicable, withdraw the molding balloon to the proximal component/distal component overlap and expand.
7. Withdraw the molding balloon to the distal covered stent and expand.
8. Open the Captor Hemostatic Valve, remove the molding balloon and replace it with an angiographic catheter to perform completion angiograms.
9. Tighten the Captor Hemostatic Valve around the angiographic catheter with gentle pressure by turning it clockwise.
10. Remove or replace all stiff wire guides to allow aorta to resume its natural position.

## Final Angiogram

1. Position angiographic catheter just above the level of the endovascular graft. Perform angiography to verify correct positioning. Verify patency of arch vessels and celiac plexus.
2. Confirm that there are no endoleaks or kinks, and verify position of proximal and distal gold radiopaque markers. Remove the sheaths, wires and catheters.  
**NOTE:** If endoleaks or other problems are observed, refer to **Section 11.2, Ancillary Devices**.
3. Repair vessels and close in standard surgical fashion.

## 10.2 Ancillary Devices

### General Use Information

Inaccuracies in device size selection or placement, changes or anomalies in patient anatomy, or procedural complications can require placement of additional endovascular grafts and extensions. Regardless of the device placed, the basic procedure(s) will be similar to the maneuvers required and described previously in this document. It is vital to maintain wire guide access.

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith TX2 TAA Endovascular Graft with Pro-Form ancillary devices.

The Zenith TX2 TAA Endovascular Graft with Pro-Form ancillary devices with the Z-Trak Plus Introduction Systems are compatible with .035 inch diameter wire guides.

### 10.2.1 Proximal Extensions

Proximal extensions are used for extending the proximal body of an *in situ* endovascular graft.

Proximal Extension Preparation/Flush

1. Remove yellow-hubbed shipping stylet. Remove cannula protector tube. Remove Peel-Away sheath from back of valve assembly. (**Fig. 7**)
2. Elevate distal tip of system and flush through the hemostatic valve until fluid emerges from the tip of the introducer sheath. (**Fig. 8**) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.  
**NOTE:** Ensure that the side-arm adapter is securely connected to the side of the valve body.  
**NOTE:** Graft flushing solution of heparinized saline is often used.
3. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal sideports and dilator tip. (**Fig. 9**)
4. Soak 4X4 gauze pads with saline and use to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

### Placement of the Proximal Extension

1. Puncture the selected artery using standard technique with an 18 gage access needle. Upon vessel entry, insert:
  - Wire guide—standard .035 inch, 260 cm, 15 mm J tip or Bentson wire guide
  - Appropriate size sheath (e.g., 5.0 French)
  - Pigtail flush catheter (often radiopaque-banded sizing catheters; i.e., Cook Centimeter Sizing CSC-20 catheter)
2. Perform angiography at the appropriate level. Using radiopaque markers, adjust position as necessary and repeat angiography.
3. Ensure introduction system has been primed with heparinized saline, and all air has been removed.
4. Give systemic heparin. Flush all catheters and wire guides with heparinized saline. This should be repeated following each exchange.
5. Replace the standard wire guide with a stiff .035 inch, 300 cm—LESDC wire guide and advance through the catheter and up to the aortic arch.
6. Remove pigtail flush catheter and sheath.  
**NOTE:** At this stage, the second femoral artery can be accessed for flush catheter placement. Alternatively, a brachial approach may be considered.
7. Introduce the freshly hydrated delivery system over the wire guide and advance until the desired graft position is reached. Ensure there is a minimum overlap of 2 stents.  
**CAUTION: To avoid twisting the endovascular graft, never rotate the delivery system during the procedure. Allow the device to conform naturally to the curves and tortuosity of the vessels.**  
**NOTE:** The dilator tip softens at body temperature.  
**NOTE:** To facilitate introduction of the wire guide into the delivery system, it may be necessary to slightly straighten the delivery system dilator tip.  
**NOTE:** The proximal extension contains barbs which should not be placed within other graft components.
8. Verify wire guide position in the aortic arch. Ensure correct graft position.  
**CAUTION: Care should be taken not to advance the sheath while the stent graft is still within it, advancing the sheath at this stage may cause the barbs to perforate the introducer sheath.**
9. Ensure that the Captor Hemostatic Valve on the Flexor Introducer Sheath is turned counter-clockwise to the open position.
10. Stabilize the grey positioner (delivery system shaft) and withdraw the sheath until the graft is fully expanded and the valve assembly docks with the control handle.  
**CAUTION: As the sheath or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.**  
**CAUTION: During sheath withdrawal, the proximal barbs are exposed and are in contact with the vessel wall. At this stage it may be possible to advance the device, but retraction may cause aortic wall damage.**  
**NOTE:** If extreme difficulty is encountered when attempting to withdraw the sheath, place the device in a less tortuous position which enables the sheath to be retracted. Very carefully withdraw the sheath until it just begins to retract, and stop instantly. Move back to original position and continue deployment.

11. Verify graft position and adjust it forward, if necessary. Recheck graft position with angiography.  
**NOTE:** If an angiographic catheter is placed parallel to the stent graft, use this to perform position angiography.
12. Loosen the safety lock from the green trigger-wire release mechanism. Withdraw the trigger-wire slowly until the proximal end of the graft opens. (**Fig. 11**) Withdrawing the trigger-wire completely will also release the distal attachment to the introducer.  
**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.
13. Remove the inner introduction system entirely, leaving the sheath and wire guide in the graft.  
**CAUTION: To avoid entangling any catheters left *in situ*, rotate the delivery system during withdrawal.**
14. Close the Captor Hemostatic Valve on the Flexor Introducer Sheath by turning it in a clockwise direction until it stops.

#### Proximal Extension Molding Balloon Insertion (Optional)

1. Prepare molding balloon as follows and/or per the manufacturer's instructions.
  - Flush wire lumen with heparinized saline
  - Remove all air from balloon
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning it counter-clockwise.
3. Advance the molding balloon over the wire guide and through the Captor Hemostatic Valve of the introduction system to the level of the proximal fixation/seal site. Maintain proper sheath positioning.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.  
**CAUTION: Do not inflate balloon in the aorta outside of graft.**
5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the proximal covered stent, starting proximally and working in the distal direction.  
**CAUTION: Confirm complete deflation of balloon prior to repositioning.**
6. Withdraw the molding balloon to the proximal extension/proximal component overlap and expand.
7. Open the Captor Hemostatic Valve, remove the molding balloon and replace it with an angiographic catheter to perform completion angiograms.
8. Tighten the Captor Hemostatic Valve around the angiographic catheter with gentle pressure by turning it clockwise.
9. Remove or replace all stiff wire guides to allow aorta to resume its natural position.

#### Final Angiogram

1. Position angiographic catheter just above the level of the endovascular graft. Perform angiography to verify correct positioning. Verify patency of arch vessels.
2. Confirm there are no endoleaks or kinks, and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.
3. Repair vessels and close in standard surgical fashion.

#### 10.2.2 Distal Extensions

Distal extensions are used for extending the distal end of an *in situ* endovascular graft or increasing the length of overlap between graft components.

##### Distal Extension Preparation/Flush

1. Remove yellow-hubbed shipping stylet. Remove cannula protector tube. Remove Peel-Away sheath from back of valve assembly. (**Fig. 7**)
2. Elevate distal tip of system and flush through the hemostatic valve until fluid emerges from the tip of the introducer sheath. (**Fig. 8**) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.  
**NOTE:** Ensure that the side-arm adapter is securely connected to the side of the valve body.  
**NOTE:** Graft flushing solution of heparinized saline is often used.
3. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal sideports and dilator tip. (**Fig. 9**)
4. Soak 4X4 gauze pads with saline and use to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

##### Placement of the Distal Extension

1. Puncture the selected artery using standard technique with an 18 gage access needle. Upon vessel entry, insert:
  - Wire guide—standard .035 inch, 260 cm, 15 mm J tip or Bentson wire guide
  - Appropriate size sheath (e.g., 5.0 French)
  - Pigtail flush catheter (often radiopaque-banded sizing catheters; i.e., Cook Centimeter Sizing CSC-20 catheter)
2. Perform angiography at the appropriate level. Using radiopaque markers, adjust position as necessary and repeat angiography.
3. Ensure graft system has been primed with heparinized saline, and all air has been removed.
4. Give systemic heparin. Flush all catheters and wire guides with heparinized saline. This should be repeated following each exchange.
5. Replace the standard wire guide with a stiff .035 inch, 300 cm—LESDC wire guide and advance through the catheter and up to the aortic arch.
6. Remove pigtail flush catheter and sheath.  
**NOTE:** At this stage, the second femoral artery can be accessed for flush catheter placement. Alternatively, a brachial approach may be considered.

7. Introduce the freshly hydrated delivery system over the wire guide and advance until the desired graft position is reached. Ensure there is a minimum overlap of two stents (plus the distal uncovered stent).

**CAUTION: To avoid twisting the endovascular graft, never rotate the delivery system during the procedure. Allow the device to conform naturally to the curves and tortuosity of the vessels.**

**NOTE:** The dilator tip softens at body temperature.

**NOTE:** To facilitate introduction of the wire guide into the delivery system, it may be necessary to slightly straighten the delivery system dilator tip.

8. Verify wire guide position in the aortic arch. Ensure correct graft position.
9. Ensure that the Captor Hemostatic Valve on the Flexor Introducer Sheath is turned counter-clockwise to the open position.
10. Stabilize the grey positioner (delivery system shaft) and withdraw the sheath until the graft is fully expanded and the valve assembly docks with the control handle.

**CAUTION: As the sheath or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.**

**NOTE:** If extreme difficulty is encountered when attempting to withdraw the sheath, place the device in a less tortuous position which enables the sheath to be retracted. Very carefully withdraw the sheath until it just begins to retract, and stop instantly. Move back to original position and continue deployment.

11. Verify graft position and adjust it forward, if necessary. Recheck graft position with angiography.  
**NOTE:** If an angiographic catheter is placed parallel to the stent graft, use this to perform position angiography.
12. Loosen the safety lock from the green trigger-wire release mechanism. Withdraw the trigger-wire slowly until the proximal end of the graft opens. (**Fig. 11**) Withdraw the trigger-wire completely to release the distal attachment to the introducer.  
**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.
13. Remove the inner introduction system entirely, leaving the sheath and wire guide in the graft.  
**CAUTION: To avoid entangling any catheters left *in situ*, rotate the delivery system during withdrawal.**
14. Close the Captor<sup>®</sup> Hemostatic Valve on the Flexor<sup>®</sup> Introducer Sheath by turning it in a clockwise direction until it stops.

#### Distal Extension Molding Balloon Insertion (Optional)

1. Prepare molding balloon as follows and/or per the manufacturer's instructions.
  - Flush wire lumen with heparinized saline.
  - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning it counter-clockwise.
3. Advance the molding balloon over the wire guide and through the Captor Hemostatic Valve of the introduction system to the level of the distal component/distal extension overlap. Maintain proper sheath positioning.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.  
**CAUTION: Do not inflate balloon in the aorta outside of graft.**
5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the overlap, starting proximally and working in the distal direction.  
**CAUTION: Confirm complete deflation of balloon prior to repositioning.**
6. Withdraw the molding balloon to the distal covered stent and expand.
7. Loosen the Captor Hemostatic Valve, remove the molding balloon and replace it with an angiographic catheter to perform completion angiograms.
8. Tighten the Captor Hemostatic Valve around the angiographic catheter with gentle pressure by turning it clockwise.
9. Remove or replace all stiff wire guides to allow aorta to resume its natural position.

#### Final Angiogram

1. Position angiographic catheter just above the level of the endovascular graft. Perform angiography to verify correct positioning. Verify patency of arch vessels.
2. Confirm there are no endoleaks or kinks, and verify position of proximal and distal gold radiopaque markers. Remove the sheaths, wires and catheters.
3. Repair vessels and close in standard surgical fashion.

## 11. IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

### 11.1 General

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive additional follow-up. Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of thoracic aortic aneurysms.
- Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in **Table 12.1**. This schedule continues to be the minimum recommendation for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.

- Annual imaging follow-up should include chest radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, chest radiographs and non-contrast CT may be used in combination with a transesophageal echocardiography for assessment of endoleak.
- The combination of contrast and non-contrast CT imaging provides information on device migration, aneurysm diameter depth change, endoleak, patency, tortuosity, progressive disease, fixation length, and other morphological changes.
- The chest radiographs provide information on device integrity (separation between components, stent fracture, and barb separation) and device migration.

Table 11.1 lists the minimum requirements for imaging follow-up for patients with the Zenith TX2 TAA Endovascular Graft with Pro-Form. Patients requiring enhanced follow-up should have interim evaluations.

**Table 11.1 Recommended Imaging Schedule for Endograft Patients**

	Angiogram	CT (contrast and non-contrast)	Chest Radiographs
Pre-procedure		X <sup>1</sup>	
Procedural	X		
1 month		X <sup>2</sup>	X
6 month		X <sup>2</sup>	X
12 month (annually thereafter)		X <sup>2</sup>	X

<sup>1</sup>Imaging should be performed within 6 months before the procedure.

<sup>2</sup>If Type I or III endoleak, prompt intervention and additional follow-up post-intervention recommended, **See Section 12.5, Additional Surveillance and Treatment.**

### 11.2 Contrast and Non-Contrast CT Recommendations

- Film sets should include all sequential images at lowest possible slice thickness ( $\leq 3$  mm). DO NOT perform large slice thickness ( $>3$  mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14" x 17" sheets if film is used.
- Both non-contrast and contrast runs are required, with matching or corresponding table positions.
- Pre-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. **Table 11.2** lists examples of acceptable imaging protocols.

**Table 11.2 Acceptable Imaging Protocols**

	Non-contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral CT or high performance MDCT capable of $>40$ seconds	Spiral CT or high performance MDCT capable of $>40$ seconds
Injection volume	n/a	Per Institutional Protocol
Injection rate	n/a	$>2.5$ cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: Smart Prep, C.A.R.E. or equivalent
Coverage-start	Neck	Subclavian aorta
Coverage-finish	Diaphragm	Profunda femoris origin
Collimation	$<3$ mm	$<3$ mm
Reconstruction	2.5 mm throughout-soft algorithm	2.5 mm throughout-soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

### 11.3 Chest Radiographs

The following views are required:

- Two films: supine-frontal (AP) and cross-table lateral
- Record the table-to-film distance and use the same distance at each subsequent examination
- Ensure entire device is captured on each single image format lengthwise
- The middle photocell should be used for all views to ensure adequate penetration of the mediastinum

If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length, including components) using 2-4X magnification visual aid.

### 11.4 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the Zenith TX2 TAA Endovascular Graft with Pro-Form is MR Conditional. It can be scanned safely under the following conditions:

#### 1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Spatial gradient field of 450 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith TX2 TAA Endovascular Graft with Pro-Form produced a temperature rise of less than 1.4 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom MR scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

### **3.0 Tesla Systems:**

- Static magnetic field of 3.0 Tesla
- Spatial gradient field of 720 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith TX2 TAA Endovascular Graft with Pro-Form produced a temperature rise of less than 1.9 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of MR scanning in a 3.0 Tesla, Excite, GE Electric Healthcare MR scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal, upper extremity, and head and neck region, depending on distance from the device to the area of interest.

Clinical information is available on six patients who received MRI scans during the course of the clinical trial. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been approximately 3,000 patients implanted with Zenith TAA Endovascular Grafts worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

### **11.5 Additional Surveillance and Treatment**

(Refer to **Section 4, WARNINGS AND PRECAUTIONS**)

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, >5 mm of maximum aneurysm diameter depth (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy, and the patient's personal choices. Patients should be counseled that subsequent reinterventions, including catheter-based and open surgical conversion, are possible following endograft placement.

## **12. REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2010