Rhapsody[®] H-30[®]

HOLMIUM LASER SYSTEM

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Table of Contents

Chapter 1: Introduction

About Cook Medical.....	1-4
Device Description.....	1-5 - 1-7
System Specifications	1-8
Labels and Symbols.....	1-9

Chapter 2: Indications

Indications	1-11
Indications for Use	1-11
Contraindications.....	1-11
Warnings and Precautions - General.....	1-12 - 1-13
Warnings and Precautions - Gynecological	1-14
Warnings and Precautions - Gastroenterological and Pulmonary	1-14

Chapter 3: Safety, Warnings, and Caution

General Warnings and Caution	1-16
General Safety Information.....	1-16
Staff Training	1-17 - 1-16
Working Area	1-17
Repair and Modifications of the Device	1-17
Eye and Skin Exposure.....	1-17 - 1-18
Hazards	1-18 - 1-19
Safety Indicators	1-20
Interference with Other Devices.....	1-20

Chapter 4: Installation

Installation	1-22
Transportation.....	1-22
Unpacking the Laser System.....	1-22 - 1-23
Cleaning.....	1-23
Contents	1-23

Ambient Operating Conditions.....	1-23
Filling the Cooling System with Water.....	1-24
Initial System Power-On.....	1-25 - 1-26
Verifying Alignment and Calibration.....	1-26 - 1-27
Electromagnetic Emission and Interference	1-27 - 1-29

Chapter 5: Operation

System Power-On.....	1-31
Checking the Optical Fibers	1-31
Laser Screens	1-32
Main Screen.....	1-32 - 1-33
Options Screen.....	1-34 - 1-35
Parameter Determination	1-35
Setting the Pulse Width.....	1-35
Adjusting the Energy per Pulse.....	1-35 - 1-36
Laser Parameters	1-36 - 1-37
Turning the System Off	1-37
Protection Against Unauthorized Use	1-37
Error Messages	1-37 - 1-38

Chapter 6: Maintenance

General Maintenance.....	1-40
Cleaning.....	1-40
Cooling Water Level.....	1-40
Power Cord	1-40
Labels.....	1-41
Replacing the Blast Shield.....	1-41 - 1-42
Troubleshooting	1-43
Error Messages	1-43 - 1-44

Chapter 7: Warranty Information

Warranty Statement	1-46
--------------------------	------

Chapter 1

Introduction

About Cook Medical

Choosing Cook Medical means choosing quality and innovation without compromise. To attain the best results with the Rhapsody H-30 Holmium Laser System and avoid the risk of dangerous faults, please thoroughly read this user manual before using the laser system.

Contact your sales representative if you have any questions regarding the operation or performance of the Rhapsody H-30 Holmium Laser System.

Cook Incorporated

750 Daniels Way
Bloomington, IN 47404 USA
Phone: 800.457.4500
www.cookmedical.com

Cook Capital Equipment Service and Repair Center

Phone: 855.207.7214
capitalservice@cookmedical.com

Device Description

Name	Laser Source
Rhapsody® H-30® Holmium Laser System	Ho:YAG

Front of Laser

Figure 1

Top Rear of Laser


Figure 2

Bottom Rear of Laser


Figure 3

Device Connections

The Rhapsody H-30 Holmium Laser System contains one laser source, emitting at 2100 nm. The laser is delivered through an optical fiber (holmium laser fiber) connected to the laser aperture. (See Figure 1, page 1-5.)

Laser Aperture

The output laser beam is delivered through an optical fiber. To connect a fiber, lift the dust shield by using the small lever located above the laser aperture. Then attach the fiber to the laser aperture. The lever is magnetized and will stay in the open position once the dust shield is lifted. (See Figure 1, page 1-5.)

Key Switch

The key switch turns the device on. There are two positions on the switch:  (on) and  (off). To power on the device, insert the key and turn it clockwise to the "on" position. To power off the device, turn the key to the "off" position and remove it. (See Figure 2, page 1-6.)

Emergency Stop Button

In case of an emergency, push the emergency stop button to immediately shut down the laser system. To reset the emergency stop button, twist the button counterclockwise and pull it out. (See Figure 1, page 1-5.)

Foot Switch

Plug the foot switch into the back panel of the laser system. When not in use, the foot switch can be stored on the back panel.


MAINS Switch

The MAINS switch controls AC power into the laser system. (See Figure 3, page 1-6.) Setting this switch to the "off" position disconnects the device from MAINS power.

System Specifications

General Specifications	
Electrical Requirements	HL-30A: 115 V ~ $\pm 10\%$ /50-60 Hz/15 A (20 W output) HL-30B: 115 V ~ $\pm 10\%$ /50-60 Hz/20 A (30 W output) HL-30E: 230 V ~ $\pm 10\%$ /50-60 Hz/10 A (30 W output)
Dimensions	20 in (W) x 23 in (D) x 51 in (H) 50 cm x 59 cm x 129 cm
Weight	225 lbs 102 kg
Temperature (operational)	50°F - 77°F, 10°C - 25°C
Temperature (storage)	41°F - 104°F, 5°C - 40°C
Relative Humidity	From 30% to 85% (noncondensing conditions only)
Cooling	Self-contained water-to-air heat exchanger
Electrical Protection	Class I
Type of Applied Part	BF (Note: Optical Fibers are Type BF Applied Parts)
Atmospheric Pressure	800-1060 hPa
Degree of Protection	IPX0 (IEC529) Not protected from solids. Not protected from liquids. IPX8 for foot pedal.
Laser Specifications	
Laser Type	Pulsed CTH:YAG (Ho:YAG)
Wavelength	2100 nm
Power to Tissue	up to 30 W
Energy/Pulse	0.5 - 3.5 J
Divergence	0.44 rad $\pm 5\%$
Application Mode	Pulsed
Pulse Duration	~150 - 800 μ s (short and long pulse - width modes)
Repetition Rate	5 - 20 Hz
Delivery Device	Wide range of flexible silica fibers
Output Deviation	Accuracy: The device contains an internal measuring system to control the actual emission of the laser energy. The device does not require user calibration. Deviation: The laser energy values can diverge from the preset values by a maximum of 20%.
Aiming Beam	Diode laser, green, 532 nm (adjustable < 5 mW)


Labels and Symbols


WARNING—Visible and invisible laser radiation.


Emergency Stop


Not for general waste


Remote InterLock

Optical Fiber
Applicator

Foot switch
laser activation device

Key Switch position

Type BF applied part

Chapter 2

Indications

Indications

This section discusses the use of the Rhapsody H-30 Holmium Laser System in clinical applications. The information is provided by specialty and includes procedural recommendations, along with specific indications and contraindications. The information provided in this section is not all-inclusive and is not intended to replace physician training or experience.

Only physicians and staff who have been appropriately trained and who are thoroughly familiar with the instructions and safety precautions provided in this manual should operate the system. A review of prior training is strongly encouraged and recommended.

Indications for Use

The Rhapsody H-30 Holmium Laser System is intended for incision, excision, resection, ablation, vaporization, coagulation, hemostasis of soft tissue, or stone fragmentation during open surgical, laparoscopic, percutaneous, or endoscopic procedures in urological, gynecological, gastroenterological, or pulmonary applications.

Contraindications

The Rhapsody H-30 Holmium Laser System should only be used in conditions where its use is appropriate and its efficacy has been proven. A qualified physician should perform clinical applications.

The use of the Rhapsody H-30 Holmium Laser System is contraindicated for patients:

- Who are inappropriate candidates for laser surgery
- Where appropriate anesthesia is contraindicated by patient history
- Whose general medical condition contraindicates surgical intervention
- Where tissue (especially tumors) is calcified
- Requiring hemostasis of vessels over approximately two millimeters in diameter
- Where laser therapy is not considered the treatment of choice
- Where any irrigant other than physiological saline or sterile water is used

The use of the Rhapsody H-30 Holmium Laser System is contraindicated for the treatment of menorrhagia in the following:

- Patients who desire future childbearing potential
- Patients who have any medical condition that contraindicates hysteroscopy, such as active pelvic inflammatory disease
- Patients who have indications of precancerous pathology

Warnings and Precautions - General

The information in this section applies to all specialties for which the Rhapsody H-30 Holmium Laser System is indicated.

The physician should become fully acquainted with the unique surgical and therapeutic effects produced with the 2100 nm wavelength before using the Rhapsody H-30 Holmium Laser System clinically. These effects may include coagulation, deep penetration, and intense cutting.

Unless otherwise stated in the specific application section, the physician should begin at the lowest power and use shorter durations. The physician should note the surgical effect and adjust the settings until the desired result is obtained.

Adverse events may result from deficiencies in operator skill level, or they may be due to a concomitantly occurring primary procedure. Careful patient selection and adherence to the manufacturer's Instructions for Use will help to mitigate risk.

- Extreme caution should be used until the physician is completely familiar with the biological interactions of the laser with various tissue types.
- Changes in the tissue texture and color are the best indicators of the laser's effect. The recommended power settings are less important than the observed effect on the tissue. Specific pulse width should depend on the tissue and is left to the surgeon's preference and best medical judgment.
- Use the lowest possible power settings that are required in order to achieve the desired effect on the tissue for treatment.
- Higher energy may be necessary to achieve the desired effect on the tissue if cooled irrigant is used. Higher power settings may increase the chance of damage to the laser fiber.
- The use of mechanical pressure on the laser fiber does not increase its cutting or vaporization effects but may induce bleeding, thermal damage, and fiber destruction.
- The risk of combustion, perforation, and laser-induced hemorrhage, any of which can cause patient death, must be fully explained to the patient prior to surgery.
- Use caution with patients who have had difficulty with previous endoscopic or laparoscopic procedures.
- The laser should be used only on a tissue target that is fully observable either by direct sight or through an endoscope.
- Steady the distal end of the fiber at the intended treatment target before depressing the foot switch.
- Vaporization or excision of anatomical structures in proximity to known critical structures, such as large arteries, large veins, bowel, ureter, or bladder should be performed carefully to avoid damage to such structures.
- Electrocautery should be easily accessible in case a bleeding artery or vein needs to be cauterized.

- The flammability of methane gas must be considered when treating the perianal area.
- Tissue perforation can occur if excessive laser energy is applied, whether through excessive pulse energy or extended laser application dwell time.
- Fiber tips are sharp; perforation may occur by mechanical means (poking).
- When risk of air embolism is present, such as when one is cutting or vaporizing organs with large vascular structures, always use fluid cooling.
- Screen tumors that are in close proximity to known arteries or veins by pretreatment workups to precisely locate the tumors.
- When radiotherapy and laser therapy are used concurrently, extra precautions should be taken, including conducting more stringent postoperative monitoring. Patients who have undergone radiotherapy present a greater risk of perforation and tissue erosion.
- Discontinue laser therapy immediately if the patient develops cardiopulmonary problems.
- Cook does not have clinical information about or experience with the use of the Rhapsody H-30 Holmium Laser System on pregnant women or nursing mothers.
- Patients who have experienced discomfort during previous laser treatment may require analgesics.
- Anesthesia considerations: If the laser is used in an oxygen-rich environment (such as in endobronchial use), oxygen concentrations should not exceed 21-30% in order to diminish the fire hazard. Endotracheal tubes should be either non-combustible or wrapped with reflective tape.
- If the physician desires to use the delivery system without a contact probe or scalpel attached, the physician must be aware of the increased risk of back-scatter and forward-scatter.
- Because of possible fluid absorption risks, patients with cardiovascular disorders should be carefully evaluated as candidates for this procedure.
- When appropriate, moistened sponges should be placed on the tissues and instrumentation surrounding the laser area to help dissipate incident laser energy and inadvertent heat production.
- Endoscopic treatment may result in adverse reactions such as stone migration, burns, hemorrhage, laceration, infection, inflammation, pain, ureteral avulsion, air embolism, pneumothorax, or bladder, ureteral, or bile duct perforation. In extreme cases, death may occur due to procedural complications, concurrent illness, or the application of laser radiation.
- Adverse events reported in clinical literature in the context of holmium laser procedures include subcapsular renal hematoma, bladder neck contracture, urethral stricture, residual adenoma, hematuria, dysuria, acute urinary retention, urinary incontinence, bladder neck stenosis, difficult micturition, urethral or bladder injury, ureteric orifice injury, hydronephrosis, steinstrasse, recatheterization, clot retention, acute epididymitis, septic shock, hypoxemia, mucosal injury, pancreatitis, cholangitis, biliary stricture, and ascites.

Warnings and Precautions - Gynecological

- During intrauterine laser surgery, do not use gas or air for cooling the laser fiber tip or for insufflation. This may cause life-threatening gas or air embolism.
- Patient monitoring: Special precautions should be taken to monitor the patient for fluid absorption (from hysteroscopic distention media). Excessive fluid absorption is worrisome if the procedure is prolonged. During the procedure venous channels seem to be open, but bleeding is controlled by the pressure of the irrigating solution. There is a large raw surface, and a large volume of fluid may be absorbed. Wider dilation of the cervix to allow irrigating fluid to flow more freely around the hysteroscope and reduce intrauterine pressure may reduce the risk. If excessive fluid is absorbed, appropriate therapy such as diuretic administration should be considered with subsequent electrolyte monitoring.
- Consideration should be given to closing the fallopian tubes prior to an intrauterine procedure to treat menorrhagia. However, if the physician, after consultation with the patient, determines that sterilization is not indicated, the patient must be advised to utilize a contraceptive technique to avoid pregnancy. Although pregnancy is highly unlikely after this type of procedure, there is no data regarding the potential risks should pregnancy occur.
- Patients are urged to have a pap smear annually. There is a possibility that in the very small percentage of women who develop uterine cancer, their symptoms of bleeding may be hidden as a result of this treatment. While other symptoms would still be expected to occur (pain, other discharges, etc.), a pap smear should be performed to reveal any possible cancerous development.
- Due to the etiology of the disease, laser treatment of endometriosis and peritoneal adhesions may result in only temporary improvement. A recurrence of the disease may require retreatment.
- Distention of the uterine cavity is at the discretion of the operating hysteroscopist. The recommended distention media for laser hysteroscopy is normal saline (0.9% sodium chloride).

Warnings and Precautions - Gastroenterological and Pulmonary

- "Flash" fire can occur during pulmonary procedures. Do not use general inhalation anesthetics. Oxygen levels in the procedure room must not be higher than 30%.
- Use caution with patients who have had previous esophageal or tracheal fistulae or episodes of aspiration.
- After esophageal procedures, swallowing may be paradoxically worsened rather than immediately improved. This is caused by secondary tissue edema. This potential problem should be explained to the patient prior to surgery.

Chapter 3

Safety, Warnings, and Caution

General Warnings and Caution

The Rhapsody H-30 is a holmium laser system, with up to 30 W, for surgical procedures. The device is manufactured in accordance with Annex II of the 93/42/EEC standard and is identified as follows:

Device	30 W Holmium Laser System
Name	Rhapsody® H-30® Holmium Laser System
Models	HL-30A (20 W), HL-30B, HL-30E
Device Category	Laser for medical use
Device Classification	II b (according to 93/42/EEC)
Laser Classification	IV
Manufactured By	Cook Incorporated 750 Daniels Way Bloomington, IN 47404, USA Phone: 800.457.4500 www.cookmedical.com

General Safety Information

- This manual contains important information about the safe use of the laser system.
- All persons who operate the laser system must understand the operating and safety instructions specified in this manual (IFU).
- Only authorized individuals with appropriate medical laser training and knowledge should operate the laser system.
- Only authorized individuals with appropriate laser training and knowledge should service the laser system.
- Only authorized service personnel should access the internal electrical or optical components.
- The laser system must be completely closed and assembled for normal operation. Only trained, authorized personnel are allowed to open the external covering panels.
- Access to the user manual should be available in the working area of the laser system.
- All of the warning labels must be kept in good condition.

DO NOT dispose of this product with residential waste. Disposal of this device should be in accordance with applicable local, regional, and national laws and regulations. This device cannot be disposed of as municipal waste, but must be separated in accordance with the WEEE (Waste Electrical and Electronic Equipment) Directive.

Staff Training

The use of the laser system is restricted to medical staff who have received proper training in laser

safety, operation of the Rhapsody H-30 Holmium Laser System, and the identification and use of appropriate personal protective equipment.

Working Area

This Rhapsody H-30 Holmium Laser System is a class 4 laser and must be used in a specific nominal hazard zone (NHZ) or working area that is defined and delimited in accordance with all applicable local and international standards (IEC60825-1:2007, ANSI Z136.3). In general clinical operation, the working area is typically defined as the operating room.

Repair and Modifications of the Device

Only authorized service personnel should perform repairs and maintenance. Constructive modifications are not permitted.

NOTE: Unauthorized repair, maintenance, and/or modifications to the laser system will void the warranty.

WARNING! Performing any procedures, operating any controls, or making any adjustments other than those specified herein may result in hazardous radiation exposure.

CAUTION: The use of controls or adjustments or the performance of procedures other than those specified herein may result in hazardous radiation exposure.

Eye and Skin Exposure

The working laser beam of the Rhapsody H-30 Holmium Laser System is invisible and, if exposed to the eyes, may cause sight loss. The skin can resist a significantly higher amount of laser energy than the eyes, but the laser beam can also seriously burn the skin.

A danger exists in any case of:

- Direct laser radiation
- Reflected laser radiation
- Diffused laser radiation

WARNING! All personnel present in the laser working area must wear appropriate laser safety eyewear.

WARNING! Incorrectly inserting the fiberoptics, sharply bending or improperly securing the fiberoptics, or failing to follow the fiberoptics installation procedure may lead to damage to the fiber or delivery system and may result in harm to the patient or operator.

Use protective eyewear with either of the following specifications:

- I 2100 L2 (EUR) @ 2100NM
- OD 4 (US) @ 2100NM

Always check the integrity and condition of protective eyewear.

Before wearing protective eyewear, be sure that the protective glass and frame are in good condition.

The worst-case nominal ocular hazard distance (NOHD) for the Rhapsody H-30 Holmium Laser System emitting laser energy from a fiber is 0.5 m, and corresponding maximum permissible exposure (MPE) for Holmium laser is $0.1 \text{ W}\cdot\text{cm}^{-2}$. During service operations, there is access to laser energy emitted directly from the laser source, which raises the NOHD to 15 m.

Hazards

Optical Hazard

The Rhapsody H-30 Holmium Laser System generates laser light at a wavelength of 2,100 nm, with a maximum energy of approximately 3.5 J, delivered from the optical fiber. Greater energies can be generated from the laser resonator, especially during service operations. At the Rhapsody H-30 Holmium Laser System's wavelength and operational energy levels, serious and permanent eye damage can occur when there is direct or indirect optical exposure.

The laser beam emitted by this laser system can cause vision loss. Laser energy transmitted by this laser that enters the eye can cause a corneal burn or can be focused on the retina. Direct absorption of laser energy by the retina can result in temporarily clouded vision, retinal lesion, long-term scotoma, or long-term photophobia.

Danger exists from any of the following:

- Direct laser radiation
- Reflected laser radiation
- Diffused laser radiation

Electrical Hazards

- Only a qualified service technician may remove the enclosure covers. Removing the enclosures exposes personnel to potential electrical hazards.
- Do not touch areas marked HIGH VOLTAGE. These and the surrounding components of the laser system contain charge voltages of up to 700 VDC.
- Do not operate the laser system if the floor is flooded.
- Service technicians should not power on the laser system unless all high-voltage areas are dry and not in danger of flooding.
- Inspect the power cord before use. Do not use a damaged power cord or a power cord with a damaged or loose power plug.

Fire/Explosion Hazard

The Rhapsody H-30 Holmium Laser System has the ability to melt, burn, or vaporize almost any material. A fire can also occur because of the nature of the laser radiation. Precautions should be taken when alcohol, gasoline, drapes, or any other ignitable materials are in or around the working area. Keep a fire extinguisher in close proximity to the laser system.

NOTE: In a clinical environment, consider using sterile water to extinguish any fire caused by laser radiation.

WARNING! High temperatures produced during the normal use of the laser system may ignite some materials (e.g., cotton wool when saturated with oxygen), or endogenous gas. Solvents of adhesives and flammable solutions used for cleaning and disinfecting should be allowed to evaporate before the laser system is used.

Fire/Explosion Hazard Precautions:

- Keep only a minimum amount of any combustible materials (e.g., alcohol) in the service room. If possible, stock combustible materials away from the working area while you are servicing the laser system.
- Always keep a fire extinguisher and/or cool water in the immediate area.
- Never direct the laser beam onto any surface except a power meter or an appropriate beam dump or heat sink.

NOTE: Water acts as an excellent heat sink for the Rhapsody H-30 Holmium Laser System's laser radiation, especially when the distal end of the fiber is immersed.

Laser-Generated Airborne Contaminants (LGAC)

The smoke plume created by laser contact with various materials can create toxic gases and vapors. If the target material is tissue, then both dead and live cellular material, including viable viruses, may be present. The potential hazardous area for LGAC exposures may often exceed the working area.

CAUTION: Laser fume and/or plume may contain viable tissue particulates.

LGAC Precautions

Use respiratory protection or ventilation (such as high-quality, well-filtered smoke evacuators) to control LGAC whenever laser radiation might generate plumes.

NOTE: Surgical masks do not protect against LGAC.

Other Hazards

- **RISK OF FLUID INGRESS**—This product is rated as IPX0 and does not protect against ingress of fluids. For cleaning, apply a cleaning agent to a cloth and wipe the external case. Do not pour or spray liquids directly onto the device. Do not operate the device if it has become wet.
- **RISK OF FIRE**—Do not point the distal end of the fiber toward flammable materials such as surgical drapes.
- **RISK OF EXPLOSION**—Do not operate the laser system in the presence of flammable substances, including flammable anesthetics.
- **RISK OF UNINTENDED MOVEMENT**—All four wheels on the laser system have non-directional locks, meaning they can be locked in any position, thus preventing unintended movement of

the laser system.

- **RISK OF FREEZING**—If any water remains in the cooling system, the laser must not be subjected to freezing temperatures.

Safety Indicators

The yellow **standby** button located on the lower left of the touch screen control panel indicates that the laser system is in the “standby” state. Depressing the foot switch will not activate the laser.

The green **ready** button indicates that the laser system is in the “ready” state. In the ready state, no system errors are present and the touch screen is displaying the user power-settings screen. Depressing the foot switch **will** activate the laser.

NOTE: When the laser system switches from “standby” to “ready,” there is a two-second delay. This delay allows time for the laser system operator to give adequate notice to personnel within the working area.

An audible warning alerts personnel that the laser is firing.

NOTE: The audible warning can be turned off in the laser options screen.

The emission-warning indicator (the blue LED atop the control panel) lights up when the laser system is in the “ready” state and blinks when the laser is firing.

The information status bar, located on the touch screen control panel, provides alerts of any internal system errors. The errors are referenced on pages 1-43 — 1-44 (see chart).

NOTE: In the event of an error, the laser system will immediately return to “standby” and the information status bar will display a corresponding error code.

Interference with Other Devices

Electromagnetic fields from other irradiants in the area surrounding the laser system can affect the laser’s performance. The laser’s electromagnetic field can also have an effect on other electrical devices.

Do not use cellular phones in the immediate vicinity of the laser system. See the section on Electromagnetic Emission and Interference (pages 1-27 — 1-29) for details on potential interference and separation distances.

Chapter 4

Installation

Installation

The Rhapsody H-30 Holmium Laser System should be installed by trained personnel upon initial receipt of the laser system from Cook Medical. Any damage that occurs during or after transportation must be corrected prior to use.

WARNING! Do not attempt to use the Rhapsody H-30 Holmium Laser System until trained personnel have performed the installation procedure.

Transportation

Transportation of the laser system must be performed in compliance with all local safety laws.

The laser system must be packaged to protect it from damage, such as in the original shipping crate. If the original shipping crate is not used for transportation, the laser system must be secured to prevent damage. The handles can be used to tie the laser system to a secure surface. The crate should be moved only by proper means of lifting and transporting, such as a pallet jack.

The laser system should not be subjected to excessive force, shock, or vibration. If excessive shock or vibration is suspected to have occurred, a trained technician should inspect the laser alignment.

If any water remains in the cooling system, the laser system must not be subjected to freezing temperatures. Freezing water within the laser cavity can damage optical components, such as the laser cavity, rod, or flash lamp, or mechanical components, such as the heat exchanger or water pump. If there is a possibility of freezing temperatures, the laser's cooling system should be drained of all water by attaching fill/drain tubes to both the air purge and water drain ports on the rear panel of the laser system. (See Figures 2 and 3, page 1-6.)

Storage and transit temperatures must not exceed 40°C (104°F). Storage humidity should not exceed 95%. If the laser is stored in a highly humid location, it must be allowed sufficient time to reach equilibrium when introduced to a cool operating environment. Failure to do this may result in damage to the optics from condensation.

Unpacking the Laser System

Your Rhapsody H-30 Holmium Laser System is shipped from the factory in a wooden crate to protect it from damage during shipment. Before removing the laser system, inspect the crate for damage. There are tip/tilt and shock indicators on both the outside and inside of the crate. If no damage is evident, continue to unpack the laser system. If there is evidence of damage, save all packing materials and notify the shipping company by filing an insurance claim report.

A Cook sales or service representative will help you with the installation of the laser system. If you have any questions, please contact Customer Relations according to your facility's location. After the Rhapsody H-30 Holmium Laser System has been unpacked, please store the shipping crate and packing materials in a safe location in case the laser system has to be transported in the future.

Not all of the latches need to be released in order for you to remove the laser system from the crate, but if you release all of the latches and unscrew the angled brackets connecting the side panels to the pallet, you can easily disassemble the crate and lay it flat for easy storage.

Cleaning

Prior to cleaning the laser system, protect the fiber connection port with the fiber applicator shutter closed. Unplug and power off the laser system prior to cleaning.

Do not use alcohol or an alcohol solution to clean the display.

Exterior enclosures may be wiped down with alcohol or another nonabrasive cleaner. Use a cloth that has been sprayed with a cleaning or disinfecting agent. Avoid spraying the disinfectant directly onto the laser system, particularly directly onto the touch screen. Do not let water enter the laser system or leak into the fiber connection port. Solutions used for cleaning and disinfecting should be allowed to evaporate before the laser system is used.

Contents

The Rhapsody H-30 Holmium Laser System is shipped with the following components and accessories. Please take inventory of the items listed below and contact your designated Customer Relations department if any items are missing or damaged.

- Laser system
- Foot switch
- Bag that contains the following:
 - Service report form with final Quality Assurance checklist
 - Two keys for the key switch
 - Remote interlock connector
- Two fill/drain tubes
- Plastic water bottle for filling the cooling system
- Vinyl cover over the laser system

Ambient Operating Conditions

The laser system should be used in a well-ventilated, air-conditioned operating room. The procedure room temperature should be kept between 18 and 24°C (65-75°F) for optimum system performance. In general, holmium lasers are very sensitive to temperature. When the laser cooling system reaches above 38°C (100°F), the safety circuits will prevent the laser system from operating. The information bar on the touch screen control panel will read "Over temperature." The laser system may overheat for the following reasons:

- There is improper ventilation.
- The ambient operating room temperature is above 25°C (77°F).
- The laser system has been running at high power for an extended period of time (more than 30 minutes).
- The input line voltage is lower than normal and has caused the cooling fans or water pump to underperform.

Filling the Cooling System with Water

1. Remove the external and internal right side panels (as you face the Rhapsody H-30 Holmium Laser System). This will expose the water reservoir and DI cartridge to view.

NOTE: 3 mm and 5 mm Allen wrenches will be needed to disassemble the side panels.

2. Insert the air purge tube into the air purge (top) port on the back of the laser system.
3. Insert the fill bottle into the water fill port on the back of the laser system.

NOTE: Use distilled water only.

4. Begin filling. Fill until water reaches the top of the reservoir.
5. Disconnect the fill bottle and remove the small tube (blue or red) with the connector from the top of the fill bottle.
6. Insert the tube from step 5 into the water drain port. Water will begin to drain. The water flow will be interrupted with burps of air. When the water flows continuously again, remove the drain tube. This will easily take place before the reservoir empties.

NOTE: The laser system has built-in safety circuitry to detect when there is too little water flow. If this circuitry engages while it fills with water, the information bar will display "Low coolant flow." Simply turn the laser system off and back on to continue.

7. Replace the tube back into the fill bottle and re-insert the fill bottle into the water fill port.
8. Continue filling until the reservoir is again full, and then stop.
9. Remove the tube again as you did in step 5.
10. Have the fill/drain tube ready near the drain port.
11. Insert the fill/drain tube into the air purge port and watch the reservoir. Turn the laser system on. When you hear the water pump turn on, and the water level suddenly drops, the pump is primed.
12. Quickly turn off the key switch.
13. Refill the reservoir to the top.
14. Turn the key switch on and watch the reservoir level drop before turning the key switch off.

Do not let the water level drop below one-third or suck air into the bottom.

15. Repeat steps 13-14 until the reservoir remains full and the system runs continuously.

NOTE: You may have to do this several times before the process is complete. The system should take about one and one-fourth liters of water when full.

16. Remove the fill bottle and air purge tube. Replace the external and internal side panels.

Initial System Power-On

Review the previous sections in this chapter and verify that all site installation requirements have been met. The laser system should only be operated in an area of restricted access. Post a laser warning sign on each access door.

1. Connect the foot switch and remote interlock connector to the laser system. (See Figure 3, page 1-6.) Find the remote interlock connector supplied separately in the bag that comes with the unit.
2. Connect the laser device to a suitable AC electrical outlet.

WARNING! To avoid the risk of electric shock, connect this laser system only to a supply MAINS with protective earth.

3. Ensure that the emergency stop button isn't engaged. (See Figure 1, page 1-5.)
4. Turn the MAINS switch on. (See Figure 3, page 1-6.)

NOTE: Position the system so both the MAINS switch and power connection to the electrical outlet can be easily reached.

5. Turn the key switch on. (See Figure 2, page 1-6.)
6. Lift the laser aperture cover and connect a 550 μm holmium laser fiber to the laser. (See Figure 1, page 1-5.)

NOTE: Do not over-tighten the laser fiber. The connector only needs to be finger tightened for proper operation.

7. After a few seconds, the information bar on the control panel screen will display "STANDBY." (See Figure 4.)

Figure 4

8. Press the **ready** button. (See Figure 6, page 1-32.) The "READY" progress bar will take about two seconds to fill.

9. Once the "READY" bar is fully lit, the information bar will display "READY to lase." This indicates that the laser system is operating normally. (See Figure 4.)
10. Follow the instructions in the next section to verify alignment and calibration.

Verifying Alignment and Calibration

The laser system has been designed with many safeguards to prevent optical misalignment during shipping and normal operation. However, it is important to verify the alignment and calibration of the laser system to ensure optimal system performance. The following test requires an energy or power meter calibrated to 2100 nm.

The following test applies to both 20 W and 30 W models. Note that some steps apply only to the 30 W models.

1. Set up the laser system for lasing.
2. Verify that the 550 μm test fiber has clean surfaces on both the proximal connector and distal ends.
3. Insert the 550 μm test fiber into the laser aperture.

WARNING! Wear appropriate safety eyewear. While the laser is in operation, all personnel in the operating room must wear appropriate laser eyewear.

4. Set up the energy or power meter in front of and parallel to the distal face of the test fiber establishing the measuring plane, as shown (see Figure 5). Be careful not to place the fiber tip too close to the energy or power meter probe, as the energy density may damage the detector's sensitive surface.


Figure 5

5. Set the power parameters to Short pulse width mode, 10 Hz and 0.5 J. (This step applies to both 20 W and 30 W models.)
6. Press the **ready** button and wait until the information bar displays "READY to lase."
7. Depress the foot switch for 5-10 seconds until you attain a stable reading on the energy or power meter. If the reading is not within 10% of the displayed settings, the laser system is either misaligned or out of tune. Please contact your sales representative or the Cook Capital Equipment Service and Repair Center.

8. Repeat steps 5 through 7 for the following settings as well:
 - a. Short pulse width mode, 10 Hz and 2.0 J. (This step applies to both 20 W and 30 W models.)
 - b. Short pulse width mode, 10 Hz and 3.0 J. (This step applies to 30 W models only)

Electromagnetic Emission and Interference

IMPORTANT! Medical electrical devices are subject to special precautions with regard to electromagnetic compatibility (EMC) according to IEC 60601-1-2. Make sure you observe the notes on EMC for installation and operation. Medical electrical devices can be influenced by mobile RF communication devices (i.e., mobile phones). If it is necessary to stack the devices or place them next to each other, and RF interference is observed, verify that the Rhapsody H-30 Holmium Laser System is operating normally. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the laser system.

IEC 60601-1-2 Table 1 Guidance and Manufacturer's Declaration - Electromagnetic Emission		
The Rhapsody H-30 Holmium Laser System is intended for use in the electromagnetic environment specified below. The customer or the end user of the laser system should ensure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment - Guidance
RF Emission CISPR 11	Group 1	The Rhapsody H-30 Holmium Laser System is not intended to use or produce RF energy. Therefore its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emission CISPR 11	Class A	The Rhapsody H-30 Holmium Laser System is not suitable for installation in all buildings, including domestic buildings and those directly connected to the public supply network in low tension. It should only be installed in buildings like hospitals with dedicated supply systems.
Harmonic Emission IEC 61000-3-2	Not Applicable	
Voltage Fluctuation / Flicker Emission IEC 61000-3-3	Not Applicable	

IEC 60601-1-2 Table 2
Guidance and Manufacturer's Declaration - Electromagnetic Emission

The Rhapsody H-30 Holmium Laser System is intended for use in the electromagnetic environment specified below. The customer or the end user of the laser system should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guide
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 KV Contact ±8 KV Air	±6 KV Contact ±8 KV Air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient / Burst EC 61000-4-4	±2 KV for power supply lines ±1 KV for I/O lines	±2 KV for power supply lines N/A for I/O lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 KV differential mode ±2 KV common mode	±1 KV differential mode ±2 KV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	<5% Ut for 0.5 cycles 40% Ut for 5 cycles 70% Ut for 25 cycles <5% Ut for 5 seconds Note: Ut is the AC mains voltage prior to the application of the test level.	<5% Ut for 0.5 cycles 40% Ut for 5 cycles 70% Ut for 25 cycles <5% Ut for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Rhapsody H-30 Holmium Laser System requires continued operation during power mains interruptions, it is recommended that the laser system be powered from an uninterruptible power supply or battery.
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

IEC 60601-1-2 Table 3
Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Rhapsody H-30 Holmium Laser System is intended for use in the electromagnetic environment specified below. The customer or the end user of the laser system should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guide
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 KHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	Portable and mobile RF communication equipment should be no closer to any part of the Rhapsody H-30 Holmium Laser System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.167 \sqrt{P}$ $d = 1.167 \sqrt{P}$ from 80 MHz to 800 MHz $d = 2.333 \sqrt{P}$ from 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcasts, and TV broadcasts cannot be predicted with accuracy. To assess the electromagnetic environment in the presence of fixed RF transmitters, consider an electromagnetic site survey. If the measured field strength in the location in which the laser system is used exceeds the applicable RF compliance level above, verify that the laser system is operating normally. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the laser system.

^b Over the 150 KHz-80 MHz frequency range, field strength should be less than 3 V/m.

Chapter 5

Operation

System Power-On

Power on the laser system and connect the optical fibers by following the steps below.

1. Verify that the laser device is connected to a suitable AC electrical outlet.
2. Ensure that the emergency stop button isn't engaged. (See Figure 1, page 1-5.)
3. Turn the MAINS switch on. (See Figure 3, page 1-6.)
4. Turn the key switch on. (See Figure 2, page 1-6.)
5. Lift the laser aperture cover and connect the holmium laser fiber to the laser. (See Figure 1, page 1-5.)

NOTE: Do not over-tighten the laser fiber. The connector only needs to be finger tightened for proper operation.

6. After a few seconds, the information bar on the control panel screen will display "STANDBY." (See Figure 4, page 1-25.)
7. Press the **ready** button. (See Figure 6, page 1-32.) The "READY" progress bar will take about two seconds to fill.
8. Once the "READY" bar is fully lit, the information bar will display "READY to lase." This indicates that the laser system is operating normally. (See Figure 4, page 1-25.)

Checking the Optical Fibers

The laser system includes a low-power green aiming diode to make it possible to identify the target during operation. This aiming diode can be used to verify the integrity of the optical fiber. Follow the instructions below to verify the integrity of the optical fiber.

1. Remove the optical fibers from their sterile packages using aseptic technique and control them to be sure that no damage occurs.
2. Plug the fiber in and check the shape of the green spot. The intensity of the aiming beam can be adjusted from the control display.
3. If the green spot from the aiming beam is well rounded, the fiber is good and working properly.
4. If the green spot from the aiming beam is misshapen, blurred, low in intensity, or missing, the fiber may be damaged. Remove and replace the fiber, or contact your service representative.

WARNING! All personnel present in the laser working area must wear appropriate laser safety eyewear.

WARNING! Incorrect insertion of the fiberoptics, sharply bending or improperly securing the fiberoptics, or failing to follow the fiberoptics installation procedure may lead to damage to the fiber or delivery system and may result in harm to the patient or operator. Refer to pages 1-18 - 1-19 for additional hazards.

Laser Screens

The control panel screen display contains two functional screens: the main screen and the options screen. Touch the tabs at the top right-hand corner of the display to toggle between these screens.

Main Screen

The main screen allows you to adjust functions for the main power settings.


Figure 6

Main Screen Items (See Figure 6)

Main Screen Tab

- Function: Displays power parameters and adjustment controls.
- Operation: Touch the tab to access the main screen.

Energy Bar

- Function: Displays the set pulse energy in joules and records total delivered energy (TOTAL) in kilojoules.
- Operation: Touch the bar to highlight orange and adjust energy per pulse with the adjustment controls.

Frequency Bar

- Function: Displays pulse frequency in Hertz.
- Operation: Touch the bar to highlight orange and adjust pulse frequency with the adjustment controls.

Preset Buttons

- Function: Stores a power parameter or jumps to a stored power parameter.
- Operation - Store a Power Parameter: Touch the S button and then touch a numbered preset button to store the power parameters that are currently displayed.
- Operation - Jump to a Power Parameter: Touch a preset button (an active preset button will be dark gray with a white number, as shown for preset #1) to jump to a stored power parameter.

Adjustment Controls

- Function: Adjust the setting that is highlighted orange: pulse energy, frequency, or visual aiming beam intensity.
- Operation: Touch the up arrow to increase the setting.
- Operation: Touch the down arrow to decrease the setting.

Information Bar

- Function: Displays simple instructions, system status, and error messages.

Pulse Width Toggle Button

- Function: Toggles between Short and Long pulse width modes.

Standby Button

- Function: Signals that the system is in "Standby" mode with a yellow bar and icon.
- Operation: Touch the button to return to "Standby" mode.

Visual Aiming Beam Button


- Function: Adjusts the intensity of the visual aiming beam with the adjustment controls when they are highlighted orange.
- Operation: Touch the button to toggle on or off. The icon will be white when on and gray when off.

Ready Button

- Function: Signals that the system is in "Ready" mode with a green bar and icon.
- Operation: Touch the button to enter "Ready" mode.

Options Screen

The Rhapsody H-30 Holmium Laser System has many options that are displayed and set via the options screen. There are six option subscreens. The general layouts of all of the subscreens are identical, but the active screen is indicated by text in the information bar as well as the white text ("GEN") in the lower buttons. Touch the desired button to display the corresponding subscreen.

Options Screen Items (See Figure 7)

Options Screen Tab

- Function: Displays subscreens of options and the controls for adjustment.
- Operation: Touch the tab to enter the options screen.

Main Text Box

- Function: Displays relevant data in the upper left of the screen.

Adjustment Controls

- Function: Controls scrolling of the text or active line (displayed in yellow text) within the main text box.
- Operation: Touch the up arrow to scroll up and the down arrow to scroll down.

Information Bar

- Function: Displays simple instructions, system status, and error messages.

Enter Button

- Function: Sets various options as directed by instructions in the information bar.
- Operation: Button is active when the button text and icon are white. The button text and icon will be gray when inactive.

Parameter Determination

Choosing the ideal laser fiber

Cook Medical's holmium laser fibers for use with the Rhapsody H-30 Holmium Laser System come in 150, 200, 273, 365, 550, and 940 μm core diameters for single-use fibers and 273, 365, 550, and 940 μm core diameters for multi-use fibers. Larger-core fibers will allow the use of higher power but are less flexible. Higher pulse energy and larger fibers result in larger ablation volumes.¹ Large-core fibers will also be consumed less during procedures because of their greater thermal mass and their ability to dissipate parasitic thermal energy faster. Users are encouraged to select the largest core size feasible considering the flexibility needed to reach the target and the size of the endoscope's working channel.

¹ Kang HW, Lee H, Teichman JMH, et al. Dependence of calculus retropulsion on pulse duration during Ho:YAG laser lithotripsy. *Lasers Surg Med.* 2006;38(8):762-772.

Setting the Pulse Width

(via the toggle button on the main screen) (See Figure 6, page 1-32.)

Lithotripsy

- Short pulse width will cause larger fragmentation particulates ("busting" the stone) but will also cause greater retropulsion.
- Long pulse width will cause smaller fragmentation particulates ("dusting" the stone) and will cause less retropulsion.

Soft Tissue

- Use the Short pulse width if you want to minimize coagulation or create a large, thick cut or vaporization area.
- Use the Long pulse width if you want to increase coagulation and decrease bleeding.

NOTE: Most soft tissue cases will benefit from using the Long pulse width mode.

NOTE: If conditions have changed, such as a piece of an impacted stone breaks off or the soft tissue is bleeding more than expected, the pulse width mode can easily be switched with the tap of a toggle button.

Adjusting the Energy per Pulse

Adjust the energy per pulse to match the desired clinical effect.

Lithotripsy

- Lower pulse energy will cause the particulate size to be smaller ("dusting" the stone) and cause less retropulsion.
- Higher pulse energy will create larger fragments ("busting" or "drilling" the stone) but will also cause more retropulsion.

Soft Tissue

- Use higher pulse energy to create larger ablation volumes than lower pulse energy.

NOTE: Higher pulse energy will also increase bleeding.

The information in the following sections provides further indications, contraindications, warnings, and precautions for the Rhapsody H-30 Holmium Laser System for specific specialties as appropriate.

Laser Parameters

The following tables display available laser parameters; possible combination cells are marked with the average power (power = energy x frequency) in watts for all available settings.

With a 150 µm fiber, the maximum power is 6 W.

With a 200 µm fiber, the maximum power is 10 W.

With a 273 µm fiber, the maximum power is 15 W.

With a 365 µm fiber, the maximum power is 20 W.

With a 550 µm fiber, the maximum power is 30 W.

With a 940 µm fiber, the maximum power is 30 W.

Combinations in Short pulse width mode:

Watts	5 Hz	8 Hz	10 Hz	12 Hz	15 Hz	20 Hz
0.5 J	2.5	4.0	5.0	6.0	7.5	10.0
0.6 J	3.0	4.8	6.0	7.2	9.0	
0.7 J	3.5	5.6	7.0	8.4	10.5	
0.8 J	4.0	6.4	8.0	9.6	12.0	
1.0 J	5.0	8.0	10.0	12.0	15.0	
1.2 J	6.0	9.6	12.0	14.4	18.0	
1.5 J	7.5	12.0	15.0	18.0	22.5	
2.0 J	10.0	16.0	20.0	24.0		
2.5 J	12.5	20.0	25.0	30.0		
3.0 J	15.0	24.0	30.0			
3.5 J	17.5					

Combinations in Long pulse width mode:

Watts	5 Hz	8 Hz	10 Hz	12 Hz	15 Hz	20 Hz
0.5 J	2.5	4.0	5.0	6.0	7.5	10.0
0.6 J	3.0	4.8	6.0	7.2	9.0	
0.7 J	3.5	5.6	7.0	8.4	10.5	

Watts	5 Hz	8 Hz	10 Hz	12 Hz	15 Hz	20 Hz
0.8 J	4.0	6.4	8.0	9.6	12.0	
1.0 J	5.0	8.0	10.0	12.0	15.0	
1.2 J	6.0	9.6	12.0	14.4	18.0	
1.5 J	7.5	12.0	15.0	18.0	22.5	
2.0 J	10.0	16.0	20.0	24.0		
2.5 J	12.5	20.0	25.0			
3.0 J	15.0					

Turning the System Off

To turn off the Rhapsody H-30 Holmium Laser System, place the system in Standby mode, turn the key switch to the off position, and remove the key.

WARNING! If the system is turned off without entering Standby mode, the parameters will not be saved.

Protection Against Unauthorized Use

When the Rhapsody H-30 Holmium Laser System is not in use, protect it against unauthorized use by removing the key from the key switch.

Error Messages

There are several error messages that can be displayed either in the information bar on the main screen or in the error log on the options screen.

The following are possible types of messages:

- Simmer failure: The flash lamp does not simmer. Call for service.
- Shutter failure: The shutter is not functioning properly. Call for service.
- Over temperature: The temperature of the cooling water is too high. Call for service.
- Low coolant flow: If water level is low, fill the reservoir with distilled water.
- No fiber attached: The fiber is not connected or is connected incorrectly. Unplug the fiber and then plug it in correctly.
- Remote open: The remote connector, near the foot switch, is not connected or is connected incorrectly. This error warning is displayed when the system is in Ready mode. Connect remote.

- Pedal not connected: The foot switch connector is not connected or is connected incorrectly. This error warning is displayed when the system is in Ready mode. Connect foot switch.
- Energy < 80%: This warning appears when the output energy is less than 80% of the set value. Call for service.
- Energy > 120%: This warning appears when the output energy is greater than 120% of the set value. Call for service.
- Energy > 200%: This warning appears when the output energy is greater than 200% of the set value. Call for service.
- Maintenance due: This warning appears when preventative maintenance is due. Call for service.
- MPS charge error: This error indicates the main power supply is not functioning properly. Call for service.
- MPS over temp: This error indicates the main power supply is not functioning properly. Call for service.
- MPS overload: This error indicates the main power supply is not functioning properly. Call for service.

Chapter 6

Maintenance

General Maintenance

If the Rhapsody H-30 Holmium Laser System is used carefully in normal operational conditions, the anticipated service life of the device is 10 years and the supplier recommends an overall system check every 6 months. Intense use, a dusty or dirty operational environment, or frequent movement necessitates more frequent maintenance. Verify the alignment and calibration using the instructions in chapter 4.

An authorized service provider is defined as: a Cook service representative, a third party representative contracted by and representing Cook, or a staff biomedical technician who has received laser technical training and received a certificate of such training by Cook.

WARNING! Only authorized personnel should provide service for the Rhapsody H-30 Holmium Laser System. Service by unauthorized personnel will void the warranty.

WARNING! Performing any procedures, operating any controls, or making any adjustments other than those specified herein may result in hazardous radiation exposure.

Cleaning

Prior to cleaning the laser system, protect the fiber connection port with the fiber applicator shutter closed. Unplug and power off the laser system prior to cleaning.

Do not use alcohol or an alcohol solution to clean the display.

Exterior enclosures may be wiped down with alcohol or another nonabrasive cleaner. Use a cloth that has been sprayed with a cleaning or disinfecting agent. Avoid spraying the disinfectant directly onto the laser system, particularly directly onto the touch screen. Do not let water enter the laser system or leak into the fiber connection port. Solutions used for cleaning and disinfecting should be allowed to evaporate before the laser system is used.

Cooling Water Level

Refill the water tank if the water level is too low.

Inspect the water level using the water level indicator on the rear panel. If the red float is not at the top of the indicator, refill the water tank with distilled water as specified in chapter 4, page 1-24.

Power Cord

The Rhapsody H-30 Holmium Laser System has a 15 foot (5 m) power cord mechanically affixed to it.

The power cord may deteriorate over time because of general use and environmental conditions. Periodically check the status of the power cord and contact your designated Customer Relations department if a problem exists.

Labels

It is up to the user to maintain the safety labels and ensure that they are in good condition.

Replacing the Blast Shield

The blast shield is a replaceable optical window that is located between the output focusing lens and the proximal surface of the fiber delivery system. This window protects the output focusing lens from damage from debris that could be ejected from a damaged fiber.


The blast shield is located on a bracket that is inserted from the top of the output focusing lens mount, just behind the laser aperture/fiber bulkhead connector.

Required Equipment: #2 Phillips-head screwdriver

Estimated Replacement Time: 1 minute

Procedure

Unplug the device from the power outlet before beginning this process.

1. Unscrew and remove the vented decorative enclosure strip at the top of the rear panel and set aside. (See two screw locations above.)

2. Completely loosen the locking screw on the bayonet latch. The locking screw is a captive screw, so it will not fall out when completely loosened.

3. Pull the bayonet latch back to disengage the top panel.

4. Lift the top panel (tabletop) and set the holding bar to hold the top panel up.

5. Unscrew the blast shield and slide it up, out of its mounting slot. The thumbscrew attached to the blast shield is a captive screw, so it will not fall off of the blast shield.

6. From both sides, look at a light source through the window at an angle to visualize the window's surfaces. Any debris, soiling, or damage to the optical coating should be evident. The blast shield can be inspected (for soiling or debris), and, if damaged, should be replaced.

7. After replacing the blast shield, close the top panel. The blast shield can only be inserted and the thumbscrew attached in one orientation.
8. Engage the bayonet latch by pushing it in.
9. Loosely tighten the bayonet latch locking screw.
10. Reattach the decorative strip at the top of the rear panel.

Troubleshooting

Trouble	Possible Cause	Solution
Key switch does nothing	- Power cord is not connected - Main line breaker is off - Emergency stop button is pressed	- Connect the power cord - Check the MAINS breaker - Release the emergency stop button
No correspondence between displayed energy and effects on skin	- Damaged fiber - Damaged blast shield	- Change the fiber - Change the blast shield

Error Messages

Error Message	Possible Cause	Solution
Low coolant flow	Not enough water	Fill the water reservoir
Over temperature	- Room temperature is too high - Airflow is obstructed - Cooling system component failure	- Lower the room temperature - Ensure nothing is obstructing the fan - Call for service
Shutter failure	Component failure	Call for service
Simmer failure	Flash lamp not functional	Call for service
No fiber attached	- Incorrect installation of the fiber - Wrong fiber is installed	- Check the proper connection of the fiber - Install suitable Rhapsody H-30 Holmium Laser System fiber
Remote open	- The remote interlock is not connected or is opened - The remote interlock is damaged	- Check the remote interlock connection - Call for service
Pedal not connected	- The foot switch is not connected - The foot switch is damaged	- Connect the foot switch properly - Call for service
No energy	Optical component damage	Call for service
Energy < 80%*	Optical component damage	If the energy error persists, call for service*
Energy > 120%*	Incorrect energy calibration	If the energy error persists, call for service*
Energy > 200%	Incorrect energy calibration	Call for service

Error Messages

Error Message	Possible Cause	Solution
Maintenance due	Preventative maintenance is due	Call for service*
MPS charge error	Main power supply malfunction	Call for service
MPS over temp	Main power supply malfunction	Call for service
MPS overload	Main power supply malfunction	Call for service

**Minor error. Schedule service, but it is okay to continue using the laser system.*

Chapter 7

Warranty Information

Warranty Statement

For warranty and service information, please refer to your authorized local representative.

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© COOK 2019, T_H-30_REV7
Bloomington, IN 47404, U.S.A