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Fig. 5


Fig. 6

TORNADO® EMBOLIZATION COILS AND MICROCOILS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Tornado Embolization Coils are available with small end exiting first (standard) or large end exiting first ("-LEF" suffix). Coil diameters of .018, .035 and .038 inch (0.46, 0.89 and 0.97 mm) are available. The embolization coil is supplied preloaded in the loading cartridge.

Longer coil length and tornado-like configuration in deployed state maximize coil exposure to cross-section of lumen for disruption of blood flow. Special platinum coil used for construction is soft, easily detected radiographically and features spaced synthetic fibers to maximize thrombogenicity.

INTENDED USE

Tornado Embolization Coils are intended for embolization of selective vessel supply to arteriovenous malformations and other vascular lesions. Tornado Embolization Coils are ideally suited for tapering vessel situations. The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known.

WARNINGS

- Positioning of Embolization Coils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- Tornado Embolization Coils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.
- If difficulties occur when deploying the embolization coil, repeat coil release procedure. If difficulties persist, withdraw the wire guide, coils and angiographic catheter simultaneously as a unit.

PRECAUTIONS

- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.

PRODUCT RECOMMENDATIONS

For .018 Inch Microcoils

- Tornado Embolization Microcoils are recommended for use through catheters that accept .018 inch (0.46 mm) maximum diameter wire guides. Catheter inner diameter (ID) must not exceed .027 inch (0.69 mm). Examples: 3.0 French nylon catheters with radiopaque tips (SCR3.0B); Cook MicroFerret™ Superselective Catheter (available outside U.S.).
- Wire guides recommended for positioning Tornado Embolization Microcoils are TFE-coated .018 inch (0.46 mm) diameter with flexible tapered tips.

For .035 and .038 Inch Coils

- .035 and .038 inch (0.89 mm and 0.97 mm) diameter, TFE-coated wire guides with flexible tapered tips are recommended for positioning Tornado Embolization Coils. Examples: TSFNA, TSFNB, TSFB and TSFBP.

Coil Size Diameter	Wire Guide Type & Size	Catheter Type & Size
.018 inch	TSF-18 STF-18	MF-3.0-18* SCR3.0B
.035 inch	TSFNB-35 TSFNA-35 TSFB-35 TSFBP-35 TSF-35	HNB(R)4.1-35 HNB(R)5.0-35 HNB6.0-35 SCBR5.0-35
.038 inch	TSFNB-38 TSFNA-38 TSFB-38 TSF-38	HNB(R)5.0-38 HNB6.0-38 SCBR5.0-38

*Available outside U.S.

INSTRUCTIONS FOR USE

For .018 Inch Microcoils:

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Push the loading cannula completely into the delivery catheter. (**Fig. 1**)
3. Advance the Luer lock connector fitting toward the catheter hub and lock into place. (**Fig. 2**)
4. Use the pusher stylet to load the Tornado Embolization Microcoil into the delivery catheter. **NOTE:** The stylet must be pushed as far as possible into the loading cannula to ensure proper loading. Remove the pusher stylet and loading cannula. (**Fig. 3**)
5. To obtain secure placement of the Tornado Embolization Microcoil, we recommend the use of the pusher stylet to push the Microcoil further into the delivery catheter. (**Fig. 4**)
6. Once coil is positioned into delivery catheter, verify position of catheter tip prior to deployment. Coil design permits delivery into the target vessel by saline/contrast flush or by push technique using an appropriately sized wire guide or pusher.
7. Perform final angiogram to confirm coil position within target vessel.

For .035 and .038 Inch Coils:

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp the loading cartridge between thumb and forefinger. Introduce the metal end of the loading cartridge into the base of the catheter hub. Lock loading cartridge onto catheter hub by turning Luer lock adapter clockwise. (**Fig. 5**)
3. Maintaining position of the cartridge, advance the stiff portion of the wire guide into the loading cannula. Push the coil into the first 20 to 30 centimeters of the angiographic catheter. (**Fig. 6**) Remove the wire guide and loading cartridge.
4. With the flexible tip of the wire guide, advance the embolization coil to the tip of the catheter. Verify position of the angiographic catheter prior to deployment.
5. Deploy the coil by advancing the wire guide past the tip of the catheter.
6. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd