EN 25	Zenith Alpha™ Abdominal Endovascular Graft Instructions for Use
CS 34	Břišní endovaskulární graft Zenith Alpha™ Návod k použití
DA 43	Zenith Alpha™ abdominal endovaskulær protese Brugsanvisning
DE 52	Zenith Alpha™ endovaskuläre Abdomenprothese Gebrauchsanweisung
EL 62	Κοιλιακό ενδαγγειακό μόσχευμα Zenith Alpha™ Οδηγίες χρήσης
ES 72	Endoprótesis vascular abdominal Zenith Alpha™ Instrucciones de uso
FR 82	Endoprothèse vasculaire abdominale Zenith Alpha™ Mode d'emploi
HU 92	Zenith Alpha™ hasi endovaszkuláris graft Használati utasítás
IT 102	Endoprotesi addominale Zenith Alpha™ Istruzioni per l'uso
NL 112	Zenith Alpha™ abdominale endovasculaire prothese Gebruiksaanwijzing
NO 122	Zenith Alpha™ abdominalt endovaskulært implantat Bruksanvisning
PL 131	Wewnątrznaczyniowy stent-graft brzuszny Zenith Alpha™ Instrukcja użycia
PT 141	Prótese endovascular abdominal Zenith Alpha™ Instruções de utilização
SV 151	Zenith Alpha™ abdominellt endovaskulärt graft Bruksanvisning

TABLE OF CONTENTS

Illustrations	11-24
1 DEVICE DESCRIPTION	25
1.1 Aortic Main Body and Iliac Leg Components	25
1.2 Main Body Delivery System	25
1.3 Iliac Leg Delivery System	25
1.4 Zenith® Low Profile AAA Endovascular Graft Ancillary Components	25
1.4.1 Main Body Extensions	25
1.4.2 Converters	25
1.4.3 Iliac Plugs	25
2 INDICATIONS FOR USE	25
3 CONTRAINDICATIONS	25
4 WARNINGS AND PRECAUTIONS	25
4.1 General	25
4.2 Patient Selection, Treatment and Follow-Up	25
4.3 Pre-Procedure Measurement Techniques and Imaging	26
4.4 Device Selection	26
4.5 Implant Procedure	26
4.6 Molding Balloon Use	27
4.7 MRI Information	27
5 POTENTIAL ADVERSE EVENTS	27
6 PATIENT SELECTION AND TREATMENT	27
6.1 Individualization of Treatment	27
7 PATIENT COUNSELING INFORMATION	28
8 HOW SUPPLIED	28
9 CLINICAL USE INFORMATION	28
9.1 Physician Training	28
9.2 Inspection Prior to Use	28
9.3 Materials Required	28
9.4 Materials Recommended	28
9.5 Device Diameter Sizing Guidelines	28
Table 9.5.1 Main Body, ZIMB (Zenith Alpha Abdominal Main Body), Graft Diameter Sizing Guide	29
Table 9.5.2 Iliac Leg, ZISL (Zenith Alpha Spiral-Z Endovascular Leg), Graft Diameter Sizing Guide	29
Table 9.5.3 Main Body Extension, ZLBE (Zenith Low Profile Body Extension), Graft Diameter Sizing Guide	29
Table 9.5.4 Converter, ZLC (Zenith Low Profile Converter), Graft Diameter Sizing Guide	29
10 DIRECTIONS FOR USE	29
Anatomical Requirements	29
General Use Information	29
Pre-Implant Determinants	29
Patient Preparation	30
10.1 Bifurcated System	30
10.1.1 Bifurcated Main Body Preparation/Flush	30
10.1.2 Iliac Leg Preparation/Flush	30
10.1.3 Vascular Access and Angiography	30
10.1.4 Main Body Placement	30
10.1.5 Contralateral Iliac Wire Guide Placement	30
10.1.6 Contralateral Iliac Leg Placement and Deployment	30
10.1.7 Main Body Distal (Bottom) Deployment	30
10.1.8 Ipsilateral Iliac Leg Placement and Deployment	31
10.1.9 Molding Balloon Insertion	31
Final Angiogram	31
10.2 Zenith® Low Profile AAA Endovascular Graft Ancillary Devices	31
General Use Information	31
10.2.1 Converter	31
10.2.2 Iliac Plugs	32
10.2.3 Main Body Extensions	32
11 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP	32
11.1 General	32
Table 11.1.1 Recommended Imaging Follow-up Schedule	32
11.2 Contrast and Non-Contrast CT Recommendations	32
Table 11.2.1 Acceptable Imaging Protocols	33
11.3 Abdominal Radiographs	33
11.4 MRI Information	33
11.5 Additional Surveillance and Treatment	33
12 MAIN BODY RELEASE TROUBLESHOOTING	33

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




45

ZENITH ALPHA™ ABDOMINAL ENDOVASCULAR GRAFT

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CAUTION: All contents of the outer pouch (including the introduction system and endovascular grafts) are supplied sterile, for single use only.

For the Zenith product line, there are several applicable Suggested Instructions for Use (IFU). This IFU describes the Suggested Instructions for Use for the Zenith Alpha Abdominal Endovascular Graft and Zenith Alpha Spiral-Z Endovascular Leg (main body, iliac legs, main body extension and converter). For information on other Zenith components, please refer to the following Suggested Instructions for Use:

- Zenith® Low Profile Endovascular Graft;
- Zenith® Branch Endovascular Graft-Iliac Bifurcation;
- Zenith® Renu™ AAA Ancillary Graft;
- Zenith Alpha™ Spiral-Z® Endovascular Leg;
- Zenith® Spiral-Z AAA Iliac Leg with the Z-Trak™ Introduction System;
- Zenith® AAA Ancillary Components;
- Zenith AAA Iliac Plug; and
- Coda® Balloon Catheter.

1 DEVICE DESCRIPTION

1.1 Aortic Main Body and Iliac Leg Components

The Zenith Alpha Abdominal Endovascular Graft is a modular system consisting of three components: a bifurcated aortic main body and two iliac legs. (Fig. 1) The graft modules are constructed of woven polyester fabric sewn to self-expanding nitinol stents with braided polyester and monofilament polypropylene suture. The modules are fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Additionally, the stents provide the necessary attachment and seal of the graft to the vessel wall.

The bare suprarenal stent at the proximal end of the graft contains barbs for additional fixation of the device. To facilitate fluoroscopic visualization of the stent graft, gold radiopaque markers are positioned on both the main body and iliac leg components. On the main body, four gold markers are positioned in a circumferential orientation within 2 mm of the most superior aspect of the graft material, a checkmark (✓) is positioned on the lateral aspect of the most distal stent of the contralateral limb, and two gold markers are positioned on the most distal stent of the ipsilateral limb. A marker also is positioned at the bifurcation of the graft.

On each iliac leg graft there are two gold markers on the proximal end, one showing minimum overlap with the main body on the contralateral side (16 mm) and one showing maximum overlap with the main body on the contralateral side (32 mm). The 16 and 32 mm gold markers are used for visualization and alignment on the ipsilateral side. A gold marker is positioned at the most distal end of the iliac leg component. Additionally, if the graft diameter increases distally, a gold marker is positioned where the graft begins to increase in diameter.

1.2 Main Body Delivery System

The Zenith Alpha Abdominal Endovascular Graft main body is shipped preloaded onto an introduction system. (Fig. 2) It has a sequential deployment method with built-in features to provide continuous control of the endovascular graft throughout the deployment procedure. The introduction system is designed for precise positioning and allows readjustment of the final graft position before deployment of the barbed suprarenal stent.

The main body graft is loaded on a 16 French (6.0 mm OD) or a 17 French (6.5 mm OD) introduction system. Dual locking mechanisms hold the endovascular graft onto the delivery system until released by the physician. All systems are compatible with a .035 inch wire guide.

The main body graft delivery system features a Flexor introducer sheath with a Captor Hemostatic Valve. For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The Flexor introducer sheath resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.3 Iliac Leg Delivery System

The Zenith Alpha Spiral-Z Endovascular Legs are shipped preloaded onto an introduction system. (Fig. 3) The introduction system is designed for ease of use with minimal preparation and uses a 12 French (4.7 mm OD) or 14 French (5.3 mm OD) introducer. All systems are compatible with a .035 inch wire guide.

For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. In addition, the delivery system features a Flexor introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.4 Zenith Low Profile AAA Endovascular Graft Ancillary Components

Zenith Low Profile ancillary components are compatible with the Zenith Low Profile AAA Endovascular Graft, as well as the Zenith Alpha Abdominal Endovascular Graft.

Additional ancillary endovascular components (main body extensions, converters and iliac plugs) are available. (Fig. 4) Ancillary components are constructed from woven polyester fabric, self-expanding nitinol stents (main body extensions and converters), self-expanding stainless steel Cook-Z® stents (iliac plugs), and braided polyester and polypropylene suture.

The aortic main body extensions can be used to provide additional length to the proximal portion of the endovascular graft. The converters and iliac plugs can be used to convert a bifurcated graft into an AUI graft, if necessary

(e.g., cases of Type III endoleak, limb occlusion or unattainable contralateral limb cannulation). Additional iliac leg components can also be used as ancillary components to provide additional length to the distal iliac portion of the graft system.

1.4.1 Main Body Extensions

Main body extensions utilize 16 French (6.0 mm OD) or 17 French (6.5 mm OD) introduction systems. (Fig. 40) The main body extension introduction system contains a single trigger-wire release mechanism. Deployment of the main body extension is achieved by sheath retraction and removing the distal trigger-wire.

To facilitate fluoroscopic visualization of the stent grafts, the main body extensions have four gold radiopaque markers positioned in a circumferential orientation within 2 mm of the most superior aspect of the graft material.

For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The main body extension delivery system features a Flexor introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.4.2 Converters

Converters utilize 16 French (6.0 mm OD) or 17 French (6.5 mm OD) introduction systems. (Fig. 30) The converter introduction system contains a single trigger-wire release mechanism. Deployment of the converter is achieved by sheath retraction and removing the distal trigger-wire.

The converter has a single gold radiopaque marker positioned within 2 mm of the most superior aspect of the graft material. For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The converter delivery system features a Flexor introducer sheath, which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.4.3 Iliac Plugs

The iliac plug is used to occlude an iliac artery, typically in conjunction with a femoral-to-femoral crossover procedure. For more information on this component, please refer to the Zenith AAA Endovascular Graft Ancillary Components Suggested Instructions for Use.

2 INDICATIONS FOR USE

The Zenith Alpha Abdominal Endovascular Graft and Zenith Low Profile ancillary components are indicated for the endovascular treatment of patients with abdominal aortic or aorto-iliac aneurysms having morphology suitable for endovascular repair, including:

- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with a 16 French (6.0 mm OD) or 17 French (6.5 mm OD) vascular introducer sheath;
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a length of at least 15 mm,
 - with a diameter measured outer wall to outer wall of no greater than 32 mm and no less than 18 mm,
 - with an angle less than 60 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal aorta.
- Iliac artery distal fixation site greater than 10 mm in length and 8 to 20 mm in diameter (measured outer wall to outer wall).

3 CONTRAINDICATIONS

The Zenith Alpha Abdominal Endovascular Graft is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, nitinol, polyester, polypropylene, urethane, PTFE, nylon or gold.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith Alpha Abdominal Endovascular Graft and introduction system should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 9.1, Physician Training**.
- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- When a 42 or 59 mm leg graft is used on the ipsilateral side, contralateral leg overlap into the contralateral main body limb should be limited to 16 mm. Failure to do so may result in occlusion of the ipsilateral limb.

4.2 Patient Selection, Treatment and Follow-Up

- The Zenith Alpha Abdominal Endovascular Graft is designed to treat aortic neck diameters no smaller than 18 mm and no larger than 32 mm. The Zenith Alpha Abdominal Endovascular Graft is designed to treat proximal aortic necks (distal to the lowest renal artery) of at least 15 mm in length. Iliac artery distal fixation site greater than 10 mm in

length and 8 to 20 mm in diameter (measured outer wall to outer wall) is required. These sizing measurements are critical to the performance of the endovascular repair.

- Key anatomical elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (>60 degrees for infrarenal neck to axis of AAA or >45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic neck (<15 mm); an inverted funnel shape (greater than 10% increase in diameter over 15 mm of proximal aortic neck length); and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. In the presence of anatomical limitations, a longer neck may be required to obtain adequate sealing and fixation. Irregular calcification and/or plaque may compromise the attachment and sealing at the fixation sites. Necks exhibiting these key anatomical elements may be more conducive to graft migration or endoleak.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 16 French (6.0 mm OD) or 17 French (6.5 mm OD) vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization, graft kinking, or thromboses. A vascular conduit technique may be necessary to achieve success in some patients.
- Pre-existing regions of stenosis/narrowing (less than approximately 20 mm ID in the aorta or 7 to 8 mm ID in the iliacs) have been shown to increase the risk of a thromboembolic event (e.g., graft limb occlusion). The potential for this increased risk in these patients may preclude placement of an endovascular graft. Dilatation of these regions with a noncompliant balloon and/or stent placement may be necessary to help assure maintained graft patency and to reduce the risk of a thromboembolic event. Additionally, the completion angiogram (with stiff wire guides removed) should be reviewed carefully to determine if further treatment in these regions is necessary (e.g., adjunctive ballooning or stenting). Failure to remove the stiff wire guide prior to the angiogram could mask any limb kinking or narrowing that might occur when the wire guide is removed.
- Follow-up imaging should be carefully reviewed for narrowing within the graft leg. Patients with a graft leg lumen of less than approximately 5 mm ID may be at increased risk of a thromboembolic event (e.g., graft limb occlusion). Reintervention (e.g., noncompliant ballooning or stenting in these regions) should be considered to help assure maintained graft patency and to reduce the risk of a thromboembolic event.
- Patients with poor outflow or a hypercoagulable state (e.g., cancer) may be at an increased risk of a thromboembolic event.
- The Zenith Alpha Abdominal Endovascular Graft is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith Alpha Abdominal Endovascular Graft is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The Zenith AAA family of grafts has not been formally tested in the following patient populations:
 - traumatic aortic injury
 - leaking, pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed endovascular grafts
 - uncorrectable coagulopathy
 - indispensable mesenteric artery
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - patients less than 18 years of age
 - morbidly obese patients
 - patients with less than 15 mm in length or greater than 60 degrees angulation of the proximal aortic neck relative to the long axis of the aneurysm
 - patients with two occluded internal iliac arteries
- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3 Pre-Procedure Measurement Techniques and Imaging**.
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in suboptimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith Alpha Abdominal Endovascular Graft. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning the x-ray C-arm during procedural

angiography such that the origins of the renal arteries, and particularly the lowest patent renal artery, are well demonstrated prior to deployment of the proximal edge of the graft material (sealing stent) of the main body. Additionally, angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are well defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg components.

Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

Lengths

Utilizing CT, length measurements should be determined to accurately assess infrarenal proximal neck length as well as planning main body sizes and leg components for the Zenith Alpha Abdominal Endovascular Graft. These reconstructions should be performed in sagittal, coronal, and 3-D.

- **The long-term performance of this endovascular graft has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm) or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- The Zenith Alpha Abdominal Endovascular Graft is not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

4.4 Device Selection

Strict adherence to the Zenith Alpha Abdominal Endovascular Graft IFU sizing guide is strongly recommended when selecting the appropriate device size (**Tables 9.5.1 thru 9.5.4**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

4.5 Implant Procedure

(Refer to **Section 10, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith Alpha Abdominal Endovascular Graft and assure accurate apposition to the aortic wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Alpha Abdominal Endovascular Graft.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- To avoid damage to the sheath, be careful to advance all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith Alpha Abdominal Endovascular Graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (exception is the inferior mesenteric artery) with the endoprosthesis. Vessel occlusion may occur.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Repositioning the stent graft distally after partial deployment of the covered proximal stent may result in damage to the stent graft and/or vessel injury.
- Inaccurate placement and/or incomplete sealing of the Zenith Alpha Abdominal Endovascular Graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Zenith Alpha Abdominal Endovascular Graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Inadequate overlap of the Zenith Alpha Spiral-Z Endovascular Leg may result in increased risk of migration of the stent graft and subsequent endoleak.
- Systemic anticoagulation should be used during the implantation procedure based on hospital- and physician-preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Fluoroscopy should be used during introduction and deployment to

confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.

- The use of the Zenith Alpha Abdominal Endovascular Graft requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise.
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- The Zenith Alpha Abdominal Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional and angiographic devices in the region of the suprarenal stent.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization or rupture of the aneurysm.
- Avoid damaging the graft or disturbing graft positioning after placement in the event reinstrumentation (secondary intervention) of the graft is necessary.
- Before deployment of the suprarenal stent, verify that the position of the access wire guide extends just distal to the aortic arch.
- Excessive overlap of 12 mm above the main body graft bifurcation may increase the risk of limb thrombosis.

4.6 Molding Balloon Use

- Do not inflate the balloon in vessel outside of the graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care when inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.
- For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

Main Body Extension and Converter

- Care should be taken not to displace the main body graft during the placement and deployment of the main body extension or converter.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.

4.7 MRI Information

NOTE: If using this device in conjunction with another endovascular graft from the Zenith family, refer to the appropriate device's IFU for additional MRI information.

Nonclinical testing has demonstrated that the Zenith Low Profile AAA Endovascular Graft (The Zenith Low Profile AAA Endovascular Graft served as a surrogate for the Zenith Alpha Abdominal Endovascular Graft) is MR Conditional. A patient with this endovascular graft can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less.
- Spatial magnetic gradient field of 1580 Gauss/cm (15.8 T/m) or less.
- The product of the spatial gradient and static magnetic field should not exceed $47.4 T^2/m$.

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 and 3.0 Tesla Systems: Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise

In nonclinical testing, the Zenith Low Profile AAA Endovascular Graft (The Zenith Low Profile AAA Endovascular Graft served as a surrogate for the Zenith Alpha Abdominal Endovascular Graft) produced a maximum temperature rise of 1.7 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Magnetom, Siemens Medical Solutions, Malvern, PA, Software Numaris/4) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise

In nonclinical testing, the Zenith Low Profile AAA Endovascular Graft (The Zenith Low Profile AAA Endovascular Graft served as a surrogate for the Zenith Alpha Abdominal Endovascular Graft) produced a maximum temperature rise of 2.0 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (Excite, GE Healthcare, Milwaukee, WI, Software G3.0-052B) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.8 W/kg).

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 5 mm of the position of the Zenith Alpha Abdominal Endovascular Graft, as found during nonclinical testing using the sequences: T1-weighted spin echo and gradient echo pulse in a 3.0 Tesla MR system (Excite, General Electric Healthcare). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US
Fax: 209-669-2450
Web: www.medicalert.org

5 POTENTIAL ADVERSE EVENTS

Adverse events associated with either the Zenith Alpha Abdominal Endovascular Graft or the implantation procedure that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion; component separation
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

Device-Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith Alpha Abdominal Endovascular Graft should be reported to Cook immediately. For customers inside the United States, to report an incident, call the Customer Relations Department at 1-800-457-4500 (24 hour) or 1-812-339-2235. For customers outside the United States, please call your distributor.

6 PATIENT SELECTION AND TREATMENT

(Refer to **Section 4, WARNINGS AND PRECAUTIONS**)

6.1 Individualization of Treatment

Cook recommends that the Zenith Alpha Abdominal Endovascular Graft component diameters be selected as described in **Tables 9.5.1, 9.5.2, 9.5.3 and 9.5.4** within **Section 9, CLINICAL USE INFORMATION**. The length of the Zenith Alpha Abdominal Endovascular Graft should extend from the lowest renal artery to just above the internal iliac (hypogastric) artery bifurcation. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy
- Comorbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Alpha Abdominal Endovascular Graft
- Patient's ability to tolerate general, regional or local anesthesia
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with a 16 French (6.0 mm OD) or 17 French (6.5 mm OD) vascular introducer sheath.
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a length of at least 15 mm,
 - with a diameter measured outer wall to outer wall of no greater than 32 mm and no less than 18 mm,
 - with an angle less than 60 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal

aorta.

- Iliac artery distal fixation site greater than 10 mm in length and 8 to 20 mm in diameter (measured outer wall to outer wall).
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

The final treatment decision is at the discretion of the physician and patient.

7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment to and compliance with postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of this endovascular graft has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.**
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements is required and should be considered a lifelong commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise every patient that it is important to seek prompt medical attention if he/she experiences signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant.
- Men who undergo endovascular or open surgical repair may experience impotence.

Physicians should refer patients to the *Patient Guide* regarding risks occurring during or after implantation of the device. Procedure-related risks include cardiac, pulmonary, neurologic, bowel and bleeding complications. Device-related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See **Section 5, POTENTIAL ADVERSE EVENTS**). The physician should complete the *Patient ID Card* and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

8 HOW SUPPLIED

- The Zenith Alpha Spiral-Z Endovascular Leg is sterilized by ethylene oxide gas, is preloaded into an introduction system, and is supplied in peel-open packages.
- The device is intended for single use only. Do not re-sterilize the device.
- The product is sterile if the package is unopened and undamaged. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- The main body grafts, main body extensions and converters are loaded in 16 French (6.0 mm OD) or 17 French (6.5 mm OD) Flexor introducer sheaths. The iliac leg devices are loaded into 12 French (4.7 mm OD) or 14 French (5.3 mm OD) Flexor introducer sheaths. The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a sterile gauze pad soaked in saline solution under sterile

conditions.

- Do not use after the expiration date printed on the label.
- Store in a cool, dry place.

9 CLINICAL USE INFORMATION

9.1 Physician Training

CAUTION: Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Alpha Abdominal Endovascular Graft should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Alpha Abdominal Endovascular Graft are outlined below.

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and comorbidities associated with AAA repair
- Knowledge of radiographic image interpretation, device selection, planning and sizing

A multidisciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Nonselective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

9.3 Materials Required

(Not included in three-piece modular system)

- Zenith Low Profile AAA Endovascular Graft Ancillary Kit
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

9.4 Materials Recommended

The following products are recommended for implantation of any component in the Zenith product line. For information on the use of these products, refer to the individual product's Suggested Instructions for Use:

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm;
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide;
 - Cook .035 inch wire guides
 - Cook Nimble™ Wire Guides
- Molding Balloons;
 - Cook Coda® Balloon Catheter
- Introducer sets;
 - Cook Check-Flo® Introducer Sets
 - Cook Extra Large Check-Flo® Introducer Sets
 - Cook Flexor® Balkin Up & Over® Contralateral Introducers
- Sizing catheter;
 - Cook Aurous® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters;
 - Cook Beacon® Tip Angiographic Catheters
 - Cook Beacon® Tip Royal Flush Catheters
- Entry needles;
 - Cook single wall entry needles

9.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Table 9.5.1 Main Body, ZIMB (Zenith Alpha Abdominal Main Body), Graft Diameter Sizing Guide*

Intended Aortic Vessel Diameter ^{1,2} (mm)	Main Body Diameter ³ (mm)	Overall Length to Contralateral Limb/ Overall Length to Ipsilateral Limb (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
18-19	22	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	16	5.3, 6.0
20-21	24	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	16	5.3, 6.0
22	26	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	16	5.3, 6.0
23-24	28	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	16	5.3, 6.0
25-26	30	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	16	5.3, 6.0
27-28	32	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	16	5.3, 6.0
28-32	36	70/94, 84/108, 98/122, 108/132, 118/142, 128/152	17	5.6, 6.5

¹ Maximum diameter along the proximal fixation site.

² Round measured aortic diameter to nearest mm.

³ Additional considerations may affect choice of diameter.

* All dimensions are nominal.

Table 9.5.2 Iliac Leg, ZISL (Zenith Alpha Spiral-Z Endovascular Leg), Graft Diameter Sizing Guide*

Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Leg Diameter ³ (mm)	Iliac Leg Label Length ⁴ (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
8	9	42, 59, 77, 93, 110, 125	12	4.0, 4.7
9-10	11	42, 59, 77, 93, 110, 125	12	4.0, 4.7
11-12	13	42, 59, 77, 93, 110, 125	12	4.0, 4.7
13-15	16	42, 59, 77, 93	12	4.0, 4.7
16-18	20	42, 59, 77, 93	14	4.7, 5.3
19-20	24	42, 59, 77, 93	14	4.7, 5.3

¹ Maximum diameter along the distal fixation site.

² Round measured iliac diameter to nearest mm.

³ Additional considerations may affect choice of diameter.

⁴ Approximate overall leg length = label length + 26 mm.

* All dimensions are nominal.

Table 9.5.3 Main Body Extension, ZLBE (Zenith Low Profile Body Extension), Graft Diameter Sizing Guide*

Intended Aortic Vessel Diameter ^{1,2} (mm)	Main Body Extension Diameter ³ (mm)	Main Body Extension Length (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
18-19	22	45, 58	16	5.5, 6.0
20-21	24	45, 58	16	5.5, 6.0
22	26	45, 58	16	5.5, 6.0
23-24	28	45, 58	16	5.5, 6.0
25-26	30	45, 58	16	5.5, 6.0
27-28	32	45, 58	16	5.5, 6.0
29-32	36	45, 58	17	5.6, 6.5

¹ Maximum diameter along the proximal fixation site.

² Round measured aortic diameter to nearest mm.

³ Additional considerations may affect choice of diameter.

* All dimensions are nominal.

Table 9.5.4 Converter, ZLC (Zenith Low Profile Converter), Graft Diameter Sizing Guide*

Main Body Diameter (mm)	Converter Diameter ¹ (mm)	Converter Length (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
22 & 24	24	66	16	5.5, 6.0
26 & 28	28	66	16	5.5, 6.0
30 & 32	32	66	16	5.5, 6.0
36	36	66	17	5.6, 6.5

¹ Additional considerations may affect choice of diameter.

* All dimensions are nominal.

10 DIRECTIONS FOR USE

Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.
- Proximal aortic neck lengths should be a minimum of 15 mm with a diameter measured outer wall to outer wall of 18 to 32 mm.
- Iliac artery distal fixation site should be greater than 10 mm in length and 8 to 20 mm in diameter (measured outer wall to outer wall).

Prior to use of the Zenith Alpha Abdominal Endovascular Graft, review this Suggested Instructions for Use booklet. The following instructions embody guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

General Use Information

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Alpha Abdominal Endovascular Graft. The Zenith Alpha Abdominal Endovascular Graft is compatible with .035 inch diameter wire guides.

Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

- Femoral artery selection for introduction of the main body system (i.e., define respective contralateral and ipsilateral iliac arteries).
- Angulation of the aortic neck, aneurysm and iliac arteries.
- Quality of the aortic neck.
- Diameters of the infrarenal aortic neck and distal iliac arteries.
- Distance from the renal arteries to the aortic bifurcation.
- Length from the aortic bifurcation to the internal iliac arteries/attachment site(s).
- Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.

8. Consider the degree of vascular calcification, stenosis and narrowing.

Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position the patient on the imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Both common femoral arteries should be prepared using standard techniques for either surgical or percutaneous access.

10.1 Bifurcated System (Fig. 2)

10.1.1 Bifurcated Main Body Preparation/Flush

1. Verify that the Captor Sleeve is inserted in the Captor Hemostatic Valve. Elevate the distal tip of the system and flush through the stopcock on the Captor Hemostatic Valve until fluid emerges through the flush groove at the proximal end of the introducer sheath. (Fig. 5a) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close the stopcock on the connecting tube.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach a syringe with heparinized saline to the clear hub on the end of the handle. Flush until fluid exits the dilator tip. (Fig. 6)

NOTE: When flushing the system, elevate the distal end of the system to facilitate the removal of air.

3. Soak sterile gauze pads in saline solution and use them to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both the sheath and dilator liberally.

10.1.2 Iliac Leg Preparation/Flush

1. Remove the Peel-Away sheath from the back of the hemostatic valve. (Fig. 8) Elevate the distal tip of the system and flush through the stopcock on the hemostatic valve until fluid emerges through the flush groove at the proximal end of the introducer sheath. (Fig. 5b) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close the stopcock on the connecting tube.

NOTE: Graft flushing solution of heparinized saline is often used.

NOTE: When removing the Peel-Away sheath from the back of the hemostatic valve, ensure that the delivery system sheath is held stationary against the dilator tip to limit possible movement.

2. Attach a syringe with heparinized saline to the black hub on the inner cannula. Flush until fluid exits the dilator tip. (Fig. 7)

NOTE: When flushing the system, elevate the distal end of the system to facilitate the removal of air.

10.1.3 Vascular Access and Angiography

1. Puncture the selected common femoral arteries using standard technique with an 18 or 19UT gage arterial needle. Upon vessel entry, insert:

- Wire guides - standard .035 inch diameter, 145 cm long
- Appropriate size sheaths (e.g., 6 or 8 French)
- Flush catheter (often radiopaque sizing catheters - e.g., Centimeter Sizing Catheter or straight flush catheter)

2. Perform angiography to identify level(s) of renals, aortic bifurcation and iliac bifurcations.

NOTE: If fluoroscopy angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

10.1.4 Main Body Placement

1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.

2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

NOTE: Monitor the patient's coagulation status throughout the procedure.

3. On the ipsilateral side, replace the J wire with a stiff wire guide (LES) .035 inch, 260 cm long, and advance through the catheter and up to the thoracic aorta. Remove the flush catheter and sheath. Maintain wire guide position.
4. Before insertion, position the main body delivery system on the patient's abdomen under fluoroscopy to determine the orientation of the contralateral limb radiopaque marker. The sidearm of the hemostatic valve may serve as an external reference to the contralateral limb radiopaque marker.
5. Insert the main body delivery system over the wire into the femoral artery, with attention to sidearm reference.

CAUTION: Maintain wire guide position during delivery system insertion.

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

6. Advance the delivery system until the four gold radiopaque markers (which are positioned 2 mm from the most proximal segment of the graft material) (Fig. 9, Illustration 1) are just inferior to the most inferior renal orifice.

7. Verify position of the wire guide in the thoracic aorta. Ensure the graft system is oriented such that the contralateral limb is positioned above and anterior to the origin of the contralateral iliac. If the contralateral limb radiopaque marker is not properly aligned, rotate the entire system until it is correctly positioned half way between a lateral and an anterior position on the contralateral side.

- A marker formation of a ✓ indicates an anterior position of the short (contralateral) limb. (Fig. 9, Illustration 4)
- A marker formation of a ∨ indicates a posterior position of the short (contralateral) limb. (Fig. 9, Illustration 5)
- A marker formation of a | line indicates a lateral position of the short (contralateral) limb. (Fig. 9, Illustration 6)

8. Repeat the angiogram to verify the four gold radiopaque markers are 2 mm or more below the most inferior renal orifice.

9. Ensure the Captor Hemostatic Valve is turned to the open position. (Fig. 10)

10. Stabilize the gray positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two covered stents by withdrawing the sheath while monitoring device location.

NOTE: The delivery system does not utilize a top cap; however, the stent graft has a suprarenal stent with barbs. The device should be accurately positioned before the outer sheath is withdrawn.

11. Without moving the table, decrease the magnification to check on the position of the contralateral limb radiopaque marker and location of renal arteries. Continue withdrawing the sheath until the contralateral limb is fully deployed. (Fig. 11) Stop withdrawing the sheath.

NOTE: Verify the contralateral limb is at least 5 mm above the aortic bifurcation and in the desired location for cannulation.

12. Repeat the angiogram and reposition if necessary.

13. While holding the black gripper, turn the black safety lock knob counterclockwise to engage the blue rotation handle. (Fig. 12)

NOTE: If the black safety lock knob is removed from the system after it has been turned counterclockwise, the blue rotation handle will remain engaged. Continue with the procedure.

14. Under fluoroscopy, turn the blue rotation handle in the direction of the arrow (clockwise) until a stop is felt. (Fig. 13)

NOTE: If the blue rotation handle stops before completing a full rotation, visually verify the position of the black safety lock knob, and if necessary, turn it to the unlock position.

NOTE: Handle system mechanism and safeties may be manually overridden; however, do not attempt to force the handle before first attempting all troubleshooting actions.

NOTE: Turning the rotation handle releases the suprarenal stent. If resistance is felt or system bowing is noticed, the device is under tension. Excessive force may cause the graft position to be altered. If excessive resistance or delivery system movement is noted, stop and assess the situation. If the stent does not fully release, see Section 12, MAIN BODY RELEASE TROUBLESHOOTING

CAUTION: During suprarenal stent deployment, verify that the position of the main body wire guide extends just distal to the aortic arch and that support of the system is maximized.

NOTE: Once the barbed suprarenal stent has been deployed, further attempts to reposition the graft are not recommended.

WARNING: The Zenith Alpha Abdominal Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

10.1.5 Contralateral Iliac Wire Guide Placement

1. Manipulate the catheter and wire guide through the open end of the contralateral limb into the body of the graft. Advance the wire guide inside the body of the graft and into the thoracic aorta. AP and oblique fluoroscopic views can aid in verification of device cannulation.
2. After cannulation, advance the angiographic catheter over the wire into the body of the endovascular graft. Remove the wire guide and perform angiography to confirm position. Reinsert the wire guide inside the body of the graft and into the thoracic aorta. Remove the angiographic catheter.

10.1.6 Contralateral Iliac Leg Placement and Deployment

NOTE: If using this device in conjunction with the Zenith Spiral-Z AAA Iliac Leg, refer to the Zenith Spiral-Z AAA Iliac Leg IFU for appropriate deployment and overlap instructions.

CAUTION: Verify that the contralateral iliac leg is selected.

NOTE: When using a 42 or 59 mm leg graft on the ipsilateral side, contralateral leg overlap into the contralateral main body limb should be limited to 16 mm.

1. Position the image intensifier to show both the contralateral internal iliac artery and contralateral common iliac artery.
2. Prior to introduction of the contralateral iliac leg delivery system, inject contrast through the contralateral femoral sheath to locate the contralateral internal iliac artery.
3. Remove the femoral sheath and introduce the contralateral iliac leg delivery system into the artery. Advance slowly until the second gold radiopaque marker on the iliac leg graft aligns with the gold checkmark of the main body graft, with 32 mm of overlap between components. (Fig. 14) If there is any tendency for the main body graft to move during this maneuver, hold it in position by stabilizing the positioner on the ipsilateral side.

NOTE: Radiopaque marker bands are positioned 16 mm from the proximal end of the iliac leg graft to identify the minimum amount of overlap and 32 mm from the proximal end of the iliac leg graft to identify the maximum amount of overlap.

NOTE: If difficulty is encountered advancing the iliac leg delivery system, exchange to a more supportive wire guide. In tortuous vessels, the anatomy may alter significantly with the introduction of the rigid wires and sheath systems.

4. Confirm position of the distal end of the contralateral iliac leg graft. Reposition the contralateral iliac leg graft if necessary to ensure both internal iliac patency and minimum overlap of 2 stents (16 mm) within the main body endovascular graft.
5. To deploy, hold the contralateral iliac leg graft in position with the gray positioner while withdrawing the sheath approximately 10 mm. (Figs. 15 and 16)
6. Check the graft position and reposition if necessary.
7. Continue to deploy the graft by withdrawing the sheath while continuously checking the position of the graft.
8. Stop withdrawing the sheath as soon as the distal end of the contralateral iliac leg graft is released.
9. Under fluoroscopy and after verification of iliac leg graft position, loosen the pin vise and retract the inner cannula to dock the tapered dilator to the positioner. Tighten the pin vise. Maintain sheath position while withdrawing the gray positioner with secured inner cannula. (Fig. 17)
10. Re-check the position of the wire guide.

10.1.7 Main Body Distal (Bottom) Deployment

1. Return to the ipsilateral side.
2. Fully deploy the ipsilateral limb of the main body by withdrawing the sheath until the most distal stent has expanded. (Fig. 18) Stop withdrawing the sheath.

NOTE: The distal stent is still secured to the delivery system.

3. While holding the black gripper, turn the gray safety lock knob counterclockwise to engage the blue rotation handle. (Fig. 19)

NOTE: If unable to release the gray safety lock knob, attempt to turn the blue rotation handle clockwise until the handle comes to a stop, then attempt to release the gray safety lock knob.

4. Turn the blue rotation handle in the direction of the arrow (clockwise) until a stop is felt. (Fig. 20)

NOTE: Handle system mechanism and safeties may be manually overridden; however, do not attempt to force the handle before first attempting all troubleshooting actions.

NOTE: Turning the rotation handle releases the distal stent of the main body graft. If the distal stent does not release, see **Section 12, MAIN BODY RELEASE TROUBLESHOOTING**.

5. Withdraw the entire inner subassembly through the sheath by holding the Captor Hemostatic Valve stationary and pulling on the black gripper. Confirm that the Captor sleeve is removed along with the inner sub assembly. (Fig. 21)

NOTE: Maintain position of sheath and wire guide.

6. Close the Captor Hemostatic Valve on the introducer sheath by turning it in a clockwise direction until it stops. (Fig. 22)

10.1.8 Ipsilateral Iliac Leg Placement and Deployment

NOTE: If using this device in conjunction with the Zenith Spiral-Z AAA Iliac Leg, refer to the Zenith Spiral-Z AAA Iliac Leg IFU for appropriate deployment and overlap instructions.

NOTE: Ensure the Captor Hemostatic Valve on the main body introducer sheath is turned to the open position. (Fig. 23)

1. Position the image intensifier to show both the ipsilateral internal iliac artery and ipsilateral common iliac artery.
2. Prior to the introduction of the ipsilateral iliac leg delivery system, inject contrast through the main body femoral sheath to locate the ipsilateral internal iliac artery.
3. Utilize the main body graft wire and sheath assembly to introduce the ipsilateral iliac leg graft. Advance the dilator and sheath assembly into the main body sheath.

NOTE: In tortuous vessels, the position of the internal iliac arteries may alter significantly with the introduction of the rigid wires and sheath systems.

4. Continue advancing slowly until the proximal edge of the ipsilateral leg graft aligns with the proximal edge of the previously-placed contralateral leg graft. (Figs. 24 and 25)
5. Confirm position of the distal end of the iliac leg graft. Using the distal gold marker as a reference point, reposition the iliac leg graft if necessary to ensure internal iliac patency.
6. To deploy, hold the iliac leg graft in position with the gray positioner while withdrawing the sheath approximately 10 mm. (Fig. 26)
7. Check the graft position and reposition if necessary.
8. Continue to deploy the graft by withdrawing the sheath while continuously checking the position of the graft. (Fig. 27)
9. Under fluoroscopy and after verification of iliac leg graft position, loosen the pin vise, and retract the inner cannula to dock the tapered dilator to the positioner. Tighten the pin vise. Maintain sheath position while withdrawing the positioner with secured inner cannula. (Fig. 28)
10. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops.
11. Re-check the position of the wire guides. Leave the sheath and wire guides in place.

10.1.9 Molding Balloon Insertion

1. Prepare the molding balloon as follows:
 - Flush the wire lumen with heparinized saline.
 - Remove all air from the balloon.
2. In preparation of the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counterclockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve of the main body introduction system to the level of the renal arteries. Maintain proper sheath positioning.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.

CAUTION: Do not inflate the balloon in the vessel outside of the graft.

5. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the most proximal covered stent and the infrarenal neck, starting proximally and working in the distal direction. (Fig. 29)

CAUTION: Confirm complete deflation of the balloon prior to repositioning.

CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.

6. Withdraw the molding balloon to the ipsilateral limb overlap region and expand.

CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.

7. Withdraw the molding balloon to the ipsilateral distal fixation site and expand.

CAUTION: Do not inflate the balloon in the vessel outside of the graft.

8. Deflate and remove the molding balloon. Transfer the molding balloon onto the contralateral wire guide and into the contralateral iliac leg introduction system. Advance the molding balloon to the contralateral limb overlap and expand.

CAUTION: Confirm complete deflation of the balloon prior to repositioning.

9. Withdraw the molding balloon to the contralateral iliac leg/vessel distal fixation site and expand. (Fig. 29)

CAUTION: Do not inflate the balloon in the vessel outside of the graft.

10. Remove the molding balloon and replace it with an angiographic catheter to perform completion angiograms.
11. Remove or replace all stiff wire guides to allow the iliac arteries to resume their natural position.

Final Angiogram

1. Position an angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of the internal iliac arteries.
 2. Confirm there are no endoleaks or kinks and verify position of the proximal gold radiopaque markers. Remove the sheaths, wires and catheters.
- NOTE:** If endoleaks or other problems are observed, refer to **Section 10.2, Zenith Low Profile AAA Endovascular Graft Ancillary Devices**.
3. Repair vessels and close in standard surgical fashion.

10.2 Zenith Low Profile AAA Endovascular Graft Ancillary Devices

General Use Information

Inaccuracies in device size selection or placement, changes or anomalies in patient anatomy or procedural complications can require placement of additional endovascular grafts, extensions, iliac plugs and converters. Regardless of the device placed, the basic procedure(s) will be similar to the maneuvers required and described previously in this document. It is vital to maintain wire guide access.

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Low Profile AAA Endovascular Graft ancillary devices. Low Profile AAA Endovascular Graft ancillary devices are compatible with .035 inch wire guides.

Low Profile AAA Endovascular Graft ancillary devices are compatible with the Zenith Low Profile AAA Endovascular Graft as well as the Zenith Alpha Abdominal Endovascular Graft.

10.2.1 Converter

Converters can be used to convert a bifurcated graft into an aortouni-iliac graft, if necessary (e.g., cases of Type III endoleak, limb occlusion or unattainable contralateral limb cannulation). (Fig. 30)

Converter Preparation/Flush

1. Remove the inner stylet (from the inner cannula), cannula protector (from the inner cannula) and dilator tip protector (from the dilator tip). Remove the Peel-Away sheath from the back of the hemostatic valve. (Fig. 31) Elevate the distal tip of the system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (Fig. 32) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach a syringe with heparinized saline to the hub of the inner cannula. Flush until fluid exits the dilator tip. (Fig. 33)

NOTE: When flushing the system, elevate the distal end of the system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use them to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both the sheath and the dilator liberally.

Converter Placement and Deployment

1. Remove the main body delivery sheath. Use the main body graft wire guide to introduce the converter into the main body.

NOTE: The converter delivery system cannot be introduced through the main body or iliac leg introducer sheath.

2. Advance slowly until the converter is at the site of the required intervention. (Fig. 34) Verify appropriate stent graft overlap to ensure proper sealing and resistance to migration. The proximal two stents should be positioned in the main body graft and the distal two stents should be positioned in the ipsilateral leg.
3. Deploy the device by withdrawing the sheath while stabilizing the gray positioner of the delivery system. (Figs. 35 and 36)
4. Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle, and then remove via the slot over the inner cannula. (Fig. 37)
5. Continue to deploy the device until the distal stent is uncovered.
6. Withdraw the tapered tip of the introducer back through the converter graft and delivery system while maintaining wire guide position. Ensure the endovascular graft is not displaced during the withdrawal of the delivery system.
7. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops. (Fig. 38)

Converter Molding Balloon Insertion

NOTE: For information on the use of recommended products, refer to the individual product's Instructions for Use.

1. Prepare the molding balloon as follows:
 - Flush the wire lumen with heparinized saline.
 - Remove all air from the balloon.

CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.

2. Advance the molding balloon over the wire guide and through the hemostatic valve to the proximal segment of the converter.
3. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.

CAUTION: Do not inflate the balloon in the vessel outside of the graft.

4. Expand the molding balloon within the proximal segment and the distal segment of the converter using diluted contrast media (as recommended by the manufacturer). (Fig. 39)

CAUTION: Confirm complete deflation of the balloon prior to repositioning.

5. Completely deflate and remove the molding balloon, replace it with an angiographic catheter and perform completion angiograms.
6. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

10.2.2 Iliac Plugs

Refer to the Zenith AAA Endovascular Graft Ancillary Components or the Zenith AAA Iliac Plug IFU.

10.2.3 Main Body Extensions

Main body extensions are used for extending the proximal body of an *in situ* endovascular graft. (Fig. 40)

Main Body Extension Preparation/Flush

1. Remove the inner stylet (from the inner cannula), cannula protector (from the inner cannula) and dilator tip protector (from the dilator tip). Remove the Peel-Away sheath from the back of the hemostatic valve. (Fig. 31) Elevate the distal tip of the system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (Fig. 32) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close the stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach a syringe with heparinized saline to the hub of the inner cannula. Flush until fluid exits the dilator tip. (Fig. 33)

NOTE: When flushing system, elevate distal tip of the system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use them to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both the sheath and the dilator liberally.

Main Body Extension Placement and Deployment

1. Remove the main body delivery sheath. Use the main body graft wire guide to introduce the main body extension into the main body.

NOTE: The main body extension delivery system cannot be introduced through the main body or iliac leg introducer sheath.

2. Advance slowly until the main body extension is at the site of the required intervention. (Fig. 41)
3. Verify the main body extension position to ensure proper sealing and resistance to migration.
4. Verify the placement with angiography to ensure that the renal arteries remain patent and that proper placement is achieved.

CAUTION: Care should be taken not to displace the main body graft during the placement and deployment of the main body extension.

5. Deploy the device by withdrawing the sheath while stabilizing the gray positioner of the delivery system. (Figs. 35 and 42) Continue to deploy the device until the most distal stent is uncovered. Stop withdrawing the sheath.
6. Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle, and then remove via the slot over the inner cannula. (Fig. 37)
7. Withdraw the tapered tip of the introducer back through the main body extension graft and delivery system while maintaining wire guide position. Ensure the main body extension and endovascular graft are not displaced during withdrawal of the delivery system.
8. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops. (Fig. 38)

Main Body Extension Molding Balloon Insertion

NOTE: For information on the use of recommended products, refer to the individual product's Instructions for Use.

1. Prepare the molding balloon as follows:
 - Flush the wire lumen with heparinized saline.
 - Remove all air from the balloon.

CAUTION: The Captor Hemostatic Valve must be open prior to repositioning of the molding balloon.

2. Advance the molding balloon over the wire guide and through the hemostatic valve of the main body introduction system to the level of the main body extension.
3. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning the Captor Hemostatic Valve clockwise.

CAUTION: Do not inflate the balloon in the vessel outside of the graft.

4. Expand the molding balloon within the proximal segment of the main body extension and then the most distal segment of the main body extension using diluted contrast media (as recommended by the manufacturer). (Fig. 43)

CAUTION: Confirm complete deflation of the balloon prior to repositioning.

5. Completely deflate and remove the molding balloon, replace it with an angiographic catheter and perform completion angiograms.
6. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

11 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

11.1 General

- **The long-term performance of this endovascular graft has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.
- Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in **Table 11.1.1**. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).

Table 11.1.1 lists the minimum imaging follow-up for patients with the Zenith Alpha Abdominal Endovascular Graft. Patients requiring enhanced follow-up should have interim evaluations.

Table 11.1.1 Recommended Imaging Follow-up Schedule

	Pre-op	Intra-op	30-Day	6-Month	12-Month ⁴
CT Scan	X ¹		X ³	X ³	X ³
Abdominal Device X-ray			X	X	X
Angiography	X ²	X			

¹ Imaging should be performed within 6 months before the procedure.

² Required only to resolve any uncertainties in anatomical measurements necessary for graft sizing.

³ Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan.

⁴ Yearly thereafter

11.2 Contrast and Non-Contrast CT Recommendations

- Film sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). DO NOT perform large slice thickness (>3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14 inch x 17 inch sheets if film is used.

It is important to follow acceptable imaging protocols during the CT exam. **Table 11.2.1** lists examples of acceptable imaging protocols

11.2.1 Acceptable Imaging Protocols

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral CT or high performance MDCT capable of >40 seconds	Spiral CT or high performance MDCT capable of >40 seconds
Injection volume	n/a	Per institutional protocol
Injection rate	n/a	>2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage – start	Diaphragm	1 cm superior to celiac axis
Coverage – finish	Proximal femur	Profunda femoris origin
Collimation	<3 mm	<3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

11.3 Abdominal Radiographs

The following views are required:

- Four films: supine-frontal (AP), lateral, 30 degree LPO and 30 degree RPO views centered on umbilicus.
- Record the table-to-film distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2-4X magnification visual aid.

11.4 MRI Information

NOTE: If using this device in conjunction with another endovascular graft from the Zenith family, refer to the appropriate device's IFU for additional MRI information.

Nonclinical testing has demonstrated that the Zenith Low Profile AAA Endovascular Graft (The Zenith Low Profile AAA Endovascular Graft served as a surrogate for the Zenith Alpha Abdominal Endovascular Graft) is MR Conditional. A patient with this endovascular graft can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less.
- Spatial magnetic gradient field of 1580 Gauss/cm (15.8 T/m) or less.
- The product of the spatial gradient and static magnetic field should not exceed 47.4 T²/m.

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 and 3.0 Tesla Systems: Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e., per scanning sequence).

1.5 Tesla Temperature Rise:

In nonclinical testing, the Zenith Low Profile AAA Endovascular Graft (The Zenith Low Profile AAA Endovascular Graft served as a surrogate for the Zenith Alpha Abdominal Endovascular Graft) produced a maximum temperature rise of 1.7 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Magnetom, Siemens Medical Solutions, Malvern, PA, Software Numaris/4) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.1 W/kg).

3.0 Tesla Temperature Rise:

In nonclinical testing, the Zenith Low Profile AAA Endovascular Graft (The Zenith Low Profile AAA Endovascular Graft served as a surrogate for the Zenith Alpha Abdominal Endovascular Graft) produced a maximum temperature rise of 2.0 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (Excite, GE Healthcare, Milwaukee, WI, Software G3.0-052B) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.8 W/kg).

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 5 mm of the position of the Zenith Alpha Abdominal Endovascular Graft, as found during nonclinical testing using the sequences: T1-weighted spin echo and gradient echo pulse in a 3.0 Tesla MR system (Excite, General Electric Healthcare). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail: MedAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

11.5 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, ≥5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's comorbidities, life expectancy and the patient's personal choices. Patients should be counseled that subsequent reinterventions, including catheter based and open surgical conversion are possible following endograft placement.

12 MAIN BODY RELEASE TROUBLESHOOTING

NOTE: Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

NOTE: In case of difficulty removing the bare stent wires from the graft during rotation of the blue handle, proceed through the normal deployment process of the ipsilateral limb release. If this still does not alleviate difficulty, proceed with the main body release troubleshooting steps below to disassemble the rotation handle.

NOTE: In case of difficulty with ipsilateral limb release, proceed directly to steps 1-4 below.

1. Place surgical forceps into the slots on the back-end clips and slide both clips out. (**Fig 44 and 45**)
 2. Remove the back-end cap from the blue rotation handle. (**Fig 46**)
 3. While holding the black gripper, remove the blue rotation handle by sliding it straight back. (**Fig. 47**) The trigger wires will be visible. (**Fig. 48**)
- NOTE:** Tension will be felt during the removal of the rotation handle.
- NOTE:** If the blue rotation handle feels as though it can not be removed, turn the blue rotation handle in the opposite direction of the white arrow (counterclockwise) and continue until the handle has been completely removed.
- NOTE:** The handle is designed such that any and all failures will not impede rotation of the blue handle and removal of the trigger wires.
4. Using forceps, grasp all the wires and withdraw them until the graft is released. (**Fig 49**)

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© COOK 2013