EN
8

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TABLE OF CONTENTS

Illustrations	1-7
1 DEVICE DESCRIPTION.....	8
1.1 Zenith Main Body and Iliac Leg Components	8
1.2 Main Body Introduction System	8
1.3 Iliac Leg Introduction System	8
1.4 Zenith AAA Endovascular Graft Ancillary Components	8
2 INTENDED USE	8
3 CONTRAINDICATIONS	8
4 WARNINGS AND PRECAUTIONS	8
4.1 General	8
4.2 Patient Selection, Treatment and Follow-Up	8
4.3 Pre-Procedure Measurement Techniques and Imaging	8
4.4 Device Selection	9
4.5 Implant Procedure	9
4.6 Molding Balloon Use	9
4.7 MRI Safety and Compatibility	9
5 ADVERSE EVENTS	9
5.1 Potential Adverse Events	9
Device-Related Adverse Event Reporting	10
6 PATIENT SELECTION AND TREATMENT	10
6.1 Individualization of Treatment	10
7 PATIENT COUNSELING INFORMATION.....	10
8 HOW SUPPLIED	10
9 CLINICAL USE INFORMATION.....	10
9.1 Physician Training	10
9.2 Inspection Prior to Use	10
9.3 Materials Required	10
9.4 Materials Recommended	10
9.5 Device Sizing Guidelines	11
Table 9.5.1 Main Body (ZAU) Graft Sizing Guide	11
Table 9.5.2 Iliac Leg (TFLE) Graft Sizing Guide	11
10 DIRECTIONS FOR USE	11
Anatomical Requirements	11
General Use Information	11
Pre-Implant Determinants	11
Patient Preparation	11
10.1 Zenith Flex AUI AAA Main Body and Zenith Flex AAA Iliac Leg	11
10.1.1 Main Body Preparation/Flush	11
10.1.2 Iliac Leg Preparation/Flush	11
10.1.3 Vascular Access and Angiography	11
10.1.4 Main Body Placement	11
10.1.5 Main Body Proximal (Top) Deployment	12
10.1.6 Main Body Distal (Bottom) Deployment	12
10.1.7 Docking of Top Cap	12
10.1.8 Iliac Leg Placement and Deployment	12
10.1.9 Iliac Plug Placement and Deployment	12
10.1.10 Molding Balloon Insertion	12
10.1.11 Final Angiogram	12
10.1.12 Femoral-to-Femoral Crossove	12
11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP	12
11.1 General	12
Table 11.1.1 Recommended Imaging Schedule for Endograft Patients	13
11.2 Contrast and Non-Contrast CT Recommendations	13
Table 11.2.1 Acceptable Imaging Protocols	13
11.3 Abdominal Radiographs	13
11.4 Ultrasound	13
11.5 MRI Safety and Compatibility	13
11.6 Additional Surveillance and Treatment	13

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ZENITH FLEX® AUI AAA ENDOVASCULAR GRAFT WITH THE Z-TRAK INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

CAUTION: All contents of the outer pouch (including the introduction system and endovascular grafts) are supplied sterile, for single use only.

For the Zenith product line there are four applicable Suggested Instructions For Use (IFU). This IFU describes the suggested Instructions for Use for the Zenith Flex AUI AAA Endovascular Graft (main body and iliac leg). For information on other Zenith components, please refer to the following Instructions for Use:

- Zenith® AAA Endovascular Graft Ancillary Components (main body extension, iliac leg extension, converter, and iliac plug);
- Zenith AAA Iliac Plug; and
- CODA® Balloon Catheter.

1 DEVICE DESCRIPTION

1.1 Zenith Main Body and Iliac Leg Components

The Zenith Flex AUI AAA Endovascular Graft consists of a tapered aortic main body. (**Fig 1**) The Zenith Flex AUI AAA Endovascular Graft may be used independently or in combination with a Zenith AAA Endovascular Graft iliac leg (TFLE). (**Fig 2**) Zenith Flex AUI AAA Endovascular Graft is constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z® stents with braided polyester and monofilament polypropylene suture. The modules are fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents provide the necessary attachment and seal of the graft to the vessel wall.

The bare suprarenal stent at the proximal end of the graft contains bars that are placed at 3 mm increments for additional fixation of the device. To facilitate fluoroscopic visualization of the stent graft, four gold radiopaque markers are positioned in a circumferential orientation within 2 mm of the most proximal aspect of the graft material.

1.2 Main Body Introduction System

The Zenith Flex AUI AAA Endovascular Graft is shipped preloaded onto the Z-Trak Introduction System. (**Fig 3**) The system has a sequential deployment method with built-in features to provide continuous control of the endovascular graft throughout the deployment procedure. The Z-Trak Introduction System is designed for precise positioning and allows readjustment of the final graft position before deployment of the barbed suprarenal stent.

The main body graft introduction system uses an 18, 20 or 22 French Z-Trak Introduction System. Dual trigger-wire release mechanisms lock the endovascular graft onto the introduction system until released by the physician. All systems are compatible with a .035 inch wire guide. For added hemostasis, the Captor® Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The introduction system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.3 Iliac Leg Introduction System

The Zenith Flex AAA Endovascular Graft iliac legs are shipped preloaded onto the Z-Trak Introduction System. (**Fig 2**) The introduction system is designed for ease of use with minimal preparation. The iliac leg introduction system uses a 14 French or 16 French Z-Trak Introduction System. All systems are compatible with a .035 inch wire guide.

For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath.

1.4 Zenith® AAA Endovascular Graft Ancillary Components

Additional ancillary endovascular components (main body extensions, iliac leg extensions and iliac plugs) are available. (**Fig. 4**) Refer to the Zenith AAA Endovascular Graft Ancillary Components Instructions for Use for more information.

2 INTENDED USE

The Zenith Flex AUI AAA Endovascular Graft is intended for treatment of patients with abdominal aortic, ruptured abdominal aortic, aortoiliac, or iliac aneurysms.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Standard techniques for placement of arterial access sheaths, angiographics catheters and wire guides should be employed.

3 CONTRAINDICATIONS

The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), polypropylene, or gold.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 9.1, Physician Training**.
- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft and/or leaks may be required to undergo secondary interventions or surgical procedures.

4.2 Patient Selection, Treatment and Follow-Up

- The Zenith Flex AUI AAA Endovascular Graft is designed to treat aortic neck diameters no smaller than 18 mm and no larger than 32 mm. The Zenith Flex AUI AAA Endovascular Graft is designed to treat proximal aortic necks (distal to the lowest renal artery) of at least 15 mm in length. Iliac artery distal fixation site greater than 10 mm in length and 7.5 – 20 mm in diameter (measured outer wall to outer wall) is required. These sizing measurements are critical to the performance of the endovascular repair.

- Key anatomical elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (>60 degrees for infrarenal neck to axis of AAA or >45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic neck (<15 mm); an inverted funnel shape (greater than 10% increase in diameter over 15 mm of proximal aortic neck length); and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. In the presence of anatomical limitations, a longer neck may be required to obtain adequate sealing and fixation. Irregular calcification and/or plaque may compromise the attachment and sealing at the fixation sites. Necks exhibiting these key anatomical elements may be more conducive to graft migration or endoleak.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and introduction systems of the profile of a 16 French to 22 French vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization. A vascular conduit technique may be necessary to achieve success in some patients
- The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System is not recommended in patients who cannot tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensible inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The safety and effectiveness of the Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System has not been evaluated in the following patient populations:

- traumatic aortic injury
- mycotic aneurysms
- pseudoaneurysms resulting from previous graft placement
- uncorrectable coagulopathy
- indispensible mesenteric artery
- genetic connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
- concomitant thoracic aortic or thoracoabdominal aneurysms
- patients with active systemic infections
- pregnant or nursing females
- morbidly obese patients
- less than 18 years of age
- patients with less than 15 mm in length or greater than 60 degrees angulation of the proximal aortic neck relative to the long axis of the aneurysm.
- patients with a neck diameter measured outer wall to outer wall greater than 32 mm and less than 18 mm
- patients with iliac distal fixation site less than 10 mm in length and less than 7.5 mm and greater than 20 mm in diameter (measured outer wall to outer wall)

- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3 Pre-Procedure Measurement Techniques and Imaging**.

- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intra-operative flexibility to achieve optimal procedural outcomes.

4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith Flex AUI AAA Endovascular Graft. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning the x-ray C-arm during procedural angiography such that the origins of the renal arteries, and particularly the lowest patent renal artery, are well demonstrated prior to deployment of the proximal edge of the graft material (sealing stent) of the main body. Additionally, angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are well defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg component.

Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

Lengths

Utilizing CT, length measurements should be determined to accurately assess infrarenal proximal neck length as well as planning main body sizes and leg components for the Zenith Flex AUI AAA Endovascular

Graft. These reconstructions should be performed in sagittal, coronal, and 3-D.

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-up.**
- The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System is not recommended in patients unable to undergo, or who will not be compliant with, the necessary pre-operative and post-operative imaging and implantation studies as described in **Section 11, Imaging Guidelines and Post-Operative Follow-up.**
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

4.4 Device Selection

- Strict adherence to the Zenith Flex AUI AAA Endovascular Graft IFU sizing guide is strongly recommended when selecting the appropriate device size (**Tables 9.5.1 and 9.5.2**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

4.5 Implant Procedure

(Refer to **Section 10, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith Flex AUI AAA Endovascular Graft and assure accurate apposition to the aortic wall.
- Do not bend or kink the introduction system. Doing so may cause damage to the introduction system and the Zenith Flex AUI AAA Endovascular Graft.
- To avoid any twist in the endovascular graft, during any rotation of the introduction system, be careful to rotate all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the introduction system if resistance is felt during advancement of the wire guide or introduction system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith Flex AUI AAA Endovascular Graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (exception is the inferior mesenteric artery) with the endoprosthesis. Vessel occlusion may occur.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Repositioning the stent graft distally after partial deployment of the covered proximal stent may result in damage to the stent graft and/or vessel injury.
- Inaccurate placement and/or incomplete sealing of the Zenith Flex AUI AAA Endovascular Graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Zenith Flex AUI AAA Endovascular Graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implant procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during introduction system insertion.
- Always use fluoroscopy for guidance, delivery and observation of any Zenith Flex AUI AAA Endovascular Graft components within the vasculature.
- The use of the Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise (e.g., adequate hydration).
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- The Zenith Flex AUI AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional and angiographic devices in the region of the suprarenal stent.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization, or rupture of the aneurysm.
- Avoid damaging the graft or disturbing graft positioning after placement in the event reinstrumentation (secondary intervention) of the graft is necessary.
- Before deployment of the suprarenal stent, verify that the position of the access wire guide extends just distal to the aortic arch.
- Verify that the predetermined occlusion device (if necessary) is selected for insertion on the contralateral side of the patient before implantation.

- Use of the Zenith Flex AUI AAA Endovascular Graft typically requires occlusion of the contralateral iliac artery to ensure exclusion of the aneurysm.

4.6 Molding Balloon Use

- Do not inflate the balloon in vessel outside of graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care in inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.
- For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

4.7 MRI Safety and Compatibility

For MRI Safety and Compatibility please refer to below description based on the Zenith AAA Endovascular Graft. Zenith Flex AUI AAA Endovascular Graft and Zenith AAA Endovascular Grafts are constructed of similar materials and therefore comparable with regards to MRI safety and compatibility.

Non-clinical testing has demonstrated that the Zenith AAA Endovascular Graft is MR Conditional. It can be scanned safely under the following conditions:

1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Spatial gradient field of 450 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris /4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Spatial gradient field of 720 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

Clinical information is available on seventeen patients who received MRI scans after stent-graft implantation. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been well over 50,000 Zenith AAA Endovascular Grafts implanted worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

5 ADVERSE EVENTS

5.1 Potential Adverse Events

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anasthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)

- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Occlusion of device or native vessel
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)

Device-Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith Flex AUI AAA Endovascular Graft should be reported to COOK immediately. Please call your distributor.

6 PATIENT SELECTION AND TREATMENT

(See **Section 4, Warnings and Precautions**)

6.1 Individualization of Treatment

Cook recommends that the Zenith Flex AUI AAA Endovascular Graft component diameters be selected as described in Tables 9.5.1 and 9.5.2. The length of the Zenith Flex AUI AAA Endovascular Graft should extend from the lowest renal artery to just above the internal iliac (hypogastric) artery bifurcation. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain.

This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The risks and benefits should be carefully considered for each patient before use of the Zenith Flex AUI AAA Endovascular Graft. Additional considerations for patient selection include but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of potential or existing aneurysm rupture compared to the risk of treatment with the Zenith Flex AUI AAA Endovascular Graft
- Ability to tolerate general, regional or local anesthesia
- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 16 French to 22 French vascular introducer sheath
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a length of at least 15 mm,
 - with a diameter measured outer wall to outer wall of no greater than 32 mm and no less than 18 mm,
 - with an angle less than 60 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal aorta
- Iliac artery distal fixation site greater than 10 mm in length and 7.5 to 20 mm in diameter (measured outer wall to outer wall)
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft
- Use of the Zenith Flex AUI AAA Endovascular Graft typically requires occlusion of the contralateral iliac artery in conjunction with a standard femoral-to-femoral-bypass procedure. For instructions for preparation and deployment of the Zenith AAA iliac plug, refer to the Zenith AAA Iliac Plug Instructions for Use.

The final treatment decision is at the discretion of the physician and patient.

7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up.**
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine post-operative follow-up requirements are required and should be considered a life-long commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise all patients that it is important to seek prompt medical attention if they experience signs of graft occlusion, aneurysm enlargement or rupture. Signs of graft occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.

- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant. Men who undergo endovascular or open surgical repair may experience impotence.

The physician should complete the Patient I.D. Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

8 HOW SUPPLIED

- The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System is supplied sterile and pre-loaded in peel-open packages.
- The device is intended for single use only. Do not re-sterilize the device.
- Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to COOK.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- The main body device is loaded onto an 18, 20 or 22 French Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a sterile gauze pad soaked in saline solution.
- Do not use after the "USE BY" (expiration) date printed on the label.
- Store in a cool, dry place.

9 CLINICAL USE INFORMATION

9.1 Physician Training

CAUTION: Always have a qualified surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System are outlined below:

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection, planning and sizing.

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to COOK. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

9.3 Materials Required

(Not included in Zenith Flex AUI AAA Endovascular Graft)

- A selection of ancillary components including main body extensions, iliac extensions and iliac plugs
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

9.4 Materials Recommended

(Not included in Zenith Flex AUI AAA Endovascular Graft)

The following products are recommended for implantation of any component in the Zenith product line. For information on the use of these products, refer to the individual product's Suggested Instructions for Use.

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide; for example:
 - Cook .035 inch wire guides
 - Cook Nimble™ Wire Guides
- Molding Balloons; for example
 - Cook CODA Balloon Catheter
- Introducer sets; for example
 - Cook Check-Flo® Introducer Set
 - Cook Extra Large Check-Flo Introducer Sets
- Sizing catheter; for example:
 - Cook Auros® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
- Entry needles; for example:
 - Cook single wall entry needles

9.5 Device Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Table 9.5.1 Main Body Graft Sizing Guide*

Intended Aortic Vessel Diameter ^{1,2} (mm)	Main Body Diameter ³ (mm)	Overall Graft Length (mm)	Introducer Sheath ID /OD
18-19	22	113, 130, 147	18 Fr/ 7.1 mm
20-21	24	113, 130, 147	18 Fr/ 7.1 mm
22	26	113, 130, 147	18 Fr/ 7.1 mm
23-24	28	113, 130, 147	20 Fr/ 7.7 mm
25-26	30	113, 130, 147	20 Fr/ 7.7 mm
27-28	32	113, 130, 147	20 Fr/ 7.7 mm
29-32	36	127, 144, 161	22 Fr/ 8.5 mm

¹Maximum diameter along the proximal fixation site.

²Round measured aortic diameter to nearest mm.

³Additional considerations may affect choice of diameter.

*All dimensions are nominal.

Table 9.5.2 Iliac Leg Graft Sizing Guide*

Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Leg Diameter ³ (mm)	Iliac Leg Working Length ⁴ (mm)	Introducer Sheath ID /OD
<8	8	37, 54, 71, 88, 105,122	14 Fr/5.3 mm
8-9	10	37, 54, 71, 88, 105,122	14 Fr/5.3 mm
10-11	12	39, 56, 73, 90, 107, 124	16 Fr/6.0mm
12-13	14	39, 56, 73, 90, 107, 124	16 Fr/6.0mm
14-15	16	39, 56, 73, 90	16 Fr/6.0mm
16-17	18	39, 56, 73, 90	16 Fr/6.0mm
18	20	39, 56, 73, 90	16 Fr/6.0mm
19	22	39, 56, 73, 90	16 Fr/6.0mm
20	24	39, 56, 73, 90	16 Fr/6.0mm

¹Maximum diameter along the proximal fixation site.

²Round measured aortic diameter to nearest mm.

³Additional considerations may affect choice of diameter.

⁴Overall leg length = working length + 22 mm stent

*All dimensions are nominal.

10 DIRECTIONS FOR USE

Illustrations

Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.
- Proximal aortic neck lengths should be a minimum of 15 mm with a diameter measured outer wall to outer wall of 18 – 32 mm
- Iliac artery distal fixation site should be greater than 10 mm in length and 7.5 – 20 mm in diameter (measured outer wall to outer wall)

Prior to use of the Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System, review this Suggested Instructions for Use booklet.

The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

General Use Information

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System. The Zenith Flex AUI AAA Endovascular Graft with the Z-Trak Introduction System is compatible with .035 inch diameter wire guides.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected.

Determinants include:

- Femoral artery selection for insertion of the introduction system
- Angulation of aortic neck, aneurysm and iliac arteries.
- Quality of the aortic neck.
- Diameters of infrarenal aortic neck and distal iliac arteries.
- Distance from the renal arteries to the internal iliac artery/attachment site.
- Length from the aortic bifurcation to the internal iliac artery/ attachment site.
- Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
- Consider the degree of vascular calcification.

Patient Preparation

- Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
- Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
- Expose both common femoral arteries using standard surgical technique.
- Establish adequate proximal and distal vascular control of both femoral vessels.

10.1 Zenith Flex AUI AAA Main Body and Zenith Flex AAA Iliac Leg (Figs. 2 and 3)

10.1.1 Main Body Preparation/Flush

- Remove grey-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away® sheath from back of hemostatic valve. (Fig. 5) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the sideport near the tip of the introducer sheath. (Fig. 6) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

NOTE: Graft flushing solution of heparinized saline is often used.

- Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (Fig. 7)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

- Soak sterile gauze pads in saline solution and use to wipe the Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

10.1.2 Iliac Leg Preparation/Flush

- Remove grey-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (Fig. 8) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the sideport near the tip of the introducer sheath. (Fig. 6) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

NOTE: Graft flushing solution of heparinized saline is often used.

- Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (Fig. 7)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

10.1.3 Vascular Access and Angiography

- Puncture the selected common femoral arteries using standard technique with an 18 or 19 UT gage arterial needle. Upon vessel entry, insert:

- Wire guides — standard .035 inch diameter, 145 cm long, J tip or Benton Wire Guide
- Appropriate size sheaths (e.g., 6 or 8 French)
- Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)

- Perform angiography to identify level(s) of renals, aortic bifurcation and iliac bifurcations.

NOTE: If fluoroscope angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

10.1.4 Main Body Placement

- Ensure the introduction system has been flushed with heparinized saline and that all air is removed from the system.
- Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

NOTE: Monitor the patient's coagulation status throughout the procedure.

- Replace J wire with stiff wire guide (LES). .035 inch, 260 cm long, and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.
- Insert main body introduction system over the wire, into the femoral artery.

CAUTION: Maintain wire guide position during introduction system insertion.

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the introduction system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

- Advance introduction system until the four gold radiopaque markers (which are positioned 2 mm from the most proximal segment of the graft material) (Fig. 9) are just inferior to the most inferior renal orifice.

- Repeat the angiogram to verify the four gold radiopaque markers are just inferior to the most inferior renal orifice, and to ensure the distal portion of the endovascular graft will not occlude the internal iliac artery.
- Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (Fig. 10)
- Use the gripper to stabilize the grey positioner (the shaft of the introduction system) while withdrawing the sheath. Continue to withdraw the sheath while verifying position until the distal segment of the Zenith AAA AUI main body is fully deployed. (Figs. 11 and 12). Stop withdrawing the sheath.

10.1.5 Main Body Proximal (Top) Deployment

- Perform angiography through an angiographic catheter to verify position of the endovascular graft with respect to the renal arteries. If necessary, carefully reposition the covered portion of the endovascular graft with respect to the renal arteries. (Repositioning can only take place over a small range of distance at this stage.)

NOTE: Ensure patency of renal arteries by confirming that the proximal graft markers are 2 mm or more below the lowest patent renal artery.

- Remove the safety lock from the top stent trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the top stent trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Fig. 13)

If unable to remove the top stent trigger-wire release mechanism from the top cap, perform the following steps under fluoroscopy:

- Remove tension on the trigger-wire by loosening the pin vise and slightly pulling the inner cannula to move the top cap down over the suprarenal stent. Avoid compressing the Zenith Flex AUI AAA Endovascular Graft main body.
 - Retighten the pin vise.
 - Remove the top stent trigger-wire release mechanism.
 - Continue with Step 3 in Section 10.1.5, Main Body Proximal (Top) Deployment.
- Loosen the pin vise. (Fig. 14) Control the position of the graft by stabilizing the grey positioner of the introducer.

CAUTION: Before deployment of suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch.

- Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the main body until the top stent is fully deployed. (Figs. 15 and 16) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise to avoid contact with the deployed suprarenal stent.

NOTE: Once the barbed suprarenal stent has been deployed, further attempts to reposition the graft are not recommended.

WARNING: The Zenith Flex AUI AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

10.1.6 Main Body Distal (Bottom) Deployment

NOTE: The distal stent is still secured by the trigger-wire.

- Remove the safety lock. Withdraw and remove the trigger-wire by sliding the trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Fig. 17)

10.1.7 Docking of Top Cap

- Loosen the pin vise. (Fig. 18)
- Secure sheath and inner cannula to avoid any movement of these components.
- Advance the grey positioner over the inner cannula until it docks with the top cap. (Figs. 19, 20 and 21)

NOTE: If resistance occurs, slightly rotate grey positioner and continue to gently advance.

- Retighten the pin vise and withdraw the entire top cap and grey positioner through the graft and through the sheath by pulling on the inner cannula. Leave the sheath and wire guide in place.

NOTE: Maintain position of sheath and wire guide.

- Close the Captor Hemostatic Valve on the main body introducer sheath by turning it in a clockwise direction until it stops. (Fig. 22)

10.1.8 Iliac Leg Placement and Deployment

NOTE: Ensure the Captor Hemostatic Valve on the main body introducer sheath is turned to the open position. (Fig. 23)

- Utilize the main body graft wire and sheath assembly to introduce the iliac leg graft. Advance dilator and sheath assembly into the main body sheath.

NOTE: In tortuous vessels, the position of the internal iliac arteries may alter significantly with the introduction of the rigid wires and sheath systems.

- Advance slowly until the iliac leg graft overlaps a minimum of two full iliac leg stents (i.e., proximal stents of iliac leg graft) inside the main body. (Fig. 24)
- Confirm position of distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure internal iliac patency.

NOTE: Ensure the Captor Hemostatic Valve on the iliac leg introducer sheath is turned to the open position. (Fig. 10)

- To deploy, stabilize the iliac leg graft with the gripper on the grey positioner while withdrawing the iliac leg sheath. (Figs. 25 and 26) If necessary, withdraw the main body sheath.
- Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to grey positioner. Tighten pin vise. Maintain main body sheath position while withdrawing iliac leg sheath and grey positioner with secured inner cannula.
- Close the Captor Hemostatic Valve on the main body introducer sheath by turning it in a clockwise direction until it stops.
- Re-check the position of the wire guides. Leave sheath and wire guides in place.

10.1.9 Iliac Plug Placement and Deployment

Use of the Zenith Flex AUI AAA Endovascular Graft typically requires occlusion of the contralateral iliac artery. (Fig. 27) For instruction on preparation and deployment of the Zenith AAA Iliac Plug, refer to the Zenith AAA Iliac Plug Endovascular Graft Instructions for Use.

10.1.10 Molding Balloon Insertion

- Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.
- In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counter-clockwise.
- Advance the molding balloon over the wire guide and through the Captor Hemostatic Valve of the main body introduction system to the level of the renal arteries. Maintain proper sheath position.
- Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.

CAUTION: Do not inflate balloon in vessel outside of graft.

- Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the most proximal covered stent and the infrarenal neck, starting proximally and working in the distal direction. (Fig. 28)

CAUTION: Confirm complete deflation of balloon prior to repositioning.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

- Withdraw the molding balloon to the main body and iliac leg overlap and expand.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

CAUTION: Do not inflate balloon in vessel outside of graft.

- Withdraw the molding balloon to the iliac leg distal sealing site and expand.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

- Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
- Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

10.1.11 Final Angiogram

- Position angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of internal iliac artery.
- Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.

NOTE: If endoleaks or other problems are observed, refer to the Suggested Instructions for Use for the Zenith AAA Endovascular Graft Ancillary Components.

- Repair vessels and close in standard surgical fashion.

10.1.12 Femoral-to-Femoral Crossover

Perform femoral-to-femoral crossover in standard surgical fashion to revascularize the contralateral limb. Repair vessels and close in standard surgical fashion.

11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

11.1 General

- The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up. Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.
- Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in Table 11.1.1. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT.

Table 11.1.1 lists the minimum imaging follow-up for patients with the Zenith Flex AUI AAA Endovascular Graft. Patients requiring enhanced follow-up should have interim evaluations.

Table 11.1.1 Recommended Imaging Schedule for Endograft Patients

	Angiogram	CT (Contrast and non-contrast)	Abdominal Radiographs
Pre-procedure	X ¹	X ¹	
Procedural	X		
Pre-discharge (within 7 days)		X ^{2,3,4}	X
1 month		X ^{2,3,4}	X
3 month		X ^{2,4,5}	
6 month		X ^{2,4}	X
12 month (annually thereafter)		X ^{2,4}	X

¹Imaging should be performed within 6 months before the procedure.

²Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan. With ultrasound, non-contrast CT is still recommended.

³Either pre-discharge or 1 month CT recommended.

⁴If Type I or III endoleak, prompt intervention and additional follow-up post-intervention recommended. See **Section 11.6, Additional Surveillance and Treatment**.

⁵Recommended if endoleak reported at pre-discharge or 1 month.

11.2 Contrast and Non-Contrast CT Recommendations

- Film sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). DO NOT perform large slice thickness (> 3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
 - All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14 inch x 17 inch (35.5 cm x 43.2 cm) sheets if film is used
 - Both non-contrast and contrast runs are required, with matching or corresponding table positions.
 - Pre-contrast and contrast run slice thickness and interval must match.
 - DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.
- Non-contrast and contrast enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. **Table 11.2.1** lists examples of acceptable imaging protocols.

Table 11.2.1 Acceptable Imaging protocols

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral CT or high performance MDCT capable of >40 seconds	Spiral CT or high performance MDCT capable of >40 seconds
Injection volume	n/a	Per institutional protocol
Injection rate	n/a	>2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage - start	Diaphragm	1 cm superior to celiac axis
Coverage - finish	Proximal femur	Profunda femoris origin
Collimation	<3 mm	<3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

11.3 Abdominal Radiographs

The following views are required:

- Four films: supine-frontal (AP), cross-table lateral, 30 degree LPO and 30 degree RPO views centered on umbilicus.
- Record the table-to-film distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2-4X magnification visual aid.

11.4 Ultrasound

Ultrasound imaging may be performed in place of contrast CT when patient factors preclude the use of image contrast media. Ultrasound may be paired with non-contrast CT. A complete aortic duplex is to be videotaped for maximum aneurysm diameter, endoleaks, stent patency and stenosis.

Included on the videotape should be the following information as outlined below:

- Transverse and longitudinal imaging should be obtained from the level of the proximal aorta demonstrating mesenteric and renal arteries to the iliac bifurcations to determine if endoleaks are present utilizing color flow and color power angiography (if accessible).
- Spectral analysis confirmation should be performed for any suspected endoleaks
- Transverse and longitudinal imaging of the maximum aneurysm should be obtained.

11.5 MRI Safety and Compatibility

For MRI Safety and Compatibility please refer to below description based on the Zenith AAA Endovascular Graft. Zenith Flex AUI AAA Endovascular Grafts and Zenith AAA Endovascular Grafts are constructed of similar materials and therefore comparable with regards to MRI safety and compatibility.

Non-clinical testing has demonstrated that the Zenith AAA Endovascular Graft is MR Conditional. It can be scanned safely under the following conditions:

1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Spatial gradient field of 450 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris /4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Spatial gradient field of 720 Gauss/cm
- Maximum whole-body-averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

Clinical information is available on seventeen patients who received MRI scans after stent-graft implantation. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been well over 50,000 Zenith AAA Endovascular Grafts implanted worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

11.6 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, ≥ 5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counseled that subsequent reinterventions including catheter based and open surgical conversion are possible following endograft placement.

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