Bakri Postpartum Balloon

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Bakri-tamponadeballon

Bakri-tamponade-ballon

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

BAKRI POSTPARTUM BALLOON

Used for control of postpartum uterine bleeding. Intended for one-time use.

CONTENTS

- Bakri Postpartum Balloon
- Syringe

CAUTION: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

WARNING: This device is intended as a temporary means of establishing hemostasis in cases indicating conservative management of postpartum uterine bleeding.

WARNING: Device should not be left indwelling for more than twenty-four (24) hours.

WARNING: Maximum inflation is 500 ml.

WARNING: Clinical data to support the safety and effectiveness of the Bakri Postpartum Balloon in the setting of uterine atony are limited. Patients in whom this device is being used should be closely monitored for signs of worsening bleeding and/or disseminated intravascular coagulation (DIC). In such cases, emergency intervention per hospital protocol should be followed.

WARNING: There are no clinical data to support use of this device in the setting of DIC.

WARNING: Patient monitoring is an integral part of managing postpartum hemorrhage. Signs of deteriorating or non-improving condition should lead to a more aggressive treatment and management of patient uterine bleeding.

IMPORTANT: Patient urine output should be monitored while Bakri Postpartum Balloon is in use.

PRECAUTIONS

- Avoid excessive force when inserting balloon into uterus.

INDICATIONS FOR USE: Use of this device is intended to provide temporary control or reduction of postpartum uterine bleeding when conservative management is warranted.

USE OF THIS DEVICE IS CONTRAINDICATED IN THE PRESENCE OF

- Arterial bleeding requiring surgical exploration or angiographic embolization
- Cases indicating hysterectomy
- Pregnancy
- Cervical cancer
- Purulent infections in the vagina, cervix, or uterus
- Untreated uterine anomaly
- Disseminated intravascular coagulation
- A surgical site which would prohibit the device from effectively controlling bleeding

IMPORTANT: In order to maximize tamponade effect, constant pressure must be affected between the balloon wall and the tissue surface. Pressure can be achieved and maintained by applying gentle traction to the balloon shaft, before securing to the patient's leg, or attaching to a weight, not to exceed 500 grams.

SUGGESTED PLACEMENT METHODS

Transvaginal Placement

1. Determine uterine volume by direct examination or ultrasound examination.
2. Using ultrasound guidance, insert the balloon portion of the catheter in the uterus making certain that the entire balloon is inserted past the cervical canal and internal ostium.

3. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.

Transabdominal Placement, Post Cesarean Section

1. Determine uterine volume by direct examination intraoperative or ultrasound examination postoperative.
2. From above via access of the Cesarean incision, pass the tamponade balloon, inflation port first, through the uterus and cervix.
3. Have an assistant pull the shaft of the balloon through the vaginal canal until the deflated balloon base comes in contact with the internal cervical ostium.
4. Close the incision per normal procedure, taking care to avoid puncturing the balloon while suturing.

SUGGESTED INSTRUCTIONS FOR BALLOON INFLATION

WARNING: Always inflate the balloon with sterile liquid. Never inflate with air, carbon dioxide or any other gas.

1. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.

WARNING: Do not overinflate the balloon. Refer to product label and instructions for use for maximum inflation volume.

IMPORTANT: To ensure that the balloon is filled to the desired volume, it is recommended that the predetermined volume of fluid be placed in a separate container, rather than relying on a syringe count to verify the amount of fluid that has been instilled into the balloon.

2. Using the enclosed syringe, begin filling the balloon to the predetermined volume through the stopcock.
3. Apply gentle traction to the balloon shaft to ensure proper contact between the balloon and tissue surface. To maintain tension, secure the balloon shaft to the patient's leg or attach to a weight, not to exceed 500 grams.

NOTE: To maximize the tamponade effect, counter pressure can be applied by packing the vaginal canal with iodine or antibiotic soaked vaginal gauze.

4. Connect the drainage port to a fluid collection bag to monitor hemostasis.

IMPORTANT: To adequately monitor hemostasis, the balloon drainage port and tubing should be flushed clear of clots with sterile isotonic saline.

5. Monitor patient continuously for signs of increased bleeding and uterine cramping.

BALLOON REMOVAL

Maximum indwell time is twenty-four (24) hours. Balloon may be removed sooner upon the physician's determination of hemostasis.

1. Remove tension from balloon shaft.
2. Remove any vaginal packing.
3. Using an appropriate syringe, aspirate the contents of the balloon until fully deflated.
4. Gently retract the balloon from the uterus and vaginal canal and discard.
5. Monitor patient for signs of bleeding.

REFERENCE

Y. N. Bakri, A. Amri, F. Abdul Jabbar: "Tamponade balloon for obstetrical bleeding," *International Journal of Gynecology and Obstetrics*, 74 (2001) 139-142.

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COOK IRELAND LTD.